# 

# ANALYTICAL METHOD VALIDATION PROTOCOL

FOR

T M TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE

HYDROCHLORIDE, CHLORPHENAMINE MALEATE

Z

PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE POWDER AND ASCORBIC ACID



Site Address: Safetab Life Science Plot No.A-67 to 72, PIPDIC Electronic Park, Thirubuvanai, Puducherry-605 107.



#### Safetab Life Science

ST/AMVAP/017 Page No. 2 of 25

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE AND ASCORBIC ACID POWDER

Revision No.:00

TITLE

19	9.7 STABILITY OF ANALYTICAL SOLUTION	
18	(iii) Intermediate Precision (Ruggedness)	
17	(ii) Method Precision	
17	(i) System precision	***************************************
17	9.6 PRECISION	
16	9.5 ACCURACY (RECOVERY)	
15	c) Oxidative degradation	
15	b) Alkali degradation	
15	a) Acid degradation	
15	9.4 INTERFERENCE FROM DEGRADANTS (forced degradation)	
ä	9.3 LINEARITY AND RANGE	
12	9.2 SPECIFICITY	
12	9.1 SYSTEM SUITABILITY	
12	9.0 VALIDATION RESULTS	10
<b>∞</b>	8.0 DESCRIPTION OF ANALYTICAL METHOD	
7	7.0 DETAILS OF INSTRUMENTS, COLUMN, SOLVENTS AND CHEMICALS TO BE USED FOR VALIDATION WORK	
o	6.0 DETAILS OF STANDARD, SAMPLES AND PLACEBO TO BE USED FOR VALIDATION WORK	
ហ	5.0 GENERAL INFORMATION  REFERENCE, TYPE OF VALIDATION, TEST TO BE VALIDATED,  COMPOSITION, BATCH NO, SPECIFICATION LIMIT, VALIDATION  STUDY, VALIDATION TEAM	1 99
И	4.0 SCOPE	
СЛ	3.0 OBJECTIVE	7.x
4	2.0 PROTOCOL APPROVAL SHEET	A 1
2	1.0 INDEX	L.A
PAGE No.		
	1.0 THE PARTY OF T	н



#### Safetab Life Science

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL,
PHENYLEPHRINE HYDROCHLORIDE,
CHLORPHENAMINE MALEATE IN PARACETAMOL,
PHENYLEPHRINE HYDROCHLORIDE,
CHLORPHENAMINE MALEATE AND ASCORBIC ACID

CHLORPHENAMINE MALEATE AND ASCORBIC ACID

POWDER

Page No. 3 of 25

9.8 FILTER PAPER STUDY 9.9 ROBUSTNESS 10.0 PARAMETERS TO BE VALIDATED 11.0 ABBREVIATION 12.0 REVISION HISTORY				INDEX PAGE No.	No.
				PAPER STUDY 20	0
				NESS 21	<b> -</b> -
11.0 ABBREVIATION 12.0 REVISION HISTORY	Ö		PARAMETERS	TO BE VALIDATED 23	ω
12.0 REVISION HISTORY	<b> </b>	Ö	<b>ABBREVIATIO</b> I	N 24	4
	<u> </u>	0	REVISION HIS	3TORY 25	U



# Safetab Life Science

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, ST/AMVAP/017

THE TEST OF ASSAY OF PARACETAMOL,
PHENYLEPHRINE HYDROCHLORIDE,
CHLORPHENAMINE MALEATE IN PARACETAMOL,
PHENYLEPHRINE HYDROCHLORIDE,
CHLORPHENAMINE MALEATE AND ASCORBIC ACID
POWDER

Page No. 4 of 25

Revision No.:00

# 2.0 PROTOCOL APPROVAL SHEET

Prepared by : Asst.Manager-QC

Name : K. SARAVAZIDZ

Signature : R

Date : 06 110 2022

Reviewed by : AGM-QC

Name: M. V. JASA FUMAR

Signature :

: 0b/10/2022

Date

Approved by : GM-QA

: A G, ICANNAZ

Name

Signature : A Common to the signature

Date : othorn

Effective Date : 10 10 2022



#### Safetab Life Science

ANALYTICAL METHOD VALIDATION PROTOCOL FOR

CHLORPHENAMINE MALEATE AND ASCORBIC ACID CHLORPHENAMINE MALEATE IN PARACETAMOL, THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, PHENYLEPHRINE HYDROCHLORIDE, POWDER

> Page No. 5 of

ST/AMVAP/

Revision No.:00

#### 3.0 OBJECTIVE:

Maleate and Ascorbic acid powder by HPLC Chlorphenamine validate the Maleate method for Ę Paracetamol, 앜 assay Phenylephrine 으 Paracetamol, Hydrochloride, Phenylephrine Chlorphenamine Hydrochloride,

#### 4.0 SCOPE:

Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Ascorbic acid powder by HPLC method. This protocol shall define the procedure, Documentation refer the acceptance criteria to be used in determination of Assay by HPLC Method. This scope of the Protocol is to evaluate the acceptability of analytical method used for the Paracetamol, Phenylephrine Hydrochloride and Chlorphenamine Maleate

#### 5.0 GENERAL INFORMATION:

REFERENCE

In-House

. .

TYPE OF VALIDATION

. Validation of non-pharmacopeial method

TEST TO BE VALIDATED

..

Phenylephrine Hydrochloride, Chlorphenamine Assay of Paracetamol, Maleate Phenylephrine Chlorphenamine Maleate Hydrochloride, Paracetamol,

and Ascorbic acid powder.

Each 4.5gm sachet contains:

COMPOSITION

Ascorbic acid BP	Chlorphenamine Maleate BP	Phenylephrine hydrochloride BP 10mg	Paracetamol BP	Content
50mg	20mg	10mg	650mg	Strength

**BATCH NO** 

ST/T/S-1322

SPECIFICATION LIMIT

.. 90.0% to 110.0% of the labeled claim

VALIDATION STUDY

QC-Laboratory, Safetab Life science, Puducherry

**VALIDATION TEAM** 

0 0 C.Albin jose

L.Parthasarathi

3. R.Vignesh



TILE

# Safetab Life Science

Protoco (No. 27 ST/AMVAP/017

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE AND ASCORBIC ACID POWDER

Revision No.:00

#### 6.0 DETAILS WORK: OF STANDARD, SAMPLES AND **PLACEBO** 70 BE USED FOR VALIDATION

NAME OF THE MATERIAL	ID NO/BATCH NO	POTENCY/PURITY
Sample	To be mentioned in report	To be mentioned in report
Plain placebo	To be mentioned in report	To be mentioned in report
Working standard Paracetamol BP	To be mentioned in report	To be mentioned in report
Phenylephrine Hydrochloride BP	To be mentioned in report	To be mentioned in report
Chlorphenamine Maleate BP	To be mentioned in report	To be mentioned in report
Ascorbic acid BP	To be mentioned in report	To be mentioned in report
API Paracetamol BP	To be mentioned in report	To be mentioned in report
Phenylephrine Hydrochloride BP	To be mentioned in report	To be mentioned in report
Chlorphenamine Maleate BP	To be mentioned in report	To be mentioned in report
Ascorbic acid BP	To be mentioned in report	To be mentioned in report



#### Safetab Life Science

CHLORPHENAMINE MALEATE AND ASCORBIC ACID ANALYTICAL METHOD VALIDATION PROTOCOL FOR Protocol Wall CHLORPHENAMINE MALEATE IN PARACETAMOL, THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, PHENYLEPHRINE HYDROCHLORIDE,

Page No. 7 of 25

ST/AMVAP/017

Revision No.:00

7.0 **VALIDATION WORK:** DETAILS 0 INSTRUMENTS, COLUMN, SOLVENTS AND CHEMICALS USED FOR

POWDER

#### Instruments:

High performance liquid chromatograph with PDA detector

Shimadzu, Model: LC-2030C 3D Prominence

performance liquid chromatograph with UV visible detector

Shimadzu, Model: LC 2030 Prominence i

### Analytical Balance

Make Sartorius, Model : Quintix-125D-10IN

#### Ţ

Make: Eutech instruments, Model No: pH 700

#### Column:

Inerstil ODS 3V, 250 mm X 4.6 mm, 5µm (or) equivalent

# Reagents, chemicals and Working standard with grade:

Paracetamol BP (Working standard)

Phenylephrine Hydrochloride BP (Working standard)

Chlorphenamine Maleate (Working standard)

1-Heptanesulphonic acid sodium salt (AR grade)

Orthophosphoric acid (AR grade)

Purified Water (Milli-Q water (or) equivalent)

Acetonitrile (HPLC grade)

Methanol (HPLC grade)



# Safetab Life Science

CHLORPHENAMINE MALEATE AND ASCORBIC ACID ANALYTICAL METHOD VALIDATION PROTOCOL FOR CHLORPHENAMINE MALEATE IN PARACETAMOL, THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, PHENYLEPHRINE HYDROCHLORIDE, POWDER

TILE

Page No. 8 of 25
Protocol Wo:
ST/AMVAP/017

Revision No.:00

# 8.0 DESCRIPTION OF ANALYTICAL METHOD:

# Preparation of Buffer solution:

Weigh filter and degas. water. And adjust pH to 3.0 $\pm$ 0.05 with Orthophosphoric acid. Filter through 0.45 $\mu$  membrane and dissolve about 2.0g of 1-heptane sulphonic acid sodium salt in 1000 mL of Milli-Q

# Preparation of Mobile Phase A:

Use buffer solution as mobile phase A.

# Preparation of Mobile Phase B:

Use acetonitrile as mobile phase B.

# Preparation of Diluent:

Prepare a degassed mixture of buffer and methanol in the ratio of 50:50 v/v.

# **Chromatographic Conditions:**

Column Inerstil ODS 3V, 250 mm X 4.6 mm, 5µm (or) equivalent

Wave length : UV at 220 nm

Column Temperature : 30°C

Flow Rate : 1.2 mL/min

Injection Volume : 50 µL

Run time : 20 Minutes

# Preparation of Blank Solution:

Use diluent as blank,

10minutes before further all the prepared using for dilution standard and filtration. sample solutions 00 bench top for

### Gradient Program:

20.0	) )	14.01	14.0		8.0		5.0		0.01	Пте	
80		80	80		50		80		80	Mobile phase A %	
20		20	50		50		20	20	30	Mobile phase B%	



Safetab

Life

Science

TITLE

CHLORPHENAMINE MALEATE AND ASCORBIC ACID ANALYTICAL METHOD VALIDATION PROTOCOL FOR CHLORPHENAMINE MALEATE IN PARACETAMOL, THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, PHENYLEPHRINE HYDROCHLORIDE, POWDER

> Protocol No. ST/AMVAP/017 Page No. 9 of 25

Revision No.:00

# Preparation of Standard Stock Solution-1:

volume with diluent and mix. Weigh and transfer accurately about 40 mg of Phenylephrine hydrochloride WS into a 200 mL clean, dry volumetric flask. Add 140 mL of diluent and sonicate to dissolve. Dilute up to the

# Preparation of Standard Stock Solution-2:

dissolve. Dilute up to the volume with diluent and mix. Paracetamol WS into a 100 mL clean dry volumetric flask. Add 70 mL of diluent and sonicate Weigh and transfer accurately about 40 mg of Chlorphenamine maleate WS and 65mg

# Preparation of Standard Solution:

volumetric flask. Dilute up to the volume with diluent and mix. Transfer each 5 mL of standard stock solution-1, standard stock solution-2 and into a 50 크

# Preparation of Sample solution-A (For Phenylephrine & Chlorphenamine maleate):

equivalent to 650mg of Paracetamol into a 500 mL volumetric flask add about 340 mL of diluent and sonicate for 20 minutes with intermittent shaking. Cool to room temperature and dilute up to the volume with diluent and mix. Filter through 0.45µm PVDF filter. Transfer and mix the contents of not less than 5 sachets. Weigh and transfer the sample

# Preparation of Sample solution-B (For Paracetamol):

volume with diluent and mix. Transfer 5 mL of above Sample solution-A in to a 100 mL volumetric flask and dilute up to the

#### Procedure:

chromatograms and measure the responses for the major peaks solution-A and Inject diluent as Sample blank solution. Inject Standard solution in five replicates, Sample solution-B in duplicates into the chromatograpl chromatograph. Inject Sample

n. Record the

purpose only. The retention times for Paracetamol, Phenylephrine and Chlorphenamine were about 4.0 minutes, 6.7 minutes and 10.4 minutes respectively and it's for information

### System suitability:

Theoretical plate

Tailing factor

Relative standard Deviation

: NLT 2000 for Paracetamol, Phenylephrine and

: NMT 2.0 for Paracetamol, Phenylephrine and

Chlorphenamine peak.

NMT 2.0% for five replicate standard injection of

Paracetamol, Phenylephrine and Chlorphenamine.

Inject 50µl of the above solution as per following sequence



Safetab

Life

TITLE

#### CHLORPHENAMINE MALEATE AND ASCORBIC ANALYTICAL METHOD VALIDATION PROTOCOL FOR CHLORPHENAMINE MALEATE IN PARACETAMOL, THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, PHENYLEPHRINE HYDROCHLORIDE, Science ACID

POWDER

Protocol No. ST/AMVAP/017 Page No. 10 of 25

Revision No.:00

## Injection sequence:

	Ú	4	ω	2	<u> </u>	S. No
	Bracketing standard	Sample solution-B	Sample solution-A	Standard solution	Diluent (Blank)	Sample Name
Home	Each after every 6 sample	2	2	И	⊢≛	No. of injections

# Calculate the assay of Paracetamol in mg/sachet as follows:

Where,

A Average area of peak due to Paracetamol in Sample solution B.

AS

SN Average area of peak due to Paracetamol in standard preparation.

Weight of Paracetamol working standard in mg.

**≶**T Weight of sample taken in mg.

AFW Average fill weight of sachet in mg.

ס Potency of Paracetamol working standard in % on as such basis.

# Calculate the assay of Paracetamol in % as follows:

5 Label claim of Paracetamol in mg/sachet.

### Calculate the assay of Phenylephrine Hydrochloride in mg/sachet as follows:



#### ANA

TITLE

# Safetab Life Science

Rage No. 11 of 25
Protocol No.:
ST/AMVAP/017

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE AND ASCORBIC ACID POWDER

Revision No.:00

Where,

A 11 solution A. Average area of peak due ç Phenylephrine Hydrochloride = Sample

SS 11 preparation. Average area of peak due to Phenylephrine Hydrochloride in standard

SNS 11 Weight of Phenylephrine Hydrochloride working standard in mg.

WT = Weight of sample taken in mg.

AFW = Average fill weight of sachet in mg.

ס 11 such basis. Potency of Phenylephrine Hydrochloride working standard 3 % 9 n as

# Calculate the assay of Phenylephrine Hydrochloride in % as follows:

5 П Label claim of Phenylephrine Hydrochloride in mg/sachet.

# Calculate the assay of Chlorphenamine maleate in mg/sachet as follows:

Where,

A solution A. Average area 으 peak due ð Chlorphenamine maleate 3 Sample

ΑS preparation. Average area 앜 peak due 6 Chlorphenamine maleate in standard

S¥ Weight of Chlorphenamine maleate working standard in mg

WT = Weight of sample taken in mg.

AFW = Average fill weight of sachet in mg.

ס 11 Potency of Chlorphenamine maleate working standard in % on as such

# Calculate the assay of Chlorphenamine maleate in % as follows:

5 Label claim of Chlorphenamine maleate in mg/sachet.



#### TILE

# Safetab Life Science

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE AND ASCORBIC ACID

POWDER

Page No. 12 of 25
Protocol No.D > Y
ST/AMVAP/017

Revision No.:00

9.0 VALIDATION RESULTS:

# 9.1 SYSTEM SUITABILITY TEST:

#### Purpose:

To establish system suitability as per methodology.

#### Study Design:

Sequence shall be in following provisional manner.

2	<b> </b>	S.No.
Standard preparation	Blank (Diluent)	Description of solution
ហ	<b>-</b>	No. of Injections

# Evaluate the following system suitability parameters:

- standard injections 1) % RSD of area 으 Paracetamol, Phenylephrine and Chlorphenamine peak in five replicate
- injection. 2 Theoretical plates for Paracetamol, Phenylephrine and Chlorphenamine peak Ξ. standard
- injection ω Tailing factor ਨ੍ਹੀ Paracetamol, Phenylephrine and Chlorphenamine peak Ξ. standard

# Acceptance Criteria:

- standard injections 1) % RSD of area for Paracetamol, should not more than 2.0%. Phenylephrine and Chlorphenamine peak in five replicate
- injection should Theoretical plates for Paracetamol, not less than 3000. Phenylephrine and Chlorphenamine peak 3 standard
- injection should not more than 2.0. Tailing factor for Paracetamol, Phenylephrine and Chlorphenamine peak ⋽. standard

### 9.2 SPECIFITY:

presence "The specificity is the ability of an analytical procedure 으 componenets that may be expected present in sample matrix". ç measure accurately an analyte 3

#### Purpose:

ОТ demonstrate that the placebo not interfering with the analyte peak.



#### Safetab Life Science

CHLORPHENAMINE MALEATE AND ASCORBIC ACID ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, PHENYLEPHRINE HYDROCHLORIDE,

POWDER

ST/AMVAP/017 Protocol No. Page No. 13 of U

Revision No.:00

#### Study design:

Sequence shall be in following provisional manner.

12 Test preparation-Soln-B
11 Test preparation-Soln-A
Plain placebo with Chlorphenamine, Phenylephrine Hcl and Paracetamol
9 Plain placebo with Chlorphenamine Maleate
8 Plain placebo with Phenylephrine Hcl
7 Plain placebo with Paracetamol
6 Chlorphenamine Maleate working Standard
5 Phenylephrine HCL working Standard
4 Paracetamol working Standard
3 Plain placebo
2 Standard preparation
1 Blank (Diluent)
S.No.

# Acceptance criteria:

- There should not be any interference due to blank, Placebo peak with analyte. Peak purity should not be less than 0.995 according to Lab solution software.

#### <u>9</u> :3 LINEARITY AND RANGE:

proportional to "The linearity of the analytical method is it's ability to elecit test resuthe concentration of the analyte in samples within give range". test results data directly

#### Purpose:

7 Establish the linearity of analyte within the specified range



#### 

# Safetab Life Science

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, CHLORPHENAMINE MALEATE AND ASCORBIC ACID CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, PHENYLEPHRINE HYDROCHLORIDE, POWDER

Page No. 14 of 25
ProtocolNo.
ST/AMVAP/017

Revision No.:00

#### Study Design:

targeted To demonstrate the linearity and range of analytical method over the range of 10% to 150% of concentration.

Sequence shall be in following provisional manner.

12		10	9	&	7	6	И	4	ω	2		S.No.
Level – 6 (150%)	Blank (Diluent)	Level - 5 (125%)	Blank (Diluent)	Level – 4 (100%)	Blank (Diluent)	Level – 3 (75%)	Blank (Diluent)	Level - 2 (50%)	Blank (Diluent)	Level - 1 (10%)	Blank (Diluent)	Description of solution
2	1	2	Д	2	1	2	₽	2	⊢	2	<b>—</b>	No. of Injections

Plot a graph Evaluate the square, slope and Y-intercept. Plot a of concentration of concentration (at X-axis) versus squared correlation coefficient  $(r^2)$ , average peak area of analyte (at Y-axis). correlation coefficient (r), residual sum of

## Acceptance criteria:

- 1) To conclude the linearity, the squared correlation coefficient should not be less than 0.995
- 2) To conclude the range. % RSD for peak area 125% and 150% should be not more than 2.0. of linearity level of 10%, 50%, 75%, 100%,



# Safetab Life Science

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, CHLORPHENAMINE MALEATE AND ASCORBIC ACID CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE PHENYLEPHRINE HYDROCHLORIDE, POWDER

Page No. 15 of 25
Protocol No.)

Revision No.:00

# **e**. INTERFERENCE FROM DEGRADANT (forced degradation)

#### Study design:

stressing the test 0 evaluate the preparation under the following maximum stress conditions interference from degradants, carry out വ forced degradation study φ

פאייממני בר מכאי מממנוסו	Oxidative degradation	Alkali degradation	Acid degradation	Degradation
80°C for 30 minutes.	Exposure to 5ml of 30% H2O2 and Heat on water hath at	Exposure to 5ml of 5N NaOH and Heat on water bath at 80°C for 30 minutes.	for 30 minutes.	Stress Condition

# Sequence shall be in following provisional manner, For forced chemical degradation:

7 Standard preparation (Bracketing)	6 Sample Solution (A&B) (Oxidative degradation)	5 Sample Solution (A&B) (Alkali degradation)	4 Sample Solution (A&B) (Acid degradation)	3 Sample Solution (A&B) (As such)	2 Standard preparation	1 Blank (Diluent)	S.No. Description of solution
1	2	2	2	2	Cī.	1	No. of Injections

samples and the degradation profiles under each stressed condition. Chromatograph the samples of chemical and physical forced degradation into HPLC system equipped with diode array detector and evaluate the peak purity for the analytes in stressed

# **Acceptance Criteria:**

- There should not be any interference due to degradants with analyte in stressed samples
- possible). The desired degradation should be 10-30% ⋽ acid, alkali and oxidative degradation, Œ



# Safetab Life Science

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, CHLORPHENAMINE MALEATE AND ASCORBIC ACID CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, PHENYLEPHRINE HYDROCHLORIDE, POWDER

TILE

Protocol NO. 7

Revision No.:00

- same shall be documented about 10% ᆼ 30% degradation is and reported. not achieved by applying above stressed condition,
- 4) Peak purity should not be less than 0.950 according to Lab solution software.

# 9.5 ACCURACY STUDY (RECOVERY STUDY)

add amount of analyte". true value. accuracy of an analytical method is the closeness of results obtained by that method to Accuracy may often be expressed as present recovery by the assay of known,

#### Purpose:

establish the accuracy of the analytical method in the specified range.

# Sequence shall be in following provisional manner

<u>1</u> 5	14	13	12	}	10	9	&	7	6	У	4	ω	2	<b> -</b>	S.No.
Standard preparation (Bkt)	Level - 3 Set - 3 (150%)	Level - 3 Set - 2 (150%)	Level - 3 Set - 1 (150%)	Blank (Diluent)	Level - 2 Set - 3 (100%)	Level – 2 Set – 2 (100%)	Level - 2 Set - 1 (100%)	Blank (Diluent)	Level - 1 Set - 3 (50%)	Level – 1 Set – 2 (50%)	Level - 1 Set - 1 (50%)	Blank (Diluent)	Standard preparation	Blank (Diluent)	Description of solution
1	1	1	Н	<b>⊢</b>	<b>⊢</b> •	<u>L</u>	<b> </b> -\$	1	1	Н	ı	1	И	<u> </u>	No. of Injections



# Safetab Life Science

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL,
PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, MALEATE POWDER AND ASCORBIC ACID

Page No. 17 of 25
Protocol No.:
ST/AMVAP/017

Revision No.:00

#### Study design:

To demonstrate the accuracy of the analytical method, prepare known quantities of drug (at level 50%, 100% and 150% c placebo. Prepare the recovery samples in triplicate for each level. prepare recovery samples by spiking 150% of targeted concentration) to

# Acceptance criteria:

The mean % recovery at each level should be 98.0 to 102.0

### 9.6 PRECISION:

repeatability and reproducibility" homogeneous factor) "The Precision of an analytical procedure express the closeness of the agreement (Degree between sample വ series under of measurements the prescribed obtained from condition. multiple Precision sampling of n may be o considered the same

## (i) System Precision

#### Purpose:

To establish the precision of the HPLC system being used for the analysis.

#### Study Design:

Sequence shall be in following provisional manner.

2	<b> </b> 4	S.No.
Standard preparation	Blank (Diluent)	Description of solution
СЛ	Ľ	No. of Injections

## Acceptance criteria:

2.0 % RSD of area of analyte peak in five replicate standard injections should not be more than

# (ii) Method Precision:

#### Purpose:

To establish the repeatability of test results obtained by the analytical method.

#### Study design:

precision by To demonstrate the method precision, analyze six sample preparations as per the methodology representing a single batch and determine the assay for the same. Evaluate the method computing the percentage and relative standard deviation of the assay results



# Safetab Life Science

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE AND ASCORBIC ACID CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE,

POWDER

_		140,	22000	- 1
	Ŋ	ש		ט
	7	3	1	Pag
	A	5		D
	AMVAP	rotoco		<b>Z</b>
	≱	No	(1000) I	27,000,000
	-	Q.	7	X
	017		)	2
	7	2000		30
		PRESE	< (	<b>У</b>

Revision No.:00

9	8	7	6	И	4	ω	2	<b> </b>	S.No.
Standard preparation (BKT)	Sample Solution (A&B) -6	Sample Solution (A&B) -5	Sample Solution (A&B) -4	Sample Solution (A&B) -3	Sample Solution (A&B) -2	Sample Solution (A&B) -1	Standard preparation	Blank (Diluent)	Description of solution
1 (after six sample injection)	2	2	2	2	2	2	σı	1	No. of Injections

## Acceptance criteria:

% RSD for assay of six preparations should not be more than 2.0.

# (iii) Intermediate Precision (Ruggedness):

#### Purpose:

standard deviation of the assay results. the variability of instrument, column (different lot no) analyst and day. Analyse six sample preparations as per the methodology representing a single batch and determine the assay for variability To demonstrate the reproducibility of test results obtained by the analytical method for the variability of instrument, column (different lot no) analyst and day. Analyse six sample same. Evaluate the intermediate precision by computing the percentage and relative



#### Safetab Life Science

Page No.

19 of

25

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL,

CHLORPI

9	8	7	6	ъ	4	ω	2	H	S.No.
Standard preparation (BKT)	Sample Solution (A&B) -6	Sample Solution (A&B) -5	Sample Solution (A&B) -4	Sample Solution (A&B) -3	Sample Solution (A&B) -2	Sample Solution (A&B) -1	Standard preparation	Blank (Diluent)	Description of solution
(after six sample injection)	2	2	2	2	2	2	ъ	1	No. of Injections

## Acceptance criteria:

- 1) % RSD for assay of six preparations should not be more than 2.0.
- 2) Cumulative % RSD for assay of twelve preparations (i.e. method precision and intermediate precision) should not be more than 2.0.

# 9.7 STABILITY OF ANALYTICAL SOLUTION:

#### Study design:

Prepare Standard and sample solution as per the methodology and store at Ambient temperature. Chromatograph this solution at regular intervals for 48 hours by using same diluent. Calculate the % difference of analyte peak area for standard and test preparations with that of initial. The study may be stopped if 2 consecutive failure of sample solution.

Sequence shall be in following provisional



# Safetab Life Science

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE AND ASCORBIC ACID CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, POWDER

TITLE

Page No. 20 of 25
Protocol No.!
ST/AMVAP/017

Revision No.:00

6	И	4	ω	2	<b>⊢</b>	S.No.
Sample solution (A&B) (Time interval)	Standard preparation (Time interval)	Sample solution (A&B) (Initial)	Standard preparation(A&B) (Initial)	Standard preparation	Blank (Diluent)	. Description of solution
Ľ	H	Н	1	ъ	<b>)</b>	No. of Injections

## Acceptance criteria:

The sample and standard solution shall be considered stable for the final period till which the area difference between initial and next periodic interval should not be more than  $\pm 2\%$ .

# 9.8 FILTER PAPER STUDY:

#### Study design:

The filter paper study of the analytical method shall perform by filtering test solution through 0.45µ PVDF filter against that of unfiltered.



# Safetab Life Science

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL,
PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE AND ASCORBIC ACID CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, POWDER

Page No. 21 of 25

ST/AMVAP/017

Revision No.:00

Sequence shall be in following provisional manner.

_		7			·····			
	<b>∞</b>	6	ហ	4	ω	2	<b>J</b> k	S.No.
	Standard preparation	Sample solution (A&B) –Filter Set 3 (0.45µ PVDF filter)	Sample solution (A&B) –Filter Set 2 (0.45µ PVDF filter)	Sample solution (A&B) –Filter Set 1 (0.45µ PVDF filter)	Sample solution (A&B) -Unfiltered (Centrifuge)	Standard preparation	Blank	Description of solution
<b> </b>	•	<b> </b>	щ	H	<b>p.k</b>	υ	⊢≛	No. of Injections

## Acceptance criteria:

The % area difference of filter solution should not differ  $\pm 2.0$  against that of unfiltered.

## 9.9 ROBUSTNESS:

#### Purpose:

To establish the robustness of the analytical method.

#### Study Design:

against deliberate changes in chromatographic conditions The robustness of the analytical method can be established by demonstrating its reliability

Sequence shall be in following provisional manner.



#### Safetab Life Science

Protocol No. ST/AMVAP/017 Page No.

22 of

25

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, CHLORPHENAMINE MALEATE AND ASCORBIC ACID CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, PHENYLEPHRINE HYDROCHLORIDE, POWDER

TITLE

Revision No.:00

н	Bracketing standard	4
2	Sample solution (A&B)	ω
VΊ	Standard preparation	2
<b>⊢</b> ≛	Blank (Diluent)	<u>н</u>
No. of Injections	Description of solution	S.No.
	According to each variable	
⊷	Bracketing standard	4
2	Sample solution (A&B)	ω
5	Standard preparation	2
<b>}</b>	Blank (Diluent)	<b>ل</b> سم
No. of Injections	Description of solution	S.No.
	As such	

parameters. Following variable shall be done according to deliberate changes in chromatographic

- a Flow rate change by  $\pm 10\%$  mean (i.e 1.1 ml/min and 1.3 ml/minute)
- b) Wave length change by  $\pm$  3nm (i.e. 217nm and 223nm)
- c) Column oven Temperature change by  $\pm$ 5.0 (i.e. 25°C and | 35°C)

## Acceptance criteria:

assay value of method precision. System suitability should comply for each variable and % of drug not differ ±2% from mean



# Safetab Life Science

Protocorno.by ST/AMVAP/017

ANALYTICAL METHOD VALIDATION PROTOCOL FOR A
THE TEST OF ASSAY OF PARACETAMOL,
PHENYLEPHRINE HYDROCHLORIDE,
CHLORPHENAMINE MALEATE IN PARACETAMOL, CHLORPHENAMINE MALEATE AND ASCORBIC ACID PHENYLEPHRINE HYDROCHLORIDE, POWDER

Revision No.:00

# 10.0 PARAMETERS TO BE VALIDATED:

œ	7	6.		5	4.	ώ			) <u></u>	No
A) Flow rate change b) Wavelength c) Temperature change	Filter paper study	Stability of Analytical solution	<ul><li>(i) System precision</li><li>(ii) Method precision</li><li>(iii) Intermediate precision</li></ul>	Precision	Accuracy (Recovery)	Linearity and range	<ul><li>(II) Interference from degradants (Forced degradation)</li><li>a) Acid degradation</li><li>b) Alkali degradation</li><li>c) Oxidative degradation</li></ul>	Specificity (Selectivity)  (i) Interference from blank and placebo	System suitability	Validation parameters

having common system suitability with bracketing standard. Note: More than one parameter can be performed at once with relevant sequence



#### Safetab Life Science

Page No. 24 of 25

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, CHLORPHENAMINE MALEATE AND ASCORBIC ACID CHLORPHENAMINE MALEATE IN PARACETAMOL, PHENYLEPHRINE HYDROCHLORIDE, POWDER

TITLE

ST/AMVAP/017

Revision No.:00

11.0 ABBREVIATION:

mg Milligram

S.No Serial Number

<u>=</u> Milliliter

% Percentage

D Identification

API Active pharmaceutical ingredient

High performance liquid chromatography

B.NO

**HPLC** 

Batch number

mm Millimeter

m Micrometer

min Minutes

ĉ Degree centigrade

nm Nanometer

RSD Relative standard deviation

三 Micro litre

HCL Hydrochloric acid

NaoH Sodium Hydroxide

H<sub>2</sub>O<sub>2</sub> Hydrogen Peroxide



#### 

Safetab Life Science

ANALYTICAL METHOD VALIDATION PROTOCOL FOR THE TEST OF ASSAY OF PARACETAMOL,
PHENYLEPHRINE HYDROCHLORIDE,
CHLORPHENAMINE MALEATE IN PARACETAMOL,
PHENYLEPHRINE HYDROCHLORIDE,
CHLORPHENAMINE MALEATE AND ASCORBIC ACID POWDER

Page No. 25 of 25
Protogol No. DY
ST/AMVAP/017

Revision No.:00

# 12.0 REVISION HISTORY:

ST/AMVAP/017	Protocol No.
10/10/2022	Effective date
New Protocol prepared.	Reason for Review

\*\* END OF THE DOCUMENT\*\*