

		Batch Filling Record (Part-A)		Page 18 of 46
PLANT-II, BADDI				
Product Name	Topotecan for Injection 4 mg (base)/vial			
Batch Number	8PK201			
Batch size	12 L	Batch Filling Code	IO720004	
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001	
Shelf life	24 Months (Provisional)	Version No.	00	

6.3 TRANSFER OF RUBBER STOPPERS:

Transfer the ready for sterilization rubber stoppers from warehouse to the unit preparation area through the static pass box (CSP/104) of equipment wash room (CSLG 014).

6.4 STERILIZATION OF RUBBER STOPPERS AS PER SOP No. SOP/CSP/037:

Cycle started at: 21:48

Cycle Completed at: 23:47

Date	Load No.	Load Pattern No.	No. of Bags	Total qty of rubber stoppers	Sterilization Hold Time		Done by ^{##}	Ckd. By ^{##}
					From	To		
14/10/08	05	6	7	7000	22:39	23:09	Shambhu	PR

(Graph attached on the reverse side of page)

6.5 STORAGE OF STERILIZED RUBBER STOPPERS (NMT 48 HOURS AFTER STERILIZATION):

Procedure	Eq. Id. No.	SOP No.	Date	Storage Time*		Done by ^{##}	Ckd. By ^{##}
				From	To		
Store the sterilized Rubber stoppers under LAF in cooling zone.	CSP/035	SOP/CSP/072	14/10/08	23:47	---	Shambhu	14/10/08
			15/10/08	---	21:45	Sayjeev	15/10/08

* Hold time starts after completion of sterilization to end of filling.

^{##} Production

[Signature]

Lead No -05


MFG.BY :MACHIN FABRIK
CUSTOMER :M/S DABUR PHARMA LTD.
EQUIPMENT :STEAM STERILIZER

BATCH DATA
STERILIZER NO :1007
BATCH NO :0
QUANTITY :1
OPERATOR CODE :1
PRINT INTERVAL :60 SEC
DATE :14/10/08
TIME :21:48:38

PROCESS : STOPPER & FLIP OFF

SET PARAMETERS
PRE VACUUM : -0.750 BAR
PRE PRESSURE : 0.500 BAR
NO OF PRE PULSE : 3 NOS
SET POINT-1 : 110.0 DEG C
SET POINT-2 : 115.0 DEG C
SET POINT-3 : 118.0 DEG C
SP-1 HOLD TIME : 2 MIN
SP-2 HOLD TIME : 3 MIN
SP-3 HOLD TIME : 5 MIN
HEAT UP CONTROL BAND : 0.2 DEG C
STER HOLD TEMP : 121.5 DEG C
STER HOLD TIME : 30 MIN
TEMP CONTROL BAND : 0.2 DEG C
OVERSHOOT TEMP : 124.0 DEG C
STER COUNT STOP TEMP : 120.9 DEG C
STER COUNT RESET TEMP : 120.8 DEG C
POST VACUUM : -0.750 BAR
POST PRESSURE : -0.100 BAR
POST VAC HOLD TIME : 30 MIN
NO OF POST PULSE : 1 NOS
EXHAUST ON TIME : 5 SEC
EXHAUST OFF TIME : 50 SEC
PROCESS END PRESSURE : -0.050 BAR
CHAM PRESSURE HIGH : 1.500 BAR

MFG.BY :MACHIN FABRIK
CUSTOMER :M/S DABUR PHARMA LTD.
DATE :14/10/08
TIME :21:49:19
EQUIPMENT :STEAM STERILIZER
PROCESS :STOPPER & FLIP OFF SEAL



* Only Rubber stopper load has been taken.

14/10/08

1/5

TIME	T1 DEG C	T2 DEG C	T3 DEG C	T4 DEG C	T5 DEG C	PRES BAR	STATUS
21:49:19	030.4	039.7	045.0	035.9	046.0	-0.021	PROCESS ON
21:49:21	030.3	039.8	045.1	036.0	046.1	-0.031	PRE VACUUM
21:50:20	029.8	037.4	043.1	033.0	044.0	-0.644	PRE VACUUM
21:50:47	029.6	036.9	042.4	032.1	043.4	-0.767	PRE PRESSURE
21:51:20	029.9	051.9	050.3	055.4	051.5	-0.586	PRE PRESSURE
21:52:20	030.0	066.1	063.0	068.0	065.2	-0.448	PRE PRESSURE
21:53:21	030.0	074.6	072.1	074.5	073.0	-0.327	PRE PRESSURE
21:54:21	030.1	080.8	078.9	080.8	079.1	-0.215	PRE PRESSURE
21:55:21	043.8	085.9	084.7	083.3	083.6	-0.104	PRE PRESSURE
21:56:20	066.1	090.4	089.4	087.4	088.2	00.008	PRE PRESSURE
21:57:22	089.4	093.9	093.7	091.0	092.1	00.085	PRE PRESSURE
21:58:21	094.8	096.5	096.4	095.4	095.0	00.151	PRE PRESSURE
21:59:20	098.1	098.9	098.9	098.7	097.8	00.218	PRE PRESSURE
22:00:21	100.9	101.2	101.3	099.6	100.3	00.284	PRE PRESSURE
22:01:22	103.2	103.2	103.4	101.2	102.4	00.348	PRE PRESSURE
22:02:21	105.3	105.0	105.2	104.2	104.4	00.411	PRE PRESSURE
22:03:22	107.1	106.8	106.9	107.4	106.1	00.473	PRE PRESSURE
22:03:49	107.8	107.5	107.6	108.2	106.9	00.499	PRE VACUUM
22:04:21	096.7	097.7	098.4	098.6	096.5	-0.201	PRE VACUUM
22:05:20	078.2	078.8	079.9	079.9	078.6	-0.593	PRE VACUUM
22:06:17	065.2	065.8	066.8	067.1	065.9	-0.759	PRE PRESSURE
22:07:23	080.0	081.7	081.6	084.5	081.8	-0.397	PRE PRESSURE
22:08:21	086.5	091.1	091.4	090.9	091.2	-0.177	PRE PRESSURE
22:09:22	093.1	097.1	097.4	096.7	097.1	-0.003	PRE PRESSURE
22:10:21	102.0	100.9	101.2	100.8	100.8	00.116	PRE PRESSURE
22:11:22	104.9	103.9	104.1	104.2	103.8	00.230	PRE PRESSURE
22:12:23	107.4	106.4	106.6	106.9	106.3	00.335	PRE PRESSURE
22:13:22	109.5	108.6	108.8	109.0	108.4	00.430	PRE PRESSURE
22:14:12	110.9	110.0	110.2	110.6	110.0	00.496	PRE VACUUM
22:14:21	109.8	109.1	109.5	109.2	108.6	00.273	PRE VACUUM
22:15:23	089.9	091.0	091.4	091.8	090.6	-0.360	PRE VACUUM

22:16:22	074.8	075.6	075.9	076.3	076.0	-0.649	PRE VACUUM
22:17:03	065.2	067.2	069.1	067.3	073.1	-0.756	PRE PRESSURE
22:17:23	069.9	070.3	070.9	074.0	075.3	-0.617	PRE PRESSURE
22:18:22	087.7	085.9	085.9	086.5	086.2	-0.295	PRE PRESSURE
22:19:23	095.9	095.3	095.7	095.3	095.4	-0.046	PRE PRESSURE
22:20:21	102.0	100.7	101.1	100.9	100.6	00.119	PRE PRESSURE
22:21:23	105.7	104.5	104.9	105.0	104.5	00.265	PRE PRESSURE
22:22:23	108.7	107.6	107.9	108.2	107.5	00.393	PRE PRESSURE
22:23:15	111.0	109.9	110.1	110.5	109.8	00.496	HEAT UP
22:24:23	109.8	109.2	109.5	109.9	109.2	00.443	HEAT UP
22:25:23	110.1	109.5	109.7	110.1	109.4	00.461	HEAT UP
22:26:22	112.6	111.6	111.9	112.3	111.7	00.589	HEAT UP
22:27:23	114.8	113.8	114.0	114.4	113.8	00.695	HEAT UP
22:28:24	115.2	114.3	114.7	115.0	114.3	00.685	HEAT UP
22:29:23	115.2	114.4	114.7	115.0	114.3	00.686	HEAT UP
22:30:23	115.0	114.3	114.6	115.0	114.3	00.674	HEAT UP
22:31:24	116.9	115.9	116.2	116.7	116.0	00.821	HEAT UP
22:32:23	118.0	117.2	117.5	117.8	117.1	00.869	HEAT UP
22:33:24	118.2	117.4	117.7	118.0	117.3	00.852	HEAT UP
22:35:24	117.9	117.2	117.5	117.9	117.2	00.855	HEAT UP
22:36:25	118.0	117.3	117.6	118.0	117.3	00.877	HEAT UP
22:37:24	118.9	118.0	118.2	118.7	118.0	00.936	HEAT UP
22:38:25	120.7	119.7	120.0	120.5	119.8	01.041	HEAT UP
22:39:24	122.1	121.1	121.5	121.8	121.1	01.128	HEAT UP
✓22:39:47	122.6	121.7	121.9	122.3	121.6	01.163	STER HOLD
22:40:25	122.5	121.7	122.0	122.3	121.6	01.152	STER HOLD
22:41:24	122.3	121.6	121.9	122.3	121.5	01.124	STER HOLD
22:42:25	122.5	121.7	122.0	122.4	121.6	01.123	STER HOLD
22:43:24	122.6	121.8	122.1	122.5	121.7	01.161	STER HOLD
22:44:26	122.4	121.6	121.9	122.3	121.6	01.145	STER HOLD
22:45:25	122.4	121.5	121.9	122.3	121.5	01.118	STER HOLD
22:46:26	122.5	121.7	122.0	122.4	121.6	01.119	STER HOLD
22:47:24	122.6	121.8	122.0	122.5	121.7	01.141	STER HOLD
22:48:25	122.5	121.7	122.0	122.5	121.6	01.161	STER HOLD
22:49:26	122.3	121.6	121.9	122.3	121.5	01.133	STER HOLD

22:50:25	122.4	121.7	122.0	122.3	121.6	01.117	STER HOLD
22:51:26	122.6	121.8	122.1	122.5	121.7	01.143	STER HOLD
22:52:25	122.5	121.8	122.0	122.5	121.7	01.166	STER HOLD
22:53:26	122.3	121.7	121.9	122.3	121.5	01.128	STER HOLD
22:54:27	122.5	121.7	122.0	122.4	121.6	01.120	STER HOLD
22:55:25	122.7	121.8	122.0	122.5	121.7	01.141	STER HOLD
22:56:26	122.6	121.8	122.0	122.5	121.7	01.166	STER HOLD
22:57:27	122.3	121.7	121.9	122.3	121.5	01.141	STER HOLD
22:58:25	122.5	121.7	122.0	122.4	121.6	01.123	STER HOLD
22:59:26	122.7	121.8	122.1	122.5	121.7	01.145	STER HOLD
23:00:27	122.6	121.7	122.0	122.5	121.7	01.166	STER HOLD
23:01:26	122.4	121.7	122.0	122.4	121.6	01.149	STER HOLD
23:02:27	122.4	121.7	122.0	122.3	121.5	01.123	STER HOLD
23:03:28	122.6	121.8	122.1	122.5	121.7	01.134	STER HOLD
23:04:26	122.7	121.8	122.1	122.5	121.7	01.160	STER HOLD
23:05:27	122.4	121.7	122.0	122.4	121.6	01.147	STER HOLD
23:06:28	122.5	121.7	122.0	122.3	121.5	01.118	STER HOLD

23:07:29	122.7	121.8	122.1	122.5	121.7	01.141	STER HOLD
23:08:28	122.5	121.8	122.0	122.5	121.6	01.164	STER HOLD
✓ 23:09:29	122.4	121.7	122.0	122.3	121.5	01.116	STER HOLD
23:09:45	122.6	121.8	122.0	122.5	121.7	01.166	POST VACUUM
23:10:28	106.5	108.7	109.5	109.2	107.7	00.138	POST VACUUM
23:11:29	087.9	090.2	090.2	090.9	091.1	-0.401	POST VACUUM
23:12:28	072.2	075.3	084.6	076.8	090.9	-0.668	POST VACUUM
23:13:00	065.3	068.8	085.0	071.3	091.1	-0.746	VACUUM HOLD
23:13:27	060.4	064.6	086.5	066.4	091.1	-0.720	VACUUM HOLD
23:14:28	052.1	062.2	089.3	060.7	091.2	-0.850	VACUUM HOLD
23:15:29	046.4	072.8	090.9	053.3	091.4	-0.882	VACUUM HOLD
23:16:28	042.5	077.9	091.8	046.6	091.7	-0.902	VACUUM HOLD
23:17:29	039.7	081.1	092.6	043.8	092.1	-0.913	VACUUM HOLD
23:18:28	037.8	082.9	093.2	043.1	092.4	-0.919	VACUUM HOLD
23:19:29	036.3	084.2	093.6	043.6	092.7	-0.924	VACUUM HOLD
23:20:29	035.3	085.2	093.9	043.5	092.9	-0.928	VACUUM HOLD

23:22:29	033.7	085.9	094.2	043.0	093.3	-0.932	VACUUM HOLD
23:23:30	033.1	086.0	094.3	042.6	093.4	-0.933	VACUUM HOLD
23:24:31	032.5	085.9	094.3	043.5	093.5	-0.935	VACUUM HOLD
23:25:30	032.0	085.8	094.3	044.2	093.6	-0.935	VACUUM HOLD
23:26:30	031.6	085.7	094.2	043.9	093.8	-0.937	VACUUM HOLD
23:27:29	031.1	085.7	094.1	047.1	093.9	-0.937	VACUUM HOLD
23:28:30	030.4	085.6	094.0	054.8	093.3	-0.938	VACUUM HOLD
23:29:30	030.5	085.5	093.9	059.5	093.0	-0.938	VACUUM HOLD
23:30:30	030.2	085.4	093.7	062.8	093.3	-0.939	VACUUM HOLD
23:31:29	029.9	085.3	093.7	065.0	093.7	-0.940	VACUUM HOLD
23:32:30	029.5	085.3	093.6	066.6	094.0	-0.940	VACUUM HOLD
23:33:31	029.2	085.2	093.5	067.8	094.2	-0.940	VACUUM HOLD
23:34:30	028.8	085.2	093.4	068.6	094.4	-0.941	VACUUM HOLD
23:35:31	028.6	085.1	093.2	069.3	094.5	-0.942	VACUUM HOLD
23:36:32	028.3	085.1	093.2	069.8	094.7	-0.942	VACUUM HOLD
23:37:32	027.9	085.1	093.1	070.3	094.8	-0.942	VACUUM HOLD

23:38:30	027.6	085.1	093.1	070.7	094.8	-0.942	VACUUM HOLD
23:39:31	027.4	085.0	093.0	071.2	094.9	-0.942	VACUUM HOLD
23:40:32	027.0	084.8	093.0	071.6	095.0	-0.942	VACUUM HOLD
23:41:33	026.8	084.7	093.0	072.0	095.1	-0.943	VACUUM HOLD
23:42:31	026.6	084.6	093.0	072.3	095.1	-0.942	VACUUM HOLD
23:43:01	026.5	084.5	092.9	072.5	095.2	-0.943	POST PRESSURE
23:43:32	028.6	084.8	092.5	073.8	095.3	-0.793	POST PRESSURE
23:44:33	029.5	086.0	091.0	073.2	096.7	-0.485	POST PRESSURE
23:45:31	030.4	085.9	090.5	071.8	096.8	-0.232	POST PRESSURE
23:46:16	030.9	084.0	090.7	071.6	096.0	-0.103	VACUUM BREAK
23:46:32	031.0	083.5	090.8	071.5	095.7	-0.072	VACUUM BREAK
23:47:31	031.4	081.1	090.7	070.2	094.4	-0.050	PROCESS END

FD VALUE:- 0030.0

Dr
14/10/08

5/5

		Batch Filling Record (Part-A)		Page 19 of 46
PLANT-II, BADDI				
Product Name	Topotecan for Injection 4 mg (base)/vial			
Batch Number	8PK201			
Batch size	12 L	Batch Filling Code	IO720004	
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001	
Shelf life	24 Months (Provisional)	Version No.	00	

6.6 PRE CLEANING DETAILS OF EQUIPMENTS USED FOR ONLINE FILTRATION AND FILLING:

Clean the equipment and filtration accessories with the WFI in the wash table lying in the equipment wash room (CSLG 014).

Cleaning started at: 22:00

Cleaning completed at: 23:00

Date: 14/10/08

S.No.	Equipment Name	Equipment Id. No.	Previous Product Detail	Done By ^{##}	Ckd. By ^{##}
1.	SS Beakers 500 mL, 1000 mL	---	Cytarabine inj. 8LA204	CS	A
2.	*Platinum Cured silicon tubes	---	---	CS	A
3.	Cartridge Filter Holder with stand	CSP/113	Cytarabine inj. 8LA204	CS	A
4.	**50 mm disc type Vent filter	---	---	CS	A
5.	Buffer Bottle	---	Cytarabine inj. 8LA204	CS	A
6.	SS forceps	---	Cytarabine inj. 8LA204	CS	A
7.	SS filling needle	---	Cytarabine inj. 8LA204	CS	A
8.	SS syringe set	---	Cytarabine inj. 8LA204	CS	A
9.	Stopper hopper, Stopper rail, chute	---	Cytarabine inj. 8LA204	CS	A
10.	Needle holding plate	---	Cytarabine inj. 8LA204	CS	A
11.	Nylon Fasteners	---	---	CS	A
12.	Stopper head set	---	Cytarabine inj. 8LA204	CS	A
13.	TC with Nozzle	---	Cytarabine inj. 8LA204	CS	A
14.	Manifold	---	Cytarabine inj. 8LA204	CS	A
15.	SS clamps, gaskets and nozzles	---	Cytarabine inj. 8LA204	CS	A

**These equipments to be used product dedicated.

* To be Used Batch Specific.

^{##}Production

	Batch Filling Record (Part-A)		Page 20 of 46
	PLANT-II, BADDI		
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

6.7 FINAL CLEANING OF EQUIPMENT AND ACCESSORIES FOR ONLINE FILTRATION AND FILLING:

A. R. No. Of WFI: 02WF08000748

Date: 14/10/08

Note: Final cleaning and rinsing of the washed equipments should be done with WFI under LAF No.: CSP/059 in unit preparation room (CSLG017).

Cleaning Started at: 22:00

Cleaning Completed at: 23:10

S.No.	Equipment Name	Equipment Id. No.	SOP No.	Qty.	Done By##	Ckd. By##
1.	SS Beakers 500 mL, 1000 mL	---	SOP/CSP/006	01 + 01	Sarungh	R
2.	*Platinum Cured silicon tubes	---	SOP/CSP/006	04 + 04 + 03	Sarungh	R
3.	Cartridge Filter Holder with stand	CSP/113	SOP/CSP/006	01	Sarungh	R
4.	**50 mm disc type Vent filter	---	SOP/CSP/006	01	Sarungh	R
5.	Buffer Bottle	---	SOP/CSP/006	01	Sarungh	R
6.	SS forceps	---	SOP/CSP/006	06	Sarungh	R
7.	SS filling needle	---	SOP/CSP/006	04	Sarungh	R
8.	SS syringe set	---	SOP/CSP/006	04 (1 set)	Sarungh	R
9.	Stopper hopper, Stopper rail, chute	---	SOP/CSP/006	01 + 01	Sarungh	R
10.	Needle holding plate	---	SOP/CSP/006	01	Sarungh	R
11.	Nylon Fasteners	---	SOP/CSP/006	15	Sarungh	R
12.	Stopper head set	---	SOP/CSP/006	02	Sarungh	R
13.	Sealing head set with plate	---	SOP/CSP/006	-	Sarungh	R
14.	TC with Nozzle	---	SOP/CSP/006	06	Sarungh	R
15.	Manifold	---	SOP/CSP/006	01	Sarungh	R
16.	SS clamps, gaskets and nozzles	---	SOP/CSP/006	06 + 06 + 06	Sarungh	R

**These equipments to be used product dedicated.

* To be Used Batch Specific

##Production

	Batch Filling Record (Part-A)		Page 21 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

6.8 FILTER HOLDER PREPARATION:

A. R. No. Of WFI: 02WFO8000748

Date: 14/10/08

Note: Preparation of Cartridge filter holder should be done under LAF No.: CSP/059

S. No.	Procedure	SOP No.	Time		Observations	Done by ^{##}	Ckd. By ^{##}
			From	To			
1.	Assemble 4 inch cartridge filter holder under LAF	SOP/CSP/022	21:26	21:30	Cartridge filter housing Eq. Id. No. <u>CSP/113</u>	Sarvesh	
2.	Arrangement of sterilization grade filter: Supplier: Millipore Type: LAGL04TP6 MOC: PVDF Porosity: 0.22 µ Place the filter in cartridge filter holder and Note down the Filter Serial No.	---	21:31	21:35	Filter Serial No.: <u>(7NN)5708</u> <u>0607</u>	Sarvesh	
3.	Carry out Pre integrity of the cartilage filter and note the observation in Step 6.9	SOP/CSP/023	23:07	23:27	53.9 PSI	Sohag	P
4.	Carry out the sterilization of the filtration equipments and accessories. Record the details in table of step 6.11.	SOP/CSP/037	15/10/08 09:20	11:25	OK	Ravinder	

^{##}Production

	Batch Filling Record (Part-A)		Page 22 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

6.9 PRE INTEGRITY TEST OF PRODUCT FILTER (BUBBLE POINT TEST):

- Perform the pre filter integrity of the Product filter as per SOP No. SOP/CSP/023.
- Properly wet the filter by flushing cool WFI before undergoing pre integrity test.
- Bubble point value is NLT 50 psi with WFI.

Date	Filter Catalogue No	Filter Serial No.	Cartridge Filter Holder No.	Observed Bubble Point	Result (Pass/Fail)	Done By ^{##}	Checked By ^{##}
14/10/08	C 8 14/10/08 LOGLO4T86	C7NN15708 0607	CSL/113	53.9	Pass	Sohan	J

6.10 PRE INTEGRITY TEST OF AIR VENT FILTER (BUBBLE POINT TEST):

- Perform the pre integrity test of the Air vent filter for buffer bottle as per SOP No. SOP/CSP/023.
- Properly wet the filter by flushing 60% IPA/WFI mixture before undergoing pre integrity test.
- Bubble point value is NLT 16 psi with 60% IPA/WFI mixture.

Date	Filter Catalogue No	Filter Serial No.	Cartridge Filter Holder No.	Observed Bubble Point	Result (Pass/Fail)	Done By ^{##}	Checked By ^{##}
14/10/08	MTA F15010	F78NO 9210 01965	—	22.48	Pass	Sohan	J 14/10/08

^{##}Production

INTEGRITEST EXACTA

Instrument Serial Number:

W20678

Date:

10/14/2008

Time:

11:07:32PM

Test Serial Number:

200810142307

Test Name:

LAGL BPT

Test Revision:

2

Test Description:

BUBBLE POINT

Bubble Point Test

Filter Catalog Number:

LAGL04TP6

Test Fluid:

Water

Filter Description:

Optiseal Cartridge Filter With
Durapore Membrane , 6/pack

Filter Size:

4"

Pore Size:

0.22

Housing Type:

OPTISEAL

Operator Inputs

Operator:

rai

Product Batch no.

TOPOTECAN INJ. 0PR201

Filter Lot no.

8PK 201

C7NN15708 ✓

User Defined Field #3 - Optional

0607

User Defined Field #4 - Optional

CSP/113 ✓

Notes:

ONLINE FILTER PRE INT.


Test Specifications

of Cartridges: 1
Minimum Bubble Point: 50.0 psig

Test Data

Gas Volume: 449 ml

Pressure psig	Flowrate ml/min
39.9	2.7
49.0	4.6
50.0	4.4
51.0	6.7
53.1	22.7
54.4	74.5
55.1	176.6



Test Results

Measured Bubble Point: 53.9 psig ✓

Immediate Retest

Test Conclusion

PASSED

Signature: _____

Sotun

14/10/08

Signature: _____

14/10/08

Palltronic Flowstar Test Report

14/Oct/2008

23:32:31

Serial Number : 15342626/3.3a
Function : BP
Test Programme : MTGR15010-BP
Production area : PROD-CSL
Filter line : BUFFER BOTTLE FILTER
Product name : TOPOTECAN FOR INJ.
Product batch no. : 8PK201 ✓
Filter part no. : MTGR
Filt. serial no. : F7PN09210 01965
Housing : ---
Wetting liquid : 60:40 IPA WATER
Test gas : AIR
BP minimum : 16.00 psi

BP value measured : 22.48 psi ✓

BP within limits

Comment : PRE-INTEGRITY

Electr. Signature : PRADEEP
Created : 14/Oct/2008 23:39:31

Sohan
14/10/08

14/10/08

14/10/08

	Batch Filling Record (Part-A)		Page 23 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

6.11 STERILIZATION OF EQUIPMENT AND ACCESSORIES FOR ONLINE FILTRATION AND FILLING AS PER SOP No. SOP/CSP/037:

Cycle started at: 09:20

Cycle end at: 11:25

Date	Load No.	Load Pattern No.	Sterilization Hold Time		Done By ^{##}	Checked By ^{##}
			From	To		
15/10/08	02	2	10:16	10:46	Ravinder	MS

(Graph attached on the reverse side of page No.: 23)

6.12 STORAGE OF STERILIZED EQUIPMENT AND ACCESSORIES FOR ONLINE FILTRATION AND FILLING (NMT 48 HOURS AFTER STERILIZATION):

Procedure	Eq. Id. No.	SOP No.	Date	Storage Time*		Done By ^{##}	Checked By ^{##}
				From	To		
Store the sterilized equipment and accessories for on line filtration and Filling under LAF in filtration area.	CSP/036	SOP/CSP/072	15/10/08	11:25	---	Narender	Shruti
			15/10/08	---	21:45	Shruti	15/10/08

* Hold time starts after completion of sterilization and completes at the end of activity.

^{##}Production

Load No-02


MFG.BY :MACHIN FABRIK
CUSTOMER :M/S DABUR PHARMA LTD.
EQUIPMENT :STEAM STERILIZER

BATCH DATA
STERILIZER NO :1007
BATCH NO :0
QUANTITY :1
OPERATOR CODE :1
PRINT INTERVAL :60 SEC
DATE :15/10/08
TIME :09:20:15

PROCESS : MACHINE PARTS

SET PARAMETERS
PRE VACUUM : -0.750 BAR
PRE PRESSURE : 0.500 BAR
NO OF PRE PULSE : 3 NOS
SET POINT-1 : 110.0 DEG C
SET POINT-2 : 115.0 DEG C
SET POINT-3 : 118.0 DEG C
SP-1 HOLD TIME : 2 MIN
SP-2 HOLD TIME : 3 MIN
SP-3 HOLD TIME : 5 MIN
HEAT UP CONTROL BAND : 0.2 DEG C
STER HOLD TEMP : 121.5 DEG C
STER HOLD TIME : 30 MIN
TEMP CONTROL BAND : 0.2 DEG C
OVERSHOOT TEMP : 124.0 DEG C
STER COUNT STOP TEMP : 120.9 DEG C
STER COUNT RESET TEMP : 120.8 DEG C
POST VACUUM : -0.750 BAR
POST PRESSURE : -0.100 BAR
POST VAC HOLD TIME : 30 MIN
NO OF POST PULSE : 1 NOS
EXHAUST ON TIME : 5 SEC
EXHAUST OFF TIME : 50 SEC
PROCESS END PRESSURE : -0.050 BAR
CHAM PRESSURE HIGH : 1.500 BAR

MFG.BY :MACHIN FABRIK
CUSTOMER :M/S DABUR PHARMA LTD.
DATE :15/10/08
TIME :09:20:43
EQUIPMENT :STEAM STERILIZER
PROCESS :MACHINE PARTS



TIME	T1 DEG C	T2 DEG C	T3 DEG C	T4 DEG C	T5 DEG C	PRES BAR	STATUS
09:20:45	028.5	032.7	035.9	029.7	036.1	-0.025	PROCESS ON
09:20:47	028.5	032.7	035.9	029.7	036.1	-0.039	PRE VACUUM
09:21:41	028.0	029.9	033.2	027.0	033.4	-0.587	PRE VACUUM
09:22:20	027.8	028.4	031.8	025.6	032.0	-0.762	PRE PRESSURE
09:22:41	027.9	037.2	036.6	041.2	038.2	-0.674	PRE PRESSURE
09:23:41	028.2	053.8	050.6	056.8	052.3	-0.581	PRE PRESSURE
09:24:42	028.3	063.1	059.6	065.2	059.6	-0.496	PRE PRESSURE
09:25:41	031.2	069.8	066.4	071.6	066.1	-0.412	PRE PRESSURE
09:26:42	044.6	075.1	072.1	076.7	072.0	-0.326	PRE PRESSURE
09:27:41	060.3	079.7	077.5	081.0	077.0	-0.238	PRE PRESSURE
09:28:42	058.9	084.3	082.8	085.4	081.5	-0.142	PRE PRESSURE
09:29:41	055.8	088.0	086.9	089.2	085.9	-0.045	PRE PRESSURE
09:30:42	064.1	091.4	091.0	092.4	089.8	00.052	PRE PRESSURE
09:31:43	087.7	094.7	094.3	094.5	092.4	00.119	PRE PRESSURE
09:32:42	094.7	096.8	096.6	096.9	094.6	00.160	PRE PRESSURE
09:33:42	097.6	098.9	098.8	099.2	097.1	00.210	PRE PRESSURE
09:34:43	100.2	100.9	100.7	099.6	099.4	00.267	PRE PRESSURE
09:35:42	102.3	102.7	102.6	100.6	101.6	00.316	PRE PRESSURE
09:36:43	104.3	104.4	104.4	103.1	103.5	00.373	PRE PRESSURE
09:37:42	106.0	105.9	106.0	106.1	105.3	00.420	PRE PRESSURE
09:38:44	107.6	107.3	107.4	108.0	106.8	00.468	PRE PRESSURE
09:39:22	108.5	108.1	108.2	108.9	107.7	00.498	PRE VACUUM
09:39:43	103.2	103.4	104.3	104.7	102.7	00.005	PRE VACUUM
09:40:44	084.5	085.4	086.3	086.4	085.0	-0.472	PRE VACUUM
09:41:44	071.1	071.5	072.9	072.6	071.5	-0.694	PRE VACUUM
09:42:13	065.3	065.8	067.2	067.1	065.9	-0.758	PRE PRESSURE
09:42:43	072.0	071.3	071.7	075.2	071.6	-0.578	PRE PRESSURE
09:43:44	086.5	086.2	086.3	087.0	086.3	-0.290	PRE PRESSURE
09:44:43	091.8	094.3	094.5	093.7	094.4	-0.085	PRE PRESSURE
09:45:44	100.3	099.5	099.7	099.2	099.4	00.065	PRE PRESSURE
09:46:44	103.7	102.7	102.9	102.9	102.6	00.184	PRE PRESSURE

09:47:43	106.3	105.4	105.6	105.8	105.3	00.222	PRE PRESSURE
09:48:44	108.7	107.8	107.9	108.2	107.6	00.394	PRE PRESSURE
09:49:44	110.5	109.7	109.9	110.2	109.6	00.477	PRE PRESSURE
09:49:58	110.9	110.0	110.2	110.6	110.0	00.500	PRE VACUUM
09:50:44	097.2	098.4	098.9	100.0	098.0	-0.149	PRE VACUUM
09:51:45	083.1	083.4	083.9	087.1	083.4	-0.506	PRE VACUUM
09:52:44	069.9	071.3	072.4	080.0	070.9	-0.707	PRE VACUUM
09:53:07	065.3	067.1	068.3	079.5	064.9	-0.756	PRE PRESSURE
09:53:45	075.7	074.7	074.5	081.4	074.8	-0.526	PRE PRESSURE
09:54:43	088.6	089.3	089.3	089.8	089.3	-0.213	PRE PRESSURE
09:55:44	097.6	097.9	098.2	097.8	097.8	00.034	PRE PRESSURE
09:56:43	103.7	102.5	102.8	102.8	102.4	00.180	PRE PRESSURE
09:57:45	106.9	105.9	106.1	106.3	105.7	00.314	PRE PRESSURE
09:58:45	109.6	108.5	108.8	109.1	108.4	00.428	PRE PRESSURE
09:59:28	111.1	110.1	110.4	110.7	110.0	00.499	HEAT UP
09:59:44	110.5	109.8	110.2	110.3	109.6	00.435	HEAT UP
10:00:43	110.1	109.3	109.5	110.0	109.3	00.459	HEAT UP
10:01:44	110.5	109.6	109.8	110.3	109.6	00.478	HEAT UP
10:02:45	112.8	111.8	112.0	112.4	111.8	00.590	HEAT UP
10:03:46	114.6	113.6	113.8	114.2	113.6	00.682	HEAT UP
10:04:45	115.2	114.3	114.5	114.9	114.2	00.682	HEAT UP
10:05:45	115.3	114.3	114.7	115.0	114.3	00.708	HEAT UP
10:06:46	114.9	114.1	114.5	114.9	114.2	00.703	HEAT UP
10:07:45	116.6	115.6	115.8	116.2	115.6	00.788	HEAT UP
10:08:45	118.1	117.1	117.4	117.7	117.1	00.870	HEAT UP
10:09:45	118.1	117.2	117.5	117.8	117.1	00.841	HEAT UP
10:10:45	118.0	117.2	117.5	117.9	117.2	00.866	HEAT UP
10:11:47	118.0	117.3	117.6	118.0	117.3	00.880	HEAT UP
10:12:47	118.3	117.4	117.7	118.1	117.3	00.869	HEAT UP
10:13:46	118.5	117.6	117.9	118.4	117.6	00.909	HEAT UP
10:14:47	120.1	119.2	119.4	119.9	119.1	01.005	HEAT UP
10:15:47	121.4	120.6	120.7	121.2	120.5	01.086	HEAT UP
✓10:16:53	122.5	121.7	121.8	122.3	121.5	01.152	STER HOLD
10:17:47	122.3	121.7	121.8	122.3	121.5	01.136	STER HOLD
10:18:47	122.1	121.5	121.7	122.1	121.3	01.113	STER HOLD

10:17:46	122.1	121.4	121.5	122.1	121.3	01.134	STER HOLD
10:20:47	122.5	121.7	122.0	122.4	121.6	01.153	STER HOLD
10:21:48	122.3	121.5	121.9	122.3	121.5	01.131	STER HOLD
10:22:47	122.6	121.7	122.1	122.5	121.7	01.138	STER HOLD
10:23:48	122.3	121.6	121.9	122.3	121.5	01.144	STER HOLD
10:24:46	122.4	121.5	122.0	122.3	121.5	01.110	STER HOLD
10:25:47	122.5	121.8	122.0	122.5	121.7	01.159	STER HOLD
10:26:48	122.3	121.7	121.8	122.3	121.5	01.144	STER HOLD
10:27:46	122.3	121.7	121.9	122.3	121.5	01.113	STER HOLD
10:28:47	122.6	121.8	122.0	122.5	121.7	01.143	STER HOLD
10:29:48	122.3	121.7	121.9	122.3	121.5	01.145	STER HOLD
10:30:47	122.3	121.6	122.0	122.3	121.5	01.115	STER HOLD
10:31:48	122.6	121.7	122.1	122.4	121.7	01.135	STER HOLD
10:32:47	122.4	121.7	122.0	122.3	121.6	01.150	STER HOLD
10:33:49	122.3	121.6	121.9	122.3	121.5	01.136	STER HOLD
10:34:47	122.5	121.7	122.0	122.4	121.6	01.124	STER HOLD
10:35:49	122.6	121.8	122.0	122.5	121.7	01.139	STER HOLD
10:36:48	122.3	121.6	121.9	122.4	121.5	01.144	STER HOLD
10:37:49	122.3	121.6	121.9	122.3	121.5	01.126	STER HOLD
10:38:49	122.6	121.8	122.0	122.5	121.7	01.138	STER HOLD
10:39:48	122.5	121.7	122.0	122.5	121.6	01.155	STER HOLD
10:40:49	122.3	121.7	121.9	122.3	121.5	01.135	STER HOLD
10:41:50	122.5	121.8	122.0	122.4	121.6	01.124	STER HOLD
10:42:49	122.6	121.8	122.0	122.5	121.7	01.167	STER HOLD
10:43:50	122.3	121.7	121.9	122.3	121.5	01.142	STER HOLD
10:44:49	122.6	121.8	122.0	122.5	121.6	01.129	STER HOLD
10:45:50	122.6	121.9	122.0	122.6	121.7	01.160	STER HOLD
10:46:49	122.3	121.7	121.9	122.3	121.5	01.145	STER HOLD
10:46:54	122.5	121.7	122.0	122.5	121.6	01.158	POST VACUUM
10:47:50	101.6	103.9	104.4	104.6	102.9	-0.038	POST VACUUM

10:48:51	086.1	087.3	087.1	088.0	087.6	-0.442	POST VACUUM
10:49:50	074.3	075.5	077.1	075.4	085.7	-0.657	POST VACUUM
10:50:31	065.3	071.7	071.0	068.2	087.5	-0.746	VACUUM HOLD

10:50:49	053.4	063.1	060.7	056.9	087.6	-0.841	VACUUM HOLD
10:51:50	047.3	069.7	056.1	050.9	090.4	-0.881	VACUUM HOLD
10:53:50	042.9	076.9	056.5	046.9	091.3	-0.901	VACUUM HOLD
10:54:51	039.8	081.2	070.3	044.3	092.3	-0.913	VACUUM HOLD
10:55:50	038.0	083.3	078.2	043.9	093.1	-0.921	VACUUM HOLD
10:56:50	037.9	084.9	083.6	047.4	093.9	-0.926	VACUUM HOLD
10:57:51	037.2	085.8	087.1	047.6	094.6	-0.930	VACUUM HOLD
10:58:52	036.2	086.6	089.4	047.5	095.1	-0.932	VACUUM HOLD
10:59:51	035.2	086.7	090.9	054.1	095.4	-0.933	VACUUM HOLD
11:00:52	034.2	087.3	092.1	060.5	095.7	-0.936	VACUUM HOLD
11:01:52	033.4	087.7	092.8	065.4	094.8	-0.937	VACUUM HOLD
11:02:51	032.7	088.0	093.4	068.8	095.2	-0.938	VACUUM HOLD
11:03:51	032.0	088.2	093.9	071.1	095.7	-0.939	VACUUM HOLD
11:04:52	031.3	088.4	094.3	072.7	096.0	-0.940	VACUUM HOLD
11:05:53	030.6	088.4	094.5	073.8	096.3	-0.940	VACUUM HOLD
11:06:51	030.1	088.5	094.7	074.4	096.5	-0.940	VACUUM HOLD
11:07:52	029.5	088.5	094.8	074.8	096.8	-0.940	VACUUM HOLD
11:08:53	029.3	088.4	094.9	075.3	096.9	-0.940	VACUUM HOLD
11:09:53	029.5	088.4	095.0	075.7	097.1	-0.940	VACUUM HOLD
11:10:52	029.7	088.4	095.0	076.1	097.2	-0.939	VACUUM HOLD
11:11:52	029.8	088.2	095.0	076.4	097.4	-0.938	VACUUM HOLD
11:12:53	029.9	088.0	095.1	076.7	097.5	-0.936	VACUUM HOLD
11:13:53	029.9	087.8	095.1	077.0	097.7	-0.934	VACUUM HOLD
11:14:54	029.9	087.6	095.1	077.2	097.8	-0.931	VACUUM HOLD
11:15:52	029.9	087.5	095.0	077.2	097.8	-0.929	VACUUM HOLD
11:16:52	030.0	087.3	095.0	077.2	097.8	-0.928	VACUUM HOLD
11:17:54	030.0	087.1	094.9	077.2	097.8	-0.927	VACUUM HOLD
11:18:54	030.0	086.9	094.8	077.2	097.8	-0.927	VACUUM HOLD
11:19:53	030.1	086.7	094.7	077.2	097.8	-0.927	VACUUM HOLD
11:20:51	030.1	086.7	094.6	077.2	097.8	-0.927	POST PRESSURE
11:20:53	030.5	084.5	094.5	077.7	097.8	-0.828	POST PRESSURE
11:21:55	031.7	087.3	093.3	077.4	098.5	-0.519	POST PRESSURE



CHECKLIST FOR FILLING ROOM
QUALITY ASSURANCE DEPARTMENT

PAGE 1 OF 1

Previous Product: *Cytidine Syg. 10mg/ml*

B. No. *SLA 204*

Product Name: *Topotecan for Inj. 1mg/ml*

B. No. *8PK 201*

Equipment No. *CSP/021*

Date: *15/10/08*

CHECKLIST FOR FILLING ROOM			
Sr. No.	Check Points	Remarks (ok / not ok)	Checked By
1.	Ensure that the filling line room for cleanliness and no remnants of previous batch/ product left out in the room.	ok	<i>Singh</i> 15/10/08
2.	Check Ensure that the filling, bunging and sealing cRAB is visibly clean and no remnants of previous batch/ product left out in the room.	ok	<i>Singh</i> 15/10/08
3.	Ensure that pressure differential of cRAB & area within the limit.	ok	<i>Singh</i> 15/10/08
4.	Ensure that Equipment usage log, room cleaning and sanitization logbook are maintained and updated.	ok	<i>Singh</i> 15/10/08
5.	Ensure the "CLEANED" label affixed on the filling, bunging, sealing is cRAB.	ok	<i>Singh</i> 15/10/08
6.	Ensure the temp. and %RH is within the limit.	ok	<i>Singh</i> 15/10/08

QA Signature

Ankur

	Batch Filling Record (Part-A)		Page 24 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

7.0 ONLINE FILTRATION, FILLING, STOPPERING AND SEALING OPERATION:

7.1 CLEANING:

Step No.	Checks	Eq. Id. No.	SOP No.	Date	Time		Done By##	Checked By##
					From	To		
1.	Cleaning of filling, Stoppering and Sealing machine	CSP/021	SOP/CSP/031	15/10/08	11:00	11:35	naresh	<i>[Signature]</i>
2.	Cleaning of cRABs 2a. cRAB1 2b. cRAB2 2c. cRAB3 2d. cRAB4 2e. RAB5	CSP/025 CSP/026 CSP/027 CSP/028 CSP/029	SOP/CSP/054	15/10/08	11:00	11:35	naresh	<i>[Signature]</i>

7.2 LINE CLEARANCE:

- Ensure the cleanliness of Vial Filling room.
- Ensure the line clearance as per the check list.
- Ensure that status label is properly written on door of the room and on holding tank.
- Ensure proper change over and setting is done as per the pack size.
- Ensure that the pressure differential in gauges of the cRAB is within the limit.
- Ensure that the pressure differential, temperature, relative humidity of filling room is within limit.

Date: 15/10/08

Step No.	Checks	Eq. Id. No.	SOP No.	Previous product		Ckd. by##	Verified by###
				Product Name	Batch No.		
1.	Line clearance for filling, stoppering and sealing machine	CSP/021	SOP/QAD/053	cytarabine in 100mg/ml	8LA204	<i>[Signature]</i>	<i>[Signature]</i> 15/10/08

##Production

###IPQA

**Batch Filling Record (Part-A)**

Page 25 of 46

PLANT-II, BADDI

Product Name Topotecan for Injection 4 mg (base)/vial

Batch Number 8PK201

Batch size 12 L

Batch Filling Code

IO720004

Department Production (CSL)

BFR No.

BFR/CSL/TOP/001

Shelf life 24 Months (Provisional)

Version No.

00

7.3 ENVIRONMENTAL MONITORING IN FILLING ROOM:

Frequency: At the start and at the end of activity

S. No.	Checks	SOP No.	Observations			
1.	Date	---	15/10/08	15/10/08		
2.	Time	---	15:20	21:55		
3.	Differential pressure reading on the magnehelic gauge: Filling room to Filling Corridor (Range: 0.5 -5.0 mm of water)	SOP/CSP/011	2.4	2.4		
4.	Filling room temperature (NMT 25°C)	SOP/CSP/011	22	23		
5.	Filling room R.H. (NMT 60%)	SOP/CSP/011	51	56		
6.	Done by ^{##}	---	Shamsher	Shamsher		
7.	Ckd. by ^{##}	---	SP	SP		

7.4 FILL WEIGHT CALCULATION DATA.Specific Gravity of the Product 1.01 as per A.R. No. 02608000084 Dated 15/10/08

Fill volume*	Target fill Volume	Min.	Max.
Range in ml	1.957	1.898	2.015
Range in Wt. (Volume X Specific Gravity of the Product)	1.976	1.916	2.035

Calculated by SP
Production 15/10/08Ckd. by Shamsher
IPQA 15/10/08

* Refer page no. 3.

^{##} ProductionShamsher

	Batch Filling Record (Part-A)		Page 26 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

7.5 ONLINE FILTRATION, FILLING AND HALF STOPPERING OPERATION:

Date: 15/10/08 First Filtration completed at: 17:10
 Filling started at: 19:50 Filling completed at: 21:45
 Duration between first filtration completion time to filling completion time 4 hrs (NMT 48 hours)
35 mins.

Step No.	Checks	SOP No.	Time		Observations	Done By##	Ckd. By##
			From	To			
1.	Note the reading of cRAB1 over infeed turntable on magnehelic gauge in 'ON' condition. (Limit: 10 to 16 Pascal)	SOP/CSP/054	---	---	At Start <u>14</u>	Sanjeev	SP
					At End <u>14</u>	Shamsher	SP
2.	Note the reading of cRAB1 over filling and Stoppering machine on magnehelic gauge in 'ON' condition. (Limit: 30 to 36 Pascal)	SOP/CSP/054	---	---	At Start <u>32</u>	Sanjeev	SP
					At End <u>32</u>	Shamsher	SP
3.	Note the reading of cRAB 2 for Sealing Machine on the magnehelic gauge in 'ON' condition (Limit: 14 to 20 Pascal)	SOP/CSP/054	---	---	At Start <u>16</u>	Sanjeev	SP
					At End <u>16</u>	Shamsher	SP
4.	Transfer the sterilized equipments; rubber stoppers and aluminum flip off overseal from the cooling zone LAF to the filling room in the mobile LAF. And note the reading of ΔP across HEPA on the magnehelic gauge in 'ON' condition (Limit: 10 to 15 mm of WC)	SOP/CSP/018	15:35	15:37	ΔP across HEPA <u>12</u> (mm of WC)	Sanjeev	SP
5.	Ensure the Stoppering unit is assembled.	SOP/CSP/031	15:40	15:43	OK	Sanjeev	SP

##Production

Form No.: SOP/QAD/038-02 Rev. 00

CONFIDENTIAL

QA Signature



	Batch Filling Record (Part-A)		Page 27 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

Step No.	Checks	SOP No.	Time		Observations	Done By ^{##}	Ckd. By ^{##}
			From	To			
6.	Check the initial weight of the tank and note the reading	---	17:36	17:38	Wt = 102.75 Kg	Sanjeev	P
7.	Connect the outlet of holding tank to the inlet of product filter through peristaltic pump and outlet of product filter to the inlet of buffer bottle and connect the outlet of buffer bottle to the manifold having connections for piston and filling needles using Platinum cured silicone tubes.	---	17:42	17:46	OK	Sanjeev	P
8.	Filter the solution and collect in buffer bottle.	SOP/CSP/031	17:49	17:53	OK	Sanjeev	P
9.	Before start of the filling of vials, purge the needles with the product removing the entrapped air from tubes and piston and collect the solution in SS container and discard.	SOP/CSP/031	17:56	18:01	OK	Sanjeev	P
10.	Adjust the volume of syringes as per the data obtained from step 7.4 and record the observation in the table 7.6.1.	SOP/CSP/031	19:37	19:45	OK	Shamshu	P
11.	Load the Rubber stoppers in hoppers.	SOP/CSP/031	19:30	19:32	OK	Shamshu	P

^{##}Production

	Batch Filling Record (Part-A)		Page 28 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

Step No.	Checks	SOP No.	Time		Observations	Done By##	Ckd. By##
			From	To			
12.	Operate the filling and Stoppering machine at validated speed and fill the solution in the vials filtering through online 0.22µ cartridge filter.	SOP/CSP/031	19:50	21:45	50 Vials per minute	Shamshu	8
13.	Check the in process Fill volume at Start and End (during process every 30 min.) of filling from each Filling head (4 vials) and record the details in table 7.6.2	---	---	---	OK	Sanjeev	8
14.	Check the final weight of the tank and note the reading (Left over solution)	---	---	---	Wt = 2.35 Kg 0.00 Kg 15/10/08	Sanjeev	37
15.	Total solution consumed = Initial Wt (Step6) -Final Wt (Step14)	---	---	---	10.4 Kg	Sanjeev	37
16.	Total rejection during filling and Stoppering.	---	---	---	Vials	24	Sanjeev 37
					Rubber stoppers	24	
17.	Send the rejected vials for deactivation and crushing.	SOP/CSP/069	---	---	24 Vials	Sanjeev	37
18.	Filled vials rejected at filling As per table 7.8	---	---	---	24 Vials	Sanjeev	37

##Production

[Signature]

	Batch Filling Record (Part-A)		Page 29 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK20)		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

7.6 INPROCESS CHECKS:

7.6.1 INITIAL FILL WEIGHT AND STOPPERING CHECK:

Specific Gravity of the Product 1.01 as per A.R. No. 020108000084 Dated 15/10/08

Fill volume*	Target fill Volume	Min.	Max.
Range in ml	^{SP} 15/10/08 1.956 1.957	1.898	2.015
Range in Wt. (Volume X Specific Gravity of the Product)	1.976	1.916	2.035

Date	Time	Filled Weight (g)				Initial Half Stoppering		Done By ^{##}	Ckd. By ^{##}
		N ₁	N ₂	N ₃	N ₄	S ₁	S ₂		
15/10/08	19:37	1.708	1.732	1.753	1.795	OK	OK	Shambhu	SP
15/10/08	19:45	1.864	1.873	1.938	1.997	OK	OK	Shambhu	SP

N₁ = nozzle 1, N₂ = nozzle 2, N₃ = nozzle 3, N₄ = nozzle 4, S₁ = Stoppering head 1, S₂ = Stoppering head 2
(Weights of filled vials sheet attached on the reverse side of page No.: 29)

* Refers page no. 3.

^{##}Production

[Signature]

PRODUCT NAME:- TOPOTECAN INJ. 4MG/2ML.
PRODUCT. BATCH.NO: 8P14201.

INITIAL FILL WEIGHT. AF CHECKS.

15.10.2008 19:37 ✓

1.708

12.432 g

1.732

12.432 g

1.753

12.419 g

1.785

12.378 g

12.216 g

12.219 g

12.273 g

15/10/08

15/10/08

15.10.2008 19:45

1.964

12.342 g

1.973

12.349 g

1.938

12.385 g

1.997

12.349 g

12.378 g

12.376 g

12.367 g

12.352 g

15/10/08

PRODUCT NAME: TOPOTECAN INJ. 4MG/2ML
 BATCH NO: 8PK201
 IN PROCESS. FILL WEIGHT CHECKS

----- Weighing -----
 15 2 15.4g 15/10/08

12.283 g
 2.000 10.283 g
 12.290 g
 2.011 10.291 g
 12.212 g
 2.027 10.212 g
 1.956 10.256 g

hmp
 15/10/08

q
 15/10/08

----- Weighing -----
 15 2 15.4g 15/10/08

12.151 g
 1.992 10.159 g
 12.335 g
 1.983 10.157 g
 12.101 g
 1.965 10.101 g
 12.248 g
 2.019 10.221 g

q
 15/10/08

hmp
 15/10/08

----- Weighing -----
 15 2 15.4g 15/10/08

12.227 g
 1.945 10.282 g
 12.855 g
 1.995 10.960 g
 12.342 g
 1.956 10.367 g
 12.102 g
 2.000 10.102 g

hmp
 15/10/08

q
 15/10/08

----- Weighing -----
 15 2 15.4g 15/10/08

12.369 g
 1.993 10.376 g
 12.862 g
 1.959 10.183 g
 12.317 g
 1.929 10.360 g
 12.181 g
 1.952 10.229 g

hmp
 15/10/08

q
 15/10/08



Batch Filling Record (Part-A)

Page 31 of 46

PLANT-II, BADDI

Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	BPK 201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

7.7 INTERVENTIONS OBSERVED DURING FILLING:

Date	Time		Interventions	Tray No.	Ckd. By ^{##}
	From	To			
15/10/08	20:30	21:13	M/C stopped for dinner break	03	J
N/A					

7.8 FILLED VIAL REJECTION DETAILS

Time	No of vials	Type of Rejection	Ckd. By ^{##}
19:31	04	Initial fill wt. check	J
19:45	04	Initial fill wt. check	J
19:50	04	In process fill wt. check	J
20:20	04	In process fill wt. check	J
21:15	04	In process fill wt. check	J
21:45	04	In process fill wt. check	J
N/A			

^{##}Production

		Batch Filling Record (Part-A)		Page 32 of 46
PLANT-II, BADDI				
Product Name	Topotecan for Injection 4 mg (base)/vial			
Batch Number	8PK 201			
Batch size	12 L	Batch Filling Code	IO720004	
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001	
Shelf life	24 Months (Provisional)	Version No.	00	

7.9 POST INTEGRITY TEST OF AIR VENT FILTER (WATER INTRUSION TEST/BUBBLE POINT TEST):

- Perform the post integrity test of the Air vent filter of holding tank as per SOP No SOP/CSP/023.
- Flow rate should be NMT 0.12mL/minute.
- Bubble point value of this filter is NLT 20psi with IPA.

Date	Filter Catalogue No.	Filter Serial No.	Cartridge Filter Holder No.	Flow Rate (mL/min.) / Bubble Point Pressure (psi)**	Result (Pass/Fail)	Done By ^{##}	Ckd. By ^{##}
15/10/08	SLK7002PFRP	FP9649 0387	CSP/016	0.10 ml/min	Pass	sohan	sp 15/10/08

7.10 POST INTEGRITY TEST OF PRODUCT FILTER (BUBBLE POINT TEST):

- Perform the post integrity test of the Product filter as per SOP No. SOP/CSP/023.
- Properly flush the filter with cool WFI before undergoing post integrity test.
- Bubble point value is NLT 50 psi with WFI.

Date	Filter Catalogue No.	Filter Serial No.	Cartridge Filter Holder No.	Observed Bubble Point	Result (Pass/Fail)	Done By ^{##}	Checked By ^{##}
15/10/08	LAGL04TP6	C7MN15708 0607	CSP/113	57.29	Pass	sohan	sp 15/10/08

7.11 POST INTEGRITY TEST OF AIR VENT FILTER (BUBBLE POINT TEST):

- Perform the post integrity test of the Air vent filter for buffer bottle as per SOP No SOP/CSP/023.
- Properly wet the filter by flushing 60% IPA/WFI mixture before undergoing post integrity test.
- Bubble point value is NLT 16 psi with 60% IPA/WFI mixture.

Date	Filter Catalogue No.	Filter Serial No.	Cartridge Filter Holder No.	Observed Bubble Point	Result (Pass/Fail)	Done By ^{##}	Checked By ^{##}
15/10/08	MTGR15010 -BP	F7PN092 1001965	-	22.48	Pass	sohan	sp 15/10/08

Note: In case of any technical problem with the integrity test apparatus, the integrity should be performed manually as per SOP No. SOP/CSP/023

** Strike off which is not applicable

(Attach the print out of the process on the backside of the page)

^{##} Production

^{###} IPQA

Palltronic Flowstar Test Report

15/Oct/2008 22:57:09
Serial Number : 15342626/3.3a
Function : BP
Test Programme : MTGR15010-BP
Production area : PROD-CSL
Filter line : BUFFER BOTTLE FILTER
Product name : TOPOTECAN FOR INJ.
Product batch no. : 8PK201
Filter part no. : MTGR
Filt. serial no. : F7PN09210 01965
Housing : -----
Wetting liquid : 60:40 IPA:WATER
Test gas : AIR
BP minimum : 16.00 psi

BP value measured : 22.48 psi

BP within limits

Comment : POST-INTEGRITY

Electr. Signature : PRADEEP
Created : 15/Oct/2008 23:04:11

Sohan 15/10/08 *****

15/10/08

15/10/08

Palltronic Flowstar Test Report

15/Oct/2008

22:30:51

Serial Number : 15342626/3.3a
Function : BP
Test Programme : LAGL04TP6
Production area : PROD-CSL
Filter line : ONLINE FILTER
Product name : TOPOTECAN FOR INJ.
Product batch no. : 8PK201
Filter part no. : LAGL04TP6 ✓
Filt. serial no. : C7MN15708 0607
Housing : CSP/113
Wetting liquid : WATER
Test gas : AIR
BP minimum : 50.00 psi

BP value measured : 57.29 psi ✓

BP within limits

Comment : POST-INTEGRITY

Electr. Signature : PRADEEP

Created : 15/Oct/2008 22:46:32

Sohan
15/10/08

Palltronic Flowstar Test Report

15/Oct/2008

23:07:12

Serial Number : 15342626/3.3a
Function : WIT
Test Programme : SLK7002PFRP-WIT ✓
Production area : PROD-CSL
Filter line : AIR VENT FILTER
Product name : TOPOTECAN FOR INJ.
Product batch no. : 8NJ201
Filter part no. : SLK7002PFRP ✓
Filt. serial no. : FP9649 0387
Housing : CSP/016
Test pressure : 36.25 psi
Test time : 236 sec (Auto)
Max. flow : 0.12 ml/min

Flow : 0.10 ml/min ✓

Flow within limits

Comment : POST-INTEGRITY
HOLDING TANK

Electr. Signature : PRADEEP
Created : 15/Oct/2008 23:19:11

Handwritten signature
15/10/08

Handwritten signature
15/10/08

Handwritten signature
15/10/08

	Batch Filling Record (Part-A)		Page 33 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	SPK 201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

7.12 POST CLEANING OF EQUIPMENTS AND ACCESSORIES FOR ONLINE FILTRATION AND FILLING:

Clean the equipments and accessories for online filtration and filling with the WFI in the wash table lying in the equipment wash room (CSLG 014).

Cleaning started at: 22:58 (15/10/08) Cleaning completed at: 23:58

S.No.	Equipment Name	Equipment Id. No.	SOP No.	Qty.	Done By ^{##}	Ckd. By ^{##}
1.	SS Beakers 500 mL, 1000 mL	---	SOP/CSP/006	01+01		
2.	Cartridge Filter Holder with stand	---	SOP/CSP/006	01		
3.	**50 mm disc type Vent filter	---	SOP/CSP/006	01		
4.	Buffer Bottle	---	SOP/CSP/006	02		
5.	SS forceps	---	SOP/CSP/006	08-		
6.	SS filling needle	---	SOP/CSP/006	04		
7.	SS syringe set	---	SOP/CSP/006	04		
8.	Stopper hopper, Stopper rail, chute	---	SOP/CSP/006	01+01 01+01		
9.	Needle holding plate	---	SOP/CSP/006	01		
10.	Nylon Fasteners	---	SOP/CSP/006	-		
11.	Stopper head set	---	SOP/CSP/006	02		
12.	TC with Nozzle	---	SOP/CSP/006	-		
13.	Manifold	---	SOP/CSP/006	01		
14.	SS clamps, gaskets and nozzles	---	SOP/CSP/006	06+06+06		

** These equipments to be used product dedicated

^{##}Production

	Batch Filling Record (Part-A)		Page 34 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK 201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

8.0 LYOPHILIZATION

8.1 Perform the vacuum leak test as per SOP No. SOP/CSP/032

Date : 14/10/08

Time : 03:05

Vacuum Achieved : 5 μ bar.

Hold Time : 60 min.

Vacuum Drop : 2 μ bar


Result : PASS / FAH

Note:- Allowable pressure rise in system should not be more than 6 μ bars.

The system should pass the leak rate not more than 0.01 mbar/litre/sec.

8.2 CIP of Lyophilizer

8.2.1 Line Clearance for CIP

S.No.	Checks	Room No.	SOP No.	Date	Previous Product		Ckd by ^{##}	Verified by ^{###}
					Product Name	Batch no.		
1.	Line Clearance of lyophilizer	CSLG/034	SOP/QAD/053	13/10/08	Topotecan Placebo	-		 13/10/08

8.2.2 Ensure the CIP of lyophilizer is complete using WFI as per SOP No. SOP/CSP/032.

8.2.3 A.R. No. of WFI 02WFE800748

8.2.4 Record the details of the operation in the table below.

Date	Details	Eq. Id. No.	Cleaning		Done by ^{##}	Checked by ^{##}
			Started at	Completed at		
14/10/08	Lyophiliser	CSL/040	06:29	08:08		

(Attach the print out of the process on the reverse side of page)

^{##}Production

Form No.: SOP/QAD/038-02 Rev. 00

CONFIDENTIAL

QA Signature



Leak Rate Test Process Report

Initiator:	SUP	Start Time:	14/10/2008 3:05:14
Batch Name:	LRT20081014030514	End Time:	14/10/2008 5:25:53

PARAMETERS:

Test vacuum: 6 ubar
Test time: 60 Min.
Drying pressure: 100 ubar
Drying time: 30 Min.
Max time to reach dry pressure: 60 Min.
Retry number: 3

RESULT:

No.: 1

Actual Pressure Rise: 2 ubar

Test Result: Success

done
14/10/08

Res
14/10/08
done
14/10/08

Leak Rate Test Process Report

Initiator:	SUP	Start Time:	14/10/2008 3:05:14
Batch Name:	LRT20081014030514	End Time:	14/10/2008 5:25:53

Date	Time	Description
Alarm: Total: 3		
14/10/2008	3:05:39	Compressor 1 cooling water low flowrate
14/10/2008	3:20:02	Vacuum pump 1 failure
14/10/2008	3:20:10	PT301 pump vacuum probe fail
Audit: Total: 5		
14/10/2008	3:05:15	Vacuum leakage test is enabled by Sup.
14/10/2008	3:20:23	SUP logged in as Application User by Sup.
14/10/2008	3:20:29	A_PUMPVACPROBEFAIL ALARM is acknowledged by Sup.
14/10/2008	3:20:29	A_VP1FAIL ALARM is acknowledged by Sup.
14/10/2008	3:20:29	A_COM1COOLWATERLOWPRES ALARM is acknowledged by Sup.

Sup
14/10/08

Ray
14/10/08
by me
14/10/08

Vacuum Leak Test

Pass

Rej 14/10/08

14/10/08

01-06ch	26.4	26.3	26.3	26.1	15	85	24	1	21	155
07-12ch	26.3	26.2	26.3	26.0	17	83	23	1	21	155
13-18ch	28.2	-75.7	40	25.4	24.5	60	28.9	80	100	
19-24ch	28.4	-48	0.06	0.05	---					
-100.0										150.0
05:00	60mm/h						CH 6			°C

01-06ch	26.1	26.0	26.4	26.0	25.9					
07-12ch	26.0	26.0	25.9	25.7	24.0					
13-18ch	28.1	-75.4	40	25.4	24.4	29.4				
19-24ch	29.1	-263	0.06	0.05	---					
-100.0										150.0
04:00	2008/10/14						CH 5			°C

2008/10/14 03:02 TREND+TPE 60mm/h / 030

2008/10/14

PART NO. 81407861-101

CIP Process Report

Initiator:	SUP	Start Time:	14/10/2008 6:29:16
Batch Name:	CIP20081014062916	End Time:	14/10/2008 8:08:23

PARAMETERS:

Chamber Wash Time:	2	Min.
Condenser Wash Time:	2	Min.
Drain time between mainfolds:	60	Sec.
Number of washes:	2	
Final drain time:	5	Min.
System drying time:	30	Min.
Flooding time:	1	Min.
Flooding Number :	1	

Batch Name: CIP20081014062916

Page 1 of 2

Arvind
14/10/08
lyn
14/10/08
[Signature]
14/10/08

CIP Process Report

Initiator:	SUP	Start Time:	14/10/2008 6:29:16
Batch Name:	CIP20081014062916	End Time:	14/10/2008 8:08:23

Date	Time	Description
Alarm: Total: 1		
14/10/2008	6:51:13	Compressor 1 cooling water low flowrate
Audit: Total: 3		
14/10/2008	6:29:17	CIP cycle is enabled by Sup.
14/10/2008	6:51:57	SUP logged in as Application User by Sup.
14/10/2008	6:52:01	A_COM1COOLWATERLOWPRES ALARM is acknowledged by Sup.

Batch Name: CIP20081014062916

Page 2 of 2

[Signature]
14/10/08

[Signature]
14/10/08
14/10/08

Manufactured By : Adam Fabriwerks Pvt. Ltd.
 Equipment : WIP Module
 Client : Dabur Pharma Ltd. Unit II, Baddi
 SIP OF WIP VESSEL SELECTED
 Operator : A
 Date : 13/10/2008 Time : 9:51

Sterilization Temperature ✓ 123.0 degC
 Stop Temperature ✓ 122.0 degC
 Reset Temperature ✓ 121.0 degC
 Sterilization Time ✓ 30 mins

TIME	STER TIME (mins:secs)	VESSEL TEMP (degC)	VES OUTLET TEMP (degC)	PROCESS STATUS
9:51	0:0	37.0	36.6	WAITING FOR TEMP.
9:52	0:0	62.6	46.8	WAITING FOR TEMP.
9:53	0:0	69.8	57.6	WAITING FOR TEMP.
9:54	0:0	74.2	66.2	WAITING FOR TEMP.
9:55	0:0	77.8	73.0	WAITING FOR TEMP.
9:56	0:0	82.4	78.6	WAITING FOR TEMP.
9:57	0:0	86.8	84.2	WAITING FOR TEMP.
9:58	0:0	91.4	89.2	WAITING FOR TEMP.
9:59	0:0	94.6	93.2	WAITING FOR TEMP.
10:00	0:0	99.2	98.2	WAITING FOR TEMP.
10:01	0:0	103.2	102.6	WAITING FOR TEMP.
10:02	0:0	107.0	106.6	WAITING FOR TEMP.
10:03	0:0	110.0	110.0	WAITING FOR TEMP.
10:04	0:0	112.4	112.6	WAITING FOR TEMP.
10:05	0:0	114.2	114.8	WAITING FOR TEMP.
10:06	0:0	115.6	116.6	WAITING FOR TEMP.
10:07	0:0	116.4	117.8	WAITING FOR TEMP.
10:08	0:0	117.6	118.6	WAITING FOR TEMP.
10:09	0:0	118.0	119.2	WAITING FOR TEMP.
10:10	0:0	118.0	119.2	WAITING FOR TEMP.
10:11	0:0	118.8	120.2	WAITING FOR TEMP.
10:12	0:0	119.0	120.2	WAITING FOR TEMP.
10:13	0:0	119.4	120.8	WAITING FOR TEMP.
10:14	0:0	119.6	121.0	WAITING FOR TEMP.
10:15	0:0	119.6	121.0	WAITING FOR TEMP.
10:16	0:0	120.2	121.8	WAITING FOR TEMP.
10:17	0:0	120.4	121.8	WAITING FOR TEMP.
10:18	0:0	120.6	122.0	WAITING FOR TEMP.
10:19	0:0	120.8	122.2	WAITING FOR TEMP.
10:20	0:0	121.0	122.4	WAITING FOR TEMP.
10:21	0:0	121.2	122.6	WAITING FOR TEMP.
10:22	0:5	121.4	123.0	STER. TIMER RUNNING
10:23	1:5	121.4	122.8	STER. TIMER RUNNING
10:24	2:5	121.4	122.8	STER. TIMER RUNNING
10:25	3:5	121.4	122.8	STER. TIMER RUNNING
10:26	4:5	121.4	123.0	STER. TIMER RUNNING
10:27	5:5	121.4	123.0	STER. TIMER RUNNING
10:28	6:5	121.4	122.8	STER. TIMER RUNNING
10:29	7:5	121.4	122.8	STER. TIMER RUNNING
10:30	8:5	121.6	123.0	STER. TIMER RUNNING
10:31	9:5	121.4	123.0	STER. TIMER RUNNING
10:32	10:5	121.4	123.0	STER. TIMER RUNNING
10:33	11:5	121.4	122.8	STER. TIMER RUNNING

10:30	11:5	121.4	122.0	STER. TIMER RUNNING
10:34	12:5	121.6	123.0	STER. TIMER RUNNING
10:35	13:5	121.4	123.0	STER. TIMER RUNNING
10:36	14:5	121.4	123.0	STER. TIMER RUNNING
10:37	15:5	121.4	123.0	STER. TIMER RUNNING
10:38	16:5	121.4	122.8	STER. TIMER RUNNING
10:39	17:5	121.4	123.0	STER. TIMER RUNNING
10:40	18:5	121.4	123.0	STER. TIMER RUNNING
10:41	19:5	121.4	123.0	STER. TIMER RUNNING
10:42	20:5	121.4	123.0	STER. TIMER RUNNING
10:43	21:5	121.4	123.0	STER. TIMER RUNNING
10:44	22:5	121.4	123.0	STER. TIMER RUNNING
10:45	23:5	121.4	123.0	STER. TIMER RUNNING
10:46	24:5	121.4	123.0	STER. TIMER RUNNING
10:47	25:5	121.4	123.0	STER. TIMER RUNNING
10:48	26:5	121.4	123.0	STER. TIMER RUNNING
10:49	27:5	121.4	123.0	STER. TIMER RUNNING
10:50	28:5	121.4	123.0	STER. TIMER RUNNING
10:51	29:5	121.4	123.0	STER. TIMER RUNNING
13/10/2008	10:52	SIP Of WIP Vessel Process Completed Successfully		

DONE BY

Arine
13/10/08

W
13/10/08

121.4
13/10/08

CHECKED BY
13/10/08

Manufactured By : Adam Fabriwerks Pvt. Ltd.
 Equipment : WIP Module
 Client : Dabur Pharma Ltd. Unit II, Baddi
 Operator : A
 P VESSEL RINSE CYCLE SELECTED
 Date : 13/10/2008 Time : 9:00

Quantity : 500 Ltrs

ME	WFI QTY.	VESSEL TEMP (degC)	PROCESS STATUS
01	0	24.6	WIP VESSEL RINSE IN PROCESS
02	80	24.4	WIP VESSEL RINSE IN PROCESS
03	140	24.2	WIP VESSEL RINSE IN PROCESS
04	196	24.2	WIP VESSEL RINSE IN PROCESS
05	240	24.0	WIP VESSEL RINSE IN PROCESS
06	284	24.0	WIP VESSEL RINSE IN PROCESS
07	320	24.0	WIP VESSEL RINSE IN PROCESS
08	364	24.0	WIP VESSEL RINSE IN PROCESS
09	400	23.8	WIP VESSEL RINSE IN PROCESS
10	436	23.8	WIP VESSEL RINSE IN PROCESS
11	468	24.0	WIP VESSEL RINSE IN PROCESS
12	500	24.0	WIP WATER LEVEL REACHED
12	500	24.0	WIP DRAIN IN PROCESS
13	276	24.0	WIP DRAIN IN PROCESS
14	180	23.8	WIP DRAIN IN PROCESS
15	84	23.8	WIP DRAIN IN PROCESS
16	0	23.6	WIP DRAIN IN PROCESS
17	0	23.6	WIP DRAIN IN PROCESS
/10/2008	9:17	Rinse Cycle Process Completed Successfully	

ONE BY

Amir
13/10/08

Amir
13/10/08

Amir
CHECKED BY
13/10/08

Manufactured By : Adam Fabriwerks Pvt. Ltd.
 Equipment : WIP Module
 Client : Dabur Pharma Ltd. Unit II, Baddi
 WIP OF LYOPHILISER SELECTED
 Operator : A
 Date : 14/10/2008 Time : 4:49

WFI Quantity : 2000 Ltrs
 Temperature : 70.0 degC
 Delay For Drain : 15 secs

TIME	WFI QTY.	VESSEL TEMP (degC)	PROCESS STATUS
4:49	0	55.8	WATER LEVEL MAKEUP IN PROCESS
4:50	0	48.4	WATER LEVEL MAKEUP IN PROCESS
4:51	0	43.4	WATER LEVEL MAKEUP IN PROCESS
4:52	0	39.6	WATER LEVEL MAKEUP IN PROCESS
4:53	116	37.0	WATER LEVEL MAKEUP IN PROCESS
4:54	204	35.6	WATER LEVEL MAKEUP IN PROCESS
4:55	260	34.0	WATER LEVEL MAKEUP IN PROCESS
4:56	312	33.0	WATER LEVEL MAKEUP IN PROCESS
4:57	364	32.6	WATER LEVEL MAKEUP IN PROCESS
4:58	404	32.8	WATER LEVEL MAKEUP IN PROCESS
4:59	448	32.8	WATER LEVEL MAKEUP IN PROCESS
5:00	492	36.2	WATER LEVEL MAKEUP IN PROCESS
5:01	532	36.2	WATER LEVEL MAKEUP IN PROCESS
5:02	560	35.8	WATER LEVEL MAKEUP IN PROCESS
5:03	604	35.4	WATER LEVEL MAKEUP IN PROCESS
5:04	644	35.0	WATER LEVEL MAKEUP IN PROCESS
5:05	676	34.6	WATER LEVEL MAKEUP IN PROCESS
5:06	716	34.2	WATER LEVEL MAKEUP IN PROCESS
5:07	748	34.0	WATER LEVEL MAKEUP IN PROCESS
5:08	792	33.6	WATER LEVEL MAKEUP IN PROCESS
5:09	828	33.4	WATER LEVEL MAKEUP IN PROCESS
5:10	864	33.2	WATER LEVEL MAKEUP IN PROCESS
5:11	900	33.0	WATER LEVEL MAKEUP IN PROCESS
5:12	944	32.8	WATER LEVEL MAKEUP IN PROCESS
5:13	980	32.6	WATER LEVEL MAKEUP IN PROCESS
5:14	1020	32.4	WATER LEVEL MAKEUP IN PROCESS
5:15	1048	32.4	WATER LEVEL MAKEUP IN PROCESS
5:16	1096	32.2	WATER LEVEL MAKEUP IN PROCESS
5:17	1128	32.0	WATER LEVEL MAKEUP IN PROCESS
5:18	1152	31.8	WATER LEVEL MAKEUP IN PROCESS
5:19	1188	31.8	WATER LEVEL MAKEUP IN PROCESS
5:20	1236	31.6	WATER LEVEL MAKEUP IN PROCESS
5:21	1272	31.4	WATER LEVEL MAKEUP IN PROCESS
5:22	1308	31.4	WATER LEVEL MAKEUP IN PROCESS
5:23	1344	31.2	WATER LEVEL MAKEUP IN PROCESS
5:24	1380	31.2	WATER LEVEL MAKEUP IN PROCESS
5:25	1416	31.2	WATER LEVEL MAKEUP IN PROCESS
5:26	1456	31.0	WATER LEVEL MAKEUP IN PROCESS
5:27	1488	31.0	WATER LEVEL MAKEUP IN PROCESS
5:28	1528	30.8	WATER LEVEL MAKEUP IN PROCESS
5:29	1560	30.8	WATER LEVEL MAKEUP IN PROCESS
5:30	1600	30.6	WATER LEVEL MAKEUP IN PROCESS
5:31	1636	30.6	WATER LEVEL MAKEUP IN PROCESS
5:32	1672	30.6	WATER LEVEL MAKEUP IN PROCESS
5:33	1712	30.6	WATER LEVEL MAKEUP IN PROCESS
5:34	1744	30.4	WATER LEVEL MAKEUP IN PROCESS
5:35	1784	30.4	WATER LEVEL MAKEUP IN PROCESS
5:36	1820	30.4	WATER LEVEL MAKEUP IN PROCESS
5:37	1856	30.2	WATER LEVEL MAKEUP IN PROCESS
5:38	1892	30.2	WATER LEVEL MAKEUP IN PROCESS
5:39	1928	30.2	WATER LEVEL MAKEUP IN PROCESS
5:40	1964	30.2	WATER LEVEL MAKEUP IN PROCESS
5:41	2000	30.2	WIP WATER LEVEL REACHED
5:41	1960	30.2	WATER LEVEL MAKEUP IN PROCESS

5:42	2240	30.2	WAITING FOR TEMP.
5:43	2116	30.8	WAITING FOR TEMP.
5:44	1988	31.8	WAITING FOR TEMP.
5:45	1960	31.6	WAITING FOR TEMP.
5:46	1968	33.0	WAITING FOR TEMP.
5:47	1956	33.0	WAITING FOR TEMP.
5:48	1972	35.0	WAITING FOR TEMP.
5:49	1984	36.2	WAITING FOR TEMP.
5:50	1960	37.0	WAITING FOR TEMP.
5:51	1996	38.0	WAITING FOR TEMP.
5:52	1968	39.8	WAITING FOR TEMP.
5:53	1980	40.8	WAITING FOR TEMP.
5:54	1984	41.4	WAITING FOR TEMP.
5:55	1980	43.2	WAITING FOR TEMP.
5:56	1968	44.2	WAITING FOR TEMP.
5:57	2000	45.0	WAITING FOR TEMP.
5:58	1964	46.4	WAITING FOR TEMP.
5:59	1996	47.0	WAITING FOR TEMP.
6:00	1972	48.8	WAITING FOR TEMP.
6:01	1984	49.0	WAITING FOR TEMP.
6:02	1984	50.2	WAITING FOR TEMP.
6:03	2024	50.6	WAITING FOR TEMP.
6:04	2004	52.6	WAITING FOR TEMP.
6:05	1960	53.8	WAITING FOR TEMP.
6:06	1984	54.4	WAITING FOR TEMP.
6:07	2016	55.8	WAITING FOR TEMP.
6:08	2060	56.2	WAITING FOR TEMP.
6:09	2040	57.4	WAITING FOR TEMP.
6:10	2004	59.0	WAITING FOR TEMP.
6:11	2024	59.6	WAITING FOR TEMP.
6:12	2016	61.0	WAITING FOR TEMP.
6:13	2072	62.0	WAITING FOR TEMP.
6:14	2036	62.8	WAITING FOR TEMP.
6:15	2036	64.0	WAITING FOR TEMP.
6:16	2028	65.0	WAITING FOR TEMP.
6:17	2056	65.8	WAITING FOR TEMP.
6:18	2080	67.2	WAITING FOR TEMP.
6:19	2040	67.8	WAITING FOR TEMP.
6:20	2000	67.6	WAITING FOR TEMP.
6:21	1896	69.2	WAITING FOR TEMP.
6:21	1896	70.0	WIP WATER TEMP. REACHED
6:22	1876	70.0	WATER RECIRCULATION
6:23	1892	71.0	WATER RECIRCULATION
6:24	1912	71.0	WATER RECIRCULATION
6:25	1944	71.4	WATER RECIRCULATION
6:26	1928	71.4	WATER RECIRCULATION
6:27	1924	71.6	WATER RECIRCULATION
6:28	1920	71.6	WATER RECIRCULATION
6:29	1920	71.6	WATER RECIRCULATION
6:30	1924	71.8	WATER RECIRCULATION
6:31	1916	71.6	WATER RECIRCULATION
6:32	1920	71.6	WATER RECIRCULATION
6:33	1920	71.6	WATER RECIRCULATION
6:34	1916	71.6	WATER RECIRCULATION
6:35	1920	71.6	WATER RECIRCULATION
6:36	1916	71.6	WATER RECIRCULATION
6:37	1920	71.8	WATER RECIRCULATION
6:38	1920	71.8	WATER RECIRCULATION
6:39	1916	71.8	WATER RECIRCULATION
6:40	1916	71.8	WATER RECIRCULATION
6:41	1912	71.8	WATER RECIRCULATION
6:41	1912	71.8	WIP PROCESS STARTED
6:42	1580	71.8	WIP IN PROCESS
6:43	1736	71.8	WATER RECIRCULATION
6:44	1724	71.8	WIP PROCESS STARTED
6:44	1544	71.8	WIP IN PROCESS
6:45	1440	71.8	WIP PROCESS STARTED
6:45	1360	71.8	WIP IN PROCESS
6:46	1228	71.6	WIP IN PROCESS
6:47	1092	71.8	WIP IN PROCESS
6:48	1284	71.6	WATER RECIRCULATION
6:49	1252	71.6	WIP PROCESS STARTED
6:49	1280	71.6	WIP PROCESS STARTED

6:48	1000	71.6	WIP IN PROCESS
6:50	872	71.6	WIP IN PROCESS
6:51	768	71.6	WIP IN PROCESS
6:52	848	71.6	WIP IN PROCESS
6:53	968	71.6	WATER RECIRCULATION
6:54	944	71.6	WIP PROCESS STARTED
6:54	716	71.6	WIP IN PROCESS
6:55	568	71.6	WIP IN PROCESS
6:56	744	71.4	WATER RECIRCULATION
6:57	720	71.4	WIP PROCESS STARTED
6:57	556	71.4	WIP IN PROCESS
6:58	556	71.4	WIP IN PROCESS
6:59	524	71.4	WIP IN PROCESS
7:00	476	71.4	WIP IN PROCESS
7:01	560	71.4	WATER RECIRCULATION
7:02	548	71.4	WIP PROCESS STARTED
7:02	308	71.4	WIP IN PROCESS
7:03	416	71.2	WIP IN PROCESS
7:04	340	70.8	WIP IN PROCESS
7:05	300	70.6	WIP IN PROCESS
7:06	400	70.8	WATER RECIRCULATION
7:07	372	70.8	WATER RECIRCULATION
7:08	360	70.8	WATER RECIRCULATION
7:09	356	70.8	WATER RECIRCULATION
7:10	348	70.8	WATER RECIRCULATION
7:11	352	70.8	WATER RECIRCULATION
7:12	348	70.8	WATER RECIRCULATION
7:13	352	70.8	WATER RECIRCULATION
7:14	356	70.8	WATER RECIRCULATION
7:15	288	70.8	WIP DRAIN IN PROCESS
7:16	224	70.4	WIP DRAIN IN PROCESS
7:17	204	70.2	WIP DRAIN IN PROCESS
7:18	168	70.0	WIP DRAIN IN PROCESS
7:19	124	70.0	WIP DRAIN IN PROCESS
7:20	76	69.6	WIP DRAIN IN PROCESS
7:21	8	69.4	WIP DRAIN IN PROCESS
7:22	0	69.2	WIP DRAIN IN PROCESS
14/10/2008	7:23		WIP Of Lyophiliser Process Completed Successfully

DONE BY

Amir
14/10/08

Amir
14/10/08

CHECKED BY

Amir
14/10/08

	Batch Filling Record (Part-A)		Page 35 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	SPK 201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

8.3 FILTER INTEGRITY TEST

8.3.1 PRE INTEGRITY TEST OF AIR VENT FILTER (HYDRO CORR TEST/BUBBLE POINT TEST)

- Prior to startup of the SIP, perform the pre filter integrity of the vent filter to be used on Lyophilizer as per SOP No. SOP/CSF/023 and record the details in the table below.
- Flow rate should be NMT 0.7mL/minute in case of 5 inch cartridge filter MDI.
- Flow rate should be NMT 0.2mL/minute in case of 4 inch cartridge filter Millipore.
- Bubble point value is NLT 16 psi with 70:30 IPA:WFI.

Date	Filter Catalogue No.	Filter Serial No.	Cartridge Filter Holder No.	Flow Rate (mL/min.) / Bubble-Point Pressure (Psi) **	Result (Pass/Fail)	Done By ^{###}	Ckd. By ^{###}
13/10/08	C6MNG7312/0189	0189	BM7AN0008K CSF/086	0.05 mL/min	Passed	Raj	Wm
14/10/08	CPTF5301ASS101	CT0198-4	M370750003	0.14 mL/min	Passed	Raj	Am
14/10/08	CPTF5301ASS101	CT0198-3	M370750015	0.13 mL/min	Passed	Raj	Wm

(Attach the print out of the process on the reverse side of page.)

** Strike off which is not applicable.

INTEGRITEST EXACTA

Instrument Serial Number:

W20678

Date:

10/13/2008

Time:

9:42:28PM

Test Serial Number:

200810132142

Test Name:

LAGR HYDROCORR

Test Revision:

3

Test Description:

LAGR 4", PTFE, 0.2micron,
6pk

HydroCorr Test

Filter Catalog Number:

LAGR04TP6

Test Fluid:

Water

Filter Description:

Optiseal Cartridge Filter With
Aervent Membrane , 6/pack

Filter Size:

4"

Pore Size:

0.2

Housing Type:

optiseal

Operator Inputs

Operator:

raj

Product Batch no.

TOPTECAN

Filter Lot no.

C6MN97312 ✓

User Defined Field #3 - Optional

0189

User Defined Field #4 - Optional

PRE INT. H.N. CSP/086

Notes:

NITROGEN LINE FILTER

Test Specifications

HydroCorr 38.0 psig

Test Pressure:

Test Time: 20.00 min

HydroCorr Specification: 0.20 ml/min

4.44 Cartridges

of Cartridges:

System Specification: 0.20 ml/min


Test Data

Gas Volume: 207 ml

Actual Test Pressure: 37.9 psig

Time	Pressure	Flowrate
min	psig	ml/min
5.29	37.8	0.12
6.08	37.7	0.15
6.92	37.7	0.12
7.75	37.7	0.11
8.58	37.7	0.10
9.42	37.7	0.09
10.25	37.6	0.08
11.08	37.6	0.07
11.92	37.6	0.07
12.75	37.6	0.07
13.58	37.6	0.06
14.42	37.6	0.06
15.25	37.5	0.06
16.08	37.5	0.06
16.92	37.5	0.06
17.75	37.5	0.06
18.58	37.5	0.06
19.42	37.5	0.05
20.00	37.5	0.05

Flowrate v. Time



Test Results

Flowrate: 0.05 ml/min

Constant Pressure 0.05 ml/min

Flowrate:

Immediate Retest

Test Conclusion

PASSED

Signature:

Raj 13/10/08

Signature:

hmm 13/10/08

hmm 13/10/08

INTEGRITEST EXACTA

Instrument Serial Number:

W20678

Date:

10/14/2008

Time:

5:30:05AM

Test Serial Number:

200810140530

Test Name:

Lyo Vent filter

Test Revision:

2

Test Description:

MDI

HydroCorr Test

Filter Catalog Number:

CPTF5301A0SS101

Test Fluid:

WFI

Filter Description:

Filter Size:

5"

Pore Size:

0.2

Housing Type:

5"

Operator Inputs

Operator:

raj

Product Batch no.

TOPOTECAN

Filter Lot no.

CPTF 5301A0SS101 ✓

User Defined Field #3 - Optional

CT0198-04

User Defined Field #4 - Optional

PRE INT. H.N. M370750003

Notes:

AIR LINE FILTER

Test Specifications

HydroCorr

38.0 psig

Test Pressure:

Test Time:

20.00 min

HydroCorr Specification: 0.70 ml/min

of Cartridges: 1

System Specification: 0.70 ml/min


Test Data

Gas Volume: 555 ml

Actual Test Pressure: 37.9 psig

Time	Pressure	Flowrate
min	psig	ml/min
5.33	37.8	-0.01
6.12	37.8	0.57
6.96	37.8	0.44
7.79	37.7	0.36
8.63	37.7	0.30
9.46	37.7	0.27
10.29	37.7	0.24
11.12	37.7	0.22
11.96	37.6	0.20
12.79	37.6	0.19
13.63	37.6	0.18
14.46	37.6	0.17
15.29	37.6	0.16
16.13	37.6	0.16
16.96	37.5	0.15
17.79	37.5	0.15
18.62	37.5	0.14
19.46	37.5	0.14
20.00	37.5	0.14

Flowrate v. Time



Test Results

Flowrate: 0.14 ml/min

Constant Pressure Flowrate: 0.14 ml/min

Test Conclusion

PASSED

Signature: *Rat* 14/10/08

Signature: *W* 14/10/08

INTEGRITEST EXACTA

Instrument Serial Number:

W20678

Date:

10/14/2008

Time:

5:59:49AM

Test Serial Number:

200810140559

Test Name:

Lyo Vent filter

Test Revision:

2

Test Description:

MDI

HydroCorr Test

Filter Catalog Number:

CPTF5301A0SS101

Test Fluid:

WFI

Filter Description:

Filter Size:

5"

Pore Size:

0.2

Housing Type:

5"

Operator Inputs

Operator:

rai

Product Batch no.

TOPOPTECAN

Filter Lot no.

CPTF5301A0SS101 ✓

User Defined Field #3 - Optional

CT 0198-03

User Defined Field #4 - Optional

PRE INT. H.N.M370750015

Notes:

AIR LINE FILTER

Test Specifications

HydroCorr 38.0 psig

Test Pressure:

Test Time: 20.00 min

HydroCorr Specification: 0.70 ml/min

of Cartridges:

1


System Specification:

0.70 ml/min

Test Data

Gas Volume: 542 ml

Actual Test Pressure: 37.9 psig

**Test Results**

Flowrate: 0.13 ml/min

Constant Pressure 0.13 ml/min

Flowrate:

Test Conclusion**PASSED**

Signature:

Red 1/11/08

Signature:

1/11/08

	Batch Filling Record (Part-A)		Page 36 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size		Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

8.4 SIP OF LYOPHILIZER

8.4.1 LINE CLEARANCE FOR SIP

S.No.	Checks	Room No.	SOP No.	Date	Previous Product		Ckd. by ^{##}	Verified by ^{###}
					Product Name	Batch no.		
1.	Line Clearance of lyophilizer.	CSLG/034	SOP/QAD/053	14/10/08	Topotecan Placebo	—		

8.4.1.1 perform the SIP as per SOP No.SOP/CSP/032 and record the details in the table below.

SIP Cycle Started at: 17:57, Completed at: 22:10

Date	Details	Eq. Id. No.	SIP Details			Done by ^{##}	Ckd. by ^{##}
			Sterilization Temperature	Ster.Hold Started at	Ster. Hold Completed at		
14/10/08	Lyophilise	CSP/040	121°C	19:17	19:47	Raj—	

(Attach the print out of the process on the reverse side of page.)

8.4.2 STERILIZED LYOPHILIZER HOLD TIME (NMT 48HOURS AFTER STERILIZATION)

Procedure	Room No.	SOP No.	Date	Time		Done by ^{##}	Ckd. by ^{##}
				From	To		
Sterilized lyophilizer hold time	CSLG/034	---	14/10/08	22:10	—	Raj—	
			15/10/08	—	17:30	Shamshu	

^{##}Production

^{###}IPQA

SIP Process Report

Initiator:	SUP	Start Time:	14/10/2008 17:57:16
Batch Name:	SIP20081014175716	End Time:	14/10/2008 22:10:23

PARAMETERS:

Sterilization temperature:	121	"C
Sterilization pressure:	2.2	BarA
Sterilization time:	30	Min.
Number of purges:	2	
Pre-evacuation pressure:	0.35	BarA
Steam discharge pressure:	1.02	BarA
Filter drying time:	15	Min.
System drying time:	45	Min.
Condenser cooling time:	10	Min.
Shelf cooling time:	15	Min.
Time to reach setpoint:	120	Min.
Time below setpoint once reached:	30	Min.
Total steaming time after setpoint reached:	240	Min.

Batch Name: SIP20081014175716

Page 1 of 4

Handwritten signatures and dates:
14/10/08
14/10/08
14/10/08

SIP Process Report

Initiator:	SUP	Start Time:	14/10/2008 17:57:16
Batch Name:	SIP20081014175716	End Time:	14/10/2008 22:10:23

Date	Time	Step No.	Step Type	Chamber Pressure (BarA)	Condenser Pressure (BarA)	Cham. Drain Temp. (°C)	Cond. Drain Temp. (°C)	Filter Drain Temp. 1 (°C)	Filter Drain Temp. 2 (°C)
14/10/2008	17:57:22	1	Retract Ram	1.03	0.92	24.6	24.9	26.6	26.3
14/10/2008	17:58:29	2	Lock Doors	0.99	1.00	24.6	24.9	26.5	26.1
14/10/2008	18:02:16			0.36	0.36	22.8	28.1	26.5	26.3
14/10/2008	18:02:26	3	Purge	0.34	0.34	22.7	28.0	26.6	26.2
14/10/2008	18:07:16			0.49	0.50	78.1	77.9	100.4	98.9
14/10/2008	18:12:16			0.59	0.61	83.6	83.6	101.6	100.6
14/10/2008	18:17:16			0.74	0.74	89.4	88.4	103.4	102.3
14/10/2008	18:22:16			0.88	0.89	93.5	93.3	105.3	104.4
14/10/2008	18:27:16			1.01	1.03	98.9	97.0	107.6	106.6
14/10/2008	18:32:16			0.41	0.44	71.6	75.4	76.6	76.6
14/10/2008	18:37:16			0.49	0.52	78.3	78.1	100.6	99.3
14/10/2008	18:42:16			0.69	0.72	88.1	87.8	102.8	101.7
14/10/2008	18:47:16			0.89	0.92	95.1	94.8	105.6	104.4
14/10/2008	18:50:19	4	Steaming	1.02	1.05	99.0	99.0	107.4	106.8
14/10/2008	18:52:16			1.09	1.11	97.6	98.6	119.4	119.3
14/10/2008	18:57:16			1.30	1.32	99.7	100.1	120.7	120.4
14/10/2008	19:02:16			1.51	1.54	110.2	110.6	122.3	122.0
14/10/2008	19:07:16			1.73	1.77	114.8	114.9	123.7	123.3
14/10/2008	19:12:16			1.95	1.99	118.6	118.8	125.4	125.1
14/10/2008	19:17:16			2.15	2.19	121.7	121.9	127.1	126.7
14/10/2008	19:22:16			2.15	2.18	121.7	121.8	126.2	125.6
14/10/2008	19:27:16			2.13	2.17	121.5	121.5	123.4	123.3
14/10/2008	19:32:16			2.14	2.18	121.6	121.8	122.3	122.8
14/10/2008	19:37:16			2.17	2.20	121.8	122.1	126.1	125.8
14/10/2008	19:42:16			2.12	2.16	121.5	121.7	122.0	122.0
14/10/2008	19:47:00	5	Steam Dischar	2.15	2.18	121.6	121.8	125.6	125.4
14/10/2008	19:47:16			2.09	2.13	121.4	121.6	124.0	124.3
14/10/2008	19:52:16			1.34	1.38	109.2	109.3	108.9	110.9
14/10/2008	19:57:16			1.06	1.09	102.4	102.9	100.6	101.1
14/10/2008	20:00:41	6	Filter Drying	1.02	1.05	99.8	101.1	98.8	98.5
14/10/2008	20:02:16			1.05	1.08	91.8	98.4	71.0	84.6
14/10/2008	20:07:16			1.07	1.10	77.9	91.3	57.2	76.6
14/10/2008	20:12:16			1.07	1.10	70.3	82.8	56.6	70.7
14/10/2008	20:15:02	7	System Drying	1.01	1.04	67.2	78.2	56.7	67.3
14/10/2008	20:17:16			0.69	0.71	89.2	96.5	88.2	86.3
14/10/2008	20:22:16			0.36	0.39	84.1	95.1	70.2	71.8
14/10/2008	20:27:16			0.19	0.21	77.9	91.6	53.1	55.3
14/10/2008	20:32:16			0.13	0.15	72.6	87.2	40.9	44.4
14/10/2008	20:37:16			0.12	0.14	67.1	80.6	35.9	37.7
14/10/2008	20:42:16			0.12	0.13	61.9	73.8	34.7	35.4
14/10/2008	20:47:16			0.12	0.13	57.8	67.8	34.1	34.4
14/10/2008	20:52:16			0.12	0.13	54.3	62.3	33.7	33.6
14/10/2008	20:57:16			0.12	0.13	51.2	57.8	33.5	33.1
14/10/2008	21:02:16			0.12	0.13	48.7	54.0	33.1	33.0
14/10/2008	21:07:08	8	Cond. Cooling	0.11	0.12	46.4	50.8	33.1	32.8
14/10/2008	21:07:16			0.12	0.13	46.3	50.7	32.9	32.8
14/10/2008	21:12:16			0.07	0.08	44.2	46.7	32.8	32.6
14/10/2008	21:17:16			0.07	0.08	42.3	43.5	32.7	32.4
14/10/2008	21:22:16			0.07	0.08	40.6	41.5	32.3	32.0
14/10/2008	21:25:03	9	Shelf Cooling	0.06	0.08	39.7	40.3	32.2	31.9
14/10/2008	21:27:16			0.07	0.08	39.0	39.7	32.0	31.8
14/10/2008	21:32:16			0.07	0.08	37.6	38.0	31.7	31.5

Batch Name: SIP20081014175716

Page 2 of 4

Handwritten:
 1000
 14/10/08
 Raj
 14/10/08
 14/10/08

SIP Process Report

Initiator:	SUP	Start Time:	14/10/2008 17:57:16
Batch Name:	SIP20081014175716	End Time:	14/10/2008 22:10:23

Date	Time	Step No.	Step Type	Chamber Pressure (BarA)	Condenser Pressure (BarA)	Cham. Drain Temp. (°C)	Cond. Drain Temp. (°C)	Filter Drain Temp. 1 (°C)	Filter Drain Temp. 2 (°C)
14/10/2008	21:37:16			0.07	0.08	36.5	37.0	31.6	31.2
14/10/2008	21:42:16			0.07	0.08	35.6	35.9	31.1	30.9
14/10/2008	21:47:16			0.07	0.08	34.6	35.2	30.8	30.7
14/10/2008	21:52:16			0.07	0.08	33.9	34.5	30.7	30.1
14/10/2008	21:57:16			0.07	0.08	33.2	33.9	30.4	30.0
14/10/2008	21:57:59	10	Unlock Doors	0.06	0.08	33.0	33.8	30.3	30.0
14/10/2008	21:58:03	11	Final Aeration	0.06	0.08	33.1	33.8	30.3	30.0
14/10/2008	22:02:16			0.49	0.49	32.5	33.3	30.2	29.8
14/10/2008	22:10:23			0.98	0.99	31.6	32.4	29.9	29.6

Handwritten:
 14/10/08
 14/10/08
 14/10/08

SIP Process Report

Initiator:	SUP	Start Time:	14/10/2008 17:57:16
Batch Name:	SIP20081014175716	End Time:	14/10/2008 22:10:23

Date	Time	Description
Alarm: Total: 3		
14/10/2008	19:14:33	PT101 chamber vacuum probe fail
14/10/2008	19:14:49	PT101 chamber vacuum probe fail
14/10/2008	21:48:52	Shelf temperature low deviation
Audit: Total: 10		
14/10/2008	18:30:32	SUP logged in as Application User by Sup.
14/10/2008	18:36:57	SUP logged in as Application User by Sup.
14/10/2008	18:42:23	SUP logged in as Application User by Sup.
14/10/2008	18:54:51	SUP logged in as Application User by Sup.
14/10/2008	19:15:57	SUP logged in as Application User by Sup.
14/10/2008	19:16:04	A_CHAMVACPROBEFAIL ALARM is acknowledged by Sup.
14/10/2008	20:47:43	SUP logged in as Application User by Sup.
14/10/2008	21:48:41	SUP logged in as Application User by Sup.
14/10/2008	21:49:02	A_SHELFTEMPPLD ALARM is acknowledged by Sup.
14/10/2008	22:10:16	SUP logged in as Application User by Sup.

Handwritten signatures and dates:
 14/10/08
 14/10/08
 14/10/08

Kas
 14/11/08
 for
 14/11/08

01-06ch	91.1	91.1	91.1	90.8	91.1
07-12ch	91.0	91.0	91.1	90.7	84.2
13-18ch	102.2	102.8	98.5	93.5	32.8
19-24ch	32.9	26	2.1	0.11	---

100.0
 21:00 60mm/h 20 40 60 80 CH 9 150-100 PC

01-06ch	100.1	99.9	100.0	100.0	100.1
07-12ch	100.1	100.1	100.1	99.9	99.9
13-18ch	114.9	100.8	101.3	98.2	64.7
19-24ch	81.4	102.1	15	1.06	1.07

100.0
 21:00 60mm/h 20 40 60 80 CH 9 150-100 PC

01-06ch	111.8	111.8	111.8	111.8	111.7
07-12ch	111.9	111.9	111.9	111.7	74.3
13-18ch	107.6	107.6	107.6	107.1	122.1
19-24ch	121.7	---	---	---	---

100.0
 19:00 60mm/h 20 40 60 80 CH 7 150-100 PC

MEMORANDUM

PART NO. S1407861-001

	Batch Filling Record (Part-A)		Page 37 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

8.5 LOADING OF LYOPHILIZER



Loading is done as per SOP No. SOP/CSP/032

Date	SOP No.	No. of Tray in Shelf	Quantity	Time		Done by ##	Ckd. by ##
				From	To		
15/10/08		1 X 6	42.40	18:10	22:20	Sharmila	 15/10/08

8.6 LYOPHILIZATION CYCLE

Start the lyophilization cycle as per SOP No. SOP/CSP/032

Recipe Name: -

Date	SOP No.	Cycle started at		Cycle end at		Done by ##	Ckd. by ##
		Date	Time	Date	Time		
16/10/08	SOP/CSP/032	16/10/08	00:41	18/10/08	18:32		 16/10/08

(Attach the print out of the process on the reverse side of page.)

Production

Form No.: SOP/QAD/038-02 Rev. 00

CONFIDENTIAL

QA Signature



	Batch Filling Record (Part-A)		Page 38 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

9.0 CLEANING OF ALUMINUM FLIP-OFF OVER SEALS:

9.1 LINE CLEARANCE:

- Ensure the cleanliness of equipment wash area.
- Ensure no remains of previous product are there.
- Ensure that proper status label is attached on machine.

Date: 17/10/08 14/10/08 17/10/08

Step No.	Checks	Eq. Id. No.	SOP No.	Previous product		Ckd. by##	Verified by###
				Product Name	Batch No.		
1.	Line clearance for washing of aluminum seals	CSP/129	SOP/QAD/053	Cytarabine inj BP/ 100mg/ml	8LA284	<u>CS</u>	<u>17/10/08</u>

9.2 CLEANING:

Cleaning of aluminum flip-off over seals done as per SOP No. SOP/CSF/053.

A.R. No. Of WFI*: 02WF08500757

Date	Quantity	Cleaning Time		Quantity		Done by##	Ckd. By##
		Started at	Completed at	No. of good seals	No. of Rejected seals		
<u>17/10/08</u> <u>14/10/08</u>	6300	15:45	16:25	6295	05	<u>CS</u>	<u>CS</u>

* Return loop A. R. No to be mentioned

Production

IPQA

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Vacuum Control Parameters:

Vacuum high limit:	3	uBarA	Vacuum preAlarm:	500	uBarA
Vacuum low limit:	3	uBarA	Vacuum alarm:	700	uBarA

Step No	Step Type	Temperature ("C)	Vacuum (ubar)	Time	
				HH:MM	MM:SS
1	Loading	10.0	--	--	--
2	Freezing	-10.0	--	0 : 20	--
3	Freezing	-10.0	--	0 : 30	--
4	Freezing	-15.0	--	0 : 5	--
5	Freezing	-15.0	--	1 : 0	--
6	Freezing	-25.0	--	0 : 10	--
7	Freezing	-25.0	--	1 : 0	--
8	Freezing	-40.0	--	0 : 15	--
9	Freezing	-40.0	--	3 : 0	--
10	Evacuation	--	200	--	--
11	Drying	-28.0	200	0 : 10	--
12	Drying	-28.0	200	1 : 30	--
13	Drying	-20.0	200	0 : 8	--
14	Drying	-20.0	200	1 : 30	--
15	Drying	-12.0	150	0 : 8	--
16	Drying	-12.0	150	4 : 0	--
17	Drying	-6.0	150	0 : 6	--
18	Drying	-6.0	150	11 : 0	--
19	Drying	0.0	100	0 : 6	--
20	Drying	0.0	100	11 : 0	--
21	Drying	7.0	80	0 : 7	--
22	Drying	7.0	80	3 : 0	--
23	Drying	15.0	80	0 : 8	--
24	Drying	15.0	80	3 : 0	--
25	Drying	30.0	60	0 : 15	--
26	Drying	30.0	60	8 : 0	--
27	Drying	45.0	60	0 : 15	--
28	Drying	45.0	60	12 : 0	--
29	Moisture Test	--	300	--	30 : 60
30	Pre-Aeration	--	--	--	--
31	Stoppering	--	--	--	0 : 0
32	Aeration	--	--	--	--
33	Storage	25.0	--	0 : 0	--

The time setpoint at the step of moisture test denotes respectively the repeat time in unit of minutes, and test time in unit of second.

Handwritten signature
18/10/08
18/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
16/10/2008	0:43:36	1	Loading	9.6	10.0	10.6	36.9	36.9	1000	--
16/10/2008	0:44:31	2	Freezing	10.4	10.0	9.1	36.9	36.9	1000	--
16/10/2008	0:46:52			6.6	7.6	8.9	36.9	37.2	1000	--
16/10/2008	0:51:52			1.5	2.6	4.3	36.9	37.2	1000	--
16/10/2008	0:56:52			-1.2	-2.4	-1.8	36.9	37.2	1000	--
16/10/2008	1:01:52			-6.8	-7.4	-6.3	36.9	36.7	1000	--
16/10/2008	1:04:26	3	Freezing	-10.7	-10.0	-7.6	36.6	36.7	1000	--
16/10/2008	1:06:52			-9.6	-9.9	-9.6	36.9	36.7	1000	--
16/10/2008	1:11:52			-9.4	-9.9	-10.7	36.9	36.7	1000	--
16/10/2008	1:16:52			-10.7	-9.9	-9.9	36.9	36.7	1000	--
16/10/2008	1:21:52			-9.6	-9.9	-11.1	36.9	36.7	1000	--
16/10/2008	1:26:52			-10.4	-9.9	-10.4	36.4	36.7	1000	--
16/10/2008	1:31:52			-9.7	-9.9	-10.8	36.4	36.7	1000	--
16/10/2008	1:34:27	4	Freezing	-9.7	-10.0	-11.3	36.3	36.3	1000	--
16/10/2008	1:36:52			-12.6	-12.4	-11.3	36.4	36.7	1000	--
16/10/2008	1:39:25	5	Freezing	-14.7	-14.9	-13.4	36.2	36.3	1000	--
16/10/2008	1:41:52			-14.7	-14.9	-15.4	36.4	36.2	1000	--
16/10/2008	1:46:52			-15.2	-14.9	-15.3	36.4	36.2	1000	--
16/10/2008	1:51:52			-15.6	-14.9	-15.3	36.4	36.2	1000	--
16/10/2008	1:56:52			-14.4	-14.9	-16.3	36.4	36.2	1000	--
16/10/2008	2:01:52			-15.1	-14.9	-15.5	36.4	36.2	1000	--
16/10/2008	2:06:52			-15.7	-14.9	-15.5	35.9	36.2	1000	--
16/10/2008	2:11:52			-14.7	-14.9	-16.2	35.9	36.2	1000	--
16/10/2008	2:16:52			-14.5	-14.9	-16.2	35.9	36.2	1000	--
16/10/2008	2:21:52			-15.8	-14.9	-15.2	35.9	35.7	1000	--
16/10/2008	2:26:52			-15.0	-14.9	-16.2	35.9	35.7	1000	--
16/10/2008	2:31:52			-14.6	-14.9	-16.1	35.9	35.7	1000	--
16/10/2008	2:36:52			-15.7	-14.9	-15.3	35.9	35.7	1000	--
16/10/2008	2:39:26	6	Freezing	-14.7	-15.0	-15.8	35.6	35.5	1000	--
16/10/2008	2:41:52			-17.2	-17.4	-17.0	35.9	35.7	1000	--
16/10/2008	2:46:52			-22.0	-22.4	-20.8	35.9	35.7	1000	--
16/10/2008	2:49:25	7	Freezing	-24.5	-24.9	-23.6	35.5	35.5	1000	--
16/10/2008	2:51:52			-25.0	-24.9	-24.5	35.4	35.7	1000	--
16/10/2008	2:56:52			-24.3	-24.9	-25.6	35.4	35.7	1000	--
16/10/2008	3:01:52			-24.9	-24.9	-25.8	35.4	35.2	1000	--
16/10/2008	3:06:52			-25.6	-24.9	-25.5	35.4	35.2	1000	--
16/10/2008	3:11:52			-24.6	-24.9	-25.1	35.4	35.2	1000	--
16/10/2008	3:16:52			-24.3	-24.9	-26.0	35.4	35.2	1000	--
16/10/2008	3:21:52			-25.1	-24.9	-25.9	35.4	35.2	1000	--
16/10/2008	3:26:52			-25.6	-24.9	-25.2	35.4	35.2	1000	--
16/10/2008	3:31:52			-25.1	-24.9	-25.1	35.4	35.2	1000	--
16/10/2008	3:36:52			-24.7	-24.9	-25.7	34.9	35.2	1000	--
16/10/2008	3:41:52			-24.4	-24.9	-25.8	34.9	35.2	1000	--
16/10/2008	3:46:52			-25.5	-24.9	-26.0	34.9	35.2	1000	--
16/10/2008	3:49:28	8	Freezing	-24.6	-25.0	-25.9	34.8	34.8	1000	--
16/10/2008	3:51:52			-26.8	-27.4	-26.4	34.9	34.7	1000	--
16/10/2008	3:56:52			-31.7	-32.4	-31.0	34.9	34.7	1000	--
16/10/2008	4:01:52			-35.1	-37.4	-34.6	34.9	34.7	1000	--
16/10/2008	4:06:52			-37.6	-40.0	-37.3	34.9	34.7	1000	--
16/10/2008	4:11:52			-38.1	-40.0	-38.5	34.9	34.7	1000	--
16/10/2008	4:13:55	9	Freezing	-39.0	-40.0	-39.1	34.6	34.5	1000	--
16/10/2008	4:16:52			-40.0	-40.0	-40.3	34.9	34.7	1000	--
16/10/2008	4:21:52			-39.6	-40.0	-40.3	34.9	34.7	1000	--
16/10/2008	4:26:52			-39.6	-40.0	-40.8	34.9	34.7	1000	--

Batch Name: AFD8PK201

Page 2 of 19

Aspin
18/10/08
by me
12/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
16/10/2008	4:31:52			-40.1	-40.0	-40.8	34.4	34.7	1000	--
16/10/2008	4:36:52			-39.9	-40.0	-40.3	34.4	34.7	1000	--
16/10/2008	4:41:52			-39.2	-40.0	-40.8	34.4	34.2	1000	--
16/10/2008	4:46:52			-40.1	-40.0	-40.1	34.4	34.2	1000	--
16/10/2008	4:51:52			-39.2	-40.0	-40.3	34.4	34.2	1000	--
16/10/2008	4:56:52			-40.3	-40.0	-40.9	34.4	34.2	1000	--
16/10/2008	5:01:52			-39.6	-40.0	-40.3	34.4	34.2	1000	--
16/10/2008	5:06:52			-39.5	-40.0	-40.9	34.4	34.2	1000	--
16/10/2008	5:11:52			-39.9	-40.0	-40.2	34.4	34.2	1000	--
16/10/2008	5:16:52			-39.2	-40.0	-41.0	34.4	34.2	1000	--
16/10/2008	5:21:52			-40.2	-40.0	-40.8	34.4	34.2	1000	--
16/10/2008	5:26:52			-39.4	-40.0	-40.2	33.9	34.2	1000	--
16/10/2008	5:31:52			-40.2	-40.0	-40.7	33.9	34.2	1000	--
16/10/2008	5:36:52			-39.6	-40.0	-40.3	33.9	34.2	1000	--
16/10/2008	5:41:52			-39.5	-40.0	-41.1	33.9	34.2	1000	--
16/10/2008	5:46:52			-40.0	-40.0	-40.3	33.9	33.7	1000	--
16/10/2008	5:51:52			-39.2	-40.0	-40.8	33.9	33.7	1000	--
16/10/2008	5:56:52			-40.1	-40.0	-40.2	33.9	33.7	1000	--
16/10/2008	6:01:52			-39.2	-40.0	-40.8	33.9	33.7	1000	--
16/10/2008	6:06:52			-40.3	-40.0	-40.8	33.9	33.7	1000	--
16/10/2008	6:11:52			-39.2	-40.0	-40.7	33.9	33.7	1000	--
16/10/2008	6:16:52			-40.1	-40.0	-40.3	33.9	33.7	1000	--
16/10/2008	6:21:52			-39.4	-40.0	-40.3	33.4	33.7	1000	--
16/10/2008	6:26:52			-40.3	-40.0	-40.7	33.4	33.7	1000	--
16/10/2008	6:31:52			-39.6	-40.0	-40.7	33.4	33.7	1000	--
16/10/2008	6:36:52			-39.7	-40.0	-40.7	33.4	33.2	1000	--
16/10/2008	6:41:52			-40.0	-40.0	-40.7	33.4	33.2	1000	--
16/10/2008	6:46:52			-39.2	-40.0	-40.7	33.4	33.2	1000	--
16/10/2008	6:51:52			-39.8	-40.0	-40.2	33.4	33.2	1000	--
16/10/2008	6:56:52			-39.4	-40.0	-40.9	33.4	33.2	1000	--
16/10/2008	7:01:52			-40.1	-40.0	-40.2	33.4	33.2	1000	--
16/10/2008	7:06:52			-39.1	-40.0	-40.5	33.4	33.2	1000	--
16/10/2008	7:11:52			-40.1	-40.0	-40.3	32.9	33.2	1000	--
16/10/2008	7:13:52	10	Evacuation	-39.4	-40.0	-40.2	33.0	32.9	1000	200
16/10/2008	7:16:52			-39.4	-40.0	-40.2	30.6	14.6	1000	200
16/10/2008	7:21:52			-39.5	-40.0	-40.2	19.3	-46.7	1000	200
16/10/2008	7:26:52			-39.3	-40.0	-40.2	-24.3	-50.9	1000	200
16/10/2008	7:31:52			-39.4	-40.0	-40.2	-33.9	-47.5	1000	200
16/10/2008	7:36:52			-39.4	-40.0	-40.7	-43.6	-56.8	1000	200
16/10/2008	7:41:52			-39.6	-40.0	-40.7	-56.1	-60.9	1000	200
16/10/2008	7:46:52			-39.7	-40.0	-40.7	-55.5	-68.6	1000	200
16/10/2008	7:51:52			-39.4	-40.0	-40.7	-60.2	-72.7	389	200
16/10/2008	7:53:23	11	Drying	-39.7	-39.5	-40.8	-60.8	-74.5	205	200
16/10/2008	7:56:52			-36.2	-35.5	-38.1	-64.4	-75.2	213	200
16/10/2008	8:01:52			-30.1	-29.7	-33.5	-68.2	-76.6	316	200
16/10/2008	8:06:52			-28.0	-28.0	-29.8	-71.5	-77.9	437	200
16/10/2008	8:11:52			-27.8	-28.0	-29.3	-70.1	-79.6	491	200
16/10/2008	8:16:52			-27.9	-28.0	-29.3	-45.8	-63.9	322	200
16/10/2008	8:19:35	12	Drying	-28.5	-28.0	-30.0	-43.2	-44.8	202	200
16/10/2008	8:21:52			-27.0	-28.0	-29.1	-55.6	-55.6	203	200
16/10/2008	8:26:52			-27.9	-28.0	-29.0	-63.8	-66.2	208	200
16/10/2008	8:31:52			-27.8	-28.0	-29.0	-72.4	-74.5	208	200
16/10/2008	8:36:52			-29.5	-28.0	-30.7	-73.1	-58.9	338	200
16/10/2008	8:41:52			-27.7	-28.0	-28.9	-74.1	-71.9	212	200

Batch Name: AFD8PK201

Page 3 of 19

Arvind
18/10/08
12/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
16/10/2008	8:46:52			-29.8	-28.0	-30.6	-71.5	-54.7	383	200
16/10/2008	8:51:52			-27.6	-28.0	-28.3	-73.4	-70.5	213	200
16/10/2008	8:56:52			-27.7	-28.0	-28.9	-73.8	-77.5	213	200
16/10/2008	9:01:52			-27.4	-28.0	-28.9	-72.6	-79.3	215	200
16/10/2008	9:06:52			-27.7	-28.0	-28.9	-67.8	-79.8	216	200
16/10/2008	9:11:52			-27.7	-28.0	-28.9	-67.9	-79.8	217	200
16/10/2008	9:16:52			-29.5	-28.0	-30.4	-67.6	-59.1	434	200
16/10/2008	9:21:52			-27.4	-28.0	-28.2	-64.6	-70.4	202	200
16/10/2008	9:26:52			-27.6	-28.0	-28.7	-63.3	-76.4	202	200
16/10/2008	9:31:52			-27.6	-28.0	-28.7	-64.4	-79.2	201	200
16/10/2008	9:36:52			-27.6	-28.0	-28.7	-64.9	-79.2	198	200
16/10/2008	9:41:52			-27.8	-28.0	-28.7	-64.3	-79.2	201	200
16/10/2008	9:46:52			-27.9	-28.0	-28.7	-65.4	-79.2	197	200
16/10/2008	9:51:52			-27.6	-28.0	-28.7	-65.9	-80.2	205	200
16/10/2008	9:56:52			-27.7	-28.0	-28.7	-65.8	-80.2	270	200
16/10/2008	10:01:52			-27.5	-28.0	-28.7	-66.4	-80.2	203	200
16/10/2008	10:06:05	13	Drying	-27.3	-28.0	-28.7	-66.1	-80.1	202	200
16/10/2008	10:06:52			-27.4	-27.2	-28.7	-65.9	-80.2	193	200
16/10/2008	10:11:52			-22.0	-22.2	-25.0	-68.1	-80.2	200	200
16/10/2008	10:14:05	14	Drying	-20.9	-20.0	-23.3	-69.0	-80.6	200	200
16/10/2008	10:16:52			-19.9	-20.0	-20.9	-68.7	-80.2	202	200
16/10/2008	10:21:52			-19.9	-20.0	-20.9	-69.2	-80.2	197	200
16/10/2008	10:26:52			-19.6	-20.0	-20.9	-72.5	-80.2	196	200
16/10/2008	10:31:52			-19.6	-20.0	-20.9	-76.2	-80.7	204	200
16/10/2008	10:36:52			-19.5	-20.0	-20.9	-73.3	-81.2	193	200
16/10/2008	10:41:52			-19.7	-20.0	-20.9	-71.5	-81.2	204	200
16/10/2008	10:46:52			-19.6	-20.0	-20.9	-69.2	-80.7	198	200
16/10/2008	10:51:52			-19.7	-20.0	-20.4	-69.8	-80.7	192	200
16/10/2008	10:56:52			-19.6	-20.0	-20.4	-68.8	-80.2	202	200
16/10/2008	11:01:52			-19.7	-20.0	-20.4	-65.4	-80.2	193	200
16/10/2008	11:06:52			-19.6	-20.0	-20.4	-67.3	-80.2	204	200
16/10/2008	11:11:52			-19.6	-20.0	-20.4	-65.8	-80.2	193	200
16/10/2008	11:16:52			-19.9	-20.0	-20.4	-67.3	-80.2	204	200
16/10/2008	11:21:52			-19.6	-20.0	-20.4	-66.8	-80.2	197	200
16/10/2008	11:26:52			-19.7	-20.0	-20.4	-64.4	-80.2	192	200
16/10/2008	11:31:52			-19.5	-20.0	-20.4	-65.7	-80.2	205	200
16/10/2008	11:36:52			-19.6	-20.0	-20.4	-66.2	-80.2	199	200
16/10/2008	11:41:52			-19.5	-20.0	-20.4	-66.0	-80.2	201	200
16/10/2008	11:44:11	15	Drying	-19.5	-20.0	-20.6	-67.3	-80.0	201	150
16/10/2008	11:46:52			-17.2	-17.3	-19.7	-67.6	-80.2	149	150
16/10/2008	11:51:52			-13.1	-12.5	-14.9	-66.2	-79.7	149	150
16/10/2008	11:52:15	16	Drying	-12.9	-12.0	-14.7	-66.7	-79.6	152	150
16/10/2008	11:56:52			-11.7	-12.0	-13.2	-67.3	-79.7	156	150
16/10/2008	12:01:52			-11.7	-12.0	-12.7	-66.0	-79.7	156	150
16/10/2008	12:06:52			-11.5	-12.0	-12.7	-68.0	-79.7	156	150
16/10/2008	12:11:52			-11.5	-12.0	-12.7	-65.5	-79.7	155	150
16/10/2008	12:16:52			-11.5	-12.0	-12.7	-67.4	-79.7	154	150
16/10/2008	12:21:52			-11.8	-12.0	-12.7	-67.9	-79.7	155	150
16/10/2008	12:26:52			-11.7	-12.0	-12.7	-65.3	-79.7	154	150
16/10/2008	12:31:52			-11.6	-12.0	-12.7	-66.4	-79.7	152	150
16/10/2008	12:36:52			-11.6	-12.0	-12.7	-66.9	-79.7	151	150
16/10/2008	12:41:52			-11.6	-12.0	-12.7	-66.9	-79.7	149	150
16/10/2008	12:46:52			-11.6	-12.0	-12.7	-66.1	-79.7	146	150
16/10/2008	12:51:52			-11.4	-12.0	-12.7	-66.6	-79.7	145	150

Batch Name: AFD8PK201

Page 4 of 19

Apurva
18/10/08
hike

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
16/10/2008	12:56:52			-11.9	-12.0	-12.7	-65.6	-79.7	145	150
16/10/2008	13:01:52			-11.6	-12.0	-12.7	-67.2	-79.7	146	150
16/10/2008	13:06:52			-11.4	-12.0	-12.7	-74.2	-80.2	145	150
16/10/2008	13:11:52			-11.6	-12.0	-12.7	-75.2	-80.2	146	150
16/10/2008	13:16:52			-11.6	-12.0	-12.7	-76.9	-80.7	151	150
16/10/2008	13:21:52			-11.4	-12.0	-12.7	-72.5	-80.7	155	150
16/10/2008	13:26:52			-11.5	-12.0	-12.7	-69.8	-80.7	155	150
16/10/2008	13:31:52			-11.7	-12.0	-12.7	-69.4	-80.7	157	150
16/10/2008	13:36:52			-11.5	-12.0	-12.7	-67.9	-80.2	157	150
16/10/2008	13:41:52			-11.5	-12.0	-12.7	-67.8	-80.2	158	150
16/10/2008	13:46:52			-11.7	-12.0	-12.7	-66.7	-79.7	154	150
16/10/2008	13:51:52			-11.5	-12.0	-12.7	-68.5	-79.7	147	150
16/10/2008	13:56:52			-11.6	-12.0	-12.7	-67.8	-79.7	144	150
16/10/2008	14:01:52			-11.5	-12.0	-12.7	-69.0	-79.7	157	150
16/10/2008	14:06:52			-11.8	-12.0	-12.7	-69.0	-79.7	157	150
16/10/2008	14:11:52			-11.6	-12.0	-12.7	-70.0	-80.2	157	150
16/10/2008	14:16:52			-11.6	-12.0	-12.7	-69.0	-80.2	156	150
16/10/2008	14:21:52			-11.7	-12.0	-12.7	-70.7	-80.2	154	150
16/10/2008	14:26:52			-11.5	-12.0	-12.7	-71.6	-80.2	154	150
16/10/2008	14:31:52			-11.7	-12.0	-12.7	-67.8	-80.2	155	150
16/10/2008	14:36:52			-11.7	-12.0	-12.7	-68.3	-80.2	155	150
16/10/2008	14:41:52			-11.6	-12.0	-12.7	-71.9	-80.2	154	150
16/10/2008	14:46:52			-11.5	-12.0	-12.7	-71.9	-80.2	153	150
16/10/2008	14:51:52			-11.6	-12.0	-12.7	-72.4	-80.2	153	150
16/10/2008	14:56:52			-11.6	-12.0	-12.7	-72.4	-80.2	153	150
16/10/2008	15:01:52			-11.7	-12.0	-12.7	-75.0	-80.2	153	150
16/10/2008	15:06:52			-11.5	-12.0	-12.7	-70.4	-80.2	152	150
16/10/2008	15:11:52			-11.6	-12.0	-12.7	-69.6	-80.2	152	150
16/10/2008	15:16:52			-11.5	-12.0	-12.7	-70.6	-80.7	155	150
16/10/2008	15:21:52			-11.7	-12.0	-12.7	-66.0	-80.7	154	150
16/10/2008	15:26:52			-11.6	-12.0	-12.7	-67.3	-80.7	154	150
16/10/2008	15:31:52			-11.7	-12.0	-12.7	-67.8	-80.2	155	150
16/10/2008	15:36:52			-11.4	-12.0	-12.7	-66.0	-80.2	154	150
16/10/2008	15:41:52			-11.7	-12.0	-12.7	-67.6	-80.2	155	150
16/10/2008	15:46:52			-11.6	-12.0	-12.7	-68.1	-80.2	156	150
16/10/2008	15:51:52			-11.6	-12.0	-12.7	-67.7	-80.2	155	150
16/10/2008	15:52:21	17	Drying	-11.7	-12.0	-12.5	-68.3	-80.3	152	150
16/10/2008	15:56:52			-8.2	-7.4	-10.0	-69.5	-80.2	156	150
16/10/2008	15:58:28	18	Drying	-5.6	-6.0	-8.2	-69.2	-80.5	153	150
16/10/2008	16:01:52			-5.2	-6.1	-6.9	-68.4	-80.2	150	150
16/10/2008	16:06:52			-5.8	-6.1	-6.9	-71.2	-80.2	145	150
16/10/2008	16:11:52			-5.8	-6.1	-6.9	-77.5	-80.7	157	150
16/10/2008	16:16:52			-5.6	-6.1	-6.4	-78.0	-80.7	148	150
16/10/2008	16:21:52			-5.6	-6.1	-6.4	-78.6	-81.2	148	150
16/10/2008	16:26:52			-5.8	-6.1	-6.4	-77.1	-81.2	154	150
16/10/2008	16:31:52			-5.6	-6.1	-6.4	-77.6	-81.2	147	150
16/10/2008	16:36:52			-5.6	-6.1	-6.4	-77.6	-81.2	154	150
16/10/2008	16:41:52			-5.7	-6.1	-6.4	-78.3	-81.2	149	150
16/10/2008	16:46:52			-5.6	-6.1	-6.4	-77.0	-81.2	150	150
16/10/2008	16:51:52			-5.4	-6.1	-6.4	-75.2	-81.2	154	150
16/10/2008	16:56:52			-5.5	-6.1	-6.4	-71.2	-81.2	145	150
16/10/2008	17:01:52			-5.6	-6.1	-6.4	-70.1	-81.2	157	150
16/10/2008	17:06:52			-5.6	-6.1	-6.4	-70.6	-81.2	144	150
16/10/2008	17:11:52			-5.5	-6.1	-6.4	-70.6	-80.1	157	150

Batch Name: AFD8PK201

Page 5 of 19

Arrived
18/10/08
18/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
16/10/2008	17:16:52			-5.8	-6.1	-6.4	-67.3	-79.6	150	150
16/10/2008	17:21:52			-5.7	-6.1	-6.4	-66.4	-79.6	152	150
16/10/2008	17:26:52			-5.7	-6.1	-6.4	-67.4	-79.6	153	150
16/10/2008	17:31:52			-5.6	-6.1	-6.4	-65.2	-79.6	145	150
16/10/2008	17:36:52			-5.4	-6.1	-6.4	-65.9	-79.6	157	150
16/10/2008	17:41:52			-5.5	-6.1	-6.4	-66.9	-79.6	146	150
16/10/2008	17:46:52			-5.6	-6.1	-6.4	-66.9	-79.6	157	150
16/10/2008	17:51:52			-5.5	-6.1	-6.4	-68.1	-80.1	149	150
16/10/2008	17:56:52			-5.4	-6.1	-6.4	-68.8	-80.1	148	150
16/10/2008	18:01:52			-5.4	-6.1	-6.4	-65.2	-80.1	154	150
16/10/2008	18:06:52			-5.7	-6.1	-6.4	-66.9	-80.1	146	150
16/10/2008	18:11:52			-5.7	-6.1	-6.4	-67.5	-80.1	156	150
16/10/2008	18:16:52			-5.5	-6.1	-6.4	-66.2	-80.1	147	150
16/10/2008	18:21:52			-5.6	-6.1	-6.4	-65.5	-80.1	156	150
16/10/2008	18:26:52			-5.6	-6.1	-6.4	-66.6	-80.1	152	150
16/10/2008	18:31:52			-5.6	-6.1	-6.4	-65.8	-80.1	146	150
16/10/2008	18:36:52			-5.6	-6.1	-6.4	-66.3	-72.8	213	150
16/10/2008	18:41:52			-5.7	-6.1	-6.4	-65.3	-74.3	146	150
16/10/2008	18:46:52			-5.6	-6.1	-6.4	-66.5	-74.5	253	150
16/10/2008	18:51:52			-5.7	-6.1	-6.4	-67.5	-74.0	150	150
16/10/2008	18:56:52			-5.6	-6.1	-6.4	-68.0	-75.4	150	150
16/10/2008	19:01:52			-5.4	-6.1	-6.4	-68.7	-75.4	153	150
16/10/2008	19:06:52			-5.6	-6.1	-6.4	-66.0	-74.7	145	150
16/10/2008	19:11:52			-5.7	-6.1	-6.4	-68.3	-75.4	157	150
16/10/2008	19:16:52			-5.5	-6.1	-6.4	-69.0	-75.4	150	150
16/10/2008	19:21:52			-5.8	-6.1	-6.4	-69.5	-75.9	146	150
16/10/2008	19:26:52			-5.6	-6.1	-6.4	-69.5	-75.4	156	150
16/10/2008	19:31:52			-5.6	-6.1	-6.4	-69.0	-76.1	150	150
16/10/2008	19:36:52			-5.5	-6.1	-6.4	-69.0	-74.7	145	150
16/10/2008	19:41:52			-5.6	-6.1	-6.4	-65.7	-74.2	158	150
16/10/2008	19:46:52			-5.6	-6.1	-6.4	-67.2	-75.7	152	150
16/10/2008	19:51:52			-5.7	-6.1	-6.4	-69.9	-75.1	149	150
16/10/2008	19:56:52			-5.7	-6.1	-6.4	-67.9	-75.2	145	150
16/10/2008	20:01:52			-5.7	-6.1	-6.4	-68.3	-75.9	152	150
16/10/2008	20:06:52			-5.6	-6.1	-6.4	-69.4	-75.8	154	150
16/10/2008	20:11:52			-5.6	-6.1	-6.4	-70.5	-76.3	157	150
16/10/2008	20:16:52			-5.7	-6.1	-6.4	-68.3	-75.7	157	150
16/10/2008	20:21:52			-5.5	-6.1	-6.4	-68.7	-76.3	154	150
16/10/2008	20:26:52			-5.5	-6.1	-6.4	-69.2	-75.7	157	150
16/10/2008	20:31:52			-5.7	-6.1	-6.4	-69.8	-75.7	157	150
16/10/2008	20:36:52			-5.5	-6.1	-6.4	-66.8	-76.7	148	150
16/10/2008	20:41:52			-5.8	-6.1	-6.4	-67.3	-75.9	146	150
16/10/2008	20:46:52			-5.4	-6.1	-6.4	-65.4	-75.4	155	150
16/10/2008	20:51:52			-5.8	-6.1	-6.4	-67.3	-75.5	146	150
16/10/2008	20:56:52			-5.5	-6.1	-6.4	-69.1	-76.5	154	150
16/10/2008	21:01:52			-5.7	-6.1	-6.4	-68.5	-75.5	144	150
16/10/2008	21:06:52			-5.5	-6.1	-6.4	-69.6	-76.6	154	150
16/10/2008	21:11:52			-5.7	-6.1	-6.4	-65.9	-77.2	156	150
16/10/2008	21:16:52			-5.6	-6.1	-6.4	-67.0	-77.2	156	150
16/10/2008	21:21:52			-5.6	-6.1	-6.4	-67.0	-79.0	147	150
16/10/2008	21:26:52			-5.4	-6.1	-6.4	-66.9	-77.1	151	150
16/10/2008	21:31:52			-5.5	-6.1	-6.4	-65.6	-76.6	147	150
16/10/2008	21:36:52			-5.5	-6.1	-6.4	-65.5	-76.6	158	150
16/10/2008	21:41:52			-5.7	-6.1	-6.4	-66.5	-76.5	154	150

Batch Name: AFD8PK201

Page 6 of 19

Arino
18/10/08
Luke
18/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
16/10/2008	21:46:52			-5.7	-6.1	-6.4	-67.7	-77.1	144	150
16/10/2008	21:51:52			-5.4	-6.1	-6.4	-66.4	-77.1	155	150
16/10/2008	21:56:52			-5.4	-6.1	-6.4	-68.3	-77.6	159	150
16/10/2008	22:01:52			-5.7	-6.1	-6.4	-68.3	-79.4	147	150
16/10/2008	22:06:52			-5.5	-6.1	-6.4	-65.2	-79.9	159	150
16/10/2008	22:11:52			-5.7	-6.1	-6.4	-67.1	-79.9	159	150
16/10/2008	22:16:52			-5.7	-6.1	-6.4	-66.9	-78.4	150	150
16/10/2008	22:21:52			-5.8	-6.1	-6.4	-66.8	-78.9	159	150
16/10/2008	22:26:52			-5.4	-6.1	-6.4	-67.4	-78.9	159	150
16/10/2008	22:31:52			-5.4	-6.1	-6.4	-66.9	-79.4	148	150
16/10/2008	22:36:52			-5.7	-6.1	-6.4	-67.0	-78.9	204	150
16/10/2008	22:41:52			-5.6	-6.1	-6.4	-68.2	-79.9	149	150
16/10/2008	22:46:52			-5.7	-6.1	-6.4	-68.7	-79.0	153	150
16/10/2008	22:51:52			-5.8	-6.1	-6.4	-68.7	-79.0	160	150
16/10/2008	22:56:52			-5.4	-6.1	-6.4	-74.7	-79.5	143	150
16/10/2008	23:01:52			-5.4	-6.1	-6.4	-77.2	-79.5	153	150
16/10/2008	23:06:52			-5.8	-6.1	-6.4	-71.6	-80.5	149	150
16/10/2008	23:11:52			-5.7	-6.1	-6.4	-71.1	-81.5	146	150
16/10/2008	23:16:52			-5.6	-6.1	-6.4	-69.8	-80.5	150	150
16/10/2008	23:21:52			-5.5	-6.1	-6.4	-70.3	-79.7	147	150
16/10/2008	23:26:52			-5.5	-6.1	-6.4	-67.7	-79.2	152	150
16/10/2008	23:31:52			-5.6	-6.1	-6.4	-68.7	-80.9	148	150
16/10/2008	23:36:52			-5.7	-6.1	-6.4	-66.9	-80.9	150	150
16/10/2008	23:41:52			-5.7	-6.1	-6.4	-67.8	-80.9	153	150
16/10/2008	23:46:52			-5.5	-6.1	-6.4	-66.8	-80.8	151	150
16/10/2008	23:51:52			-5.5	-6.1	-6.4	-68.3	-80.2	148	150
16/10/2008	23:56:52			-5.4	-6.1	-6.4	-67.3	-81.3	148	150
17/10/2008	0:01:52			-5.6	-6.1	-6.4	-68.4	-82.3	148	150
17/10/2008	0:06:52			-5.4	-6.1	-6.4	-68.9	-82.3	148	150
17/10/2008	0:11:52			-5.5	-6.1	-6.4	-66.8	-82.3	149	150
17/10/2008	0:16:52			-5.5	-6.1	-6.4	-68.5	-82.3	150	150
17/10/2008	0:21:52			-5.7	-6.1	-6.4	-67.1	-82.3	150	150
17/10/2008	0:26:52			-5.7	-6.1	-6.4	-68.9	-82.3	150	150
17/10/2008	0:31:52			-5.6	-6.1	-6.4	-68.9	-82.8	152	150
17/10/2008	0:36:52			-5.5	-6.1	-6.4	-68.9	-82.8	152	150
17/10/2008	0:41:52			-5.5	-6.1	-6.4	-67.9	-83.3	149	150
17/10/2008	0:46:52			-5.5	-6.1	-6.4	-69.1	-82.8	148	150
17/10/2008	0:51:52			-5.6	-6.1	-6.4	-68.4	-82.8	150	150
17/10/2008	0:56:52			-5.6	-6.1	-6.4	-69.5	-82.8	150	150
17/10/2008	1:01:52			-5.5	-6.1	-6.4	-69.5	-82.8	152	150
17/10/2008	1:06:52			-5.5	-6.1	-6.4	-67.7	-82.8	153	150
17/10/2008	1:11:52			-5.6	-6.1	-6.4	-67.6	-82.3	146	150
17/10/2008	1:16:52			-5.6	-6.1	-6.4	-68.8	-82.3	146	150
17/10/2008	1:21:52			-5.6	-6.1	-6.4	-68.2	-82.3	146	150
17/10/2008	1:26:52			-5.6	-6.1	-6.4	-68.0	-82.3	147	150
17/10/2008	1:31:52			-5.5	-6.1	-6.4	-66.2	-82.3	147	150
17/10/2008	1:36:52			-5.6	-6.1	-6.4	-68.1	-82.3	150	150
17/10/2008	1:41:52			-5.5	-6.1	-6.4	-68.7	-82.3	151	150
17/10/2008	1:46:52			-5.6	-6.1	-6.4	-66.4	-81.8	150	150
17/10/2008	1:51:52			-5.6	-6.1	-6.4	-67.9	-81.8	149	150
17/10/2008	1:56:52			-5.7	-6.1	-6.4	-68.4	-81.8	150	150
17/10/2008	2:01:52			-5.7	-6.1	-6.4	-69.1	-81.8	151	150
17/10/2008	2:06:52			-5.6	-6.1	-6.4	-67.3	-81.8	151	150
17/10/2008	2:11:52			-5.5	-6.1	-6.4	-68.3	-81.8	152	150

Batch Name: AFD8PK201

Page 7 of 19

Approved
18/10/08

Lydon
18/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
17/10/2008	2:16:52			-5.4	-6.1	-6.4	-68.2	-81.8	153	150
17/10/2008	2:21:52			-5.5	-6.1	-6.4	-69.0	-81.8	154	150
17/10/2008	2:26:52			-5.5	-6.1	-6.4	-69.5	-81.3	153	150
17/10/2008	2:31:52			-5.5	-6.1	-6.4	-66.6	-81.3	152	150
17/10/2008	2:36:52			-5.4	-6.1	-6.4	-68.7	-81.3	152	150
17/10/2008	2:41:52			-5.5	-6.1	-6.4	-69.2	-80.3	150	150
17/10/2008	2:46:52			-5.4	-6.1	-6.4	-66.1	-80.3	152	150
17/10/2008	2:51:52			-5.5	-6.1	-6.4	-68.0	-80.8	154	150
17/10/2008	2:56:52			-5.6	-6.1	-6.4	-68.5	-80.8	152	150
17/10/2008	3:01:52			-5.4	-6.1	-6.4	-66.5	-80.8	152	150
17/10/2008	3:05:23	19	Drying	-5.8	-6.0	-6.2	-68.4	-80.9	150	150
17/10/2008	3:08:52			-5.4	-4.7	-6.4	-68.8	-80.8	121	100
17/10/2008	3:11:36	20	Drying	0.3	0.0	-2.4	-69.1	-80.9	102	100
17/10/2008	3:11:52			0.3	-0.1	-1.5	-69.4	-80.8	101	100
17/10/2008	3:16:52			0.2	-0.1	-0.9	-71.4	-80.8	104	100
17/10/2008	3:21:52			0.2	-0.1	-0.3	-72.8	-80.8	97	100
17/10/2008	3:26:52			0.4	-0.1	-0.3	-72.9	-80.8	102	100
17/10/2008	3:31:52			0.4	-0.1	-0.3	-73.5	-80.8	105	100
17/10/2008	3:36:52			0.5	-0.1	-0.3	-71.0	-80.8	97	100
17/10/2008	3:41:52			0.3	-0.1	-0.3	-69.7	-80.8	102	100
17/10/2008	3:46:52			0.3	-0.1	-0.3	-70.7	-80.8	104	100
17/10/2008	3:51:52			0.3	-0.1	-0.3	-70.7	-80.8	98	100
17/10/2008	3:56:52			0.4	-0.1	-0.3	-70.7	-80.8	101	100
17/10/2008	4:01:52			0.4	-0.1	-0.3	-68.1	-80.8	102	100
17/10/2008	4:06:52			0.5	-0.1	-0.3	-70.0	-80.8	103	100
17/10/2008	4:11:52			0.5	-0.1	-0.3	-69.2	-80.8	104	100
17/10/2008	4:16:52			0.5	-0.1	-0.3	-69.2	-80.3	103	100
17/10/2008	4:21:52			0.3	-0.1	-0.3	-70.2	-80.3	102	100
17/10/2008	4:26:52			0.5	-0.1	-0.3	-71.3	-80.3	97	100
17/10/2008	4:31:52			0.3	-0.1	-0.3	-71.8	-80.3	96	100
17/10/2008	4:36:52			0.5	-0.1	-0.3	-71.8	-80.3	98	100
17/10/2008	4:41:52			0.5	-0.1	-0.3	-72.4	-80.3	97	100
17/10/2008	4:46:52			0.3	-0.1	-0.3	-71.0	-80.8	96	100
17/10/2008	4:51:52			0.6	-0.1	-0.3	-72.1	-80.8	99	100
17/10/2008	4:56:52			0.3	-0.1	-0.3	-68.7	-80.8	101	100
17/10/2008	5:01:52			0.2	-0.1	-0.3	-69.8	-81.3	104	100
17/10/2008	5:06:52			0.3	-0.1	-0.3	-70.3	-81.3	103	100
17/10/2008	5:11:52			0.4	-0.1	-0.3	-69.2	-81.3	101	100
17/10/2008	5:16:52			0.5	-0.1	-0.3	-68.9	-81.3	100	100
17/10/2008	5:21:52			0.3	-0.1	-0.3	-69.5	-81.3	97	100
17/10/2008	5:26:52			0.5	-0.1	-0.3	-68.5	-81.3	102	100
17/10/2008	5:31:52			0.5	-0.1	-0.3	-70.0	-81.3	104	100
17/10/2008	5:36:52			0.3	-0.1	-0.3	-69.3	-81.3	104	100
17/10/2008	5:41:52			0.5	-0.1	-0.3	-70.2	-81.3	100	100
17/10/2008	5:46:52			0.5	-0.1	-0.3	-69.1	-81.3	97	100
17/10/2008	5:51:52			0.2	-0.1	-0.3	-70.3	-81.3	103	100
17/10/2008	5:56:52			0.3	-0.1	-0.3	-70.3	-81.3	103	100
17/10/2008	6:01:52			0.5	-0.1	-0.3	-68.2	-81.3	102	100
17/10/2008	6:06:52			0.4	-0.1	-0.3	-69.9	-81.3	99	100
17/10/2008	6:11:52			0.5	-0.1	-0.3	-69.9	-81.3	101	100
17/10/2008	6:16:52			0.4	-0.1	-0.3	-69.6	-80.8	104	100
17/10/2008	6:21:52			0.3	-0.1	-0.3	-70.6	-80.8	101	100
17/10/2008	6:26:52			0.5	-0.1	-0.3	-68.6	-81.3	96	100
17/10/2008	6:31:52			0.5	-0.1	-0.3	-70.5	-81.3	104	100

Batch Name: AFD8PK201

Page 8 of 19

Arvin
18/10/08

Lyons
18/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
17/10/2008	6:36:52			0.5	-0.1	-0.3	-69.1	-81.3	103	100
17/10/2008	6:41:52			0.4	-0.1	-0.3	-69.0	-81.3	100	100
17/10/2008	6:46:52			0.3	-0.1	-0.3	-69.5	-81.3	102	100
17/10/2008	6:51:52			0.4	-0.1	-0.3	-70.0	-81.3	104	100
17/10/2008	6:56:52			0.4	-0.1	-0.3	-68.1	-81.3	100	100
17/10/2008	7:01:52			0.6	-0.1	-0.3	-69.4	-81.3	101	100
17/10/2008	7:06:52			0.5	-0.1	-0.3	-69.9	-81.3	104	100
17/10/2008	7:11:52			0.4	-0.1	-0.3	-68.3	-80.8	101	100
17/10/2008	7:16:52			0.3	-0.1	-0.3	-70.5	-80.8	102	100
17/10/2008	7:21:52			0.5	-0.1	-0.3	-69.6	-81.3	102	100
17/10/2008	7:26:52			0.5	-0.1	-0.3	-71.0	-81.3	98	100
17/10/2008	7:31:52			0.2	-0.1	-0.3	-70.5	-81.8	103	100
17/10/2008	7:36:52			0.4	-0.1	-0.3	-71.2	-81.8	99	100
17/10/2008	7:41:52			0.5	-0.1	-0.3	-70.2	-81.8	104	100
17/10/2008	7:46:52			0.3	-0.1	-0.3	-72.2	-81.8	101	100
17/10/2008	7:51:52			0.5	-0.1	-0.3	-72.7	-81.8	102	100
17/10/2008	7:56:52			0.5	-0.1	-0.3	-73.3	-81.8	102	100
17/10/2008	8:01:52			0.3	-0.1	-0.3	-73.8	-81.3	101	100
17/10/2008	8:06:52			0.7	-0.1	-0.3	-69.3	-81.3	102	100
17/10/2008	8:11:52			0.3	-0.1	-0.3	-70.5	-81.3	102	100
17/10/2008	8:16:52			0.3	-0.1	-0.3	-69.4	-81.3	102	100
17/10/2008	8:21:52			0.3	-0.1	-0.3	-69.2	-81.3	101	100
17/10/2008	8:26:52			0.5	-0.1	-0.3	-70.3	-81.3	102	100
17/10/2008	8:31:52			0.5	-0.1	-0.3	-70.3	-81.3	103	100
17/10/2008	8:36:52			0.3	-0.1	-0.3	-68.7	-81.3	101	100
17/10/2008	8:41:52			0.4	-0.1	-0.3	-70.6	-81.3	104	100
17/10/2008	8:46:52			0.5	-0.1	-0.3	-71.4	-81.3	101	100
17/10/2008	8:51:52			0.5	-0.1	-0.3	-72.0	-81.3	103	100
17/10/2008	8:56:52			0.5	-0.1	-0.3	-72.0	-81.3	101	100
17/10/2008	9:01:52			0.3	-0.1	-0.3	-72.5	-81.3	104	100
17/10/2008	9:06:52			0.5	-0.1	-0.3	-73.0	-80.8	98	100
17/10/2008	9:11:52			0.5	-0.1	-0.3	-73.7	-80.8	103	100
17/10/2008	9:16:52			0.6	-0.1	-0.3	-70.1	-80.8	97	100
17/10/2008	9:21:52			0.5	-0.1	-0.3	-71.3	-80.3	103	100
17/10/2008	9:26:52			0.5	-0.1	-0.3	-71.8	-80.8	101	100
17/10/2008	9:31:52			0.5	-0.1	-0.3	-69.5	-80.8	102	100
17/10/2008	9:36:52			0.5	-0.1	-0.3	-70.8	-81.3	104	100
17/10/2008	9:41:52			0.5	-0.1	-0.3	-70.3	-81.3	99	100
17/10/2008	9:46:52			0.5	-0.1	-0.3	-70.4	-81.3	104	100
17/10/2008	9:51:52			0.5	-0.1	-0.3	-70.9	-81.3	98	100
17/10/2008	9:56:52			0.3	-0.1	-0.3	-68.6	-81.3	102	100
17/10/2008	10:01:52			0.2	-0.1	-0.3	-70.1	-81.3	104	100
17/10/2008	10:06:52			0.3	-0.1	-0.3	-70.6	-81.3	98	100
17/10/2008	10:11:52			0.3	-0.1	-0.3	-68.1	-81.3	104	100
17/10/2008	10:16:52			0.5	-0.1	-0.3	-69.7	-80.8	96	100
17/10/2008	10:21:52			0.6	-0.1	-0.3	-70.2	-79.8	102	100
17/10/2008	10:26:52			0.5	-0.1	-0.3	-70.3	-80.9	101	100
17/10/2008	10:31:52			0.5	-0.1	-0.3	-71.9	-80.9	101	100
17/10/2008	10:36:52			0.5	-0.1	-0.3	-70.6	-80.9	103	100
17/10/2008	10:41:52			0.6	-0.1	-0.3	-70.4	-80.9	96	100
17/10/2008	10:46:52			0.6	-0.1	-0.3	-70.9	-81.4	102	100
17/10/2008	10:51:52			0.5	-0.1	-0.3	-68.9	-81.4	104	100
17/10/2008	10:56:52			0.3	-0.1	-0.3	-70.1	-81.4	100	100
17/10/2008	11:01:52			0.3	-0.1	-0.3	-68.6	-81.4	101	100

Batch Name: AFD8PK201

Page 9 of 19

Alvin
K 18/10/08
Lupin
18/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A (°C)	Control Temp. (°C)	TE202 (°C)	TE203 (°C)	TE204 (°C)	PT101 (ubar)	Control Vacuum (ubar)
17/10/2008	11:06:52			0.5	-0.1	-0.3	-69.1	-81.4	104	100
17/10/2008	11:11:52			0.3	-0.1	-0.3	-69.6	-80.9	98	100
17/10/2008	11:16:52			0.6	-0.1	-0.3	-69.6	-80.9	103	100
17/10/2008	11:21:52			0.6	-0.1	-0.3	-68.8	-80.9	98	100
17/10/2008	11:26:52			0.2	-0.1	-0.3	-70.0	-80.9	102	100
17/10/2008	11:31:52			0.3	-0.1	-0.3	-67.8	-80.9	103	100
17/10/2008	11:36:52			0.5	-0.1	-0.3	-68.8	-80.9	100	100
17/10/2008	11:41:52			0.4	-0.1	-0.3	-69.3	-80.9	104	100
17/10/2008	11:46:52			0.4	-0.1	-0.3	-68.6	-80.9	96	100
17/10/2008	11:51:52			0.5	-0.1	-0.3	-70.4	-80.9	102	100
17/10/2008	11:56:52			0.5	-0.1	-0.3	-70.4	-80.9	102	100
17/10/2008	12:01:52			0.5	-0.1	-0.3	-68.2	-80.9	100	100
17/10/2008	12:06:52			0.5	-0.1	-0.3	-70.1	-81.4	104	100
17/10/2008	12:11:52			0.6	-0.1	-0.3	-70.6	-81.4	96	100
17/10/2008	12:16:52			0.6	-0.1	-0.3	-72.3	-81.4	102	100
17/10/2008	12:21:52			0.5	-0.1	-0.3	-72.3	-81.4	103	100
17/10/2008	12:26:52			0.5	-0.1	-0.3	-72.6	-81.9	99	100
17/10/2008	12:31:52			0.5	-0.1	-0.3	-71.7	-81.9	104	100
17/10/2008	12:36:52			0.5	-0.1	-0.3	-71.0	-81.9	96	100
17/10/2008	12:41:52			0.3	-0.1	-0.3	-69.7	-81.9	102	100
17/10/2008	12:46:52			0.3	-0.1	-0.3	-71.2	-81.9	103	100
17/10/2008	12:51:52			0.5	-0.1	-0.3	-68.6	-81.9	100	100
17/10/2008	12:56:52			0.4	-0.1	-0.3	-70.2	-81.9	104	100
17/10/2008	13:01:52			0.6	-0.1	-0.3	-70.2	-81.9	96	100
17/10/2008	13:06:52			0.5	-0.1	-0.3	-69.8	-81.4	102	100
17/10/2008	13:11:52			0.5	-0.1	-0.3	-70.2	-80.9	103	100
17/10/2008	13:16:52			0.5	-0.1	-0.3	-70.7	-80.9	101	100
17/10/2008	13:21:52			0.3	-0.1	-0.3	-70.7	-80.9	104	100
17/10/2008	13:26:52			0.3	-0.1	-0.3	-70.4	-81.4	98	100
17/10/2008	13:31:52			0.3	-0.1	-0.3	-70.9	-81.4	102	100
17/10/2008	13:36:52			0.5	-0.1	-0.3	-68.5	-81.4	104	100
17/10/2008	13:41:52			0.3	-0.1	-0.3	-69.8	-81.4	99	100
17/10/2008	13:46:52			0.5	-0.1	-0.3	-71.6	-81.4	102	100
17/10/2008	13:51:52			0.5	-0.1	-0.3	-73.7	-81.9	104	100
17/10/2008	13:56:52			0.3	-0.1	-0.3	-76.4	-81.9	99	100
17/10/2008	14:01:52			0.3	-0.1	-0.3	-73.5	-81.9	104	100
17/10/2008	14:06:52			0.4	-0.1	-0.3	-73.0	-81.9	96	100
17/10/2008	14:11:45	21	Drying	0.6	0.1	-0.1	-70.1	-81.4	102	80
17/10/2008	14:11:52			0.5	0.2	-0.3	-70.0	-81.4	103	100
17/10/2008	14:16:52			4.9	5.2	2.9	-69.9	-81.4	79	80
17/10/2008	14:18:39	22	Drying	6.4	7.0	5.4	-69.6	-81.4	84	80
17/10/2008	14:21:52			7.2	6.9	6.4	-69.9	-81.4	79	80
17/10/2008	14:26:52			7.2	6.9	6.4	-69.6	-80.9	85	80
17/10/2008	14:31:52			7.4	6.9	6.9	-71.2	-80.9	82	80
17/10/2008	14:36:52			7.5	6.9	6.9	-69.4	-81.4	78	80
17/10/2008	14:41:52			7.3	6.9	6.9	-70.5	-81.4	81	80
17/10/2008	14:46:52			7.4	6.9	6.9	-70.5	-81.4	85	80
17/10/2008	14:51:52			7.3	6.9	6.9	-69.8	-81.4	81	80
17/10/2008	14:56:52			7.4	6.9	6.9	-69.3	-81.4	77	80
17/10/2008	15:01:52			7.5	6.9	6.9	-70.3	-81.4	83	80
17/10/2008	15:06:52			7.5	6.9	6.9	-68.9	-81.4	85	80
17/10/2008	15:11:52			7.3	6.9	6.9	-70.1	-81.4	81	80
17/10/2008	15:16:52			7.5	6.9	6.9	-71.1	-80.9	79	80
17/10/2008	15:21:52			7.4	6.9	6.9	-68.2	-80.9	78	80

Batch Name: AFD8PK201

Page 10 of 19

Arino
K 18/10/08
by
12/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
17/10/2008	15:26:52			7.2	6.9	6.9	-70.1	-80.9	85	80
17/10/2008	15:31:52			7.6	6.9	6.9	-70.1	-80.9	82	80
17/10/2008	15:36:52			7.2	6.9	6.9	-69.1	-80.9	80	80
17/10/2008	15:41:52			7.4	6.9	6.9	-70.2	-80.9	77	80
17/10/2008	15:46:52			7.4	6.9	6.9	-68.4	-80.9	78	80
17/10/2008	15:51:52			7.4	6.9	6.9	-70.1	-80.9	86	80
17/10/2008	15:56:52			7.4	6.9	6.9	-70.1	-80.9	85	80
17/10/2008	16:01:52			7.3	6.9	6.9	-71.3	-80.9	81	80
17/10/2008	16:06:52			7.2	6.9	6.9	-69.0	-80.9	78	80
17/10/2008	16:11:52			7.5	6.9	6.9	-70.3	-80.9	80	80
17/10/2008	16:16:52			7.3	6.9	6.9	-69.4	-80.4	85	80
17/10/2008	16:21:52			7.5	6.9	6.9	-71.0	-80.9	82	80
17/10/2008	16:26:52			7.3	6.9	6.9	-71.5	-80.9	80	80
17/10/2008	16:31:52			7.4	6.9	6.9	-69.4	-80.9	76	80
17/10/2008	16:36:52			7.5	6.9	6.9	-70.5	-80.9	86	80
17/10/2008	16:41:52			7.3	6.9	6.9	-71.0	-80.9	85	80
17/10/2008	16:46:52			7.3	6.9	6.9	-71.0	-80.9	82	80
17/10/2008	16:51:52			7.3	6.9	6.9	-70.4	-80.9	80	80
17/10/2008	16:56:52			7.6	6.9	6.9	-71.5	-80.9	77	80
17/10/2008	17:01:52			7.2	6.9	6.9	-70.1	-80.9	84	80
17/10/2008	17:06:52			7.4	6.9	6.9	-70.6	-80.9	84	80
17/10/2008	17:11:52			7.4	6.9	6.9	-70.0	-80.9	82	80
17/10/2008	17:16:52			7.5	6.9	6.9	-70.5	-80.9	78	80
17/10/2008	17:18:52	23	Drying	7.4	7.0	6.9	-69.8	-80.8	83	80
17/10/2008	17:21:52			10.3	10.0	9.0	-71.1	-80.9	77	80
17/10/2008	17:26:52			14.6	15.0	12.4	-72.4	-80.9	86	80
17/10/2008	17:26:58	24	Drying	14.7	15.0	12.8	-72.5	-80.9	82	80
17/10/2008	17:31:52			15.4	15.0	14.8	-73.6	-80.9	83	80
17/10/2008	17:36:52			15.4	15.0	14.8	-74.1	-80.9	80	80
17/10/2008	17:41:52			15.4	15.0	15.3	-76.0	-81.4	77	80
17/10/2008	17:46:52			15.5	15.0	15.3	-75.4	-81.4	81	80
17/10/2008	17:51:52			15.5	15.0	15.3	-75.4	-81.4	85	80
17/10/2008	17:56:52			15.2	15.0	15.3	-75.4	-81.4	85	80
17/10/2008	18:01:52			15.2	15.0	15.3	-70.5	-81.4	84	80
17/10/2008	18:06:52			15.4	15.0	15.3	-71.7	-81.4	82	80
17/10/2008	18:11:52			15.4	15.0	15.3	-72.2	-81.4	79	80
17/10/2008	18:16:52			15.4	15.0	15.3	-71.7	-81.4	77	80
17/10/2008	18:21:52			15.4	15.0	15.3	-70.6	-81.4	76	80
17/10/2008	18:26:52			15.3	15.0	15.3	-71.1	-81.4	81	80
17/10/2008	18:31:52			15.3	15.0	15.3	-70.5	-81.4	86	80
17/10/2008	18:36:52			15.3	15.0	15.3	-70.0	-81.4	86	80
17/10/2008	18:41:52			15.5	15.0	15.3	-71.1	-81.4	183	80
17/10/2008	18:46:52			15.4	15.0	15.3	-69.3	-80.9	87	80
17/10/2008	18:51:52			15.5	15.0	15.3	-71.8	-80.9	84	80
17/10/2008	18:56:52			15.3	15.0	15.3	-70.9	-81.4	84	80
17/10/2008	19:01:52			15.3	15.0	15.3	-70.2	-82.0	84	80
17/10/2008	19:06:52			15.5	15.0	15.3	-72.2	-82.5	82	80
17/10/2008	19:11:52			15.3	15.0	15.3	-72.2	-82.5	81	80
17/10/2008	19:16:52			15.3	15.0	15.3	-73.9	-82.5	80	80
17/10/2008	19:21:52			15.4	15.0	15.3	-70.2	-82.5	79	80
17/10/2008	19:26:52			15.5	15.0	15.3	-70.8	-82.5	78	80
17/10/2008	19:31:52			15.5	15.0	15.3	-70.2	-82.5	78	80
17/10/2008	19:36:52			15.4	15.0	15.3	-71.2	-82.5	77	80
17/10/2008	19:41:52			15.7	15.0	15.3	-71.2	-82.5	76	80

Batch Name: AFD8PK201

Page 11 of 19

Handwritten:
 18/10/08
 12/10/18

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
17/10/2008	19:46:52			15.3	15.0	15.3	-69.6	-82.0	76	80
17/10/2008	19:51:52			15.5	15.0	15.3	-71.9	-82.0	76	80
17/10/2008	19:56:52			15.5	15.0	15.3	-71.1	-82.0	80	80
17/10/2008	20:01:52			15.4	15.0	15.3	-70.1	-82.0	79	80
17/10/2008	20:06:52			15.5	15.0	15.3	-71.4	-82.0	82	80
17/10/2008	20:11:52			15.3	15.0	15.3	-72.4	-82.0	86	80
17/10/2008	20:16:52			15.3	15.0	15.3	-72.4	-82.0	86	80
17/10/2008	20:21:52			15.4	15.0	15.3	-71.9	-82.0	85	80
17/10/2008	20:26:52			15.3	15.0	15.3	-71.4	-81.5	84	80
17/10/2008	20:31:49	25	Drying	15.2	15.0	15.4	-70.1	-81.8	82	60
17/10/2008	20:31:52			15.4	15.1	15.3	-70.2	-81.5	83	80
17/10/2008	20:36:52			19.7	20.1	17.9	-72.3	-81.5	58	60
17/10/2008	20:41:52			24.1	25.1	23.5	-73.6	-81.5	66	60
17/10/2008	20:46:45	26	Drying	29.5	30.0	28.4	-75.2	-81.7	59	60
17/10/2008	20:46:52			29.6	29.9	28.4	-74.9	-82.0	59	60
17/10/2008	20:51:52			30.2	29.9	31.1	-74.9	-82.0	63	60
17/10/2008	20:56:52			30.1	29.9	30.6	-74.4	-82.0	65	60
17/10/2008	21:01:52			30.1	29.9	30.6	-72.3	-82.0	62	60
17/10/2008	21:06:52			30.1	29.9	30.6	-72.8	-82.0	60	60
17/10/2008	21:11:52			30.4	29.9	30.6	-72.8	-82.0	58	60
17/10/2008	21:16:52			30.4	29.9	31.1	-73.8	-82.0	58	60
17/10/2008	21:21:52			30.2	29.9	31.1	-71.5	-82.0	57	60
17/10/2008	21:26:52			30.5	29.9	31.1	-72.3	-82.0	56	60
17/10/2008	21:31:52			30.3	29.9	31.1	-72.3	-82.0	57	60
17/10/2008	21:36:52			30.5	29.9	31.1	-71.8	-82.0	57	60
17/10/2008	21:41:52			30.4	29.9	31.1	-73.4	-82.0	59	60
17/10/2008	21:46:52			30.5	29.9	31.1	-70.4	-82.0	61	60
17/10/2008	21:51:52			30.4	29.9	31.1	-71.6	-81.5	63	60
17/10/2008	21:56:52			30.6	29.9	31.1	-70.5	-81.5	65	60
17/10/2008	22:01:52			30.4	29.9	31.1	-71.7	-81.5	66	60
17/10/2008	22:06:52			30.4	29.9	31.1	-72.7	-81.5	58	60
17/10/2008	22:11:52			30.2	29.9	31.1	-70.7	-81.0	58	60
17/10/2008	22:16:52			30.4	29.9	31.1	-72.0	-81.5	59	60
17/10/2008	22:21:52			30.2	29.9	31.1	-72.5	-81.5	62	60
17/10/2008	22:26:52			30.1	29.9	31.1	-70.8	-81.5	63	60
17/10/2008	22:31:52			30.4	29.9	31.1	-70.8	-81.5	66	60
17/10/2008	22:36:52			30.4	29.9	31.1	-69.8	-81.5	162	60
17/10/2008	22:41:52			30.2	29.9	31.1	-71.5	-81.0	81	60
17/10/2008	22:46:52			30.4	29.9	31.1	-72.0	-81.5	61	60
17/10/2008	22:51:52			30.4	29.9	31.1	-74.2	-82.0	62	60
17/10/2008	22:56:52			30.2	29.9	31.1	-74.7	-82.0	62	60
17/10/2008	23:01:52			30.4	29.9	31.1	-79.5	-82.0	64	60
17/10/2008	23:06:52			30.4	29.9	31.1	-79.0	-82.0	64	60
17/10/2008	23:11:52			30.6	29.9	31.1	-80.3	-82.0	64	60
17/10/2008	23:16:52			30.5	29.9	31.1	-80.8	-82.0	65	60
17/10/2008	23:21:52			30.4	29.9	31.1	-80.8	-82.0	65	60
17/10/2008	23:26:52			30.5	29.9	31.1	-81.4	-82.0	65	60
17/10/2008	23:31:52			30.2	29.9	31.1	-82.5	-82.0	65	60
17/10/2008	23:36:52			30.1	29.9	31.1	-82.4	-82.0	64	60
17/10/2008	23:41:52			30.6	29.9	31.1	-81.9	-82.0	64	60
17/10/2008	23:46:52			30.4	29.9	31.1	-81.4	-82.0	64	60
17/10/2008	23:51:52			30.3	29.9	31.1	-81.8	-82.0	63	60
17/10/2008	23:56:52			30.2	29.9	31.1	-82.3	-81.5	63	60
18/10/2008	0:01:52			30.6	29.9	31.1	-81.8	-82.0	62	60

Batch Name: AFD8PK201

Page 12 of 19

Avila
18/10/08
hmc
18/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
18/10/2008	0:06:52			30.4	29.9	31.1	-82.3	-82.0	60	60
18/10/2008	0:11:52			30.4	29.9	31.1	-81.8	-82.0	60	60
18/10/2008	0:16:52			30.4	29.9	31.1	-81.3	-82.0	59	60
18/10/2008	0:21:52			30.6	29.9	31.1	-82.5	-82.0	58	60
18/10/2008	0:26:52			30.6	29.9	31.1	-81.3	-82.0	61	60
18/10/2008	0:31:52			30.4	29.9	31.1	-81.4	-82.0	66	60
18/10/2008	0:36:52			30.2	29.9	31.1	-81.4	-82.5	65	60
18/10/2008	0:41:52			30.2	29.9	31.1	-81.4	-83.0	64	60
18/10/2008	0:46:52			30.4	29.9	31.1	-82.0	-83.0	61	60
18/10/2008	0:51:52			30.5	29.9	31.1	-82.0	-82.5	57	60
18/10/2008	0:56:52			30.2	29.9	31.1	-82.0	-82.5	64	60
18/10/2008	1:01:52			30.1	29.9	31.1	-82.6	-81.9	64	60
18/10/2008	1:06:52			30.6	29.9	31.1	-82.6	-81.8	61	60
18/10/2008	1:11:52			30.4	29.9	31.1	-81.8	-82.3	58	60
18/10/2008	1:16:52			30.2	29.9	31.1	-82.3	-82.3	67	60
18/10/2008	1:21:52			30.4	29.9	31.1	-81.7	-82.8	65	60
18/10/2008	1:26:52			30.4	29.9	31.1	-81.7	-82.8	60	60
18/10/2008	1:31:52			30.4	29.9	31.1	-81.8	-82.8	57	60
18/10/2008	1:36:52			30.4	29.9	31.1	-81.8	-82.8	66	60
18/10/2008	1:41:52			30.1	29.9	31.1	-82.4	-82.8	63	60
18/10/2008	1:46:52			30.6	29.9	31.1	-82.9	-82.8	58	60
18/10/2008	1:51:52			30.4	29.9	31.1	-81.9	-82.8	63	60
18/10/2008	1:56:52			30.3	29.9	31.1	-81.6	-82.8	64	60
18/10/2008	2:01:52			30.3	29.9	31.1	-81.1	-82.8	61	60
18/10/2008	2:06:52			30.4	29.9	31.1	-81.1	-82.8	56	60
18/10/2008	2:11:52			30.1	29.9	31.1	-80.8	-82.8	66	60
18/10/2008	2:16:52			30.4	29.9	31.1	-81.4	-82.8	62	60
18/10/2008	2:21:52			30.4	29.9	31.1	-81.4	-82.8	57	60
18/10/2008	2:26:52			30.4	29.9	31.1	-80.5	-82.3	66	60
18/10/2008	2:31:52			30.2	29.9	31.1	-81.1	-82.3	63	60
18/10/2008	2:36:52			30.3	29.9	31.1	-81.7	-82.3	58	60
18/10/2008	2:41:52			30.4	29.9	31.1	-80.9	-82.3	66	60
18/10/2008	2:46:52			30.4	29.9	31.1	-81.4	-82.3	63	60
18/10/2008	2:51:52			30.2	29.9	31.1	-81.9	-82.3	59	60
18/10/2008	2:56:52			30.1	29.9	31.1	-81.4	-82.3	67	60
18/10/2008	3:01:52			30.4	29.9	31.1	-81.3	-82.3	63	60
18/10/2008	3:06:52			30.4	29.9	31.1	-81.3	-82.3	57	60
18/10/2008	3:11:52			30.3	29.9	31.1	-81.8	-82.3	66	60
18/10/2008	3:16:52			30.5	29.9	31.1	-82.4	-82.3	62	60
18/10/2008	3:21:52			30.4	29.9	31.1	-82.3	-82.3	57	60
18/10/2008	3:26:52			30.2	29.9	31.1	-82.3	-82.3	65	60
18/10/2008	3:31:52			30.1	29.9	31.1	-82.3	-82.3	61	60
18/10/2008	3:36:52			30.6	29.9	31.1	-82.3	-82.3	59	60
18/10/2008	3:41:52			30.4	29.9	31.1	-82.3	-82.3	64	60
18/10/2008	3:46:52			30.1	29.9	31.1	-82.3	-82.8	59	60
18/10/2008	3:51:52			30.5	29.9	31.1	-82.3	-82.8	67	60
18/10/2008	3:56:52			30.3	29.9	31.1	-81.8	-82.8	62	60
18/10/2008	4:01:52			30.2	29.9	31.1	-82.3	-82.8	57	60
18/10/2008	4:06:52			30.4	29.9	31.1	-82.3	-82.8	65	60
18/10/2008	4:11:52			30.4	29.9	31.1	-82.3	-82.8	59	60
18/10/2008	4:16:52			30.1	29.9	31.1	-82.3	-82.8	63	60
18/10/2008	4:21:52			30.4	29.9	31.1	-78.8	-82.8	63	60
18/10/2008	4:26:52			30.4	29.9	31.1	-80.0	-82.8	58	60
18/10/2008	4:31:52			30.4	29.9	31.1	-80.1	-82.8	66	60

Batch Name: AFD8PK201

Page 13 of 19

Archive
18/10/08
18/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
18/10/2008	4:36:52			30.2	29.9	31.1	-81.3	-83.3	60	60
18/10/2008	4:41:52			30.5	29.9	31.1	-81.4	-82.8	61	60
18/10/2008	4:46:52			30.4	29.9	31.1	-80.9	-82.8	64	60
18/10/2008	4:51:52			30.4	29.9	31.1	-81.6	-82.8	58	60
18/10/2008	4:52:38	27	Drying	30.4	30.0	31.0	-81.7	-82.7	61	60
18/10/2008	4:56:52			33.1	34.1	33.4	-82.7	-82.8	67	60
18/10/2008	5:01:52			39.1	39.2	38.6	-83.2	-82.8	62	60
18/10/2008	5:06:52			43.3	44.2	42.7	-83.7	-82.8	59	60
18/10/2008	5:07:40	28	Drying	44.8	45.0	44.2	-84.0	-82.9	62	60
18/10/2008	5:11:52			45.0	45.0	46.4	-83.1	-82.8	58	60
18/10/2008	5:16:52			45.3	45.0	46.4	-81.2	-82.8	59	60
18/10/2008	5:21:52			45.1	45.0	46.4	-81.8	-82.8	64	60
18/10/2008	5:26:52			45.1	45.0	46.4	-82.8	-82.8	62	60
18/10/2008	5:31:52			45.1	45.0	46.4	-83.3	-82.8	63	60
18/10/2008	5:36:52			45.3	45.0	46.4	-83.3	-82.3	59	60
18/10/2008	5:41:52			45.2	45.0	46.4	-83.3	-82.3	65	60
18/10/2008	5:46:52			45.3	45.0	46.4	-83.3	-82.3	66	60
18/10/2008	5:51:52			45.4	45.0	46.4	-83.3	-82.3	66	60
18/10/2008	5:56:52			45.4	45.0	46.4	-83.8	-82.3	65	60
18/10/2008	6:01:52			45.4	45.0	46.4	-83.8	-82.3	66	60
18/10/2008	6:06:52			44.9	45.0	45.5	-83.8	-82.3	66	60
18/10/2008	6:11:52			44.4	45.0	46.5	-82.0	-82.3	66	60
18/10/2008	6:16:52			44.6	45.0	46.6	-78.7	-82.3	66	60
18/10/2008	6:21:52			45.9	45.0	46.2	-77.1	-82.3	66	60
18/10/2008	6:26:52			45.4	45.0	46.1	-77.1	-82.3	66	60
18/10/2008	6:31:52			44.3	45.0	46.7	-77.6	-82.3	66	60
18/10/2008	6:36:52			45.3	45.0	46.6	-77.1	-82.3	66	60
18/10/2008	6:41:52			45.8	45.0	45.6	-76.6	-81.8	66	60
18/10/2008	6:46:52			44.8	45.0	45.7	-77.1	-81.8	66	60
18/10/2008	6:51:52			44.5	45.0	46.8	-77.0	-81.8	64	60
18/10/2008	6:56:52			45.9	45.0	46.2	-76.1	-81.8	63	60
18/10/2008	7:01:52			44.9	45.0	45.5	-76.6	-81.8	62	60
18/10/2008	7:06:52			45.1	45.0	46.8	-75.5	-81.8	60	60
18/10/2008	7:11:52			45.3	45.0	46.8	-77.0	-81.8	57	60
18/10/2008	7:16:52			45.3	45.0	46.8	-78.1	-81.3	67	60
18/10/2008	7:21:52			45.3	45.0	46.8	-78.1	-81.3	66	60
18/10/2008	7:26:52			45.5	45.0	46.8	-79.8	-81.3	63	60
18/10/2008	7:31:52			45.3	45.0	46.8	-79.8	-81.3	60	60
18/10/2008	7:36:52			45.2	45.0	46.8	-81.3	-82.3	64	60
18/10/2008	7:41:52			45.3	45.0	46.8	-81.8	-82.8	65	60
18/10/2008	7:46:52			45.4	45.0	46.8	-80.7	-82.8	60	60
18/10/2008	7:51:52			45.9	45.0	46.6	-80.7	-82.8	60	60
18/10/2008	7:56:52			45.3	45.0	46.7	-80.2	-82.8	64	60
18/10/2008	8:01:52			45.4	45.0	46.7	-81.3	-82.8	59	60
18/10/2008	8:06:52			45.4	45.0	46.7	-81.8	-82.8	67	60
18/10/2008	8:11:52			45.1	45.0	46.7	-81.8	-82.8	62	60
18/10/2008	8:16:52			45.4	45.0	46.7	-81.7	-82.8	61	60
18/10/2008	8:21:52			45.6	45.0	46.7	-81.2	-82.8	63	60
18/10/2008	8:26:52			45.5	45.0	46.7	-81.2	-82.8	56	60
18/10/2008	8:31:52			45.3	45.0	46.7	-80.7	-82.1	158	60
18/10/2008	8:36:52			44.9	45.0	46.7	-80.4	-81.3	102	60
18/10/2008	8:41:52			45.3	45.0	46.8	-81.4	-81.3	57	60
18/10/2008	8:46:52			45.4	45.0	46.8	-81.9	-81.8	62	60
18/10/2008	8:51:52			45.4	45.0	46.8	-82.4	-81.6	66	60

Batch Name: AFD8PK201

Page 14 of 19

Abir
18/10/08
hishu
18/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A (°C)	Control Temp. (°C)	TE202 (°C)	TE203 (°C)	TE204 (°C)	PT101 (ubar)	Control Vacuum (ubar)
18/10/2008	8:56:52			45.4	45.0	46.8	-81.4	-81.6	57	60
18/10/2008	9:01:52			45.4	45.0	46.8	-81.3	-81.6	60	60
18/10/2008	9:06:52			45.5	45.0	46.8	-81.8	-81.2	63	60
18/10/2008	9:11:52			45.4	45.0	46.8	-81.8	-81.2	67	60
18/10/2008	9:16:52			45.2	45.0	46.8	-81.2	-81.2	56	60
18/10/2008	9:21:52			45.4	45.0	46.8	-82.2	-82.2	59	60
18/10/2008	9:26:52			45.3	45.0	46.8	-82.7	-82.2	61	60
18/10/2008	9:31:52			45.5	45.0	46.8	-82.2	-82.7	64	60
18/10/2008	9:36:52			45.4	45.0	46.8	-82.2	-82.7	66	60
18/10/2008	9:41:52			45.6	45.0	46.8	-82.7	-83.2	66	60
18/10/2008	9:46:52			45.4	45.0	46.8	-81.7	-82.6	56	60
18/10/2008	9:51:52			45.4	45.0	46.8	-81.7	-82.6	58	60
18/10/2008	9:56:52			45.5	45.0	46.8	-82.2	-82.5	58	60
18/10/2008	10:01:52			45.7	45.0	46.8	-82.2	-82.5	60	60
18/10/2008	10:06:52			45.8	45.0	46.8	-82.7	-82.5	61	60
18/10/2008	10:11:52			45.9	45.0	46.8	-82.7	-82.5	61	60
18/10/2008	10:16:52			46.0	45.0	47.3	-82.7	-83.5	62	60
18/10/2008	10:21:52			46.1	45.0	47.3	-82.7	-84.0	63	60
18/10/2008	10:26:52			46.3	45.0	47.3	-82.7	-84.0	65	60
18/10/2008	10:31:52			46.4	45.0	47.3	-82.7	-84.0	65	60
18/10/2008	10:36:52			46.5	45.0	47.8	-82.2	-84.0	66	60
18/10/2008	10:41:52			46.6	45.0	47.8	-82.2	-84.0	66	60
18/10/2008	10:46:52			46.8	45.0	47.8	-81.7	-84.0	64	60
18/10/2008	10:51:52			47.1	45.0	47.8	-81.7	-83.5	57	60
18/10/2008	10:56:52			47.0	45.0	48.3	-81.7	-83.0	57	60
18/10/2008	11:01:52			47.1	45.0	48.3	-81.2	-83.0	58	60
18/10/2008	11:06:52			47.3	45.0	48.3	-81.2	-83.0	59	60
18/10/2008	11:11:52			47.4	45.0	48.3	-81.2	-83.0	60	60
18/10/2008	11:16:52			47.4	45.0	48.8	-80.7	-83.0	61	60
18/10/2008	11:21:52			47.7	45.0	48.8	-80.7	-82.5	62	60
18/10/2008	11:26:52			47.8	45.0	48.8	-80.7	-82.5	62	60
18/10/2008	11:31:52			47.8	45.0	48.8	-80.2	-83.0	63	60
18/10/2008	11:36:52			47.8	45.0	48.8	-80.2	-83.0	64	60
18/10/2008	11:41:52			47.8	45.0	49.3	-80.2	-83.0	65	60
18/10/2008	11:46:52			48.1	45.0	49.3	-79.7	-83.0	65	60
18/10/2008	11:51:52			48.1	45.0	49.3	-79.7	-83.0	66	60
18/10/2008	11:56:52			48.2	45.0	49.3	-79.7	-82.5	66	60
18/10/2008	12:01:52			48.4	45.0	49.3	-79.7	-82.5	66	60
18/10/2008	12:06:52			48.4	45.0	49.3	-79.7	-82.5	61	60
18/10/2008	12:11:52			48.4	45.0	49.8	-79.2	-82.5	56	60
18/10/2008	12:16:52			48.6	45.0	49.8	-79.2	-82.5	56	60
18/10/2008	12:21:52			48.6	45.0	49.8	-79.2	-82.0	57	60
18/10/2008	12:26:52			48.8	45.0	49.8	-78.7	-81.0	58	60
18/10/2008	12:31:52			48.9	45.0	49.8	-78.7	-82.1	57	60
18/10/2008	12:36:52			48.9	45.0	50.3	-78.7	-82.1	58	60
18/10/2008	12:41:52			48.9	45.0	50.3	-78.7	-82.1	58	60
18/10/2008	12:46:52			49.0	45.0	50.3	-78.2	-82.6	58	60
18/10/2008	12:51:52			49.2	45.0	50.3	-78.2	-82.6	60	60
18/10/2008	12:56:52			49.3	45.0	50.8	-78.2	-82.6	58	60
18/10/2008	13:01:52			49.3	45.0	50.8	-78.2	-82.6	59	60
18/10/2008	13:06:52			49.5	45.0	50.8	-78.2	-82.6	58	60
18/10/2008	13:11:52			49.5	45.0	50.8	-78.2	-82.6	58	60
18/10/2008	13:16:52			49.7	45.0	50.8	-77.7	-82.6	57	60
18/10/2008	13:21:52			49.8	45.0	51.3	-77.7	-82.6	57	60

Batch Name: AFD8PK201

Page 15 of 19

Arvin
18/10/08
Arvin
18/10/08


Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
18/10/2008	13:26:52			50.0	45.0	51.3	-77.7	-82.6	57	60
18/10/2008	13:31:52			50.0	45.0	51.3	-77.7	-82.1	56	60
18/10/2008	13:36:52			50.2	45.0	51.3	-77.7	-82.1	59	60
18/10/2008	13:41:52			50.2	45.0	51.3	-77.7	-82.1	63	60
18/10/2008	13:46:52			50.3	45.0	51.8	-77.7	-82.1	66	60
18/10/2008	13:51:52			50.4	45.0	51.8	-77.2	-82.1	66	60
18/10/2008	13:56:52			50.5	45.0	51.8	-77.2	-82.1	65	60
18/10/2008	14:01:52			50.6	45.0	51.8	-77.2	-82.1	63	60
18/10/2008	14:06:52			49.6	45.0	51.8	-77.2	-82.1	117	60
18/10/2008	14:11:52			48.0	45.0	50.0	-77.2	-81.6	66	60
18/10/2008	14:16:52			46.9	45.0	48.8	-77.2	-81.6	64	60
18/10/2008	14:21:52			45.7	45.0	47.3	-77.2	-81.6	61	60
18/10/2008	14:26:52			45.4	45.0	46.8	-67.1	-81.6	59	60
18/10/2008	14:31:52			45.5	45.0	46.8	-72.1	-82.1	65	60
18/10/2008	14:36:52			45.7	45.0	46.8	-75.4	-82.1	64	60
18/10/2008	14:41:52			45.9	45.0	46.8	-75.9	-82.1	58	60
18/10/2008	14:46:52			46.0	45.0	47.3	-75.9	-82.1	66	60
18/10/2008	14:51:52			46.2	45.0	47.3	-75.9	-82.1	62	60
18/10/2008	14:56:52			46.4	45.0	47.3	-75.9	-82.1	57	60
18/10/2008	15:01:52			46.6	45.0	47.8	-75.9	-82.1	65	60
18/10/2008	15:06:52			46.6	45.0	47.8	-75.9	-82.1	59	60
18/10/2008	15:11:52			46.8	45.0	47.8	-75.9	-82.1	66	60
18/10/2008	15:16:52			46.9	45.0	48.3	-75.9	-81.6	61	60
18/10/2008	15:21:52			47.1	45.0	48.3	-75.9	-81.6	62	60
18/10/2008	15:26:52			47.2	45.0	48.3	-75.9	-81.6	63	60
18/10/2008	15:31:52			47.4	45.0	48.3	-75.9	-81.6	58	60
18/10/2008	15:36:52			47.4	45.0	48.8	-75.9	-82.1	63	60
18/10/2008	15:41:52			47.6	45.0	48.8	-75.9	-82.1	57	60
18/10/2008	15:46:52			47.8	45.0	48.8	-75.9	-82.1	63	60
18/10/2008	15:51:52			48.0	45.0	49.3	-75.9	-82.1	57	60
18/10/2008	15:56:52			48.1	45.0	49.3	-75.9	-82.6	65	60
18/10/2008	16:01:52			48.1	45.0	49.3	-75.9	-82.1	57	60
18/10/2008	16:06:52			48.3	45.0	49.3	-75.9	-82.1	65	60
18/10/2008	16:11:52			48.4	45.0	49.8	-75.9	-82.1	58	60
18/10/2008	16:16:52			48.5	45.0	49.8	-75.9	-82.1	66	60
18/10/2008	16:21:52			48.7	45.0	49.8	-75.9	-82.1	58	60
18/10/2008	16:26:52			48.8	45.0	49.8	-75.9	-82.1	66	60
18/10/2008	16:31:52			49.0	45.0	49.8	-75.9	-82.1	60	60
18/10/2008	16:36:52			49.0	45.0	50.3	-75.9	-82.1	66	60
18/10/2008	16:41:52			49.0	45.0	50.3	-75.9	-82.1	60	60
18/10/2008	16:46:52			49.2	45.0	50.3	-75.9	-82.1	66	60
18/10/2008	16:51:52			49.3	45.0	50.8	-75.9	-82.1	61	60
18/10/2008	16:56:52			49.5	45.0	50.8	-75.9	-82.1	66	60
18/10/2008	17:01:52			49.6	45.0	50.8	-75.9	-82.1	61	60
18/10/2008	17:06:52			49.6	45.0	50.8	-75.9	-81.6	61	60
18/10/2008	17:11:52			49.7	45.0	50.8	-75.9	-81.6	63	60
18/10/2008	17:16:52			49.9	45.0	51.3	-75.9	-81.6	58	60
18/10/2008	17:21:52			50.2	45.0	51.3	-75.9	-81.6	63	60
18/10/2008	17:26:52			49.6	45.0	51.3	-75.9	-81.6	65	60
18/10/2008	17:31:52			48.9	45.0	50.8	-75.9	-81.6	91	60
18/10/2008	17:36:52			48.0	45.0	49.8	-75.9	-81.6	60	60
18/10/2008	17:41:52			46.7	45.0	48.6	-75.9	-82.1	66	60
18/10/2008	17:44:17	29	Moisture Test	46.0	45.0	47.6	-75.9	-82.1	63	60
18/10/2008	17:45:33	30	Pre-Aeration	45.6	45.0	47.4	-75.9	-82.1	183	--

Batch Name: AFD8PK201

Page 16 of 19


 18/10/08
 18/10/08

Auto Freeze Drying Process Report

Start Time:	18/10/2008 0:41:52	Product Quantity:	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Step No.	Step Type	TE201A ("C)	Control Temp. ("C)	TE202 ("C)	TE203 ("C)	TE204 ("C)	PT101 (ubar)	Control Vacuum (ubar)
18/10/2008	17:46:52			45.1	45.0	45.6	-75.1	-66.4	1000	--
18/10/2008	17:51:52			44.9	45.0	46.4	-70.1	-60.3	1000	--
18/10/2008	17:56:52			45.1	45.0	46.4	-64.3	-57.9	1000	--
18/10/2008	18:01:52			45.0	45.0	46.4	-60.2	-56.1	1000	--
18/10/2008	18:02:58	31	Stoppering	45.2	45.0	46.3	-59.7	-56.9	1000	--
18/10/2008	18:06:16	32	Aeration	44.9	45.0	46.3	-58.0	-55.9	1000	--
18/10/2008	18:06:52			44.9	45.0	46.4	-57.7	-55.3	1000	--
18/10/2008	18:11:52			45.1	45.0	46.4	-49.2	-51.3	1000	--
18/10/2008	18:16:00	33	Storage	44.9	45.0	46.3	-38.4	-45.3	1000	--
18/10/2008	18:16:52			44.1	25.0	46.4	-38.0	-44.1	1000	--
18/10/2008	18:21:52			37.6	25.0	40.0	-35.8	-40.9	1000	--
18/10/2008	18:32:52			25.7	25.0	27.5	-31.0	-34.0	1000	--

Arrived
 18/10/08
 18/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity	4615
End Time:	18/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Description
Alarm: Total: 91		
16/10/2008	4:09:19	Circulation pump 1 failure
16/10/2008	4:09:23	Circulation pump 2 failure
16/10/2008	4:09:30	Compressor 1 failure
16/10/2008	4:09:40	Compressor 2 failure
16/10/2008	4:10:37	Compressor 1 cooling water low flowrate
16/10/2008	4:10:37	Compressor 2 cooling water low flowrate
16/10/2008	6:07:53	Circulation pump 2 failure
16/10/2008	6:07:53	Circulation pump 1 failure
16/10/2008	6:08:00	Compressor 1 failure
16/10/2008	6:08:17	Compressor 1 cooling water low flowrate
16/10/2008	6:08:17	Compressor 2 cooling water low flowrate
16/10/2008	7:30:01	Vacuum pump 1 failure
16/10/2008	7:30:11	PT301 pump vacuum probe fail
16/10/2008	7:54:13	Vacuum less than low limit
16/10/2008	8:13:33	Vacuum failure pre-alarm
16/10/2008	8:20:28	Vacuum less than low limit
16/10/2008	8:25:23	Vacuum less than low limit
16/10/2008	8:32:14	FV5 gas bleed bypass valve fail
16/10/2008	8:32:25	FV5 gas bleed bypass valve fail
16/10/2008	8:33:39	FV5 gas bleed bypass valve fail
16/10/2008	8:34:07	Vacuum failure pre-alarm
16/10/2008	8:34:10	Vacuum failure alarm
16/10/2008	8:38:36	Vacuum less than low limit
16/10/2008	8:44:11	Vacuum failure pre-alarm
16/10/2008	8:44:18	Vacuum failure alarm
16/10/2008	8:49:05	Vacuum less than low limit
16/10/2008	9:14:25	FV5 gas bleed bypass valve fail
16/10/2008	9:14:45	Vacuum failure pre-alarm
16/10/2008	9:14:45	Vacuum failure alarm
16/10/2008	9:46:43	Vacuum less than low limit
16/10/2008	9:49:34	Vacuum failure pre-alarm
16/10/2008	9:52:29	Vacuum less than low limit
16/10/2008	9:54:51	Vacuum failure pre-alarm
16/10/2008	18:33:57	FV25 filter isolation valve fail
16/10/2008	18:33:57	Circulation pump 2 failure
16/10/2008	18:33:57	Circulation pump 1 failure
16/10/2008	18:33:57	FV4 aerate isolation valve fail
16/10/2008	18:33:57	FV29 chamber inlet valve 1 fail
16/10/2008	18:34:07	Compressor 2 failure
16/10/2008	18:34:14	Compressor 1 failure
16/10/2008	18:35:29	Vacuum failure pre-alarm
16/10/2008	18:41:57	Vacuum less than low limit
16/10/2008	18:45:29	Vacuum failure pre-alarm
16/10/2008	22:34:52	FV25 filter isolation valve fail
16/10/2008	22:34:52	Circulation pump 1 failure
16/10/2008	22:34:52	FV4 aerate isolation valve fail
16/10/2008	22:34:52	FV29 chamber inlet valve 1 fail
16/10/2008	22:34:56	Circulation pump 2 failure
16/10/2008	22:35:03	Compressor 1 failure
16/10/2008	22:35:11	Compressor 2 failure
16/10/2008	22:35:40	Compressor 1 cooling water low flowrate
16/10/2008	22:35:40	Compressor 2 cooling water low flowrate
16/10/2008	22:37:07	Vacuum failure pre-alarm
17/10/2008	18:39:24	Circulation pump 1 failure
17/10/2008	18:39:28	FV25 filter isolation valve fail
17/10/2008	18:39:28	Circulation pump 2 failure

Batch Name: AFD8PK201

Page 18 of 19

Armed
K 18/10/08
by
18/10/08

Auto Freeze Drying Process Report

Start Time:	16/10/2008 0:41:52	Product Quantity:	4615
End Time:	16/10/2008 18:32:52	Product Name:	Topotecan For Injection
Initiator:	SUP	Product Model:	
Batch Name:	AFD8PK201	Recipe Number:	25

Date	Time	Description
17/10/2008	18:39:28	FV4 aerate isolation valve fail
17/10/2008	18:39:28	FV29 chamber inlet valve 1 fail
17/10/2008	18:39:37	Compressor 1 failure
17/10/2008	18:39:46	Compressor 2 failure
17/10/2008	18:40:04	Compressor 1 cooling water low flowrate
17/10/2008	18:40:04	Compressor 2 cooling water low flowrate
17/10/2008	18:42:08	Vacuum failure pre-alarm
17/10/2008	22:34:14	FV25 filter isolation valve fail
17/10/2008	22:34:14	Circulation pump 1 failure
17/10/2008	22:34:14	FV4 aerate isolation valve fail
17/10/2008	22:34:14	FV29 chamber inlet valve 1 fail
17/10/2008	22:34:19	Circulation pump 2 failure
17/10/2008	22:34:23	Compressor 2 failure
17/10/2008	22:34:36	Compressor 1 failure
17/10/2008	22:34:59	Compressor 1 cooling water low flowrate
17/10/2008	22:34:59	Compressor 2 cooling water low flowrate
17/10/2008	22:36:54	Vacuum failure pre-alarm
18/10/2008	8:28:43	FV25 filter isolation valve fail
18/10/2008	8:28:43	Circulation pump 2 failure
18/10/2008	8:28:43	Circulation pump 1 failure
18/10/2008	8:28:43	FV4 aerate isolation valve fail
18/10/2008	8:28:43	FV29 chamber inlet valve 1 fail
18/10/2008	8:28:52	Compressor 1 failure
18/10/2008	8:29:05	Compressor 2 failure
18/10/2008	8:31:06	Compressor 1 cooling water low flowrate
18/10/2008	8:31:06	Compressor 2 cooling water low flowrate
18/10/2008	8:32:35	Vacuum failure pre-alarm
18/10/2008	13:25:34	Shelf temperature high deviation
18/10/2008	14:02:28	Circulation pump 1 failure
18/10/2008	14:02:32	FV25 filter isolation valve fail
18/10/2008	14:02:32	FV29 chamber inlet valve 1 fail
18/10/2008	14:03:44	Vacuum failure pre-alarm
18/10/2008	17:25:04	Circulation pump 1 failure
18/10/2008	17:25:04	FV29 chamber inlet valve 1 fail
18/10/2008	17:27:04	Vacuum failure pre-alarm
Events: Total: 1		
18/10/2008	17:45:42	Pressure rise test passed. The final chamber vacuum rise in this test is :9ubar



 18/10/08

 18/10/08

	Batch Filling Record (Part-A)		Page 39 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

9.3 STERILIZATION OF ALUMINUM FLIP-OFF OVER SEALS AS PER SOP No. SOP/CSP/037:

Cycle started at: 16:59

Cycle Completed at: 18:49

Date	Load No.	Load Pattern No.	No. of Bags	Total qty of aluminum flip off Over seals.	Sterilization Hold Time		Done by ^{##}	Ckd. By ^{##}
					From	To		
17/10/08	04	6	04	6295	17:42	18:12	Sohan	17/10/08

(Graph attached on the reverse side of page No.: 29)

9.4 STORAGE OF STERILIZED ALUMINUM FLIP-OFF OVER SEALS: (NMT 48 HOURS AFTER STERILIZATION):

Procedure	Eq. Id. No.	SOP No.	Date	Storage Time*		Done by ^{##}	Ckd. By ^{##}
				From	To		
Store the sterilized aluminum flip-off over seals under LAF in cooling zone.	CSP/035	SOP/CSP/072	17/10/08	18:50	---	Sohan	17/10/08
			18/10/08	---	16:20	Sanjeev	18/10/08

* Hold time starts after completion of sterilization and completes at the end of sealing.

9.5 TRANSFER OF STERILIZED ALUMINUM FLIP-OFF OVER SEALS:

Procedure	Eq. Id. No.	Date	Time		Done by ^{##}	Ckd. By ^{##}
			From	To		
Transfer the sterilized aluminum flip-off over seals from cool zone to vial washing area through the corridor in closed material transfer trolley. Load the sterilized aluminum flip-off over seals to the hopper of the sealing machine inside the unloading cRAB	CSP/030	18/10/08	16:20	16:32	Ranjit	18/10/08

^{##} Production

	Batch Filling Record (Part-A)		Page 40 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

10.0 CLEANING OF UNLOADING cRAB AND SEALING MACHINE:

10.1 Cleaning of sealing Unloading cRAB is done as per SOP No. SOP/CSP/111.

Cleaning started at: 17:00


Cleaning completed at: 17:30

Date: 18/10/08

S.No.	Equipment Name	Equipment Id. No.	Previous Product Name	Batch No.	Time		Done By ^{##}	Ckd. By ^{##}
					From	To		
1.	Unloading cRAB	CSP/030	R & D Trial Batch	-	17:00	17:15	Ranjit	
2.	Sealing machine	CSP/171	R & D Trial Batch	-	17:16	17:30	Ranjit	

10.2 UNLOADING AND SEALING OF LYOPHILIZED VIALS:

10.3 LINE CLEARANCE:

Step No.	Checks	Equipment ID	SOP No.	Date	Previous Product		Ckd. by ^{##}	Verified by ^{###}
					Product Name	Batch no.		
1.	Line Clearance of Sealing Unloading cRAB	CSP/030	SOP/QAD/053	18/10/08	R & D Trial Batch	-		
2.	Line Clearance of Sealing machine	CSP/171	SOP/QAD/053	18/10/08	R & D Trial Batch	-		


^{##} Production

^{###} IPQA

Form No.: SOP/QAD/038-02 Rev. 00

CONFIDENTIAL

QA Signature



LOAD No- 04


MFG.BY : MACHIN FABRIK
 CUSTOMER : M/S DABUR PHARMA LTD.
 EQUIPMENT : STEAM STERILIZER

BATCH DATA
 STERILIZER NO : 1007
 BATCH NO : 0
 QUANTITY : 1
 OPERATOR CODE : 1
 PRINT INTERVAL : 60 SEC
 DATE : 17/10/08
 TIME : 16:59:06

PROCESS : STOPPER & FLIP

SET PARAMETERS
 PRE VACUUM : -0.750 BAR
 PRE PRESSURE : 0.500 BAR
 NO OF PRE PULSE : 3 NOS
 SET POINT-1 : 110.0 DEG C
 SET POINT-2 : 115.0 DEG C
 SET POINT-3 : 118.0 DEG C
 SP-1 HOLD TIME : 2 MIN
 SP-2 HOLD TIME : 3 MIN
 SP-3 HOLD TIME : 5 MIN
 HEAT UP CONTROL BAND : 0.2 DEG C
 STER HOLD TEMP : 121.5 DEG C
 STER HOLD TIME : 30 MIN
 TEMP CONTROL BAND : 0.2 DEG C
 OVERSHOOT TEMP : 124.0 DEG C
 STER COUNT STOP TEMP : 120.9 DEG C
 STER COUNT RESET TEMP : 120.8 DEG C
 POST VACUUM : -0.750 BAR
 POST PRESSURE : -0.100 BAR
 POST VAC HOLD TIME : 30 MIN
 NO OF POST PULSE : 1 NOS
 EXHAUST ON TIME : 5 SEC
 EXHAUST OFF TIME : 50 SEC
 PROCESS END PRESSURE : -0.050 BAR
 CHAM PRESSURE HIGH : 1.500 BAR

MFG.BY : MACHIN FABRIK
 CUSTOMER : M/S DABUR PHARMA LTD.
 DATE : 17/10/08
 TIME : 16:59:23
 EQUIPMENT : STEAM STERILIZER
 PROCESS : STOPPER & FLIP OFF SEAL



* Only Flip off Seal
 has been taken

1/5

TIME	T1 DEG C	T2 DEG C	T3 DEG C	T4 DEG C	T5 DEG C	PRES BAR	STATUS
16:59:25	030.6	046.4	048.1	039.6	052.7	-0.020	PROCESS ON
16:59:27	030.6	046.3	048.3	039.7	052.9	-0.036	PRE VACUUM
17:00:03	030.3	046.8	049.7	039.2	053.7	-0.479	PRE VACUUM
17:00:53	029.9	046.1	049.8	038.3	053.5	-0.765	PRE PRESSURE
17:01:02	030.0	047.5	050.4	045.1	054.4	-0.643	PRE PRESSURE

17:02:01	030.2	070.6	068.9	070.7	070.3	-0.385	PRE PRESSURE
17:03:02	030.2	081.5	080.3	081.4	081.4	-0.192	PRE PRESSURE
17:04:03	030.5	088.6	087.0	086.1	088.8	-0.021	PRE PRESSURE
17:05:01	090.4	093.5	092.7	091.9	093.4	00.088	PRE PRESSURE
17:06:03	096.8	097.6	097.2	096.9	097.0	00.186	PRE PRESSURE
17:07:01	100.4	100.9	100.8	101.2	100.3	00.287	PRE PRESSURE
17:08:02	103.9	103.9	103.9	104.7	103.3	00.385	PRE PRESSURE
17:09:03	106.7	106.5	106.5	107.4	105.8	00.478	PRE PRESSURE
17:09:15	107.2	107.0	107.0	107.8	106.3	00.498	PRE VACUUM
17:10:02	092.2	092.5	093.2	093.4	092.1	-0.316	PRE VACUUM
17:11:03	073.8	073.9	074.5	074.7	074.6	-0.672	PRE VACUUM
17:11:33	065.4	066.3	067.3	067.1	071.9	-0.757	PRE PRESSURE
17:12:02	074.1	072.8	073.0	074.9	075.5	-0.550	PRE PRESSURE
17:13:03	090.5	089.0	089.0	089.5	088.9	-0.220	PRE PRESSURE
17:14:02	097.9	097.1	097.3	097.6	096.8	-0.003	PRE PRESSURE
17:15:02	102.6	101.6	101.8	102.1	101.3	00.143	PRE PRESSURE
17:16:03	106.1	105.0	105.3	105.6	104.9	00.280	PRE PRESSURE
17:17:04	109.0	108.0	108.3	108.6	107.8	00.408	PRE PRESSURE
17:17:49	110.8	109.7	110.0	110.5	109.7	00.496	PRE VACUUM
17:18:03	108.3	107.9	108.4	107.9	107.5	00.201	PRE VACUUM
17:19:04	086.6	088.1	088.7	088.7	089.7	-0.442	PRE VACUUM
17:20:04	068.0	071.3	073.5	069.8	085.9	-0.733	PRE VACUUM
17:20:14	065.6	070.6	073.5	067.6	085.6	-0.757	PRE PRESSURE
17:21:04	082.1	080.5	080.8	081.5	087.2	-0.409	PRE PRESSURE
17:22:03	095.4	093.5	093.7	094.2	093.8	-0.084	PRE PRESSURE
17:23:04	102.1	100.7	101.0	101.3	100.6	00.126	PRE PRESSURE
17:24:05	106.1	104.9	105.1	105.5	104.7	00.281	PRE PRESSURE
17:25:03	109.0	108.0	108.2	108.6	107.8	00.410	PRE PRESSURE
17:25:46	110.9	109.9	110.1	110.6	109.8	00.496	HEAT UP
17:26:04	110.5	109.7	110.0	110.3	109.6	00.437	HEAT UP
17:27:03	110.1	109.4	109.6	110.1	109.2	00.459	HEAT UP
17:28:04	110.4	109.6	109.9	110.3	109.6	00.482	HEAT UP
17:29:03	113.0	111.9	112.2	112.6	111.9	00.604	HEAT UP

215

17:30:05	115.0	114.0	114.3	114.8	114.0	00.715	HEAT UP
17:31:04	114.9	114.2	114.5	114.9	114.2	00.695	HEAT UP
17:32:05	114.9	114.2	114.5	115.0	114.2	00.697	HEAT UP
17:33:04	114.9	114.3	114.5	115.0	114.3	00.708	HEAT UP
17:34:05	117.2	116.1	116.5	116.9	116.2	00.827	HEAT UP
17:35:04	117.9	117.1	117.4	117.8	117.1	00.856	HEAT UP
17:36:05	117.9	117.2	117.5	117.9	117.2	00.872	HEAT UP
17:37:06	117.9	117.2	117.5	118.0	117.2	00.877	HEAT UP
17:38:04	118.0	117.3	117.6	118.0	117.3	00.881	HEAT UP
17:39:06	118.2	117.4	117.7	118.1	117.4	00.871	HEAT UP
17:40:05	118.9	117.9	118.1	118.7	117.9	00.933	HEAT UP
17:41:06	120.8	119.7	120.0	120.5	119.8	01.050	HEAT UP
17:42:05	122.1	121.1	121.3	121.8	121.1	01.125	HEAT UP
17:42:32	122.5	121.6	121.9	122.3	121.5	01.152	STER HOLD
17:43:06	122.3	121.5	121.9	122.2	121.5	01.140	STER HOLD
17:44:07	122.5	121.7	121.9	122.3	121.6	01.116	STER HOLD
17:45:06	122.5	121.7	121.9	122.5	121.6	01.162	STER HOLD
17:46:07	122.3	121.5	121.9	122.3	121.5	01.133	STER HOLD
17:47:06	122.6	121.8	122.0	122.5	121.7	01.132	STER HOLD
17:48:06	122.6	121.8	122.0	122.5	121.7	01.166	STER HOLD
17:49:08	122.4	121.6	121.9	122.3	121.5	01.148	STER HOLD
17:50:06	122.4	121.6	122.0	122.3	121.5	01.111	STER HOLD
17:51:06	122.6	121.8	122.0	122.5	121.7	01.149	STER HOLD
17:52:06	122.5	121.7	121.9	122.4	121.6	01.155	STER HOLD
17:53:07	122.3	121.6	121.9	122.3	121.5	01.133	STER HOLD
17:54:08	122.5	121.7	122.0	122.4	121.6	01.126	STER HOLD
17:55:06	122.6	121.8	122.0	122.5	121.7	01.162	STER HOLD
17:56:07	122.3	121.6	121.9	122.3	121.5	01.145	STER HOLD
17:57:08	122.5	121.7	122.0	122.3	121.6	01.121	STER HOLD
17:58:08	122.6	121.8	122.1	122.5	121.7	01.151	STER HOLD
17:59:07	122.4	121.7	121.9	122.3	121.6	01.148	STER HOLD
18:00:07	122.5	121.7	122.0	122.3	121.5	01.117	STER HOLD
18:01:08	122.6	121.8	122.0	122.4	121.6	01.131	STER HOLD

18:02:07	122.6	121.8	122.0	122.5	121.7	01.167	STER HOLD
18:03:08	122.4	121.7	121.9	122.3	121.6	01.152	STER HOLD
18:04:08	122.3	121.7	121.9	122.3	121.5	01.127	STER HOLD
18:05:07	122.6	121.8	122.0	122.4	121.7	01.134	STER HOLD
18:06:08	122.6	121.8	122.1	122.5	121.7	01.147	STER HOLD
18:07:07	122.4	121.7	122.0	122.4	121.6	01.157	STER HOLD
18:08:09	122.3	121.7	121.9	122.3	121.5	01.134	STER HOLD
18:09:08	122.6	121.8	122.0	122.5	121.6	01.129	STER HOLD
18:10:09	122.6	121.8	122.0	122.6	121.7	01.166	STER HOLD
18:11:08	122.3	121.6	121.9	122.4	121.5	01.142	STER HOLD
18:12:09	122.5	121.7	122.0	122.5	121.5	01.118	STER HOLD
18:12:32	122.7	121.8	122.1	122.7	121.7	01.155	POST VACUUM
18:13:08	108.0	110.2	111.2	111.7	109.5	00.173	POST VACUUM
18:14:09	085.6	087.3	090.0	090.8	095.9	-0.476	POST VACUUM
18:15:08	070.5	079.8	090.3	073.2	094.8	-0.735	POST VACUUM
18:15:13	070.8	080.1	090.4	072.3	094.7	-0.745	VACUUM HOLD
18:16:09	072.0	083.7	091.4	063.4	094.3	-0.834	VACUUM HOLD
18:17:08	068.7	086.3	092.7	055.3	094.4	-0.878	VACUUM HOLD
18:18:09	064.5	088.3	093.7	051.6	094.8	-0.905	VACUUM HOLD
18:19:11	060.1	090.0	094.9	050.4	095.5	-0.918	VACUUM HOLD
18:20:10	056.4	090.9	095.5	057.8	096.2	-0.924	VACUUM HOLD
18:21:10	053.1	091.8	096.1	063.0	096.7	-0.929	VACUUM HOLD
18:22:11	050.3	092.3	096.5	063.9	097.1	-0.932	VACUUM HOLD
18:23:09	048.1	092.7	096.9	067.4	097.5	-0.935	VACUUM HOLD
18:24:11	046.2	093.0	097.0	070.5	097.7	-0.935	VACUUM HOLD
18:25:11	044.5	093.0	097.1	071.8	097.8	-0.937	VACUUM HOLD
18:26:11	043.0	093.1	097.2	073.3	098.0	-0.938	VACUUM HOLD
18:27:10	041.6	093.2	097.2	074.5	098.1	-0.939	VACUUM HOLD
18:28:11	040.3	093.3	097.2	075.2	098.2	-0.940	VACUUM HOLD
18:29:10	039.2	093.3	097.2	075.6	098.3	-0.941	VACUUM HOLD
18:30:11	038.1	093.4	097.1	076.1	098.4	-0.942	VACUUM HOLD
18:31:12	037.1	093.4	097.0	076.4	098.4	-0.942	VACUUM HOLD
18:32:11	036.2	093.4	097.0	076.6	098.4	-0.942	VACUUM HOLD

18:33:12	093.3	093.5	096.9	076.8	098.4	-0.944	VACUUM HOLD
18:34:12	094.5	093.5	096.8	077.0	098.5	-0.944	VACUUM HOLD
18:35:12	093.8	093.4	096.6	077.2	098.4	-0.943	VACUUM HOLD
18:36:11	093.1	093.4	096.6	077.4	098.4	-0.944	VACUUM HOLD
18:37:11	092.5	093.2	096.4	077.5	098.4	-0.944	VACUUM HOLD
18:38:12	091.9	093.2	096.4	077.7	098.4	-0.944	VACUUM HOLD
18:39:13	091.3	093.1	096.4	077.7	098.4	-0.945	VACUUM HOLD
18:41:13	090.4	092.9	096.3	077.8	098.3	-0.942	VACUUM HOLD
18:42:12	090.1	092.8	096.2	077.7	098.3	-0.942	VACUUM HOLD
18:43:13	029.7	092.6	096.0	077.5	098.2	-0.945	VACUUM HOLD
18:44:14	029.5	092.5	095.9	077.4	098.2	-0.946	VACUUM HOLD
18:45:12	029.2	092.3	095.8	077.4	098.1	-0.946	VACUUM HOLD
18:45:15	029.2	092.3	095.7	077.4	098.1	-0.946	POST PRESSURE
18:46:13	029.3	092.5	094.6	077.7	098.9	-0.648	POST PRESSURE
18:47:14	029.8	093.0	094.6	075.3	099.7	-0.351	POST PRESSURE
18:48:13	030.2	091.3	095.2	073.1	099.2	-0.140	POST PRESSURE
18:48:29	030.2	090.9	095.4	072.7	099.0	-0.101	VACUUM BREAK
18:49:05	030.4	090.4	095.5	072.0	098.3	-0.050	PROCESS END

S/S

FD VALUE:- 0030.0

		Batch Filling Record (Part-A)		Page 41 of 46
PLANT-II, BADDI				
Product Name	Topotecan for Injection 4 mg (base)/vial			
Batch Number	8 PK 201			
Batch size	12 L	Batch Filling Code	IO720004	
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001	
Shelf life	24 Months (Provisional)	Version No.	00	

10.4 UNLOADING AND SEALING OF LYOPHILIZED VIALS:

Sealing started at: 19:45.

Sealing completed at: 21:05.

Step No.	Checks	SOP No.	Time		Observations	Done By##	Ckd. By##
			From	To			
1	Note the reading of magnehelic gauge over unloading cRAB.	SOP/CSP/111	19:13	19:14	ΔP across HEPA 5 (mm of WC)	Ranjit	18/10/08
2	Check the in process sealing quality of sealed vials from each sealing head (2 vials) at Start and End of sealing (during process every 30 min.) and record the details in table 10.5.	SOP/CSP/111	19:45	21:05	Vial No. 1. OK / NOT OK 2. OK / NOT OK	Ranjit	18/10/08

10.5 IN PROCESS SEALING QUALITY CHECK

Frequency: Start and End of filling (during process every 30 min.)

Date: 18/10/08

Checks		Time							
		19:45	20:15	20:45					
Sealing quality	*Crimping	OK	OK	OK	N/A				
	Loose seals	NIL	NIL	NIL					
Ckd. By##		18/10/08	18/10/08	18/10/08					

Checks		Time							
Sealing quality	*Crimping	N/A							
	Loose seals								
Ckd. By##									

* To be checked for Every Head

Production

	Batch Filling Record (Part-A)		Page 42 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	10720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00


11.0 EXTERNAL DECONTAMINATION OF FILLED AND SEALED VIALS:

11.1 CLEANING:

Step No.	Checks	Eq. Id. No.	SOP No.	Date	Time		Done By##	Ckd. By##
					From	To		
1.	Cleaning of External Decontaminating machine	CSP/023	SOP/CSP/033	18/10/08	17:40	18:05	Rangect	 18/10/08

11.2 LINE CLEARANCE:

- Ensure the cleanliness of Vial washing area.
- Ensure no remains of previous product are there.
- Ensure that status label is properly attached with the machine.
- Ensure proper change over and setting is done as per the pack size.

Step No.	Checks	SOP No.	Date	Previous product		Ckd. by##	Verified by###
				Product Name	Batch No.		
1.	Line clearance of External Decontaminating machine	SOP/QAD/053	18/10/08	R & D Trial Batch	—	 18/10/08	 18/10/08

11.3 EXTERNAL DECONTAMINATING MACHINE OPERATION:

Date: 18/10/08

External decontamination started at: 19:53

External decontamination completed at: 21:13

Step No.	Checks	SOP No.	Time		Observations	Done By##	Ckd. By##
			From	To			
1.	Pressure of Purified Water (NLT 0.5 bar)	SOP/CSP/033	19:30	19:32	0.6 bar	Rangect	
2.	Compressed Air Pressure (NLT 3.0 bar)	SOP/CSP/033	19:33	19:34	3.0 bar	Rangect	
3.	No. of vials washed	SOP/CSP/033	19:53	21:13	Qty. 4866 vials	Rangect	
4.	No. of vials broken	---	---	---	Qty. — vials	Rangect	
5.	No. of vials rejected	---	---	---	Qty. — vials	Rangect	
6.	No. of good vials	---	---	---	Qty. 4840 vials	Rangect	

Production

20/10/08 ### IPQA



	Batch Filling Record (Part-A)		Page 43 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

12.0 CODING OF SEALED VIALS ON ALUMINUM FLIP-OFF OVER SEALS FOR BATCH IDENTIFICATION:

12.1 CLEANING:

Step No.	Checks	Eq. Id. No.	SOP No.	Date	Time		Done By ^{##}	Ckd. By ^{##}
					From	To		
1.	Cleaning of conveyor belt	CSP/164	SOP/CSP/058	18/10/08	18:10	18:12	Ranjit	 18/10/08
2.	Cleaning of Ink-jet Printer	CSP/054	SOP/CSP/057	18/10/08	18:20	18:30	Ranjit	 18/10/08


12.2 LINE CLEARANCE:

- Ensure the cleanliness of Vial washing area.
- Ensure no remains of previous product are there.
- Ensure that status label is properly attached with the machine.
- Ensure proper instructions are feeded in the printer and previous details have been omitted.

Date: 18/10/08

Step No.	Checks	Eq. Id. No.	SOP No.	Previous product		Ckd. by ^{##}	Verified by ^{###}
				Product Name	Batch No.		
1.	Line clearance for conveyor belt.	CSP/164	SOP/QAD/053	Gemcitabine Inj	8KA205	 18/10/08	 18/10/08
2.	Line clearance for Ink-jet Printer.	CSP/054	SOP/QAD/053	Cytarabine Inj	8LH206	 18/10/08	 18/10/08

12.3 BATCH OVER PRINTING INSTRUCTIONS:

Batch No	Ckd. by ^{##}	Verified by ^{###}
8 PK 201	 18/10/08	 18/10/08

Date: 18/10/08

Over-Printing started at: 19:57



Over-Printing completed at: 21:22

^{##} Production

^{###} IPQA

	Batch Filling Record (Part-A)		Page 44 of 46
	PLANT-II, BADDI		
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

12.4 INITIAL VIALS CODING DETAILS:



Date	Time	Over-Printing Details	Ckd by ^{##}	Verified by ^{###}
18/10/08	19:57	8PK201		

12.5 INPROCESS CHECKS OF OVER PRINTING:

Frequency: Every 15 min. during operation.



Date	Time	Batch No.	Observation	Done By ^{##}	Ckd. By ^{##}
18/10/08	19:57	8PK201	OK	Rangset	
18/10/08	20:12	8PK201	OK	Rangset	
18/10/08	20:27	8PK201	OK	Rangset	
18/10/08	20:42	8PK201	OK	Rangset	
18/10/08	20:57	8PK201	OK	Rangset	
18/10/08	21:12	8PK201	OK	Rangset	

13.0 QUALITY CONTROL SAMPLING:

Date	Quantity of sample		Stage	Sampled by ^{###}	Remarks
	No. Of Vials	Solution (mL)			
15/10/08	3+3+3=9 vials	170 ml + 40 ml ^{***}	After 2nd filtration		-
18/10/08	53+62+83=198 vials	-	After lyophilization		-

^{##} Production

^{###} IPQA

^{**} 40 ml (25 ml + 15 ml) taken after end of filling: - 
^{*} 9 vial taken during filling 
^{***} 83 (80+3) vial taken for stability and BET validation 

	Batch Filling Record (Part-A)		Page 45 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

14.0 TRANSFER OF LYOPHILIZED AND SEALED VIALS FOR QUARANTINE:

Date	Qty. transferred for quarantine	Handed over by	Received by
20/10/08	4668 vials		Sub. hdu
—	—	—	—

15.0 YIELD AFTER FILLING AND SEALING STAGE:

- Effective Batch size (Page 2) 10.297 L
- Actual batch size step 15 (Page 28) 10.297 L
- Target Fill Volume 1.957 mL
- Quality control Bulk sample for analysis (Step 13) 210 mL
- Sample taken for fill volume check 46.968 mL
- Total No. of good vials after external washing (Step 6 of 11.3 page 42) 4866 Nos.
- Yield at filling stages after accounting for quality control samples:


% Yield with respect to Actual batch size =

$$[(d + e) + (f \times c) / (b \times 1000)] \times 100 = \underline{94.97} \%$$

(Limit: Tentative yield should be in between 90 - 100% and standard yield to be finalized after observing the trend of 10 batches. The yield less than 90 % and more than 100% should be investigated as per SOP No. SOP/QAD/030)

Variance with explanation: _____

Calculated by: 

Ckd. by: 
(Production Head)

Date: 20/10/08

Date: 20/10/2008

	Batch Filling Record (Part-A)		Page 46 of 46
PLANT-II, BADDI			
Product Name	Topotecan for Injection 4 mg (base)/vial		
Batch Number	8PK201		
Batch size	12 L	Batch Filling Code	IO720004
Department	Production (CSL)	BFR No.	BFR/CSL/TOP/001
Shelf life	24 Months (Provisional)	Version No.	00

16.0 ABBREVIATIONS:

%	Percent
°C	Degree Celsius
A.R. No.	Analytical Reference Number
Approx.	Approximately
B.No.	Batch Number
NMT	Not More Than
Ckd. by	Checked By
CSL	Cyto Sterile Liquid
Eq.Id.No.	Equipment Identification number
Exp.	Expiry
HDPE	High Density Polyethylene
IPQA	In process quality assurance
Issd.	Issued
TC	Triclover
L	Liters
LAF	Laminar Air Flow
Max.	Maximum
mbar	millibar
MFG./Mfg.	Manufacturing
Min./mins	Minutes
mL	Milliliter
mm	Millimeter
psi	Pounds per square inch
CIP	Cleaning in place
IP	Indian Pharmacopoeia
Kg	Kilogram
IPA	Isopropyl alcohol
HEPA	High efficiency particulate air

NLT	Not Less Than
No.	Number
QC	Quality Control
q.s.	Quantity Sufficient
Qty.	Quantity
Q.A.	Quality assurance
R.H.	Relative Humidity
CSP	Cyto Sterile Production
Mpa	Mega Pascal
QAD	Quality Assurance department
BFR	Batch Filling Record
g/gm	Grams
SOP	Standard Operating Procedure
S. No.	Serial Number
WFI	Water For Injection
Std. Qty.	Standard Quantity
Temp.	Temperature
Dept.	Department
USP	United States Pharmacopoeia
μ	Micron
mg	Milligram
cRAB	Closed Restricted Access Barrier
MOC	Material of Construction
SIP	Sterilization in place
Ph.Eur.	Pharma Europa
NA	Not applicable
ΔP	Differential pressure
--	---

[Signature]

01-06ch	27.4	27.5	27.4	27.6	27.1
07-12ch	27.5	27.3	27.4	27.1	24.9
13-18ch	28.3	-26.3	24.8	21.5	27.0
19-24ch	26.7	117	1.83	1.83	---

10:00 2003/10/28 1802

Z1

01-06ch	26.6	26.1	26.1	23.2	26.7
07-12ch	26.4	23.1	27.3	26.6	25.6
13-18ch	28.4	24.8	24.7	22.3	27.8
19-24ch	27.7	115	1.83	1.81	---

15:00 60mm/h CH21

20

01-06ch	36.8	36.2	36.8	37.6	35.8
07-12ch	37.3	37.3	36.9	35.8	46.2
13-18ch	46.6	-55.9	24.6	27.0	28.3
19-24ch	28.1	112	8.80	8.85	---

15:00 60mm/h CH21

1

01-06ch	38.1	38.3	37.8	38.5	36.2
07-12ch	38.5	38.4	38.2	37.3	58.8
13-18ch	51.1	-81.7	24.4	22.9	28.2
19-24ch	28.0	-53	8.86	8.85	---

15:00 60mm/h CH19

18

01-06ch	37.3	37.5	37.8	37.7	35.5
07-12ch	37.7	37.4	37.4	37.6	49.4
13-18ch	49.6	-82.1	22.6	22.6	27.8
19-24ch	27.6	-35	8.86	8.85	---

15:00 2003/10/28 CH18

17

00 08 09 0 02 0

01-06ch	36.6	36.1	36.8	36.4	37.8	34.9
07-12ch	37.0	36.7	36.9	36.9	36.9	47.8
13-18ch	48.8	-79.0	-82.1	23.2	21.6	27.7
19-24ch	27.4	41	-46	0.86	0.05	---

100.0
11:00 60mm/h CH17 150.0 °C

01-06ch	38.2	38.5	38.4	37.9	38.6	36.2
07-12ch	38.7	38.5	38.6	38.4	37.4	51.9
13-18ch	52.1	-77.1	-81.9	23.7	21.7	27.7
19-24ch	27.5	52	141	0.86	0.05	---

100.0
12:00 60mm/h CH18 150.0 °C

01-06ch	37.4	37.6	37.6	37.8	35.5
07-12ch	37.9	37.8	37.6	36.7	50.7
13-18ch	51.0	-82.5	23.2	21.7	27.6
19-24ch	27.3	-51	0.86	0.05	---

100.0
13:00 60mm/h CH19 150.0 °C

01-06ch	36.6	36.1	36.8	36.4	37.0	34.8
07-12ch	37.1	37.0	36.8	36.9	49.5	
13-18ch	49.7	-82.9	23.1	21.8	27.5	
19-24ch	27.2	-43	0.86	0.05	---	

100.0
14:00 60mm/h CH20 150.0 °C

01-06ch	35.8	36.0	35.6	37.2	34.1
07-12ch	36.3	36.2	36.3	37.2	48.4
13-18ch	48.6	-82.9	23.0	21.7	27.4
19-24ch	27.1	-51	0.86	0.05	---

100.0
15:00 60mm/h CH21 150.0 °C

12

001 09 09 09 09 02 0

01-06ch	35.0	34.0	35.2	34.4	33.4
07-12ch	35.4	35.1	35.3	34.4	46.8
13-18ch	47.1	-82.7	-82.4	22.7	27.3
19-24ch	27.1	34	-57	0.05	---

100.0 150.0
10:00 2008/10/18 CH12 °C

11

01-06ch	34.8	34.9	34.7	33.1	33.1
07-12ch	35.2	35.1	34.9	34.1	46.6
13-18ch	46.8	-81.7	-82.3	22.6	27.4
19-24ch	27.1	-58	0.05	---	---

100.0 150.0
09:00 60mm/h CH11 °C

10

01-06ch	34.4	34.6	34.1	32.8	32.8
07-12ch	34.8	34.7	34.8	33.7	46.5
13-18ch	46.7	-82.9	22.3	22.7	27.4
19-24ch	27.1	-61	0.05	---	---

100.0 150.0
08:00 2008/10/18 CH10 °C

09

01-06ch	33.8	33.6	34.0	32.6	32.3
07-12ch	34.3	34.1	34.0	33.2	44.8
13-18ch	46.7	-76.6	-81.9	22.4	27.6
19-24ch	27.3	44	-58	0.05	---

100.0 150.0
07:00 60mm/h CH 9 °C

08

01-06ch	33.1	33.3	33.0	31.6	31.6
07-12ch	33.5	33.2	33.5	32.3	46.6
13-18ch	46.8	-83.8	-82.2	22.8	27.8
19-24ch	27.5	32	-52	0.05	---

100.0 150.0
06:00 2008/10/18 CH 8 °C

07

001 00 09 07 02 0

01-06ch	27.2	27.2	27.4	27.6	26.4
07-12ch	27.4	27.5	27.3	26.7	40.2
13-18ch	39.4	-83.0	22.0	22.6	27.8
19-24ch	27.4	-50	0.05	0.05	---

100.0
05:00 60mm/h CH 7 150.0 °C

06

01-06ch	26.1	26.1	26.0	26.2	25.6
07-12ch	26.2	26.1	26.2	25.8	31.0
13-18ch	31.0	-83.0	22.7	22.6	27.8
19-24ch	27.4	-37	0.05	0.05	---

100.0
04:00 2000/10/18 CH 6 150.0 °C

05

01-06ch	25.9	26.0	25.8	26.0	25.4
07-12ch	26.1	26.0	26.3	25.6	31.1
13-18ch	31.3	-82.1	22.7	22.6	27.7
19-24ch	27.4	-50	0.05	0.05	---

100.0
03:00 60mm/h CH 5 150.0 °C

04

01-06ch	25.7	25.5	25.9	25.7	25.3
07-12ch	26.0	25.8	25.8	25.5	31.3
13-18ch	31.2	-80.7	-82.7	22.5	22.6
19-24ch	27.5	-54	0.05	0.05	---

100.0
02:00 2000/10/18 CH 4 150.0 °C

03

01-06ch	25.5	25.7	25.6	25.4	25.0
07-12ch	25.7	25.6	25.6	25.2	31.0
13-18ch	31.0	-81.1	-81.1	22.1	22.6
19-24ch	27.4	-48	0.05	0.05	---

100.0
01:00 60mm/h CH 3 150.0 °C

PART NO. 840701-001

MEMORANDUM

22:00 2008/10/17

CH23

01-06ch	21.8	21.8	22.1	21.6
07-12ch	21.9	22.8	21.5	38.9
13-18ch	30.8	-82.1	22.7	27.9
19-24ch	27.4	-52	0.05	---

01-06ch	17.3	17.2	17.2	17.9
07-12ch	17.2	17.1	17.4	15.8
13-18ch	15.6	-81.7	22.8	27.9
19-24ch	27.6	-17	0.05	---

20:00 2008/10/17

CH21

01-06ch	16.9	16.9	16.9	17.6
07-12ch	16.9	16.8	17.1	15.8
13-18ch	15.5	-82.7	23.0	27.8
19-24ch	27.4	-1	0.05	---

01-06ch	16.1	16.1	16.1	16.8
07-12ch	16.1	16.0	16.3	15.7
13-18ch	15.5	-81.7	23.0	27.8
19-24ch	27.5	-1	0.05	---

18:00 2008/10/17

CH19

01-06ch	13.2	13.2	13.1	14.5
07-12ch	13.1	13.0	13.6	7.6
13-18ch	7.2	-81.0	22.8	27.7
19-24ch	27.5	-13	0.05	---

100.0

150.0

17:00 6000/h

CH18

01-06ch	13.0	13.0	12.8	14.2
07-12ch	12.9	12.8	13.4	7.5
13-18ch	7.1	-81.1	22.8	27.5
19-24ch	27.2	-3	0.85	---

139.8

14:00 2000/10/17

CH17

159.8

01-06ch	12.5	12.5	12.4	13.8
07-12ch	12.3	12.3	12.9	7.5
13-18ch	7.1	-81.5	22.9	27.5
19-24ch	27.1	-13	0.85	---

139.8

14:00 2000/10/17

CH16

159.8

15

01-06ch	10.8	10.8	9.8	11.8
07-12ch	9.8	9.7	10.7	0.4
13-18ch	0.1	-81.5	22.9	27.5
19-24ch	27.2	34	0.85	---

139.8

14:00 2000/10/17

CH15

159.8

01-06ch	9.8	9.9	9.8	11.5
07-12ch	9.6	9.9	10.5	0.3
13-18ch	0.0	-81.7	23.8	27.5
19-24ch	27.2	19	0.85	---

139.8

13:00 6000/h

CH14

159.8

01-06ch	9.6	9.6	9.4	11.3
07-12ch	9.4	9.2	10.5	0.2
13-18ch	0.0	-81.2	23.1	27.6
19-24ch	27.2	24	0.85	---

139.8

159.8

10:20:00 2008/10/17						CH13	°C
12							
01-06ch	9.5		9.4		9.2	11.2	
07-12ch	9.5		9.1		10.0	0.3	
13-18ch	0.0		-80.2		23.4	27.6	
19-24ch	27.2		26		0.05	---	
-100.0							150.0
11:00 60mm/h						CH12	°C
11							
01-06ch	9.3		9.3		9.1	11.0	
07-12ch	9.1		9.0		10.0	0.4	
13-18ch	0.0		-81.3		22.8	27.3	
19-24ch	27.0		36		0.05	---	
-100.0							150.0
12:00 2008/10/17						CH11	°C
0							
01-06ch	9.1		9.1		8.9	10.8	
07-12ch	8.9		8.8		9.8	0.4	
13-18ch	0.1		-81.0		22.5	27.0	
19-24ch	26.7		36		0.05	---	
-100.0							150.0
14:00 50mm/h						CH10	°C
09							
01-06ch	9.0		9.0		8.8	10.7	
07-12ch	8.8		8.7		9.7	0.3	
13-18ch	0.1		-81.0		22.2	26.9	
19-24ch	26.7		22		0.05	---	
-100.0							150.0
16:00 2008/10/17						CH 9	°C
0							
01-06ch	9.0		8.9		8.7	10.6	
07-12ch	8.8		8.6		9.6	0.2	
13-18ch	-0.2		-81.4		22.1	27.0	
19-24ch	26.8		10		0.05	---	
-100.0							150.0

07

01-06ch	8.9		8.3		8.6	10.5
07-12ch	8.6		8.5		9.5	0.3
13-18ch	9.1		-81.2		22.1	27.1
19-24ch	26.9		23		0.05	---

-100.0

06:00 2000/10/17

CH 7 150.0
°C

06

01-06ch	8.7		8.6		8.5	10.4
07-12ch	8.5		8.3		9.3	0.3
13-18ch	8.0		-80.4		22.1	27.1
19-24ch	26.9		12		0.05	---

-100.0

05:00 60mm/h

CH 6 150.0
°C

05

01-06ch	8.3	8.4	8.2	8.5	8.1	10.0
07-12ch	8.1	8.4	8.0	8.3	8.9	0.1
13-18ch	8.0	-82.4	-80.5	8.8	22.1	27.1
19-24ch	26.9	11	21	0.6	0.05	---

-100.0

04:00 2000/10/17

CH 5 150.0
°C

04

01-06ch	6.5	6.4	6.4	6.8	6.2	8.6
07-12ch	6.2	6.5	6.1	6.5	7.5	-6.0
13-18ch	-6.3	-82.2	-80.8	8.0	22.2	27.2
19-24ch	26.9	13	131	0.6	0.05	---

-100.0

03:00 60mm/h

CH 4 150.0
°C

03

01-06ch	6.4		6.3	6.7	6.1	8.5
07-12ch	6.1		6.8	6.4	7.4	-5.7
13-18ch	-6.1		-82.4	8.1	22.4	27.3
19-24ch	26.9		12	0.1	0.05	---

-100.0

02:00 2000/10/17

CH 3 150.0
°C

2008-10-17

01-06ch	6.3	6.4	5.8	8.4
07-12ch	5.8	5.3	7.3	-5.9
13-18ch	-6.1	-82.3	22.5	27.3
19-24ch	26.9	12.3	0.05	---

1

01-06ch	6.1	6.0	5.8	8.2
07-12ch	5.8	5.3	7.1	-6.1
13-18ch	-6.2	-82.5	22.6	27.2
19-24ch	26.8	11.3	0.05	---

-180.0

2008-10-17

CH 1 150.0

23

01-06ch	5.4	5.7	5.4	5.2	7.7	
07-12ch	5.2	5.5	5.8	5.4	6.4	-5.7
13-18ch	-6.1	-80.3	22.8	22.6	27.3	
19-24ch	26.8	100	0.05	0.05	---	

0.00

23:00 60mm/h

CH23

23

01-06ch	4.7	4.4	5.4	4.5	7.1
07-12ch	4.5	4.4	4.8	5.7	-5.8
13-18ch	-6.1	-80.3	22.8	22.8	27.3
19-24ch	26.9	52	0.05	0.05	---

0.00

2008-10-14

CH22

21

01-06ch	4.2	4.3	5.6	4.1	6.7
07-12ch	4.8	4.8	4.7	5.2	-5.9
13-18ch	-6.3	-76.3	22.6	22.8	27.4
19-24ch	26.9	54	0.05	0.05	---

000

09

09

07

02

10000

01-06ch 3.9 0.9 3.8 6.3
 07-12ch 3.6 2.9 4.9 -5.9
 13-18ch -6.5 -25.3 12.8 27.4
 19-24ch 26.9 33 0.05 ---

19

01-06ch 3.6 0.9 3.5 6.0
 07-12ch 3.3 2.3 4.6 -5.8
 13-18ch -6.1 -25.3 23.8 27.3
 19-24ch 26.8 33 0.05 ---

-100.0

150.0

19:00 60mm/h

CH19 °C

18

01-06ch 3.2 3.3 5.5 3.1 5.7
 07-12ch 2.9 2.9 4.2 -6.0
 13-18ch -6.4 -79.9 23.8 27.5
 19-24ch 27.8 67 21.05 0.05 ---

-100.0

150.0

18:00 2000/10/16

CH18 °C

07

01-06ch 2.6 2.6 3.8 2.5 5.1
 07-12ch 2.4 2.4 2.4 3.7 -5.8
 13-18ch -6.2 -81.3 23.3 27.4
 19-24ch 26.8 56 0.05 0.05 ---

-100.0

150.0

17:00 60mm/h

CH17 °C

01-06ch 0.8 0.8 0.8 0.8 3.5
 07-12ch 0.6 0.6 0.6 1.8 -6.2
 13-18ch -6.8 -88.6 23.3 27.3
 19-24ch 26.8 65 0.05 0.05 ---

-100.0

150.0

01-06ch	-0.2	-0.3	-0.4	2.5
07-12ch	-0.6	-0.6	0.9	-12.0
13-18ch	-12.4	-00.5	23.2	27.1
19-24ch	26.6	0.2	0.06	---

133.0 150.0
13:00 60mm/h CH15

01-06ch	-0.7	-0.0	-0.9	2.0
07-12ch	-1.1	-1.1	2.4	-11.8
13-18ch	-12.2	-00.0	22.0	27.0
19-24ch	26.5	0.0	0.05	---

01-06ch	-1.3	-1.4	-0.9	-1.5	1.4
07-12ch	-1.7	-1.6	-0.3	-0.1	-11.9
13-18ch	-12.3	-00.1	22.0	22.2	27.0
19-24ch	26.5	0.0	0.05	0.06	---

100.0 150.0
13:00 60mm/h CH13 °C

01-06ch	-2.9	-3.0	-2.6	-3.0	-0.1
07-12ch	-3.3	-3.2	-2.0	-1.8	-12.0
13-18ch	-12.5	-00.6	22.0	21.8	26.9
19-24ch	26.4	0.0	0.05	0.06	---

100.0 150.0
12:00 2000/10/16 CH12 °C

01-06ch	-5.0	-5.1	-5.0	-5.3	-1.8
07-12ch	-5.5	-5.5	-5.0	-3.6	-20.1
13-18ch	-20.5	-00.9	23.0	21.8	26.6
19-24ch	26.2	0.0	0.05	0.06	---

00:00 0 20 150.0

00100 60mm/h

CH 1

01-06ch	-8.0	-8.0	-8.0	-8.4	-4.4
07-12ch	-8.5	-8.5	-8.5	-8.4	-28.2
13-18ch	-28.6	-28.6	-28.6	22.1	26.8
19-24ch	25.8	25.8	25.8	0.06	---

100.0

00100 60mm/h

CH 3

01-06ch	-10.0	-10.0	-10.0	-10.3	-6.4
07-12ch	-10.6	-10.6	-10.6	-8.5	-28.3
13-18ch	-28.8	-28.8	-28.8	22.3	25.5
19-24ch	25.4	25.4	25.4	0.06	---

100.0

00100 60mm/h

CH 9

150.0

150.0

01-06ch	-16.9	-16.9	-16.9	-17.3	-12.2
07-12ch	-17.5	-17.5	-17.5	-14.6	-28.9
13-18ch	-31.0	-31.0	-31.0	23.0	25.4
19-24ch	25.3	25.3	25.3	0.05	---

100.0

00100 60mm/h

CH 9

150.0

150.0

01-06ch	-18.2	-18.2	-18.2	-20.4	-13.3
07-12ch	-18.5	-18.5	-18.5	-16.4	-40.2
13-18ch	-40.9	-40.9	-40.9	24.8	25.6
19-24ch	25.5	25.5	25.5	0.05	---

100.0

00100 60mm/h

CH 7

150.0

150.0

86

01-06ch	-16.8	-16.8	-16.8	-19.0	-12.8
07-12ch	-17.2	-17.2	-17.2	-15.1	-40.7
13-18ch	-40.5	-40.5	-40.5	22.1	25.9
19-24ch	25.9	25.9	25.9	0.06	---

100.0

00100 60mm/h

CH 7

150.0

150.0

06:00 2008/10/16

CH 6 °C

01-06ch	-14.9	-17.0	-18.0
07-12ch	-13.3	-31.1	-41.1
13-18ch	-40.5	25.3	26.1
19-24ch	25.1	0.88	---
00+100.0			
03:00 50mm/h			

CH 5 °C

01-06ch	-9.4	-11.9	-4.8
07-12ch	-9.9	-7.7	-36.5
13-18ch	-35.7	25.6	26.3
19-24ch	26.4	0.86	---

01-06ch	-3.5	-5.1	0.5
07-12ch	-3.8	-2.1	-26.0
13-18ch	-25.4	25.8	26.4
19-24ch	26.6	0.87	---
00+100.0			
03:00 50mm/h			

CH 3 °C

01-06ch	2.4	1.3	6.0
07-12ch	2.2	3.8	-14.8
13-18ch	-15.5	25.8	26.5
19-24ch	26.7	0.87	---

CH 2 °C

01-06ch	11.4	9.7	14.2
07-12ch	11.1	12.6	-9.9
13-18ch	-8.2	25.5	26.4
19-24ch	26.6	0.88	---

00+100.0

150.0



DISPENSED MATERIAL

		Date : 14/10/08
Name of Material	:	Aluminium Flip off Seal
A.R. No.	:	02PM08000062 Purple 20mm
Product Name	:	Topotecan for Inj 4mg
Batch No.	:	8PK201
Batch Size	:	12 L
Container No.	:	1 of 1
Gross weight	:	-
Tare weight	:	-
Net weight / Qty.	:	6300 Nbs
Dispensed By (Warehouse)	:	14/10/08
Checked & Verified By (Production)	:	14/10/08
Form No. : SOP/QAD/033-01 Rev.01		



DISPENSED MATERIAL

		Date : 14/10/08
Name of Material	:	Clonidine Tablets P-1
A.R. No.	:	02PM08000042
Product Name	:	Topotecan for Inj 4mg
Batch No.	:	8PK201
Batch Size	:	12 L
Container No.	:	1 of 1
Gross weight	:	-
Tare weight	:	-
Net weight / Qty.	:	6696 Nbs
Dispensed By (Warehouse)	:	14/10/08
Checked & Verified By (Production)	:	14/10/08
Form No. : SOP/QAD/033-01 Rev.01		



DISPENSED MATERIAL

		Date : 14/10/08
Name of Material	:	R-Plug. Pluritec 20mm
A.R. No.	:	02PM08000038
Product Name	:	Topotecan for Inj 4mg
Batch No.	:	8PK201
Batch Size	:	12 L
Container No.	:	1 of 1
Gross weight	:	-
Tare weight	:	-
Net weight / Qty.	:	7000 Nbs
Dispensed By (Warehouse)	:	14/10/08
Checked & Verified By (Production)	:	14/10/08
Form No. : SOP/QAD/033-01 Rev.01		