




 Safetab Life Science BATCH PACKING RECORD		Page 1 of 109			
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

BATCH PACKING RECORD COVER SHEET:

PRODUCT NAME	LITACOLD FLU TABLETS		
LABEL CLAIM	Each Uncoated tablet contains: Chlorphenamine Maleate BP 2mg Phenylephrine Hydrochloride BP 5mg Caffeine Anhydrous BP 30mg Paracetamol BP 500mg Colour: Approved colour used		
PACK SIZE	TP: 1X4's strip pack		
STORAGE CONDITION	Store in cool and dry place Protect from light and moisture.		
SHELF LIFE	36 MONTHS		
CUSTOMER NAME	GENERIC HEALTHCARE PVT.LTD.		
EFFECTIVE BATCH NUMBER	G17240415	EFFECTIVE SIGN AND DATE:	 11/05/24
MANUFACTURING LICENCE No.	21 14 4378		
MFG. DATE:	01/2025	BPR. No	STBPR:434/02
EXP. DATE:	12/2027	Supersedes BPR.No	STBPR:434/01
BPR ISSUED BY/DATE (QA)	K. King 01/01/25		
BPR CHECKED & RECEIVED BY/DATE (PRODUCTION)	B. Sai 01/01/25		
Allotted By: QA Executive Sign/Date:  01/01/25			
DATE OF COMMENCEMENT:	01/01/25	DATE OF COMPLETION:	17/02/25
This BPR contains 109 Pages + Annex - 10 Additional Pages			
Checked by Competent Staff Sign/Date:  17/02/25			

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature		 11/05/24	 11/05/24

UNCONTROLLED COPY

CONTROLLED COPY

Safetab Life Science
BATCH PACKING RECORD

Page 2 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	078	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

1.0 LIST OF EQUIPMENTS:

S.No	Name of Equipments	Equipment ID No. to be used
1	Balance -31 Kg Capacity	ST/WB/196
2	Balance -15 Kg Capacity	ST/WB/184
3	Balance - 600 g Capacity	ST/WB/192
4	Balance - 300g Capacity	ST/WB/191
5	Leak test apparatus	ST/IPL/003-02 ,ST/QALT/001
6	Strip packing machine	ST/PPST/001orST/PRSR/001 or ST/PRSR/002
7	Carton coding machine	ST/PPCC/002

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature				
UNCONTROLLED COPY				

Safetab Life Science BATCH PACKING RECORD

Page 3 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	078	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

2.0 DISPENSING:

2.1 LINE CLEARANCE CHECKS:

S.No	CHECKS	Observation store	IPQA
1	Ensure that all empty containers, labels, used Polybags, hand gloves of previous dispensed Product/ batch are removed from the area.	✓	✓
2	Ensure that all dispensed material of previous product / batch removed from the area.	✓	✓
3	Ensure that the previous batch BMR / documents are removed from the area.	✓	✓
4	Check the housekeeping general cleanliness of the floor, walls, ceiling, pass box, doors and door view panels of man and material entry of dispensing area is clean.	✓	✓
5	Check the cleanliness of the working table, cupboard, washing area and waste bin.	✓	✓
6	Check the cleanliness of the drain points.	✓	✓
7	Ensure that the dispensing booth is cleaned and it's status.	✓	✓
8	Check the cleanliness of the supply and return grill.	✓	✓
9	Check the cleanliness of side wall of the dispensing booth.	✓	✓
10	Check the differential pressure of pre filter, intermediate filter and HEPA filter are within the limit	✓	✓
11	Check the cleanliness of the utility line (compressed air, First RO water and purified water).	✓	✓
12	Check the cleanliness of the SS side panel and SS Stool.	✓	✓
13	Ensure that the balances are cleaned and it's calibration status.	✓	✓
14	Check the status board of the area is updated.	✓	✓
15	Check the temperature & relative humidity of the area meet as per BMR requirement and ensure that the records are updated.	✓	✓
16	Check the differential pressure of the area should be within the limit.	✓	✓
17	Ensure that dispensing tools or accessories like scoops, spoons are cleaned and labeled with "Clean" label.	✓	✓
18	Check proper gowning procedure followed as per the SOP	✓	✓
19	Ensure the concurrent entries of log book.	✓	✓
20	Ensure that the Laminar of dispensing booth run Not less than 15 minutes after switch on.	✓	✓

Note: Tick mark (✓) if Complies and (X) if Does not comply

Checked by (store) sign/date: *[Signature]* 08/02/25

Verified by (IPQA) sign/date: *[Signature]* 08/02/25

Equipment	Dispensing Booth	Weighing Balance	AREA
Equipment. ID No.		<i>[Signature]</i>	
Previous Product & Batch No.		<i>AB-PW-SR-3 25/1701/01A</i>	<i>F-118</i>
Cleaned By		<i>Ray's</i>	
Checked by Stores Sign/Date		<i>[Signature]</i>	
Counter checked by IPQA/Sign/Date		<i>[Signature]</i>	

2.2 ENVIRONMENTAL CONDITIONS:

Area: Dispensing	Temperature: 23°C to 27°C	Observation	Stores Sign
	Relative Humidity: 45% to 55%	<i>23.5</i>	<i>[Signature]</i>
	Differential pressure NLT 0.5mm of water or NLT 5 pascal	<i>5</i>	
	Pre Filter 1-5mm of water	<i>NA</i>	
	Fine Filter 4-10mm of water	<i>NA</i>	
	HEAP Filter 8-20mm of water	<i>NA</i>	

2.3 APPROVAL FOR DISPENSING BY IPQA: Sign *[Signature]* Date: *08/02/25*

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature	<i>[Signature]</i>	<i>[Signature]</i>	<i>[Signature]</i>

CONTROLLED COPY

UNCONTROLLED COPY

Safetab Life Science
BATCH PACKING RECORD

Page 4 of 109

Product	LITACOLD FLU TABLETS	Product ID No.	978	Batch Size	1250000 Tablets
Generic	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

2.4 RECORD FOR PACKING MATERIALS: (12 TRACK)

Material Name	Unit	Material code	Standard Qty.	A/R No.	Issued Qty	Issued by store	Checked by production	Verified by IQA
Printed strip foil:302mm	Kg	PML00081	162.500	alp243479	162.000	A	Pi	(S)
				NA				
Plain strip back foil: 302mm	Kg	PMST0044	162.500	alp243480	162.000	A	Pi	(S)
				NA				
Catch cover	No's	PML00084	312,500	alp243864	312.000	A	Pi	(S)
				NA				
Printed cartons 1 X 4's	No's	PML00082	15625	alp250012	15625	A	Pi	(S)
				NA				
LITACOLD FLU leaflet	No's	PML00083	15625	alp243539	15625	A	Pi	(S)
				NA				
Shippers (75X20X1X4's)	No's	PTMT0027	213	alp250009	44	A	Pi	(S)
				alp250093	169	A	Pi	(S)
				NA				
Transparent GHPL printed Bopp Tape 2"	No's	PMG00071	13	alp250191	13	A	Pi	(S)
				NA				
Shipper label	No's		213		205		Pi	(S)

NOTE: If any discrepancies observed in received material inform to QA and take necessary corrective action.

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature			
UNCONTROLLED COPY		CONTROLLED COPY	

 Safetab Life Science BATCH PACKING RECORD		Page 5 of 108	
Product:	LITACOLD FLU TABLETS	Product ID No.:	976
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3
Batch No.:	G17250103	Mfg. Date:	01/2025
		Batch Size	1250000 Tablets
		Supersedes MPR No:	ST/MPR/004/R3
		Exp. Date:	12/2027

2.4.1 RECORD FOR ADDITIONAL PACKING MATERIALS

Reason for Additional Quantity Request

Approved by Production Head/Designee Sign/date:

Authorized by QA Head/Designee Sign/date:

Material Name	Unit	Material code	A.R. No.	Issued Qty	Issued by store	Checked by production	Verified by IQA
Printed strip foil: 302mm	Kg	PML00081					
Plain strip back foil: 302mm	Kg	PMST0044					
catch cover	No's	PML00084					
Printed cartons 1 X 4's	No's	PML00082					
LITACOLD FLU leaflet	No's	PML00083					
Shippers (75X20X1X4's)	No's	PTMT0027					
Transparent GHPL printed Bopp Tape 2"	No's	PMG00071					
Shipper label	No's	*****					

12/02/25

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature			
UNCONTROLLED COPY		CONTROLLED COPY	



Safetab Life Science
BATCH PACKING RECORD

Page 6 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

2.5 RECORD FOR PACKING MATERIALS: (6 TRACK)

Material Name	Unit	Material code	Standard Qty.	A.R. No.	Issued Qty	Issued by store	Checked by production	Verified by IPQA
Printed strip foil: 180mm	Kg	PML00081	162.500					
Plain strip back foil: 180mm	Kg	PMST0035	162.500					
catch cover	No's	PMST0044	312.500					
Printed cartons 1 X 4's	No's	PML00082	15625					
LITACOLD FLU leaflet	No's	PML00083	15625					
Shippers (75X20X1X4's)	No's	PTMT0027	213					
Transparent GHPL printed Bopp Tape 2"	No's	PMG00071	13					
Shipper label	No's	-----	213					

8/11/25
14/02/25

NOTE: If any discrepancies observed in received material inform to QA and take necessary corrective action.

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature				
UNCONTROLLED COPY				CONTROLLED COPY

Safetab Life Science BATCH PACKING RECORD				Page 7 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

2.5.1 RECORD FOR ADDITIONAL PACKING MATERIALS

Reason for Additional Quantity Request

Approved by Production Head/Designee Sign/date:

Authorized by QA Head/Designee Sign/date:

Material Name	Unit	Material code	A.R. No.	Issued Qty	Issued by store	Checked by production	Verified by IPQA
Printed strip foil:180mm	Kg	PML00081					
Plain strip back foil: 180mm	Kg	PMST0035					
catch cover	No's	PML00084					
Printed cartons 1 X 4's	No's	PML00082					
LITACOLD FLU leaflet	No's	PML00083					
Shippers (75X20X1X4's)	No's	PTMT0027					
Transparent GHPL printed Bopp Tape 2"	No's	PMG00071					
Shipper label	No's						

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature	<i>[Signature]</i>	<i>[Signature]</i>	<i>[Signature]</i>	
UNCONTROLLED COPY				



Safetab Life Science
BATCH PACKING RECORD

Page 8 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

2.6 MATERIAL RETURN NOTE:

Material Name	Unit	size	A.R. No.	Quantity returned	Returned by production	Verified by IPQA	Received by store
Printed strip foil	Kg						
Plain strip back foil	Kg						
catch cover	No's						
Printed cartons 1 X 4's	No's						
LITACOLD FLU leaflet	No's						
Transparent GHPL printed Bopp Tape 2"	No's						
1/2 inch cello tape	No's						
Shipper label	No's						

14/02/25

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY				

Safetab Life Science BATCH PACKING RECORD				Page 9 of 109
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size: 1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No: ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date: 12/2027

3.0 AFFIX DISPENSING CARD FOR PRIMARY PACKING MATERIAL:

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PML00081/01 | Category: NA

Item Name: 302MM LITACOLD FLU SALES PRINTED FOIL

Operation: Packing

AR No.: QCP243479 | Retest Date: 21/11/25

Mfr: DAGA POLY LAMINATORS (V0197)

Lot No.: 1 | Pack 1 of 1 | Bal. Id.: ST/SRWB/049
 Req. Qty.: 162.5 Kg | Issd. Qty.: 16.0 Kg | Conv. Fa.: 1
 Gr. Wt.: 16.3 Kg | Tr. Wt.: 0.3 Kg | Nt. Wt.: 16.0 Kg



Dispensed By
(Store)

LBST:001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0064 - 1

Verified By
(Prod.)

Sat Feb 08 25 17:40:07

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PML00081/01 | Category: NA

Item Name: 302MM LITACOLD FLU SALES PRINTED FOIL

Operation: Packing

AR No.: QCP243479 | Retest Date: 21/11/25

Mfr: DAGA POLY LAMINATORS (V0197)

Lot No.: 1 | Pack 2 of 1 | Bal. Id.: ST/SRWB/049
 Req. Qty.: 162.5 Kg | Issd. Qty.: 32.0 Kg | Conv. Fa.: 1
 Gr. Wt.: 16.3 Kg | Tr. Wt.: 0.3 Kg | Nt. Wt.: 16.0 Kg



Dispensed By
(Store)

LBST:001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0064 - 3

Verified By
(Prod.)

Sat Feb 08 25 17:49:22

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PML00081/01 | Category: NA

Item Name: 302MM LITACOLD FLU SALES PRINTED FOIL

Operation: Packing

AR No.: QCP243479 | Retest Date: 21/11/25

Mfr: DAGA POLY LAMINATORS (V0197)

Lot No.: 1 | Pack 3 of 1 | Bal. Id.: ST/SRWB/049
 Req. Qty.: 162.5 Kg | Issd. Qty.: 48.0 Kg | Conv. Fa.: 1
 Gr. Wt.: 16.3 Kg | Tr. Wt.: 0.3 Kg | Nt. Wt.: 16.0 Kg



Dispensed By
(Store)

LBST:001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0064 - 4

Verified By
(Prod.)

Sat Feb 08 25 17:49:38

Dispensed materials label pasted by (production sign/date): 11/02/25 verified by (IPQA) sign/date: 11/02/25

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	



Safetab Life Science
BATCH PACKING RECORD

Page 10 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

3.OA AFFIX DISPENSING CARD FOR PRIMARY PACKING MATERIAL:

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PML00081/01 | Category: NA

Item Name: 302MM LITACOLD FLU SALES PRINTED FOIL

Operation: Packing

AR No.: QCP243479 | Retest Date: 21/11/25

Mfr: DAGA POLY LAMINATORS (V0197)

Lot No.: 1 | Pack 4 of 1 | Bal. Id.: ST/SRWB/049
Req. Qty.: 162.5 Kg | Issd. Qty.: 64.0 Kg | Conv. Fa.: 1
Gr. Wt.: 16.3 Kg | Tr. Wt.: 0.3 Kg | Nt. Wt.: 16.0 Kg



Dispensed By
(Store)

LBST:001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0004 - 5

Verified By
(Prod.)

Sat Feb 08 25 17:49:53

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PML00081/01 | Category: NA

Item Name: 302MM LITACOLD FLU SALES PRINTED FOIL

Operation: Packing

AR No.: QCP243479 | Retest Date: 21/11/25

Mfr: DAGA POLY LAMINATORS (V0197)

Lot No.: 1 | Pack 5 of 1 | Bal. Id.: ST/SRWB/049
Req. Qty.: 162.5 Kg | Issd. Qty.: 81.0 Kg | Conv. Fa.: 1
Gr. Wt.: 17.3 Kg | Tr. Wt.: 0.3 Kg | Nt. Wt.: 17.0 Kg



Dispensed By
(Store)

LBST:001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0004 - 6

Verified By
(Prod.)

Sat Feb 08 25 17:50:14

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets

Item Code: PML00081/01 | Category: NA

Item Name: 302MM LITACOLD FLU SALES PRINTED FOIL

Operation: Packing

AR No.: QCP243479 | Retest Date: 21/11/25

Mfr: DAGA POLY LAMINATORS (V0197)

Lot No.: 1 | Pack 6 of 1 | Bal. Id.: ST/SRWB/049
Req. Qty.: 162.5 Kg | Issd. Qty.: 98.0 Kg | Conv. Fa.: 1
Gr. Wt.: 17.3 Kg | Tr. Wt.: 0.3 Kg | Nt. Wt.: 17.0 Kg



Dispensed By
(Store)

LBST:001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0004 - 7

Verified By
(Prod.)

Sat Feb 08 25 17:50:35

Dispensed materials label pasted by(production sign/date: 11/02/25 verified by(IPQA) sign/date: 01/02/25

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY			CONTROLLED COPY

Safetab Life Science BATCH PACKING RECORD					
Page 11 of 109					
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/003/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

3.0B AFFIX DISPENSING CARD FOR PRIMARY PACKING MATERIAL:

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PML00081/01 | Category: NA

Item Name: 302MM LITACOLD FLU SALES PRINTED FOIL

Operation: Packing

AR No.: QCP243479 | Retest Date: 21/11/25

Mir: DAGA POLY LAMINATORS (V0197)

Lot No.: 1 | Pack 7 of 1 | Bal. Id.: ST/SRWB/049

Req. Qty.: 162.5 Kg | Issd. Qty.: 114.5 Kg | Conv. Fa.: 1

Gr. Wt.: 16.8 Kg | Tr. Wt.: 0.3 Kg | NI. Wt.: 16.5 Kg



DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets

Item Code: PML00081/01 | Category: NA

Item Name: 302MM LITACOLD FLU SALES PRINTED FOIL

Operation: Packing

AR No.: QCP243479 | Retest Date: 21/11/25

Mir: DAGA POLY LAMINATORS (V0197)

Lot No.: 1 | Pack 9 of 1 | Bal. Id.: ST/SRWB/049

Req. Qty.: 162.5 Kg | Issd. Qty.: 147.5 Kg | Conv. Fa.: 1

Gr. Wt.: 16.8 Kg | Tr. Wt.: 0.3 Kg | NI. Wt.: 16.5 Kg

Verified By (Prod.)

Checked By (QA)

Dispensed By (Store)

Sat Feb 08 25 17:51:38

ID: DE-REC-2025-0084 - 10

LBST:001 Revision No: R1

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PML00081/01 | Category: NA

Item Name: 302MM LITACOLD FLU SALES PRINTED FOIL

Operation: Packing

AR No.: QCP243479 | Retest Date: 21/11/25

Mir: DAGA POLY LAMINATORS (V0197)

Lot No.: 1 | Pack 8 of 1 | Bal. Id.: ST/SRWB/049

Req. Qty.: 162.5 Kg | Issd. Qty.: 131.0 Kg | Conv. Fa.: 1

Gr. Wt.: 16.8 Kg | Tr. Wt.: 0.3 Kg | NI. Wt.: 16.5 Kg



Dispensed By (Store)

Checked By (QA)

Verified By (Prod.)

LBST:001 Revision No: R1

ID: DE-REC-2025-0084 - 9

Sat Feb 08 25 17:51:17

Dispensed materials label pasted by (production sign/date):

verified by (IPQA) sign/date:

11/02/25

11/02/25

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	



Safetab Life Science BATCH PACKING RECORD

Page 12 of 100

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

3.0C AFFIX DISPENSING CARD FOR PRIMARY PACKING MATERIAL:

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)
Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PML00081/01 | Category: NA
Item Name: 302MM LITACOLD FLU SALES PRINTED FOIL
Operation: Packing



AR No.: QCP243479 | Retest Date: 21/11/25
Mfr: DAGA POLY LAMINATORS (V0197)

Lot No.: 1 | Pack 10 of 1 | Bal. Id.: ST/SRWB/049
Req. Qty.: 162.5 Kg | Issd. Qty.: 163.0 Kg | Conv. Fa.: 1
Gr. Wt.: 15.8 Kg | Tr. Wt.: 0.3 Kg | Nt. Wt.: 15.5 Kg

Dispensed By
(Store)

LBST-001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0004 - 11

Verified By
(Prod.)

Sat Feb 08 25 17:52:01

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)
Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PMST0044 | Category: NA
Item Name: 302MM STRIP FOIL BSO
Operation: Packing



AR No.: QCP243480 | Retest Date: 21/11/25
Mfr: DAGA POLY LAMINATORS (V0197)

Lot No.: 1 | Pack 1 of 1 | Bal. Id.: ST/SRWB/049
Req. Qty.: 162.5 Kg | Issd. Qty.: 17.5 Kg | Conv. Fa.: 1
Gr. Wt.: 17.8 Kg | Tr. Wt.: 0.3 Kg | Nt. Wt.: 17.5 Kg

Dispensed By
(Store)

LBST-001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0004 - 2

Verified By
(Prod.)

Sat Feb 08 25 17:52:28



DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)
Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PMST0044 | Category: NA

Item Name: 302MM STRIP FOIL BSO

Operation: Packing

AR No.: QCP243480 | Retest Date: 21/11/25

Mfr: DAGA POLY LAMINATORS (V0197)

Lot No.: 1 | Pack 2 of 1 | Bal. Id.: ST/SRWB/049
Req. Qty.: 162.5 Kg | Issd. Qty.: 35.0 Kg | Conv. Fa.: 1
Gr. Wt.: 17.8 Kg | Tr. Wt.: 0.3 Kg | Nt. Wt.: 17.5 Kg

Dispensed By
(Store)

LBST-001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0004 - 12

Verified By
(Prod.)

Sat Feb 08 25 17:52:41

Dispensed materials label pasted by (production sign/date):

H1102125

verified by (IPQA) sign/date:

a11/02/25

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY			CONTROLLED COPY

Safetab Life Science
BATCH PACKING RECORD

Page 13 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

4.0 AFFIX DISPENSING CARD FOR SECONDARY PACKING MATERIAL:

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PMQ00071 | Category: NA

Item Name: GHPL PRINTED TRANSPERENT BOPP TAPE 2"

Operation: Packing

AR No.: QCP260191 | Retest Date: 21/01/27

Mfr: SRI VIGNESHWARA PACKAGE (V0220)

Lot No.: 1 | Pack 1 of 1 | Bal. Id.: NA

Req. Qty.: 13.0 Nos | Issd. Qty.: 13.0 Nos | Conv. Fa.: 1

Gr. Wt.: NA | Tr. Wt.: NA | Nt. Wt.: 13.0 Nos



Dispensed By
(Store)

LBST:001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0003 - 0

Verified By
(Prod.)

Sat Feb 08 25 17:58:45

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PML00082 | Category: NA

Item Name: LITACOLD FLU SALES PRINTED CARTON

Operation: Packing

AR No.: QCP260012 | Retest Date: 03/01/27

Mfr: ANUP PRINTERS PVT. LTD. (V0285)

Lot No.: 1 | Pack 1 of 1 | Bal. Id.: NA

Req. Qty.: 15625.0 Nos | Issd. Qty.: 15625.0 Nos | Conv. Fa.: 1

Gr. Wt.: NA | Tr. Wt.: NA | Nt. Wt.: 15625.0 Nos



DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets

Item Code: PML00083 | Category: NA

Item Name: LITACOLD FLU LEAFLET

Operation: Packing

AR No.: QCP243539 | Retest Date: 26/11/26

Mfr: ANUP PRINTERS PVT. LTD. (V0285)

Lot No.: 1 | Pack 1 of 1 | Bal. Id.: NA

Req. Qty.: 15625.0 Nos | Issd. Qty.: 15625.0 Nos | Conv. Fa.: 1

Gr. Wt.: NA | Tr. Wt.: NA | Nt. Wt.: 15625.0 Nos

Checked By
(QA)

ID: DE-REC-2025-0003 - 3

Dispensed By
(Store)

LBST:001 Revision No: R1

Verified By
(Prod.)

Sat Feb 08 25 17:57:05

Dispensed materials label pasted by (production sign/date):

[Signature] 11/02/25

verified by (IPQA) sign/date:

[Signature] 11/02/25

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with	<i>[Signature]</i>	<i>[Signature]</i> 11/05/24	<i>[Signature]</i>

UNCONTROLLED COPY

CONTROLLED COPY

Safetab Life Science BATCH PACKING RECORD				Page 14 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

4.0A AFFIX DISPENSING CARD FOR SECONDARY PACKING MATERIAL:

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PML00084 | Category: NA

Item Name: LITACOLD FLU PRINTED CATCH COVER

Operation: Packing

AR No.: QCP243864 | Retest Date: 04/01/27

Mfr: ANUP PRINTERS PVT. LTD. (V0285)



Lot No.: 1 | Pack 1 of 1 | Bal. Id.: NA
 Req. Qty.: 312500.0 Nos | Issd. Qty.: 312500.0 Nos | Conv. Fa.: 1
 Gr. Wt.: NA | Tr. Wt.: NA | Nt. Wt.: 312500.0 Nos

Dispensed By
(Store)

LBST:001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0003 - 1

Verified By
(Prod.)

Sat Feb 08 25 17:57:15

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | Batch Size: 1250000 Tablets



Item Code: PTMT0027 | Category: NA

Item Name: 390X345X395MM GHPL LITACOLD FLU PRINTED CATCH COVER

Operation: Packing

AR No.: QCP250092 | Retest Date: 11/01/27

Mfr: JJS PACKAGING (V0219)



Lot No.: 1 | Pack 1 of 1 | Bal. Id.: NA
 Req. Qty.: 213.0 Nos | Issd. Qty.: 44.0 Nos | Conv. Fa.: 1
 Gr. Wt.: NA | Tr. Wt.: NA | Nt. Wt.: 44.0 Nos

Dispensed By
(Store)

LBST:001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0003 - 4

Verified By
(Prod.)

Sat Feb 08 25 17:57:25

Dispensed materials label pasted by(production sign/date:

11/02/25

verified by(IPQA) sign/date:

11/02/25

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY				

 Safetab Life Science BATCH PACKING RECORD						Page 15 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets		
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2		
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027		

4.08 AFFIX DISPENSING CARD FOR SECONDARY PACKING MATERIAL:

DISPENSED MATERIAL

Product: LITACOLD FLU(4's) (BPR-978-01)

Batch No.: G17250103 | **Batch Size:** 1250000 Tablets



Item Code: PTMT0027 | **Category:** NA

Item Name: 390X345X395MM GHPL LITACOLD FLU PRINTED 4's

Operation: Packing

AR No.: QCP250098 | **Retest Date:** 11/01/27

Mfr: JJS PACKAGING (V0219)

Lot No.: 1 | **Pack 1 of** 1 | **Bal. Id.:** NA

Req. Qty.: 213.0 Nos | **Issd. Qty.:** 169.0 Nos | **Conv. Fa.:** 1

Gr. Wt.: NA | **Tr. Wt.:** NA | **Nt. Wt.:** 169.0 Nos



Dispensed By
(Store)

LB67:001 Revision No: R1

Checked By
(QA)

ID: DE-REC-2025-0003 - 5

Verified By
(Prod.)

8th Feb 08 25 17:57:30

Dispensed materials label pasted by (production sign/date):

[Signature] 11/02/25

verified by (IPQA) sign/date:

[Signature] 11/02/25

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with	<i>[Signature]</i>	<i>[Signature]</i>	<i>[Signature]</i>	

UNCONTROLLED COPY

 Safetab Life Science BATCH PACKING RECORD		Page 16 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No:	978
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3
Batch No.:	G17250103	Mfg.Date:	01/2025
Batch Size:	1250000 Tablets	Supersedes MPR No:	ST/MPR/005/R2
		Exp.Date:	12/2027

5.0 CARTON CHECK PRIOR TO OVER CODING:

AREA NAME	Inkjet Printing	AREA CODE	F124		
5.1.0 CARTON PRE-PRINTED LINE CLEARANCE:					
1	Check the cleanliness of the checking area.				
2	Ensure the Left over carton of previous batch/product if any is removed from area.				
3	Check the waste bin for the absence of previous product cartons.				
Checked by / date (Packing):	[Signature] 11/02/25				
Verified by /date (IPQA):	[Signature] 11/02/25				
5.1.1 CARTON CHECKING DETAILS : (ENTER THE DETAILS FROM CARTON)					
Product Name :	Litacold flu				
Manufacturing License Number :	21144378				
Packing Code:	1511326				
Carton Received By/On:	[Signature] 11/02/25				
5.1.2 CARTON CHECKS					
Total carton received quantity in number's	15625	Total number of box's	15		
Received carton A.R.Number with quantity	QCP 25 0012		15625		
NA [Signature] 11/02/25					
5.1.3 CHECKING PROCESS					
Date	Shift	Time (Hrs)		Cartons checked by	Cartons verified by
		From	To		
11/02/25	II	15:50	15:57	Viji	sneha
NA [Signature] 11/02/25					

Job title:	Prepared by:	Checked by:	Approved by:
Designation:	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with	[Signature]	[Signature]	[Signature]
UNCONTROLLED COPY		CONTROLLED COPY	

 Safetab Life Science BATCH PACKING RECORD		Page 17 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3
Batch No.:	G17250103	Mfg. Date:	01/2025
		Batch Size	1250000 Tablets
		Supersedes MPR No:	ST/MPR/005/R2
		Exp. Date:	12/2027

5.1.4 PRE-PRINTED CARTON RECONCILIATION DETAILS IN NUMBERS


Date	Quantity of carton received	Quantity of carton Inspected	Rejected carton due to brand name / generic name	Other Rejected carton	Good carton issued for packing	Checked by Packing	Verified by IPQA
11/02/25	15625	15625	-	-	15625	ff	Q
NA ff 11/02/25							
Total	15625	15625	-	-	15625	ff	Q

5.2 CARTON OVER PRINTING RECORD:

5.2.1 CARTON OVER PRINTING DETAILS	To be entered by referring the Job Order	Checked by Production	Verified by IPQA
PRODUCT NAME :	Litacold flu	ff	Q
BATCH NO :	G17250103		
MFG DT :	01/2025		
EXP DT :	12/2027		
MRP :	-NA-		

5.2.2 CARTON OVER PRINTING LINE CLEARANCE:

1	Check the cleanliness of the printing area.
2	Ensure the Left over carton of previous batch/product if any is removed from area.
3	Check the waste bin for the absence of previous product cartons.
4	Check the machine log book are updated and signed.
5	Ensure free from previous product /batch details in inject printing machine.

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

 Safetab Life Science BATCH PACKING RECORD		Page 18 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3
Batch No.:	G17250103	Mfg. Date:	01/2025
		Exp. Date:	12/2027

5.2.3 INITIAL CHECKS

	DATE: 11/02/25		DATE: 12/02/25		DATE: 13/02/25	
	HEAD 1	HEAD 2	HEAD 1	HEAD 2	HEAD 1	HEAD 2
PRODUCT NAME :	Litacold flu	Litacold flu	Litacold flu		Litacold flu	Litacold flu
BATCH NO :	G17250103	G17250103	G17250103		G17250103	G17250103
MFG.DT :	01/2025	01/2025	01/2025		01/2025	01/2025
EXP.DT :	12/2027	12/2027	12/2027		12/2027	12/2027
MRP :	-NA-	-NA-	-NA-		-NA-	-NA-
CHECKED BY PRODUCTION:	ft	ft	ft	ft	ft	ft
CHECKED BY IPQA :	Q	Q	Q	12/02/25	Q	Q
REMARKS :	-NIL-	-Nil-	-NIL-		Nil	Nil

	DATE: 13/02/25	
	HEAD 1	HEAD 2
PRODUCT NAME :	Litacold flu	
BATCH NO :	G17250103	
MFG.DT :	01/2025	
EXP.DT :	12/2027	
MRP :	-	
CHECKED BY PRODUCTION:	Q	
CHECKED BY IPQA :	Q	
REMARKS :	Nil	

Note: If any deviation observed inform to QA head and production head for further course of action. Attach a specimen of each head duly signed by production and IPQA with the BPR.

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

QME ft 11/02/25



Safetab Life Science

BATCH PACKING RECORD


Page 19 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.3A INITIAL CHECKS

	DATE:		DATE:		DATE:	
	HEAD 1	HEAD 2	HEAD 1	HEAD 2	HEAD 1	HEAD 2
PRODUCT NAME :						
BATCH NO :						
MFG.DT :						
EXP.DT :						
MRP :						
CHECKED BY PRODUCTION:						
CHECKED BY IPQA :						
REMARKS :						
	DATE:		DATE:		DATE:	
	HEAD 1	HEAD 2	HEAD 1	HEAD 2	HEAD 1	HEAD 2
PRODUCT NAME :						
BATCH NO :						
MFG.DT :						
EXP.DT :						
MRP :						
CHECKED BY PRODUCTION:						
CHECKED BY IPQA :						
REMARKS :						

Note: If any deviation observed inform to QA head and production head for further course of action. Attach a specimen of each head duly signed by production and IPQA with the BPR.

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science BATCH PACKING RECORD

Page 20 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.4 CARTON PRINTING INPROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR

DATE :	11/02/25	11/02/25	11/02/25	11/02/25
TIME :	16:05	16:11	17:10	17:16
PRODUCT NAME :	Litacold flu	Litacold flu	Litacold flu	Litacold flu
BATCH NO :	G17250103	G17250103	G17250103	G17250103
MFG.DT :	01/2025	01/2025	01/2025	01/2025
EXP.DT :	12/2027	12/2027	12/2027	12/2027
MRP :	-NA-	-NA-	-NA-	-NA-
PRINTING CLARITY :	OK	OK	OK	OK
CARTON PRINTED BY	Viji, selvakumar	Viji, selvakumar	Viji, selvakumar	Viji, selvakumar
PRINTING CHECKED BY	Srija, sneha	Srija, sneha	Srija, sneha	Srija, sneha
CHECKED BY PRODUCTION/IPQA SIGN/DATE :	11/02/25	11/02/25	11/02/25	11/02/25
REMARKS :	-NIL-	-NIL-	-NIL-	-NIL-
DATE :	11/02/25	11/02/25	12/02/25	12/02/25
TIME :	18:13	18:20	18:00	18:05
PRODUCT NAME :	Litacold flu	Litacold flu	Litacold flu	Litacold flu
BATCH NO :	G17250103	G17250103	G17250103	G17250103
MFG.DT :	01/2025	01/2025	01/2025	01/2025
EXP.DT :	12/2027	12/2027	12/2027	12/2027
MRP :	-NA-	-NA-	-NA-	-NA-
PRINTING CLARITY :	OK	OK	OK	OK
CARTON PRINTED BY	Viji, selvakumar	Viji, selvakumar	Vaitheswaran	Vaitheswaran
PRINTING CHECKED BY	Srija, sneha	Srija, sneha	Srija	Srija
CHECKED BY PRODUCTION/IPQA SIGN/DATE :	11/02/25	11/02/25	12/02/25	12/02/25
REMARKS :	-NIL-	-NIL-	-NIL-	-NIL-
TIME :	07:05 (13/02/25)	07:10 (13/02/25)	08:10	08:15
PRODUCT NAME :	Litacold flu	Litacold flu	Litacold flu	Litacold flu
BATCH NO :	G17250103	G17250103	G17250103	G17250103
MFG.DT :	01/2025	01/2025	01/2025	01/2025
EXP.DT :	12/2027	12/2027	12/2027	12/2027
MRP :	-NA-	-NA-	-NA-	-NA-
PRINTING CLARITY :	OK	OK	OK	OK
CARTON PRINTED BY	Poorkalai, Gnan	Poorkalai, Gnan	Poorkalai, Gnan	Poorkalai, Gnan
PRINTING CHECKED BY	Santhya, Maha	Santhya, Maha	Santhya, Maha	Santhya, Maha
CHECKED BY PRODUCTION/IPQA SIGN/DATE :	13/02/25	13/02/25	13/02/25	13/02/25
REMARKS :	-NIL-	-NIL-	-NIL-	-NIL-

Job title:	Prepared by:	Checked by:	Approved by:
Production	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science

BATCH PACKING RECORD

Page 21 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.4 (A) CARTON PRINTING INPROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR

DATE :	13/02/25	13/02/2025	13/02/25	13/02/2025
TIME :	09:20	09:25	12:30	12:35
PRODUCT NAME :	Litacold Flu	Litacold Flu	Litacold Flu	Litacold Flu
BATCH NO :	G17250103	G17250103	G17250103	G17250103
MFG.DT :	01/2025	01/2025	01/2025	01/2025
EXP.DT :	12/2027	12/2027	12/2027	12/2027
MRP :	-NA-	-NA-	-NA-	-NA-
PRINTING CLARITY :	OK	OK	OK	OK
CARTON PRINTED BY	Perkalai, Gouri	Perkalai, Gouri	Rathin	Rathin, Sanchya
PRINTING CHECKED BY	Sanchya, Hahn	Sanchya, Hahn	Karun	Karun, Perkalai
CHECKED BY PRODUCTION/IPQA SIGN/DATE :	P. 13/02/25	P. 13/02/2025	P. 13/02/25	P. 13/02/2025
REMARKS :	NIL	NIL	NIL	NIL
DATE :	13/02/25	13/02/25	13/02/25	13/02/25
TIME :	13:35	13:40	14:40	14:45
PRODUCT NAME :	Litacold Flu	Litacold Flu	Litacold Flu	Litacold Flu
BATCH NO :	G17250103	G17250103	G17250103	G17250103
MFG.DT :	01/2025	01/2025	01/2025	01/2025
EXP.DT :	12/2027	12/2027	12/2027	12/2027
MRP :	-NA-	-NA-	-NA-	-NA-
PRINTING CLARITY :	OK	OK	OK	OK
CARTON PRINTED BY	Perkalai, Gouri	Perkalai, Gouri	Viji, Sanchya	Viji, Sanchya
PRINTING CHECKED BY	Sanchya, Hahn	Sanchya, Hahn	Sanchya, Hahn	Sanchya, Hahn
CHECKED BY PRODUCTION/IPQA SIGN/DATE :	P. 13/02/25	P. 13/02/25	P. 13/02/25	P. 13/02/25
REMARKS :	NIL	NIL	NIL	NIL
DATE :	13/02/25	13/02/25		
TIME :	15:45	15:50		
PRODUCT NAME :	Litacold Flu	Litacold Flu		
BATCH NO :	G17250103	G17250103		
MFG.DT :	01/2025	01/2025		
EXP.DT :	12/2027	12/2027		
MRP :	-NA-	-NA-		
PRINTING CLARITY :	OK	OK		
CARTON PRINTED BY	Viji, Sanchya	Viji, Sanchya		
PRINTING CHECKED BY	Sanchya, Hahn	Sanchya, Hahn		
CHECKED BY PRODUCTION/IPQA SIGN/DATE :	P. 13/02/25	P. 13/02/25		
REMARKS :	NIL	NIL		

W.F. 13/02/25

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	



Safetab Life Science
BATCH PACKING RECORD

Page 22 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR No:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.4 (B) CARTON PRINTING INPROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR

DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				CONTROLLED COPY
UNCONTROLLED COPY				

Safetab Life Science BATCH PACKING RECORD

Page 23 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.4 (C) CARTON PRINTING IN PROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR

DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				


Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY 11/05/24 11/05/24				

		Safetab Life Science BATCH PACKING RECORD		Page 24 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.4 (D) CARTON PRINTING IN PROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR


DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY				

 Safetab Life Science BATCH PACKING RECORD		Page 25 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3
Batch No.:	G17250103	Mfg.Date:	01/2025
		Exp.Date:	12/2027
		Batch Size	1250000 Tablets
		Supersedes MPR No:	ST/MPR/005/R2

5.2.4 (E) CARTON PRINTING INPROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR

DATE :			
TIME :			
PRODUCT NAME :			
BATCH NO :			
MFG.DT :			
EXP.DT :			
MRP :			
PRINTING CLARITY :			
CARTON PRINTED BY			
PRINTING CHECKED BY			
CHECKED BY PRODUCTION/IPQA			
SIGN/DATE :			
REMARKS :			
DATE :			
TIME :			
PRODUCT NAME :			
BATCH NO :			
MFG.DT :			
EXP.DT :			
MRP :			
PRINTING CLARITY :			
CARTON PRINTED BY			
PRINTING CHECKED BY			
CHECKED BY PRODUCTION/IPQA			
SIGN/DATE :			
REMARKS :			
DATE :			
TIME :			
PRODUCT NAME :			
BATCH NO :			
MFG.DT :			
EXP.DT :			
MRP :			
PRINTING CLARITY :			
CARTON PRINTED BY			
PRINTING CHECKED BY			
CHECKED BY PRODUCTION/IPQA			
SIGN/DATE :			
REMARKS :			

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	





Safetab Life Science
BATCH PACKING RECORD

Page 26 of 100

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.4 (F) CARTON PRINTING IN PROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR

DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science

BATCH PACKING RECORD

Page 27 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.4 (G) CARTON PRINTING INPROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR


DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				
UNCONTROLLED COPY			CONTROLLED COPY	

 Safetab Life Science BATCH PACKING RECORD						Page 28 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets		
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2		
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027		

5.2.4 (H) CARTON PRINTING INPROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR

DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				

UNCONTROLLED COPY

Safetab Life Science
BATCH PACKING RECORD

Page 29 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.4 (I) CARTON PRINTING IN PROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR


DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY :				
PRINTING CHECKED BY :				
CHECKED BY PRODUCTION/IPQA :				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY :				
PRINTING CHECKED BY :				
CHECKED BY PRODUCTION/IPQA :				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY :				
PRINTING CHECKED BY :				
CHECKED BY PRODUCTION/IPQA :				
SIGN/DATE :				
REMARKS :				

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY				

Safetab Life Science BATCH PACKING RECORD						Page 30 of 109
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets	
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2	
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027	

5.2.4 (J) CARTON PRINTING INPROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR

DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY				



Safetab Life Science
BATCH PACKING RECORD

Page 31 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/P3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.4 (K) CARTON PRINTING INPROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR

DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				

14/02/25

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY				

Safetab Life Science
BATCH PACKING RECORD

Page 32 of 109

Product	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.4 (L) CARTON PRINTING INPROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR

DATE :	
TIME :	
PRODUCT NAME :	
BATCH NO :	
MFG.DT :	
EXP.DT :	
MRP :	
PRINTING CLARITY :	
CARTON PRINTED BY	
PRINTING CHECKED BY	
CHECKED BY PRODUCTION/IPQA SIGN/DATE :	
REMARKS :	
DATE :	
TIME :	
PRODUCT NAME :	
BATCH NO :	
MFG.DT :	
EXP.DT :	
MRP :	
PRINTING CLARITY :	
CARTON PRINTED BY	
PRINTING CHECKED BY	
CHECKED BY PRODUCTION/IPQA SIGN/DATE :	
REMARKS :	

NOTE: Signed specimen carton shall be attached during start and end of the day of printing process and batch end

5.2.5 PRINTED CARTON RECONCILIATION DETAILS IN NUMBERS

Date	Quantity Of over Printed Carton	Quantity of carton Inspected	Quantity Of Printing rejected Carton	Good carton issued for packing	Checked by Packing	Verified by IPQA
14/12/25	15625	15625	301	15324	✓	✓
						14/12/25
Total	15625	15625	301	15324	✓	✓

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			

UNCONTROLLED COPY

CONTROLLED COPY



Safetab Life Science
BATCH PACKING RECORD

MA

Page 33 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.4 (M) CARTON PRINTING INPROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR

DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY				
PRINTING CHECKED BY				
CHECKED BY PRODUCTION/IPQA				
SIGN/DATE :				
REMARKS :				

NOTE: Signed specimen carton shall be attached during start and end of the day of printing process and batch end

5.2.5 PRINTED CARTON RECONCILIATION DETAILS IN NUMBERS

Date	Quantity Of over Printed Carton	Quantity of carton Inspected	Quantity Of Printing rejected Carton	Good carton issued for packing	Checked by Packing	Verified by IPQA
Total					8/1/25 14/02/25	

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY				

Safetab Life Science

BATCH PACKING RECORD

Page 34 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tab. etc
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR No:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

5.2.4 (N) CARTON PRINTING INPROCESS CHECKS PRODUCTION AND IPQA: EVERY ONE HOUR

DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY :				
PRINTING CHECKED BY :				
CHECKED BY PRODUCTION/IPQA :				
SIGN/DATE :				
REMARKS :				
DATE :				
TIME :				
PRODUCT NAME :				
BATCH NO :				
MFG.DT :				
EXP.DT :				
MRP :				
PRINTING CLARITY :				
CARTON PRINTED BY :				
PRINTING CHECKED BY :				
CHECKED BY PRODUCTION/IPQA :				
SIGN/DATE :				
REMARKS :				


NOTE: Signed specimen carton shall be attached during start and end of the day of printing process and batch end

5.2.5 PRINTED CARTON RECONCILIATION DETAILS IN NUMBERS

Date	Quantity Of over Printed Carton	Quantity of carton Inspected	Quantity Of Printing rejected Carton	Good carton issued for packing	Checked by Packing	Verified by IPQA
Total						

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				

UNCONTROLLED COPY

		Safetab Life Science BATCH PACKING RECORD		 Page 35 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125,000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

6.0 PACKING PROCESS INSTRUCTIONS:

S.NO	INSTRUCTION
1	Wear hand gloves during operation.
2	Transfer approved bulk product and packing material to packing area after verification, record the verification details.
3	Set and operate the machine as per current version of SOP.
4	Check the initial machine setting parameters and record the observations.
5	Check and attach specimen of printed foil, carton and overprinted shipper label before starting the packing activity.
6	Check and attach specimen of overprinted /embossed foil, carton and overprinted shipper label before starting the packing activity.
7	Check the initially run blister for all the parameters recorded in the process checks. It should be counter checked by IPQA.
8	Carry out in-process checks and record observations in respective tables as per frequency.
9	If any of the parameters does not comply, stop the machine, rectify and recheck the parameters for compliance. Fill process interruption sheet.
10	Record the shipper weight for individual shipper.
11	After completing the batch packing, fill Finished Goods Transfer Note and transfer the batch to Finished goods store.

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				

UNCONTROLLED COPY

		Safetab Life Science BATCH PACKING RECORD		Page 36 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

6.1 LINE CLEARANCE CHECKS - PRIMARY PACKING:

S. No	Checks	DAY	Observation			
			Prod.		IPQA	
			I	II	I	II
1	Check whether area status board is updated for "CLEANED" status.		✓	✓	✓	✓
2	Left over Printed blister foil and Base foil/PVC film rolls of previous batch/product if any is removed from machine and cubicle.		✓	✓	✓	✓
3	Scrap (empty packs, trims) from rejection box below machine is removed.		✓	✓	✓	✓
4	Scrap of defoiled packs is removed from the defoiling machine.		✓	✓	✓	✓
5	Check whether empty containers of previous product/ batch are removed from the cubicle.		✓	✓	✓	✓
6	Check for the absence of Capsules/tablets of previous batch/ product from hopper, below vibrator and feeding box/channel.		✓	✓	✓	✓
7	Check for the absence of Capsules/tablets of previous batch/ product at punching tool and area below machine and machine conveyor belt.		✓	✓	✓	✓
8	Check for the absence of Capsules/tablets of previous batch/ product from collection box of defoiling machine.		✓	✓	✓	✓
9	Check for the absence of Capsules/tablets of previous batch/product from rejected tablets box.		✓	✓	✓	✓
10	Check for cleaning of hopper, feeding channel, vibrator plate		✓	✓	✓	✓
11	Check for cleaning of area below machine, backside of machine by opening the machine guards.		✓	✓	✓	✓
12	Check for removal of previous product stereos from the roller.		✓	✓	✓	✓
13	Check for cleaning of defoiling machine and its collection box.		✓	✓	✓	✓
14	Check the working desk for absence of container label, dispensing labels, other labels and document of previous batch/product.		✓	✓	✓	✓
15	Check whether the cubicle (floor, wall, ceiling, door, windows) is clean.		✓	✓	✓	✓
16	Check whether the daily records are updated and signed.		✓	✓	✓	✓
17	In case of materials / Intermediate products transferred between Block A and Block B, ensure that the containers received are clean and free from contamination.		✓	✓	✓	✓
18	Ensure that for previous product stereos removed from the embossing unit.		✓	✓	✓	✓

Note: Tick mark (✓) if Complies and (X) if Does not comply

Checked by (production) sign/date: *P. S. 12/2/25* (DAY I) Verified by (IPQA) sign/date: *12/2/25* (DAY I)

Checked by (production) sign/date: *12/2/25* (DAY II) Verified by (IPQA) sign/date: *12/2/25* (DAY II)

Equipment	Strip packing machine		AREA
Equipment ID. No.	ST/PRSR/002		
Previous Product	Litacold Flu		
Batch No.	G17250103		
Cleaned by	<i>Navin</i>		
Checked by (Production):	<i>P. S.</i>		
Verified by (IPQA):	<i>12/2/25</i>		

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with	<i>P. S.</i>	<i>11/2/25</i>	<i>12/2/25</i>

CONTROLLED COPY

UNCONTROLLED COPY

Safetab Life Science					
BATCH PACKING RECORD					
Page 37 of 109					
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	125000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

6.1A LINE CLEARANCE CHECKS - PRIMARY PACKING:

S. No	Checks	DAY	Observation			
			Prod		IPQA	
			III	IV	III	IV
1	Check whether area status board is updated for "CLEANED" status.		✓	✓	✓	✓
2	Left over Printed blister foil and Base foil/PVC film rolls of previous batch/product if any is removed from machine and cubicle.		✓	✓	✓	✓
3	Scrap (empty packs, trims) from rejection box below machine is removed.		✓	✓	✓	✓
4	Scrap of defoiled packs is removed from the defoiling machine.		✓	✓	✓	✓
5	Check whether empty containers of previous product/ batch are removed from the cubicle.		✓	✓	✓	✓
6	Check for the absence of Capsules/tablets of previous batch/ product from hopper, below vibrator and feeding box/channel.		✓	✓	✓	✓
7	Check for the absence of Capsules/tablets of previous batch/ product at punching tool and area below machine and machine conveyor belt.		✓	✓	✓	✓
8	Check for the absence of Capsules/tablets of previous batch/ product from collection box of defoiling machine.		✓	✓	✓	✓
9	Check for the absence of Capsules/tablets of previous batch/product from rejected tablets box.		✓	✓	✓	✓
10	Check for cleaning of hopper, feeding channel, vibrator plate		✓	✓	✓	✓
11	Check for cleaning of area below machine, backside of machine by opening the machine guards.		✓	✓	✓	✓
12	Check for removal of previous product stereos from the roller.		✓	✓	✓	✓
13	Check for cleaning of defoiling machine and its collection box.		✓	✓	✓	✓
14	Check the working desk for absence of container label, dispensing labels, other labels and document of previous batch/product.		✓	✓	✓	✓
15	Check whether the cubicle (floor, wall, ceiling, door, windows) is clean.		✓	✓	✓	✓
16	Check whether the daily records are updated and signed.		✓	✓	✓	✓
17	In case of materials / Intermediate products transferred between Block A and Block B, ensure that the containers received are clean and free from contamination.		✓	✓	✓	✓
18	Ensure that for previous product stereos removed from the embossing unit.		✓	✓	✓	✓

Note: Tick mark (✓) if Complies and (X) if Does not comply

Checked by (production) sign/date: 14/02/25 (DAY I) Verified by (IPQA) sign/date: 14/02/25 (DAY I)

Checked by (production) sign/date: 14/02/25 (DAY II) Verified by (IPQA) sign/date: 14/02/25 (DAY II)

Equipment	Strip packing machine	AREA
Equipment ID No.		
Previous Product		
Batch No.		
Cleaned by		
Checked by (Production):		
Verified by (IPQA):		

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with	<u>11/05/24</u>	<u>11/05/24</u>	<u>11/05/24</u>	
UNCONTROLLED COPY				



Safetab Life Science
BATCH PACKING RECORD

Page 38 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Pseudoephedrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

6.1B LINE CLEARANCE CHECKS - PRIMARY PACKING:

S. No.	Checks	Observation				
		DAY	Prod.	IPQA		
		V	VI	V	VI	
1	Check whether area status board is updated for "CLEANED" status.					
2	Left over Printed blister foil and Base foil/PVC film rolls of previous batch/product if any is removed from machine and cubicle.					
3	Scrap (empty packs, trims) from rejection box below machine is removed.					
4	Scrap of defoiled packs is removed from the defoiling machine.					
5	Check whether empty containers of previous product/ batch are removed from the cubicle.					
6	Check for the absence of Capsules/tablets of previous batch/ product from hopper, below vibrator and feeding box/channel.					
7	Check for the absence of Capsules/tablets of previous batch/ product at punching tool and area below machine and machine conveyor belt.					
8	Check for the absence of Capsules/tablets of previous batch/ product from collection box of defoiling machine.					
9	Check for the absence of Capsules/tablets of previous batch/product from rejected tablets box.					
10	Check for cleaning of hopper, feeding channel, vibrator plate					
11	Check for cleaning of area below machine, backside of machine by opening the machine guards.					
12	Check for removal of previous product stereotypes from the roller.					
13	Check for cleaning of defoiling machine and its collection box.					
14	Check the working desk for absence of container label, dispensing labels, other labels and document of previous batch/product.					
15	Check whether the cubicle (floor, wall, ceiling, door, windows) is clean.					
16	Check whether the daily records are updated and signed.					
17	In case of materials / Intermediate products transferred between Block A and Block B, ensure that the containers received are clean and free from contamination.					
18	Ensure that for previous product stereotypes removed from the embossing unit.					

Note: Tick mark (√) if Complies and (X) if Does not comply

Checked by (production) sign/date:.....(DAY I) Verified by (IPQA) sign/date:.....(DAY I)

Checked by (production) sign/date :.....(DAY II) Verified by (IPQA) sign/date:..... (DAY II)

Equipment	Strip packing machine	AREA 14/02/25
Equipment ID. No.		
Previous Product		
Batch No.		
Cleaned by		
Checked by (Production):		
Verified by (IPQA):		

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science

BATCH PACKING RECORD

Page 39 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	9/8	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

6.2 ENVIRONMENTAL CONDITIONS:	Observation					
	Date:	Date:	Date:	Date:	Date:	Date:
Temperature: 23°C to 27°C	12/02/25	13/02/25	14/02/25			
Relative Humidity: 45% to 55 %	24.2°C	24.2°C	24.3°C			
Differential Pressure NLT 5pascal or NLT 0.5 mm of water	53.1	53.1	53.1			
Production Sign	0.6m	0.6m	8			
IPQA Sign	Pri	Pri	Pri			
Primary Packing started at: 07:15		Completed at: 15:00				
6.3 APPROVAL FOR PRIMARY PACKING:		IPQA Sign/Date: 14/02/25				
Packing Sign/Date: Pri 12/02/25						


6.4 RECEIPT OF PRODUCT AND MATERIALS:

Check the following before start of packing activity and record the observation

Check Points	Date: 12/02/25		
	Observation	Done by (Packing)	Verified By (IPQA)
Temperature observed on display (Limit: 23°C - 27°C).	24.3°C	Pri	NO1
Relative Humidity observed on display (Limit: 45% - 55%).	53.1	Pri	NO1
Bulk product containers are of the required product /final analytical request in present.	✓	Pri	NO2
Status label of machine / line is for same product.	✓	Pri	NO1
All secondary packing materials are of the required product.	✓	Pri	NO1

6.5 VERIFICATION OF PRODUCT:

	Observations	Checked By/date (Packing)	Observations	Verified By/date (IPQA)
1. Average Weight	635 mg	Pri 12/02/25	635 mg	NO1 10/02/25
2. Quantity	779.260	Pri 12/02/25	779.260	NO1 10/02/25
3. Number of containers	23 no's	Pri 12/02/25	23 no's	NO1 10/02/25
If any Remarks:	NIL			

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				

UNCONTROLLED COPY

CONTROLLED COPY

Safetab Life Science
BATCH PACKING RECORD

Page 40 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125000
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

7.0 RECEIPT WEIGHING RECORD OF BULK PRODUCT:

Container No	Gross weight in kg	Tare weight in kg	Net weight in kg	Container No	Gross weight in kg	Tare weight in kg	Net weight in kg
1	37.70	2.60	35.10	36			
2	37.10	2.60	34.50	37			
3	35.50	2.50	33.00	38			
4	34.80	2.60	32.20	39			
5	34.00	2.60	31.40	40			
6	33.20	2.60	30.60	41			
7	32.70	2.60	30.10	42			
8	32.30	2.60	29.70	43			
9	33.50	2.60	30.90	44			
10	36.90	2.60	34.30	45			
11	37.70	2.20	35.50	46			
12	42.10	2.00	40.10	47			
13	42.15	2.40	39.75	48			
14	34.48	2.20	32.28	49			
15	41.75	2.60	39.15	50			
16	36.62	2.50	34.12	51			
17	33.68	2.50	31.18	52			
18	36.02	2.50	33.52	53			
19	42.72	2.50	40.22	54			
20	31.88	2.50	29.38	55			
21	33.54	2.50	31.04	56			
22	45.85	2.50	43.35	57			
23	29.45	2.50	26.95	58			
24				59			
25				60			
26				61			
27				62			
28				63			
29				64			
30				65			
31				66			
32				67			
33				68			
34				69			
35				70			
Total (kg)				Total (kg)			
				Grand total (kg)			
				779.260			
Checked By/date (production):				Checked by/date (IPQA):			
B. 12/02/25				N. 12/02/25			

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science



BATCH PACKING RECORD

Page 41 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

7.0A RECEIPT WEIGHING RECORD OF BULK PRODUCT:

Container No	Gross weight in kg	Tare weight in kg	Net weight in kg	Container No	Gross weight in kg	Tare weight in kg	Net weight in kg
1				36			
2				37			
3				38			
4				39			
5				40			
6				41			
7				42			
8				43			
9				44			
10				45			
11				46			
12				47			
13				48			
14				49			
15				50			
16				51			
17				52			
18				53			
19				54			
20				55			
21				56			
22				57			
23				58			
24				59			
25				60			
26				61			
27				62			
28				63			
29				64			
30				65			
31				66			
32				67			
33				68			
34				69			
35				70			
Total (kg)				Total (kg)			
				Grand total (kg)			
Checked By/date (production):..				Verified by/date (IPQA):..			

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	



Safetab Life Science BATCH PACKING RECORD

Page 42 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

7.1 STRIP PACKING MACHINE PARAMETERS:

Parameters	Standard and Limits	Observation			
		Date: 12/02/25	Date: 12/02/25	Date: 14/02/25	Date:
Temperature of Sealing Roller	100° C to 120° C	118° C	118° C	118° C	
Stereo impression quality	OK / NOT OK (Should be legible)	OK	OK	OK	
Machine code	10 to 45 cycle/minute	38	38	38	
Overprinting details:	Batch No.:	G17250103	G17250103	G17250103	
	MFG. Date:	01/2025	01/2025	01/2025	
	EXP. Date:	12/2027	12/2027	12/2027	
	Code number:	Pow Dose 21144378	Pow Dose 21144378	Pow Dose 21144378	
Leak test	OK / Not OK (No Strip should fail in leak test)	OK	OK	OK	
Checked by Packing Sign / Date		Pi 12/02/25	Pi 13/02/25	8 14/02/25	
Checked by IPQA Sign / Date		Pi 12/02/25	Pi 12/02/25	SP 14/02/25	Pi 14/02/25

NOTE: Check the over printed foil for correction of batch coding details. Attach a specimen foil with BPR duly signed by production and IPQA after checking.

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	


		Safetab Life Science BATCH PACKING RECORD		Page 43 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	973	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

8.0 PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:	12/02/25															
Machine Operator Name:	Sudha Lakshmi															
EVERY ONE HOUR																
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q
Time of inspection:	07:20	08:25	08:30	08:35	09:10	09:15	10:15	10:20	11:20	11:25	12:15	12:20	13:20	13:25	14:10	14:15
Sealing temperature. 100°C to 120°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C
Leak test on strip pack : (9 strips per check)	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Knurling quality:	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Horizontal cutting quality:	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Vertical cutting quality:	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Quality of over printing:	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Code no:	Pow Drugs 21144348	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Batch No.	G17250103	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Mfg. Date	01/2025	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Exp. Date	12/2027	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Machine speed (10-45 Cycles per minute)	38	38	38	38	38	38	38	38	38	38	38	38	38	38	38	38
Packing Supervisor / IPQA Signature:	Pj	Pj	Pj	Pj	Pj	Pj	Pj	Pj	Pj	Pj	Pj	Pj	Pj	Pj	Pj	Pj

• If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.

- OK / Not OK mention to particular process.
- P - Production
- Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	



Safetab Life Science BATCH PACKING RECORD

Page 44 of 109


Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

8.0A PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:		12/02/25								12/02/25								12/02/25	
Machine Operator Name:		Lalitha								Lalitha								Susha	
EVERY ONE HOUR																			
Checked by:		P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q		
Time of inspection:		15:16	15:21	16:25	16:30	17:28	17:33	18:36	18:41	20:15	20:21	21:20	21:25	06:30	06:35	08:00	08:05		
Sealing temperature, 100°C to 120°C		118°C	118°C	117°C	117°C	117°C	117°C	117°C	117°C	118°C	117°C	118°C	118°C	118°C	118°C	118°C	118°C		
Leak test on strip pack : (9 strips per check)		ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok		
Knurling quality:		ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok		
Horizontal cutting quality:		ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok		
Vertical cutting quality:		ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok		
Quality of over printing:		ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok		
Code no:	PON DRUMS 21144378	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok		
Batch No.	G17250103	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok		
Mfg. Date	01/2025	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok		
Exp. Date	12/2027	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok		
Machine speed (10-45 Cycles per minute)		38	38	38	38	38	38	38	38	38	38	38	38	38	38	38	38		
Packing Supervisor / IPQA Signature:		S	S	S	S	S	S	S	S	S	S	S	S	S	S	S	S		

• If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.

- OK / Not OK mention to particular process.
- P - Production
- Q - Quality Assurance.

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			
UNCONTROLLED COPY		CONTROLLED COPY	



Safetab Life Science BATCH PACKING RECORD

Page 45 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Pseudoephedrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

8.0B PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:		13/02/25								13/02/25							
Machine Operator Name:		Sudha								Lalitha							
EVERY ONE HOUR																	
Checked by:		P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q
Time of inspection:		09:10 11:00	09:15	10:08	10:10	11:15	11:20	12:10	12:15	13:15	13:20	14:20	14:20	15:25	15:30	16:30	16:30
Sealing temperature, 100°C to 120°C		118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C
Leak test on strip pack : (9 strips per check)		OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Knurling quality:		OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Horizontal cutting quality:		OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Vertical cutting quality:		OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Quality of over printing:		OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Code no:	Pond Drugs 21144378	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Batch No.	G17250103	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Mfg. Date	01/2025	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Exp. Date	12/2027	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
Machine speed (10-45 Cycles per minute)		38	30	30	30	30	30	30	30	30	30	38	38	38	38	38	38
Packing Supervisor / IPQA Signature:		81	81	81	81	81	81	81	81	81	81	81	81	81	81	81	81

* If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.

- OK / Not OK mention to particular process.
- P - Production
- Q - Quality Assurance

* 13/02/25

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

		Safetab Life Science BATCH PACKING RECORD		Page 46 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

8.0C PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:	13/02/25								14/02/25							
Machine Operator Name:	Lalitha								Sudha							
EVERY ONE HOUR																
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q
Time of inspection:	17:35	17:45	18:40	18:45	20:10	20:15	21:15	21:22	06:30	6:35	07:35	7:40	08:25	8:30	09:30	09:35
Sealing temperature, 100°C to 120°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C
Leak test on strip pack : (9 strips per check)	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok
Knurling quality:	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok
Horizontal cutting quality:	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok
Vertical cutting quality:	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok
Quality of over printing:	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok
Code no:	PON DRMS 21144878	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok
Batch No.	G17250103	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok
Mfg. Date	01/2025	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok
Exp. Date	12/2027	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok	ok
Machine speed (10-45 Cycles per minute)	38	38	38	38	38	38	38	38	38	38	38	38	38	38	38	38
Packing Supervisor / IPQA Signature:																

• If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.

- OK / Not OK mention to particular process.
- P - Production
- Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY				



Safetab Life Science

BATCH PACKING RECORD

Page 47 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

8.0D PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:	14/02/25	14/02/25														
Machine Operator Name:	Sudha	Lakshmi														
EVERY ONE HOUR																
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q
Time of inspection:	10:25	10:30	11:30	11:35	12:35	12:40	13:35	13:40	14:35	14:40						
Sealing temperature, 100°C to 120°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C	118°C						
Leak test on strip pack : (9 strips per check)	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK						
Knurling quality:	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK						
Horizontal cutting quality:	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK						
Vertical cutting quality:	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK						
Quality of over printing:	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK						
Code no:	PN/DRMS 21144378	OK	OK	OK	OK	OK	OK	OK	OK	OK						
Batch No.	G17250103	OK	OK	OK	OK	OK	OK	OK	OK	OK						
Mfg. Date	01/2025	OK	OK	OK	OK	OK	OK	OK	OK	OK						
Exp. Date	12/2027	OK	OK	OK	OK	OK	OK	OK	OK	OK						
Machine speed (10-45 Cycles per minute)	38	38	38	38	38	38	38	38	38	38						
Packing Supervisor / IPQA Signature:	S	S	S	S	S	S	S	S	S	S						

If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.

- OK / Not OK mention to particular process.
- P - Production
- Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY				


		Safetab Life Science BATCH PACKING RECORD		Page 48 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

8.0E PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:																
Machine Operator Name:																
EVERY ONE HOUR																
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q
Time of inspection:																
Sealing temperature, 100°C to 120°C																
Leak test on strip pack : (9 strips per check)																
Knurling quality:																
Horizontal cutting quality:																
Vertical cutting quality:																
Quality of over printing:																
Code no:																
Batch No.																
Mfg. Date																
Exp.Date																
Machine speed (10-45 Cycles per minute)																
Packing Supervisor / IPQA Signature:																

• If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.

- OK / Not OK mention to particular process.
- P - Production
- Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	



Safetab Life Science

BATCH PACKING RECORD

Page 49 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1,250,000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

8.0F PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:																	
Machine Operator Name:																	
EVERY ONE HOUR																	
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	
Time of inspection:																	
Sealing temperature: 100°C to 120°C																	
Leak test on strip pack : (9 strips per check)																	
Knurling quality:																	
Horizontal cutting quality:																	
Vertical cutting quality:																	
Quality of over printing:																	
Code no:																	
Batch No.																	
Mfg. Date																	
Exp. Date																	
Machine speed (10-45 Cycles per minute)																	
Packing Supervisor / IPQA Signature:																	

• If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.

- OK / Not OK mention to particular process,
P - Production
Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science
BATCH PACKING RECORD

Page 50 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

8.00 PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:																
Machine Operator Name:																
EVERY ONE HOUR																
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q
Time of inspection:																
Sealing temperature. 100°C to 120°C																
Leak test on strip pack : (9 strips per check)																
Knurling quality:																
Horizontal cutting quality:																
Vertical cutting quality:																
Quality of over printing:																
Code no:																
Batch No:																
Mfg. Date																
Exp.Date																
Machine speed (10-45 Cycles per minute)																
Packing Supervisor / IPQA Signature:																

- If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.
- OK / Not OK mention to particular process.
- P - Production
- Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

 Safetab Life Science BATCH PACKING RECORD		Page 51 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	078
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3
Batch No.:	G17250103	Mfg. Date:	01/2025
		Exp. Date:	12/2027
		Batch Size:	175000 TABLETS
		Supersedes MPR No:	ST/MPR/005/R2

8.0H PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:																
Machine Operator Name:																
EVERY ONE HOUR																
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q
Time of inspection:																
Sealing temperature: 100°C to 120°C																
Leak test on strip pack : (9 strips per check)																
Knurling quality:																
Horizontal cutting quality:																
Vertical cutting quality:																
Quality of over printing:																
Code no:																
Batch No.																
Mfg. Date																
Exp. Date																
Machine speed (10-45 Cycles per minute)																
Packing Supervisor / IPQA Signature:																

- If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.
- OK / Not OK mention to particular process.
- P – Production
- Q – Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			
UNCONTROLLED COPY		CONTROLLED COPY	



Safetab Life Science

BATCH PACKING RECORD

Page 52 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

8.01 PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:																
Machine Operator Name:																
EVERY ONE HOUR																
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q
Time of inspection:																
Sealing temperature, 100°C to 120°C																
Leak test on strip pack : (9 strips per check)																
Knurling quality:																
Horizontal cutting quality:																
Vertical cutting quality:																
Quality of over printing:																
Code no:																
Batch No.																
Mfg. Date																
Exp.Date																
Machine speed (10-45 Cycles per minute)																
Packing Supervisor / IPQA Signature:																

• If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.

• OK / Not OK mention to particular process.

P = Production

Q = Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science BATCH PACKING RECORD

Page 53 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

8.03 PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:																
Machine Operator Name:																
EVERY ONE HOUR																
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q
Time of inspection:																
Sealing temperature, 100°C to 120°C																
Leak test on strip pack : (9 strips per check)																
Knurling quality:																
Horizontal cutting quality:																
Vertical cutting quality:																
Quality of over printing:																
Code no:																
Batch No.																
Mfg. Date																
Exp. Date																
Machine speed (10-45 Cycles per minute)																
Packing Supervisor / IPQA Signature:																

- If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.
- OK / Not OK mention to particular process.
- P - Production
- Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	CONTROLLED COPY
Signature with				
UNCONTROLLED COPY				



Safetab Life Science
BATCH PACKING RECORD

Page 54 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

8.0K PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:		
Machine Operator Name:		
EVERY ONE HOUR		
Checked by:	P	Q
Time of inspection:		
Sealing temperature: 100°C to 120°C		
Leak test on strip pack : (9 strips per check)		
Knurling quality:		
Horizontal cutting quality:		
Vertical cutting quality:		
Quality of over printing:		
Code no:		
Batch No.		
Mfg. Date		
Exp.Date		
Machine speed (10-45 Cycles per minute)		
Packing Supervisor / IPQA Signature:		

- If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.

- OK / Not OK mention to particular process.

P – Production

Q – Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date	 11/01/24	 11/01/24	 11/01/24

UNCONTROLLED COPY


CONTROLLED COPY

 Safetab Life Science BATCH PACKING RECORD						Page 55 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets		
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2		
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027		

8.01 PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:																
Machine Operator Name:																
EVERY ONE HOUR																
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q
Time of inspection:																
Sealing temperature, 100°C to 120°C																
Leak test on strip pack : (9 strips per check)																
Knurling quality:																
Horizontal cutting quality:																
Vertical cutting quality:																
Quality of over printing:																
Code no:																
Batch No.																
Mfg. Date																
Exp. Date																
Machine speed (10-45 Cycles per minute)																
Packing Supervisor / IPQA Signature:																

- If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.
- OK / Not OK mention to particular process.
- P - Production
- Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				

UNCONTROLLED COPY

11/05/24

 Safetab Life Science BATCH PACKING RECORD		Page 56 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3
Batch No.:	G17250103	Mfg. Date:	01/2025
		Batch Size	1250000 Tablets
		Supersedes MPR No:	ST/MPR/005/R2
		Exp. Date:	12/2027

8.0M PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:																	
Machine Operator Name:																	
EVERY ONE HOUR																	
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	
Time of inspection:																	
Sealing temperature. 100°C to 120°C																	
Leak test on strip pack : (9 strips per check)																	
Knurling quality:																	
Horizontal cutting quality:																	
Vertical cutting quality:																	
Quality of over printing:																	
Code no:																	
Batch No.																	
Mfg. Date																	
Exp. Date																	
Machine speed (10-45 Cycles per minute)																	
Packing Supervisor / IPQA Signature:																	

- If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.
- OK / Not OK mention to particular process.
 - P – Production
 - Q – Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	


Safetab Life Science

BATCH PACKING RECORD

Page 57 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

8. ON PRIMARY PACKING LINE IN-PROCESS CHECKS:

Date:																
Machine Operator Name:																
EVERY ONE HOUR																
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q
Time of inspection:																
Sealing temperature. 100°C to 120°C																
Leak test on strip pack : (9 strips per check)																
Knurling quality:																
Horizontal cutting quality:																
Vertical cutting quality:																
Quality of over printing:																
Code no:																
Batch No.																
Mfg. Date																
Exp. Date																
Machine speed (10-45 Cycles per minute)																
Packing Supervisor / IPQA Signature:																

- If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.
- OK / Not OK mention to particular process.
- P - Production
- Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				

UNCONTROLLED COPY

Safetab Life Science BATCH PACKING RECORD

Page 58 of 100

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	100000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

8.0- O P R I M A R Y P A C K I N G L I N E I N - P R O C E S S C H E C K S :

Date:																
Machine Operator Name:																
EVERY ONE HOUR																
Checked by:	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q	P	Q
Time of inspection:																
Sealing temperature: 100°C to 120°C																
Leak test on strip pack : (9 strips per check)																
Knurling quality:																
Horizontal cutting quality:																
Vertical cutting quality:																
Quality of over printing:																
Code no:																
Batch No.																
Mfg. Date																
Exp. Date																
Machine speed (10-45 Cycles per minute)	P-1 14/02/25 ✓															
Packing Supervisor / IPQA Signature:																

- If vacuum test fails perform every 15 minutes for next one hour and all the shippers to be checked for leak test from the last Successful leak test. All the shippers packed from the time of last successful leak test shall be segregated and random checked by QA and will be released after approval from IPQA.
- OK / Not OK mention to particular process.
- P - Production
- Q - Quality Assurance

9.0 DEFOILING: Inspection of De-foiled tablets:

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				
	11/02/25		11/02/25	

UNCONTROLLED COPY



Safetab Life Science

BATCH PACKING RECORD

Page 59 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

Date	Inspection started at	Inspection completed at	Weight of tablets taken for inspection	Inspected by	Quantity of good tablets		Packed in shipper No.	Weight of rejected tablets	Checked by
					Kg	No's			
14/02/25	13:30	13:50	0.900	Amulya	0.700	1102	204	0.200	By
Total			0.900	-	0.700	1102	204	0.200	By
Total quantity of rejected tablets destroyed: 314 Nos.				Done by/Date (Packing)			Checked By/Date (IPQA)		
				By 14/02/25			(Su) 14/02/25		

10.0 LINE CLEARANCE CHECKS - SECONDARY PACKING:

SECONDARY PACKING		Observation						
S.No	Checks	Production						
		DAY	I	II	III	I	II	III
1	Inspect the packaging area and the packing line for cleanliness. The area and packing line should be visibly clean.		✓	✓	✓	✓	✓	✓
2	Batch Packing record checked.		✓	✓	✓	✓	✓	✓
3	Check the status labels of the area/packing line and ensure their correctness with respect to product name, batch number, date, signature of production personnel.		✓	✓	✓	✓	✓	✓
4	Removal of Previous Secondary Packing Material.		✓	✓	✓	✓	✓	✓
5	Removal of Previous Documents.		✓	✓	✓	✓	✓	✓
6	Ensure that the balances are cleaned and it's calibration status.		✓	✓	✓	✓	✓	✓
7	Check whether the balance to be used are within calibration due date.		✓	✓	✓	✓	✓	✓
CONVEYOR BELT								
8	Cleanliness of the Conveyor Belt Checked		✓	✓	✓	✓	✓	✓
9	Cleanliness of guard checked		✓	✓	✓	✓	✓	✓
10	Check & ensure that the area is free from dosage units (strips, tablets etc.)/packaging components (labels, leaflets, cartons etc.) & documents of the previous batch/product.		✓	✓	✓	✓	✓	✓
11	Ensure that any tools/materials required for equipment maintenance/operation are kept properly under control.		✓	✓	✓	✓	✓	✓
12	Check and ensure that labeled bins are available for disposal of rejects of packing material.		✓	✓	✓	✓	✓	✓
Note: Tick mark (✓) if Complies and (X) if Does not comply.								
Checked by packing sign /date <u>8/12/25</u> Verified by IPQA sign /Date <u>9/12/2025</u>								

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date	By 11/05/24	By 11/05/24	By 11/05/24	
UNCONTROLLED COPY				

		Safetab Life Science BATCH PACKING RECORD		Page 60 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

10.0A LINE CLEARANCE CHECKS - SECONDARY PACKING:

S.No	Checks	Observation					
		Production			IPQA		
	DAY	IV	V	VI	IV	V	VI
1	Inspect the packaging area and the packing line for cleanliness. The area and packing line should be visibly clean.						
2	Batch Packing record checked.						
3	Check the status labels of the area/packing line and ensure their correctness with respect to product name, batch number, date, signature of production personnel.						
4	Removal of Previous Secondary Packing Material.						
5	Removal of Previous Documents.						
6	Ensure that the balances are cleaned and it's calibration status.						
7	Check whether the balance to be used are within calibration due date.						
CONVEYOR BELT							
8	Cleanliness of the Conveyor Belt Checked						
9	Cleanliness of guard checked						
10	Check & ensure that the area is free from dosage units (strips, tablets etc.)/packaging components (labels, leaflets, cartons etc.) & documents of the previous batch/product.						
11	Ensure that any tools/materials required for equipment maintenance/operation are kept properly under control.						
12	Check and ensure that labeled bins are available for disposal of rejects of packing material.						

Note: Tick mark (✓) if Complies and (X) if Does not comply.

Checked by packing sign /date: _____ Verified by IPQA sign /Date: _____

11.0 PACKING INSTRUCTIONS :(SALES)

Pack 4 tablets in one Strip and check the quality of sealing.

At the end of packing operation remove rejected tablets and destroy in presence of IPQA.

Pack 1 strip into the 1 catch cover. Pack such 20 catch cover and 1 leaflet in a carton. Pack 75 cartons into one 5-ply shipper. 3 cartons per row and 5 cartons per column and pack 5 carton layers per shipper. Seal the 5 ply shipper using tape and carry out the in process checks as listed in the batch record. Weight and record the weight of individual shipper.

Destroy the rejected printed packing materials at the end of the shift in presence of IPQA.

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				

UNCONTROLLED COPY



Safetab Life Science
BATCH PACKING RECORD

Page 61 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/P2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

12.0 WEIGHT CHECKS ON CARTON:

12.1. STRIP WEIGHT:

S.No	1	2	3	4	5	6	7	8	9	10
Weight in gm	3.42	3.44	3.38	3.39	3.45	3.44	3.46	3.37	3.38	3.49

Average of strips = $\frac{34.2}{10}$ = A $\frac{3.42}{10}$ (g)

12.2 LEAFLET WEIGHT:

S.No	1	2	3	4	5	6	7	8	9	10
Weight in gm	1.87	1.87	1.83	1.84	1.81	1.90	1.91	1.93	1.94	1.84

Average of empty leaflet = $\frac{18.6}{10}$ = B $\frac{1.86}{10}$ (g)

12.3 EMPTY CARTON WEIGHT:

S.No	1	2	3	4	5	6	7	8	9	10
Weight in gm	18.10	18.04	18.06	18.04	18.07	17.97	18.38	18.23	18.34	18.37

Average weight of empty carton = $\frac{181.6}{10}$ = C $\frac{18.16}{10}$ (g)

12.4 AVERAGE FILLED CARTON WEIGHT:

D = A + B + C

Average of filled carton = $\frac{1412.3}{10}$ = D $\frac{141.23}{10}$ (g)

I = $\pm 50\%$ average weight of one blister (A/2) = $\frac{1.71}{10}$ (g)

II = average weight of one leaflet (B) = $\frac{1.86}{10}$ (g)

The limit of carton weight (E) shall be fixed based on I or II which is less in weight.

LOWER LIMIT E = D - I or D - II = $\frac{139.52}{10}$ (g)

HIGHER LIMIT F = D + I or D + II = $\frac{142.94}{10}$ (g)

Checked By Production:

Sign & Date:

[Signature]
12/02/25

verified By IPQA:

sign & Date:

[Signature]
12/02/25

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date	<i>[Signature]</i> 12/02/25	<i>[Signature]</i> 11/02/24	<i>[Signature]</i> 12/02/25

UNCONTROLLED COPY

CONTROLLED COPY

Safetab Life Science
BATCH PACKING RECORD

Page 52 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

13.0 SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date/ Time	ACTIVITIES	DONE BY
12/02/25 07:30	Strip printing checking	Kowsalya D, Chithra, Nirmala, Valarmathi	12/02/25 14:06	Strip printing checking	Ganathi, Nivetha, Uimala, Ganavalli
	Carton folding & printing checking	Rajakumari, Kala, Sumala		Carton folding & printing checking	mmu, kala
	Checking of Carton (Brand Name)	Devi, Ravathi		Checking of Carton (Brand Name)	Hirsha, Sundhar
	Packing of 1 strips per catch cover	Ranjani, Nisasha, Senbagavalli		Packing of 1 strips per catch cover	Kabirya, Naitthis, Vani JA
	Packing of 20 catch cover per carton	Sathiyakala, Grestha		Packing of 20 catch cover per carton	Vijayalakshmi, Arisami
	Packing of leaflet per carton	Thaashini, Aadhilakshmi		Packing of leaflet per carton	Ambika, Parmani
	Cartons weighing	V. Kowsalya		Cartons weighing	A. Jayabari
	Carton put into Sply shipper	Kowsalya V		Carton put into Sply shipper	Jayabari
	Shipper weighing	Kowsalya V		Shipper weighing	Jayabari
	Shipper label pasting	Kowsalya V		Shipper label pasting	Jayabari
	Checked by Packing Supervisor sign/date	Praveen 12/02/25		Checked by Packing Supervisor sign/date	H 12/02/25

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				
UNCONTROLLED COPY				



Safetab Life Science BATCH PACKING RECORD

Page 63 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125 Right Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

13.0A SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date / Time	ACTIVITIES	DONE BY
13/02/25 01:30	Strip printing checking	Anitha.v Chithra Vakrathri, Sarasu	13/02/25 14:05	Strip printing checking	Ammu . Ambika Gomathi, Ganavalli
	Carton folding & printing checking	Nirmala . Thenmozhi Kowsalya		Carton folding & printing checking	A.Vaithishwari Laila . vijalakshmi
	Checking of Carton (Brand Name)	Ravathi, Magashani Krishnaveni		Checking of Carton (Brand Name)	Uma . Abirami Usha
	Packing of 1 strips per catch cover	Anuma . N Nivedha		Packing of 1 strips per catch cover	Sundhara
	Packing of 20 catch cover per carton	Tanasi Gowri . P		Packing of 20 catch cover per carton	Nivetha . nirmala Ponmani Sowmiya
	Packing of leaflet per carton	pavithra		Packing of leaflet per carton	Latha
	Cartons weighing	Anitha . R		Cartons weighing	A . Jayasri
	Carton put into Sply shipper	Anitha . R		Carton put into Sply shipper	A . Jayasri
	Shipper weighing	Anitha . R		Shipper weighing	A . Jayasri
	Shipper label pasting	Anitha . R		Shipper label pasting	Sayasri
	Checked by Packing Supervisor sign/date	P 13/02/25		Checked by Packing Supervisor sign/date	P 13/02/25

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY		CONTROLLED COPY		

Safetab Life Science
BATCH PACKING RECORD

Page 64 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

13.08 SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date / Time	ACTIVITIES	DONE BY
14/02/25 6:40	Strip printing checking	Chitra, Azuma, Soundarya, Selamathi	14/02/25 14:10	Strip printing checking	Ammu, Ambika, Granavalli, Nivasha
	Carton folding & printing checking	Abinaya		Carton folding & printing checking	Kaila, Vijalakshmi
	Checking of Carton (Brand Name)	Thennaschi, Nivasha		Checking of Carton (Brand Name)	Uma, Abiricani, Usha
	Packing of 1 strips per catch cover	Nivasha, Anita, Sujatha, Gowri, Ravathi		Packing of 1 strips per catch cover	Nimla, Ponnani, Sowmya, Kavasei
	Packing of 20 catch cover per carton	Hakshwari		Packing of 20 catch cover per carton	Sundharini
	Packing of leaflet per carton	Asha Bharathi		Packing of leaflet per carton	Latha
	Cartons weighing	B. Selvi		Cartons weighing	A. Jayasri
	Carton put into Sply shipper	B. Selvi		Carton put into Sply shipper	A. Jayasri
	Shipper weighing	B. Selvi		Shipper weighing	A. Jayasri
	Shipper label pasting			Shipper label pasting	A. Jayasri
	Checked by Packing Supervisor sign/date	14/02/25		Checked by Packing Supervisor sign/date	14/02/25

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				CONTROLLED COPY
				UNCONTROLLED COPY

Safetab Life Science
BATCH PACKING RECORD

Page 65 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

13.0C SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date/ Time	ACTIVITIES	DONE BY
	Strip printing checking			Strip printing checking	
	Carton folding & printing checking			Carton folding & printing checking	
	Checking of Carton (Brand Name)			Checking of Carton (Brand Name)	
	Packing of 1 strips per catch cover			Packing of 1 strips per catch cover	
	Packing of 20 catch cover per carton			Packing of 20 catch cover per carton	
	Packing of leaflet per carton			Packing of leaflet per carton	
	Cartons weighing			Cartons weighing	
	Carton put into Sply shipper			Carton put into Sply shipper	
	Shipper weighing			Shipper weighing	
	Shipper label pasting			Shipper label pasting	
	Checked by Packing Supervisor sign/date			Checked by Packing Supervisor sign/date	

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				CONTROLLED COPY
				UNCONTROLLED COPY



Safetab Life Science
BATCH PACKING RECORD

Page 66 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

13.00 SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date / Time	ACTIVITIES	DONE BY
	Strip printing checking			Strip printing checking	
	Carton folding & printing checking			Carton folding & printing checking	
	Checking of Carton (Brand Name)			Checking of Carton (Brand Name)	
	Packing of 1 strips per catch cover			Packing of 1 strips per catch cover	
	Packing of 20 catch cover per carton			Packing of 20 catch cover per carton	
	Packing of leaflet per carton			Packing of leaflet per carton	
	Cartons weighing			Cartons weighing	
	Carton put into 5ply shipper			Carton put into 5ply shipper	
	Shipper weighing			Shipper weighing	
	Shipper label pasting			Shipper label pasting	
	Checked by Packing Supervisor sign/date			Checked by Packing Supervisor sign/date	

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				

UNCONTROLLED COPY

 Safetab Life Science BATCH PACKING RECORD		Page 67 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3
Batch No.:	G17250103	Mfg.Date:	01/2025
		Batch Size:	1250000 Tablets
		Supersedes MPR No:	ST/MPR/005/R2
		Exp.Date:	12/2027

13.0E SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date / Time	ACTIVITIES	DONE BY
	Strip printing checking			Strip printing checking	
	Carton folding & printing checking			Carton folding & printing checking	
	Checking of Carton (Brand Name)			Checking of Carton (Brand Name)	
	Packing of 1 strips per catch cover			Packing of 1 strips per catch cover	
	Packing of 20 catch cover per carton			Packing of 20 catch cover per carton	
	Packing of leaflet per carton			Packing of leaflet per carton	
	Cartons weighing			Cartons weighing	
	Carton put into Sply shipper			Carton put into Sply shipper	
	Shipper weighing			Shipper weighing	14/12/25
	Shipper label pasting			Shipper label pasting	
	Checked by Packing Supervisor sign/date			Checked by Packing Supervisor sign/date	

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			
UNCONTROLLED COPY		CONTROLLED COPY	




Safetab Life Science
BATCH PACKING RECORD

Page 68 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

13. OF SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date/ Time	ACTIVITIES	DONE BY
	Strip printing checking			Strip printing checking	
	Carton folding & printing checking			Carton folding & printing checking	
	Checking of Carton (Brand Name)			Checking of Carton (Brand Name)	
	Packing of 1 strips per catch cover			Packing of 1 strips per catch cover	
	Packing of 20 catch cover per carton			Packing of 20 catch cover per carton	
	Packing of leaflet per carton			Packing of leaflet per carton	
	Cartons weighing			Cartons weighing	
	Carton put into 5ply shipper			Carton put into 5ply shipper	
	Shipper weighing			Shipper weighing	
	Shipper label pasting			Shipper label pasting	
	Checked by Packing Supervisor sign/date			Checked by Packing Supervisor sign/date	

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			

UNCONTROLLED COPY

CONTROLLED COPY

Safetab Life Science


BATCH PACKING RECORD

Page 69 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

13.0G SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date / Time	ACTIVITIES	DONE BY
	Strip printing checking			Strip printing checking	
	Carton folding & printing checking			Carton folding & printing checking	
	Checking of Carton (Brand Name)			Checking of Carton (Brand Name)	
	Packing of 1 strips per catch cover			Packing of 1 strips per catch cover	
	Packing of 20 catch cover per carton			Packing of 20 catch cover per carton	
	Packing of leaflet per carton			Packing of leaflet per carton	
	Cartons weighing			Cartons weighing	
	Carton put into Sply shipper			Carton put into Sply shipper	
	Shipper weighing			Shipper weighing	14/12/25
	Shipper label pasting			Shipper label pasting	
	Checked by Packing Supervisor sign/date			Checked by Packing Supervisor sign/date	

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			
UNCONTROLLED COPY		CONTROLLED COPY	



Safetab Life Science
BATCH PACKING RECORD

Page 70 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

13.0H SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date / Time	ACTIVITIES	DONE BY
	Strip printing checking			Strip printing checking	
	Carton folding & printing checking			Carton folding & printing checking	
	Checking of Carton (Brand Name)			Checking of Carton (Brand Name)	
	Packing of 1 strips per catch cover			Packing of 1 strips per catch cover	
	Packing of 20 catch cover per carton			Packing of 20 catch cover per carton	
	Packing of leaflet per carton			Packing of leaflet per carton	
	Cartons weighing			Cartons weighing	
	Carton put into Sply shipper			Carton put into Sply shipper	
	Shipper weighing			Shipper weighing	14/02/25
	Shipper label pasting			Shipper label pasting	
	Checked by Packing Supervisor sign/date			Checked by Packing Supervisor sign/date	

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				CONTROLLED COPY
				UNCONTROLLED COPY

Safetab Life Science					
BATCH PACKING RECORD					
Page 71 of 109					
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

13.01 SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date / Time	ACTIVITIES	DONE BY
	Strip printing checking			Strip printing checking	
	Carton folding & printing checking			Carton folding & printing checking	
	Checking of Carton (Brand Name)			Checking of Carton (Brand Name)	
	Packing of 1 strips per catch cover			Packing of 1 strips per catch cover	
	Packing of 20 catch cover per carton			Packing of 20 catch cover per carton	
	Packing of leaflet per carton			Packing of leaflet per carton	
	Cartons weighing			Cartons weighing	
	Carton put into 5ply shipper			Carton put into 5ply shipper	
	Shipper weighing			Shipper weighing	14/02/25
	Shipper label pasting			Shipper label pasting	
	Checked by Packing Supervisor sign/date			Checked by Packing Supervisor sign/date	

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			
UNCONTROLLED COPY CONTROLLED COPY			

Safetab Life Science


BATCH PACKING RECORD

Page 72 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

13.03 SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date/ Time	ACTIVITIES	DONE BY
	Strip printing checking			Strip printing checking	
	Carton folding & printing checking			Carton folding & printing checking	
	Checking of Carton (Brand Name)			Checking of Carton (Brand Name)	
	Packing of 1 strips per catch cover			Packing of 1 strips per catch cover	
	Packing of 20 catch cover per carton			Packing of 20 catch cover per carton	
	Packing of leaflet per carton			Packing of leaflet per carton	
	Cartons weighing			Cartons weighing	
	Carton put into 5ply shipper			Carton put into 5ply shipper	
	Shipper weighing			Shipper weighing	11/05/24
	Shipper label pasting			Shipper label pasting	
	Checked by Packing Supervisor sign/date			Checked by Packing Supervisor sign/date	

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			

UNCONTROLLED COPY

CONTROLLED COPY



Safetab Life Science
BATCH PACKING RECORD

Page 73 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

13.0K SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date/ Time	ACTIVITIES	DONE BY
	Strip printing checking			Strip printing checking	
	Carton folding & printing checking			Carton folding & printing checking	
	Checking of Carton (Brand Name)			Checking of Carton (Brand Name)	
	Packing of 1 strips per catch cover			Packing of 1 strips per catch cover	
	Packing of 20 catch cover per carton			Packing of 20 catch cover per carton	
	Packing of leaflet per carton			Packing of leaflet per carton	
	Cartons weighing			Cartons weighing	
	Carton put into Sply shipper			Carton put into Sply shipper	
	Shipper weighing			Shipper weighing	me 19/10/25
	Shipper label pasting			Shipper label pasting	
	Checked by Packing Supervisor sign/date			Checked by Packing Supervisor sign/date	

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	CONTROLLED COPY
Signature with Date				
	UNCONTROLLED COPY			



Safetab Life Science
BATCH PACKING RECORD

Page 74 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

13.0L SECONDARY PACKING LINE CHECKS:

Date / Time	ACTIVITIES	DONE BY	Date / Time	ACTIVITIES	DONE BY
	Strip printing checking			Strip printing checking	
	Carton folding & printing checking			Carton folding & printing checking	
	Checking of Carton (Brand Name)			Checking of Carton (Brand Name)	
	Packing of 1 strips per catch cover			Packing of 1 strips per catch cover	
	Packing of 20 catch cover per carton			Packing of 20 catch cover per carton	
	Packing of leaflet per carton			Packing of leaflet per carton	
	Cartons weighing			Cartons weighing	
	Carton put into 5ply shipper			Carton put into 5ply shipper	
	Shipper weighing			Shipper weighing	
	Shipper label pasting			Shipper label pasting	
	Checked by Packing Supervisor sign/date			Checked by Packing Supervisor sign/date	

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				CONTROLLED COPY
				UNCONTROLLED COPY

Safetab Life Science

BATCH PACKING RECORD

Page 75 of 109


Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

14.0 SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One-Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date		12/02/25	12/02/25	12/02/25	12/02/25	12/02/25	12/02/25	12/02/25	12/02/25
Time		07:40	07:45	08:00	08:53	9:15	09:20	10:20	10:35
Parameter	Standards and Limits								
strip in catch cover	1X4's	OK	OK	OK	OK	OK	OK	OK	OK
No. of catch cover and leaf per carton	20X1X4's and 1 Leaf per carton	OK	OK	OK	OK	OK	OK	OK	OK
Carton checking (Brand Name)		OK	OK	OK	OK	OK	OK	OK	OK
Carton weight check (3 carton) Limit:		OK	OK	OK	OK	OK	OK	OK	OK
Print Quality	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No: G17250103								
	Mfg.Dt: 01/2025	OK	OK	OK	OK	OK	OK	OK	OK
	Exp.Dt: 12/2027								
Packing Supervisor/ IPQA Signature		Pi	W	Pi	W	Pi	W	Pi	W
SHIPPERS									
Shipper No.		01	01	02	3	05	06	07	07
Print Quality of Shipper Label	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No: G17250103								
	Mfg.Dt: 01/2025	OK	OK	OK	OK	OK	OK	OK	OK
	Exp.Dt: 12/2027								
	Qty: (75x20x1x4'S)	OK	OK	OK	OK	OK	OK	OK	OK
Shipper sealing quality	Uniform	OK	OK	OK	OK	OK	OK	OK	OK
Packing Supervisor/ IPQA Signature		Pi	W	Pi	W	Pi	W	Pi	W

• OK / Not OK mention to particular process.

P - Production Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

		Safetab Life Science BATCH PACKING RECORD		Page 76 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO.:	ST/MPR/005/R3	Supersedes MPR No.:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027




14.0A SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One-Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date		12/02/25	12/02/25	12/02/25	12/02/25	12/02/25	12/02/25	12/02/25	12/02/25
Time		11:25	11:30	12:20	12:30	13:15	13:28	14:15	14:20
Parameter	Standards and Limits								
strip in catch cover	1X4's	OK	OK	OK	OK	OK	OK	OK	OK
No. of catch cover per carton	20X1X4's and 1 Leaf per carton	OK	OK	OK	OK	OK	OK	OK	OK
Print Quality	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Carton checking (Brand Name)		OK	OK	OK	OK	OK	OK	OK	OK
Carton weight check (3 carton) Limit:		OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No: G17250103	OK							
	Mfg. Dt: 01/2025		OK	OK	OK	OK	OK	OK	OK
	Exp. Dt: 12/2027								
Packing Supervisor/ IPQA Signature		OK	OK	OK	OK	OK	OK	OK	OK
SHIPPERS									
Shipper No.		12	13	17	18	21	22	25	25
Print Quality of Shipper Label	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No: G17250103	OK							
	Mfg. Dt: 01/2025		OK	OK	OK	OK	OK	OK	OK
	Exp. Dt: 12/2027								
	Qty: (75x20x1x4'S)	OK	OK	OK	OK	OK	OK	OK	OK
Shipper sealing quality	Uniform	OK	OK	OK	OK	OK	OK	OK	OK
Packing Supervisor/ IPQA Signature		OK	OK	OK	OK	OK	OK	OK	OK

• OK / Not OK mention to particular process.

P - Production

Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				
UNCONTROLLED COPY				



Safetab Life Science

BATCH PACKING RECORD

Page 77 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	12x40x10 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

14.0B SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One-Hour

Checked By		P	Q	P	Q	P	Q	P	Q
Date		12/02/25	12/02/25	12/02/25	12/02/25	12/02/25	12/02/25	12/02/25	12/02/25
Time		15:20	15:26	16:26	16:31	17:32	17:38	18:40	18:43
Parameter	Standards and Limits								
strip in catch cover	1X4's	OK	OK	OK	OK	OK	OK	OK	OK
No. of catch cover per carton	20X1X4's and 1 Leaf per carton	OK	OK	OK	OK	OK	OK	OK	OK
Print Quality	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Carton checking (Brand Name)		OK	OK	OK	OK	OK	OK	OK	OK
Carton weight check (3 carton) Limit:		OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No:								
	Mfg.Dt:								
	Exp.Dt:								
Packing Supervisor/ IPQA Signature		ft	ft	ft	ft	ft	ft	ft	ft
SHIPPERS									
Shipper No.		26	27	28	29	30	31	32	33
Print Quality of Shipper Label	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No:								
	Mfg.Dt:								
	Exp.Dt:								
Shipper sealing quality	Uniform	OK	OK	OK	OK	OK	OK	OK	OK
Packing Supervisor/ IPQA Signature		ft	ft	ft	ft	ft	ft	ft	ft

• OK / Not OK mention to particular process.

P - Production

Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				
UNCONTROLLED COPY				



Safetab Life Science
BATCH PACKING RECORD

Page 78 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027


14.00 SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One-Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date		12/02/25	12/02/25	12/02/25	12/02/25	13/02/25	13/02/25	13/02/25	13/02/25
Time		20:20	20:25	21:25	21:30	06:25	06:40	07:40	07:45
Parameter	Standards and Limits								
strip in catch cover	1X4's	OK	OK	OK	OK	OK	OK	OK	OK
No. of catch cover per carton	20X1X4's and 1 Leaf per carton	OK	OK	OK	OK	OK	OK	OK	OK
Print Quality	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Carton checking (Brand Name)		OK	OK	OK	OK	OK	OK	OK	OK
Carton weight check (3 carton) Limit:		OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No:								
	G17250103								
	Mfg.Dt:								
Exp.Dt:	01/2025	OK	OK	OK	OK	OK		OK	OK
	12/2027								
Packing Supervisor/ IPQA Signature		ff	ff	ff	ff	ff	ff	ff	ff
SHIPPERS									
Shipper No.		34	38	43	49	44	45	48	48
Print Quality of Shipper Label	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No:								
	G17250103								
	Mfg.Dt:								
Exp.Dt:	01/2025	OK	OK	OK	OK	OK		OK	OK
	12/2027								
Qty: (75x20x1x4'S)		OK	OK	OK	OK	OK	OK	OK	OK
Shipper sealing quality	Uniform	OK	OK	OK	OK	OK	OK	OK	OK
Packing Supervisor/ IPQA Signature		ff	ff	ff	ff	ff	ff	ff	ff

• OK / Not OK mention to particular process.

P - Production

Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			
UNCONTROLLED COPY		CONTROLLED COPY	



Safetab Life Science
BATCH PACKING RECORD

Page 79 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	50000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027


14.00 SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One Hour

Checked By		P	Q	P	Q	P	Q	P	Q
Date		12/2/25	10/10	12/2/25	13/01/25	13/02/25	13/02/25	13/02/25	13/02/25
Time		07:05	9:10	10:10	10:12	11:05	11:19	12:20	12:25
Parameter	Standards and Limits								
strip in catch cover	1X4's	OK	OK	OK	OK	OK	OK	OK	OK
No. of catch cover per carton	20X1X4's and 1 Leaf per carton	OK	OK	OK	OK	OK	OK	OK	OK
Print Quality	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Carton checking (Brand Name)		OK	OK	OK	OK	OK	OK	OK	OK
Carton weight check (3 carton) Limit:		OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No:	OK							
	Mfg. Dt:			OK		OK	OK	OK	
	Exp. Dt:		OK		OK	OK	OK	OK	OK
Packing Supervisor/ IPQA Signature		81	81	81	81	81	81	81	81
SHIPPERS									
Shipper No.		53	54	60	61	66	67	75	75
Print Quality of Shipper Label	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No:								
	Mfg. Dt:	OK	OK	OK	OK	OK	OK	OK	OK
	Exp. Dt:	OK	OK	OK	OK	OK	OK	OK	OK
Qty: (75x20x1x4'S)		OK	OK	OK	OK	OK	OK	OK	OK
Shipper sealing quality	Uniform	OK	OK	OK	OK	OK	OK	OK	OK
Packing Supervisor/ IPQA Signature		81	81	81	81	81	81	81	81

• OK / Not OK mention to particular process.

P - Production

Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science

BATCH PACKING RECORD

Page 80 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

14.0E SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date		13/02/25	13/02/25	13/02/25	13/02/25	13/02/25	13/02/25	13/02/25	13/02/25
Time		13:15	13:20	14:20	14:26	15:25	15:31	16:30	16:36
Parameter	Standards and Limits								
strip in catch cover	1X4's	OK	OK	OK	OK	OK	OK	OK	OK
No. of catch cover per carton	20X1X4's and 1 Leaf per carton	OK	OK	OK	OK	OK	OK	OK	OK
Print Quality	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Carton checking (Brand Name)		OK	OK	OK	OK	OK	OK	OK	OK
Carton weight check (3 carton) Limit:		OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No:								
	G17250103								
	Mfg. Dt:								
	01/2025	OK	OK	OK	OK	OK	OK	OK	OK
	Exp. Dt:	OK	OK						
	12/2027								
Packing Supervisor/ IPQA Signature		A	W	h	Q	h	Q	h	Q
SHIPPERS									
Shipper No.		90	91	105	106	112	113	119	119
Print Quality of Shipper Label	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No:								
	G17250103								
	Mfg. Dt:								
	01/2025	OK	OK	OK	OK	OK	OK	OK	OK
	Exp. Dt:	OK	OK						
	12/2027								
	Qty: (75x20x1x4'S)	OK	OK	OK	OK	OK	OK	OK	OK
Shipper sealing quality	Uniform	OK	OK	OK	OK	OK	OK	OK	OK
Packing Supervisor/ IPQA Signature		A	W	h	Q	h	Q	h	Q

- OK / Not OK mention to particular process.

P - Production

Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				CONTROLLED COPY
UNCONTROLLED COPY				

Safetab Life Science

BATCH PACKING RECORD

Page 34 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125/000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

14.0F SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One-Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date		13/02/25	13/02/25	13/02/25	13/02/25	13/02/25	13/02/25	13/02/25	13/02/25
Time		17:34	17:40	18:40	18:46	20:35	20:40	21:40	21:47
Parameter	Standards and Limits								
strip in catch cover	1X4's	OK	OK	OK	OK	OK	OK	OK	OK
No. of catch cover per carton	20X1X4's and 1 Leaf per carton	OK	OK	OK	OK	OK	OK	OK	OK
Print Quality	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Carton checking (Brand Name)		OK	OK	OK	OK	OK	OK	OK	OK
Carton weight check (3 carton) Limit:		OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No: G17250103 Mfg.Dt: 01/2025 Exp.Dt: 12/2027	OK	OK	OK	OK	OK	OK	OK	OK
Packing Supervisor/ IPQA Signature		h	h	h	h	h	h	h	h
SHIPPERS									
Shipper No.		126	127	134	135	147	148	159	159
Print Quality of Shipper Label	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No: G17250103 Mfg.Dt: 01/2025 Exp.Dt: 12/2027 Qty: (75x20x1x4'S)	OK	OK	OK	OK	OK	OK	OK	OK
Shipper sealing quality	Uniform	OK	OK	OK	OK	OK	OK	OK	OK
Packing Supervisor/ IPQA Signature		h	h	h	h	h	h	h	h

• OK / Not OK mention to particular process.

P - Production Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				
UNCONTROLLED COPY				



Safetab Life Science

BATCH PACKING RECORD

Page 82 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	928	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

14.06 SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One-Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date		14/02/25	14/02/25	14/02/25	14/02/25	14/02/25	14/02/25	14/02/25	14/02/25
Time		6:40	6:45	7:45	7:50	9:10	9:15	10:15	10:20
Parameter	Standards and Limits								
strip in catch cover	1X4's	OK	OK	OK	OK	OK	OK	OK	OK
No. of catch cover per carton	20X1X4's and 1 Leaf per carton	OK	OK	OK	OK	OK	OK	OK	OK
Print Quality	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Carton checking (Brand Name)		OK	OK	OK	OK	OK	OK	OK	OK
Carton weight check (3 carton) Limit:		OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No: G17250103 Mfg.Dt: 01/2025 Exp.Dt: 12/2027	OK	OK	OK	OK	OK	OK	OK	OK
Packing Supervisor/ IPQA Signature		A	(S)	A	(S)	A	(S)	A	(S)
SHIPPERS									
Shipper No.		163	164	167	169	171	173	175	175
Print Quality of Shipper Label	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No: G17250103 Mfg.Dt: 01/2025 Exp.Dt: 12/2027 Qty: (75x20x1x4'S)	OK	OK	OK	OK	OK	OK	OK	OK
Shipper sealing quality	Uniform	OK	OK	OK	OK	OK	OK	OK	OK
Packing Supervisor/ IPQA Signature		A	(S)	A	(S)	A	(S)	A	(S)

• OK / Not OK mention to particular process.

P - Production

Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with	(S)	(S)	(S)
UNCONTROLLED COPY		CONTROLLED COPY	

		Safetab Life Science BATCH PACKING RECORD		Page 82 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

14.0H SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date		14/02/25	14/02/25	14/02/25	14/02/25	14/02/25	14/02/25	14/02/25	14/02/25
Time		11:20	11:25	12:12	12:17	13:15	13:20	14:20	14:25
Parameter	Standards and Limits								
strip in catch cover	1X4's	OK	OK	OK	OK	OK	OK	OK	OK
No. of catch cover per carton	20X1X4's and 1 Leaf per carton	OK	OK	OK	OK	OK	OK	OK	OK
Print Quality	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Carton checking (Brand Name)		OK	OK	OK	OK	OK	OK	OK	OK
Carton weight check (3 carton) Limit:		OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No:	OK	OK	OK	OK	OK	OK	OK	OK
	G17250103								
	Mfg.Dt:								
	01/2025								
	Exp.Dt:		OK	OK	OK	OK	OK	OK	OK
	12/2027								
Packing Supervisor/ IPQA Signature		R	(SIC)	R	(SIC)	R	(SIC)	R	(SIC)
SHIPPERS									
Shipper No.		180	181	186	187	194	195	196	196
Print Quality of Shipper Label	Clear & Legible	OK	OK	OK	OK	OK	OK	OK	OK
Overprinting Details	Batch No:	OK	OK	OK	OK	OK	OK	OK	OK
	G17250103								
	Mfg.Dt:								
	Exp.Dt:								
	01/2025								
	12/2027		OK	OK	OK	OK	OK	OK	OK
	Qty: (75x20x1x4'S)		OK	OK	OK	OK	OK	OK	OK
Shipper sealing quality	Uniform	OK	OK	OK	OK	OK	OK	OK	OK
Packing Supervisor/ IPQA Signature		R	(SIC)	R	(SIC)	R	(SIC)	R	(SIC)

• OK / Not OK mention to particular process.

P -- Production Q -- Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				
UNCONTROLLED COPY				

		Safetab Life Science BATCH PACKING RECORD		Page 84 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1230000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

14.01 SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One-Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date		14/02/25	14/02/25						
Time		15:25	15:30						
Parameter	Standards and Limits								
strip in catch cover	1X4's	OK	OK						
No. of catch cover per carton	20X1X4's and 1 Leaf per carton	OK	OK						
Print Quality	Clear & Legible	OK	OK						
Carton checking (Brand Name)		OK	OK						
Carton weight check (3 carton) Limit:		OK	OK						
Overprinting Details	Batch No: G17250103								
	Mfg. Dt: 01/2025	OK	OK						
	Exp. Dt: 12/2027								
Packing Supervisor/ IPQA Signature		✓	✓						
SHIPPERS									
Shipper No.		205	205						
Print Quality of Shipper Label	Clear & Legible	OK	OK						
Overprinting Details	Batch No: G17250103								
	Mfg. Dt: 01/2025								
	Exp. Dt: 12/2027	OK	OK						
Qty: (75x20x1x4's)		OK	OK						
Shipper sealing quality	Uniform	OK	OK						
Packing Supervisor/ IPQA Signature		✓	✓						

• OK / Not OK mention to particular process.

P - Production Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science

BATCH PACKING RECORD

Page 85 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO.:	ST/MPR/005/R3	Supersedes MPR No.:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

14.01 SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date									
Time									
Parameter	Standards and Limits								
strip in catch cover	1X4's								
No. of catch cover per carton	20X1X4's and 1 Leaf per carton								
Print Quality	Clear & Legible								
Carton checking (Brand Name)									
Carton weight check (3 carton) Limit:									
Overprinting Details	Batch No:								
	Mfg.Dt:								
	Exp.Dt:								
Packing Supervisor/ IPQA Signature									
SHIPPERS									
Shipper No.									
Print Quality of Shipper Label	Clear & Legible								
Overprinting Details	Batch No:								
	Mfg.Dt:								
	Exp.Dt:								
	Qty: (75x20x1x4'S)								
Shipper sealing quality	Uniform								
Packing Supervisor/ IPQA Signature									

• OK / Not OK mention to particular process.

P - Production

Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science
BATCH PACKING RECORD

Page 86 of 169

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	20X1X4's
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

14.0K SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One-Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date									
Time									
Parameter	Standards and Limits								
strip in catch cover	1X4's								
No. of catch cover per carton	20X1X4's and 1 Leaf per carton								
Print Quality	Clear & Legible								
Carton checking (Brand Name)									
Carton weight check (3 carton) Limit:									
Overprinting Details	Batch No:								
	Mfg.Dt:								
	Exp.Dt:								
Packing Supervisor/ IPQA Signature									
SHIPPERS									
Shipper No.									
Print Quality of Shipper Label	Clear & Legible								
Overprinting Details	Batch No:								
	Mfg.Dt:								
	Exp.Dt:								
	Qty: (75x20x1x4'S)								
Shipper sealing quality	Uniform								
Packing Supervisor/ IPQA Signature									

• OK / Not OK mention to particular process.

P - Production Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science

BATCH PACKING RECORD

Page 82 of 109


Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R.3	Supersedes MPR No:	ST/MPR/005/R.2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

14.01 SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One-Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date:									
Time:									
Parameter	Standards and Limits								
strip in catch cover	1X4's								
No. of catch cover per carton	20X1X4's and 1 Leaf per carton								
Print Quality	Clear & Legible								
Carton checking (Brand Name)									
Carton weight check (3 carton) Limit:									
Overprinting Details	Batch No:								
	Mfg. Dt:								
	Exp. Dt:								
Packing Supervisory/ IPQA Signature									
SHIPPERS									
Shipper No.									
Print Quality of Shipper Label	Clear & Legible								
Overprinting Details	Batch No:								
	Mfg. Dt:								
	Exp. Dt:								
	Qty: (75x20x1x4'S)								
Shipper sealing quality	Uniform								
Packing Supervisor/ IPQA Signature									

• OK / Not OK mention to particular process.

P - Production Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				

UNCONTROLLED COPY



Safetab Life Science

BATCH PACKING RECORD

Page 88 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

14.0M SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One-Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date									
Time									
Parameter	Standards and Limits								
strip in catch cover	1X4's								
No. of catch cover per carton	20X1X4's and 1 Leaf per carton								
Print Quality	Clear & Legible								
Carton checking (Brand Name)									
Carton weight check (3 carton) Limit:									
Overprinting Details	Batch No:								
	Mfg.Dt:								
	Exp.Dt:								
Packing Supervisor/ IPQA Signature									
SHIPPERS									
Shipper No.									
Print Quality of Shipper Label	Clear & Legible								
Overprinting Details	Batch No:								
	Mfg.Dt:								
	Exp.Dt:								
	Qty: (75x20x1x4'S)								
Shipper sealing quality	Uniform								
Packing Supervisor/ IPQA Signature									

• OK / Not OK mention to particular process.

P - Production

Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				
UNCONTROLLED COPY				



Safetab Life Science

BATCH PACKING RECORD

Page 89 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

14.0N SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One-Hour

Checked by		P	Q	P	Q	P	Q	P	Q
Date									
Time									
Parameter	Standards and Limits								
strip in catch cover	1X4's								
No. of catch cover per carton	20X1X4's and 1 Leaf per carton								
Print Quality	Clear & Legible								
Carton checking (Brand Name)									
Carton weight check (3 carton) Limit:									
Overprinting Details	Batch No:								
	Mfg.Dt:								
	Exp.Dt:								
Packing Supervisor/ IPQA Signature:									
SHIPPERS									
Shipper No.									
Print Quality of Shipper Label	Clear & Legible								
Overprinting Details	Batch No:								
	Mfg.Dt:								
	Exp.Dt:								
	Qty: (75x20x1x4'S)								
Shipper sealing quality	Uniform								
Packing Supervisor/ IPQA Signature									

• OK / Not OK mention to particular process.

P - Production

Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY				

Safetab Life Science

BATCH PACKING RECORD

Page 90 of 109

Product:	LITACOLD FLU TABLETS	Product ID No:	978	Batch Size:	250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

14.0-0 SECONDARY PACKING LINE IN-PROCESS CHECKS: Frequency: Every One-Hour

Checked by	P	Q	P	Q	P	Q	P	Q
Date								
Time								
Parameter	Standards and Limits							
strip in catch cover	1X4'S							
No. of catch cover per carton	20X1X4'S and 1 Leaf per carton							
Print Quality	Clear & Legible							
Carton checking (Brand Name)								
Carton weight check (3 carton) Limit:								
Overprinting Details	Batch No:							
	Mfg.Dt:							
	Exp.Dt:							
Packing Supervisor/ IPQA Signature								
SHIPPERS								
Shipper No.								
Print Quality of Shipper Label	Clear & Legible							
Overprinting Details	Batch No:							
	Mfg.Dt:							
	Exp.Dt:							
Qty: (75x20x1x4'S)								
Shipper sealing quality	Uniform							
Packing Supervisor/ IPQA Signature								

• OK / Not OK mention to particular process.

P - Production

Q - Quality Assurance

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY CONTROLLED COPY			



Safetab Life Science
BATCH PACKING RECORD

Page 91 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

15.0 WEIGHT CHECKS ON SHIPPERS:

15.1 LIMITS:

Weight of Ten empty shipper

S.No	1	2	3	4	5	6	7	8	9	10
Weight in gm	1.142	1.170	1.166	1.154	1.212	1.168	1.166	1.170	1.172	1.152

Average weight of empty shipper = $\frac{11.66}{10} = (A) 1.166$

Filled carton weight (Refer 12.4) (B) * 141.23

Weight of one filled shipper (C) = A + (B X 75) + Tape Weight
= 11.766

SHIPPER LIMIT (C±B/2) = (11.696 TO 11.836)

15.1.1 PACKING PROFILE:

Pack 4 tablet in 1 Strip → 1 strip in 1 catch cover → 20 catch cover into 1 carton → 75 carton in 1 shipper (75x20x1x4'S)

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	CONTROLLED COPY
Signature with				

UNCONTROLLED COPY

Safetab Life Science BATCH PACKING RECORD

Page 92 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	125000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

16.0 WEIGHT CHECKS:

Date: 12/02/25			Date: 12/02/25			Date: 13/02/25			Date: 13/02/25		
Shipper No	Gross weight (kg)	Done by packing	Shipper No	Gross weight (kg)	Done by packing	Shipper No	Gross weight (kg)	Done by packing	Shipper No	Gross weight (kg)	Done by packing
1	11.800	Sayari	23	11.768	Sayari	45	11.794	Anita	67	11.768	Anita
2	11.768	Sayari	24	11.756	Sayari	46	11.740	Anita	68	11.740	Anita
3	11.756	Sayari	25	11.764	Sayari	47	11.742	Anita	69	11.780	Anita
4	11.768	Sayari	26	11.770	Sayari	48	11.764	Anita	70	11.776	Anita
5	11.764	Sayari	27	11.742	Sayari	49	11.756	Anita	71	11.768	Anita
6	11.794	Sayari	28	11.756	Sayari	50	11.784	Anita	72	11.730	Anita
7	11.770	Sayari	29	11.800	Sayari	51	11.752	Anita	73	11.748	Anita
8	11.742	Sayari	30	11.814	Sayari	52	11.750	Anita	74	11.760	Anita
9	11.764	Sayari	31	11.742	Sayari	53	11.760	Anita	75	11.758	Anita
10	11.756	Sayari	32	11.756	Sayari	54	11.766	Anita	76	11.756	Anita
11	11.784	Sayari	33	11.764	Sayari	55	11.770	Anita	77	11.740	Anita
12	11.756	Sayari	34	11.771	Sayari	56	11.742	Anita	78	11.780	Anita
13	11.752	Sayari	35	11.730	Sayari	57	11.730	Anita	79	11.788	Anita
14	11.760	Sayari	36	11.742	Sayari	58	11.728	Anita	80	11.760	Anita
15	11.771	Sayari	37	11.780	Sayari	59	11.768	Anita	81	11.768	Anita
16	11.730	Sayari	38	11.800	Elavari	60	11.756	Anita	82	11.770	Anita
17	11.742	Sayari	39	11.791	Elavari	61	11.748	Anita	83	11.760	Elavari
18	11.756	Sayari	40	11.761	Elavari	62	11.558	Anita	84	11.788	Elavari
19	11.770	Sayari	41	11.770	Elavari	63	11.750	Anita	85	11.770	Elavari
20	11.768	Sayari	42	11.764	Elavari	64	11.780	Anita	86	11.770	Elavari
21	11.756	Sayari	43	11.768	Anita	65	11.786	Anita	87	11.756	Elavari
22	11.756	Sayari	44	11.768	Anita	66	11.760	Anita	88	11.760	Elavari
Date		Packed in Shipper No		Weighed by		Checked by		Packed quantity		Loose pack details	
		From	To								
12/02/25		01	42	Sayari		m		3150 X 802		Shipper No:	
13/02/25		43	88	m				3450 X 803		Quantity:	
										Total packed quantity:	528000
Checked by production: 										Checked by production	
Verified by IPQA: 										Verified by IPQA	
Remarks (If any): N/A											

Note - If any shipper weight found out of limit, the same shall be checked physically by opening the shipper.

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science

BATCH PACKING RECORD

Page 93 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	12000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

16.0A WEIGHT CHECKS:

Shipper No	Gross weight (kg)	Done by packing	Shipper No	Gross weight (kg)	Done by packing	Shipper No	Gross weight (kg)	Done by packing	Shipper No	Gross weight (kg)	Done by packing
89	11.770	Ekanesi	111	11.720	Rajavalli	133	11.898	Jayasee	155	11.722	R. Santhi
90	11.730	Ekanesi	112	11.732	Rajavalli	134	11.742	Jayasee	156	11.708	R. Santhi
91	11.732	Jayasee	113	11.740	Rajavalli	135	11.786	Jayasee	157	11.706	R. Santhi
92	11.764	Jayasee	114	11.762	Rajavalli	136	11.800	Jayasee	158	11.730	R. Santhi
93	11.768	Jayasee	115	11.768	Rajavalli	137	11.810	Jayasee	159	11.710	R. Santhi
94	11.770	Jayasee	116	11.756	Rajavalli	138	11.760	Jayasee	160	11.746	Jaya
95	11.758	Jayasee	117	11.740	Rajavalli	139	11.780	Jayasee	161	11.700	Jaya
96	11.768	Jayasee	118	11.738	Rajavalli	140	11.764	Jayasee	162	11.732	B. Selvi
97	11.800	Jayasee	119	11.740	Rajavalli	141	11.788	Jayasee	163	11.720	B. Selvi
98	11.810	Jayasee	120	11.756	Rajavalli	142	11.756	Jayasee	164	11.718	B. Selvi
99	11.794	Jayasee	121	11.750	Rajavalli	143	11.746	Jayasee	165	11.716	B. Selvi
100	11.794	Jayasee	122	11.748	Rajavalli	144	11.742	R. Santhi	166	11.742	B. Selvi
101	11.776	Jayasee	123	11.762	Rajavalli	145	11.710	R. Santhi	167	11.780	B. Selvi
102	11.764	Jayasee	124	11.736	Rajavalli	146	11.740	R. Santhi	168	11.732	B. Selvi
103	11.784	Jayasee	125	11.756	Jayasee	147	11.716	R. Santhi	169	11.742	B. Selvi
104	11.768	Jayasee	126	11.740	Jayasee	148	11.720	R. Santhi	170	11.736	B. Selvi
105	11.770	Jayasee	127	11.732	Jayasee	149	11.724	R. Santhi	171	11.732	B. Selvi
106	11.740	Jayasee	128	11.748	Jayasee	150	11.716	R. Santhi	172	11.740	B. Selvi
107	11.748	Jayasee	129	11.720	Jayasee	151	11.736	R. Santhi	173	11.742	B. Selvi
108	11.762	Jayasee	130	11.734	Jayasee	152	11.708	R. Santhi	174	11.740	B. Selvi
109	11.736	Jayasee	131	11.780	Jayasee	153	11.710	R. Santhi	175	11.748	B. Selvi
110	11.720	Jayasee	132	11.746	Jayasee	154	11.720	R. Santhi	176	11.750	B. Selvi

Date	Packed in Shipper No		Weighed by	Checked by	Packed quantity	Loose pack details	
	From	To					
13/02/25	89	159	Jayasee	✓	5325 X 805	Shipper No:	-
14/02/25	160	176	B. Selvi	✓	1275 X 805	Quantity:	-
				mm		Total packed quantity:	5,28,000
Checked by production: ✓						Checked by production	R.
Verified by IPQA: ✓						Verified by IPQA	✓
Remarks (If any): Nil							

Note - If any shipper weight found out of limit, the same shall be checked physically by opening the shipper.

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			

UNCONTROLLED COPY

CONTROLLED COPY

Safetab Life Science

BATCH PACKING RECORD

Page 94 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 TABLETS
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO.:	ST/MPR/005/R3	Supersedes MPR No.:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

16.0B WEIGHT CHECKS:

Date:	Shipper No	Gross weight (kg)	Done by packing	Date:	Shipper No	Gross weight (kg)	Done by packing	Date:	Shipper No	Gross weight (kg)	Done by packing
14/02/25	177	11.760	B. Selvi	14/02/25	149	11.736	Jayashri				
	178	11.772	B. Selvi		200	11.740	Jayashri				
	179	11.732	B. Selvi		201	11.732	Jayashri				
	180	11.738	B. Selvi		202	11.742	Jayashri				
	181	11.740	B. Selvi		203	11.800	Jayashri				
	182	11.722	B. Selvi		204	11.764	Jayashri				
	183	11.718	B. Selvi		205	11.456	Jayashri				
	184	11.710	B. Selvi								
	185	11.740	B. Selvi								
	186	11.784	B. Selvi								
	187	11.802	B. Selvi								
	188	11.748	B. Selvi								
	189	11.754	B. Selvi								
	190	11.766	B. Selvi								
	191	11.774	B. Selvi								
	192	11.776	B. Selvi								
	193	11.780	B. Selvi								
	194	11.760	B. Selvi								
	195	11.768	Jayashri								
	196	11.784	Jayashri								
	197	11.800	Jayashri								
	198	11.772	Jayashri								
Date	Packed in Shipper No		Weighed by	Checked by	Packed quantity	Loose pack details					
	From	To				Shipper No:					
14/02/25	177	205	Selvi Jay	✓	2102 x 80's	Shipper No:	205 ✓				
						Quantity:	2 x 20 x 4's ✓				
						Total packed quantity:	168160				
Checked by production: ✓						Checked by production	Beng ✓				
Verified by IPQA: ✓						Verified by IPQA	✓				
Remarks (If any) nil											

Note - If any shipper weight found out of limit, the same shall be checked physically by opening the shipper.

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date			
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science

BATCH PACKING RECORD

Page 95 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	125000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

16.0C WEIGHT CHECKS:

[illegible]

Note - If any shipper weight found out of limit, the same shall be checked physically by opening the shipper.

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date	 11/05/24	 11/05/24	 11/05/24
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science
BATCH PACKING RECORD

Page 95 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

16.0D WEIGHT CHECKS:

[illegible]

Note - If any shipper weight found out of limit, the same shall be checked physically by opening the shipper.

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date	 11/05/24	 11/05/24	 11/05/24	
UNCONTROLLED COPY				

[illegible]

16.0E WEIGHT CHECKS:


Note - If any shipper weight found out of limit, the same shall be checked physically by opening the shipper.

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date	 14/05/24	 14/05/24	 14/05/24

CONTROLLED COPY

UNCONTROLLED COPY



Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date	 11/07/24	 11/07/24	 11/07/24

CONTROLLED COPY

UNCONTROLLED COPY

Page 99 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

16.0G WEIGHT CHECKS:

[illegible]

Note - If any shipper weight found out of limit, the same shall be checked physically by opening the shipper.

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date	 24	 11/05/24	 11/05/24

CONTROLLED COPY

UNCONTROLLED COPY



Safetab Life Science
BATCH PACKING RECORD

Page 100 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	150000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027


17.0 PACKING MATERIALS RECONCILIATION:

S. No.	Material	Sales							
		Printed strip Foil (Kg)	Plain strip Foil (Kg)	Catch cover (no's)	Carton (No's)	5 Ply Shipper (No's)	Transparent GHPL printed Bopp Tape 2"	Leaflets	Labels (No's)
A	Quantity received for the batch	163.00	162.65	312500	15625	213	13	15625	206
B	Quantity packed	161.00	161.00	306040	15302	204	13	15302	205
C	Quantity rejected (Online) & Destroyed	2.00	1.65	6259	301	9	-	312	-
D	Quantity Returned materials	-	-	-	-	-	-	-	-
E	Sample affixed on BPR	-	-	1	12	-	-	1	1
F	QC samples	-	-	200	10	-	-	10	-
G	Other samples (If any)	-	-	-	-	-	-	-	-
H	Total accountable Quantity (B + C + D + E + F + G)	163.00	162.65	312500	15625	213	13	15625	206
I	Variation (A - I)	-	-	-	-	-	-	-	-
J	% Variation ($\frac{I}{A} \times 100$)	0%	0%	0%	0%	0%	0%	0%	0%
K	Acceptable limit	$\pm 2\%$			$\pm 1\%$				

If % Variation is out of the acceptable limit, investigate & specify the reason:

17.1 DESTRUCTION OF REJECTED PACKING MATERIALS AND STEREOS:

Date	Stereos (No's)	Printed Foil (Kg)	Plain BSO foil (Kg)	Catch cover (No's)	Printed Cartons (No's)	Shipper Labels (No's)	Shipper (No's)	Leaflets	Destroyed by Packing	Verified by IPQA
14/10/25	30	2.00	1.65	6259	301	-	9	312	A	100
Total Qty. Destroyed	30	2.00	1.65	6259	301	-	9	312		

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date	 11/05/24	 11/05/24	 11/05/24
UNCONTROLLED COPY		CONTROLLED COPY	

Safetab Life Science
BATCH PACKING RECORD

Page 101 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

18.0 SHIPPER LABEL:

18.1 Reference Shipper label format:


SHIPPING MARK		
Product Name		
Batch No.		
Mfg.Date		
Exp.Date		
Quantity		
Gross Weight		kg
Net Weight		kg
Case No.		
Supplier. Generic Healthcare Pvt. Ltd.		
Barcode details shall be printed		


18.2 GHPL Printed Transparent Bopp Tape H type sealing for both side and White color Patti H type.

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				
UNCONTROLLED COPY				

		Safetab Life Science BATCH PACKING RECORD		Page 102 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR No:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

19.0 PACKING LABEL PASTE:

SHIPPING MARK	
Product Name :	Litacold CPR B/80
Batch No.	G17250103
Mfg. Date	01/2025
Exp. Date	12/2027
Quantity	75 X 20 X 1 X 4T
Gross Weight	11.94 KGS
Net Weight	10.76 KGS
Case No.	1
 	
Supplier: Generic Healthcare Pvt Ltd. [FOR EXPORT ONLY]	
 (01)58906067080252(17)271231(10)G17250103	
 (00)089060670805959598	

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	CONTROLLED COPY
Signature with Date				
UNCONTROLLED COPY				

Safetex Life Science
BATCH PACKING RECORD

Page 1 of 1

Product: LITACOLD FLU TABLET	Product ID No.: 378	Batch Size:	Batch Label:
Generic: Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate	MPR No.: S1/MPR/005/R3	Supersedes MPR No.:	MPR/005/R2
Batch No.: G17250103	Mfg. Date: 01/2025	Exp. Date:	12/2027

BMR COPY/QUALITY CONTROL COPY

20.0 TEST REQUEST FORM FOR FINISHED PRODUCT:

From: Manufacturing Department	To: Q.A. Department
To be filled by production	To be filled by IPQA
Date of Mfg.: 01/2025 Date of exp.: 12/2027	Quantity Sampled: 1 X 13 X 4'S
Test Request for: <u>Content</u>	Sampled by: <u>A. B.</u> on date: 12/02/25
TRF raised by: <u>ft</u> on date: 12/02/25	
To be filled by QC	
TRF sample received sign & date: <u>N. J. J.</u> 12/02/25	
Sample tested on: 12/02/25	
Test Results	Limits
IMPACTS parameters complies as per specification	
Analytical Reference no.: <u>SFP20185</u>	Space for QC Approval (Stamping and Signing)
Analysed By: <u>Km</u> on: 12/02/25	

[Signature] 12/02/25

Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with Date	<u>[Signature]</u> 12/02/25	<u>[Signature]</u> 12/02/25	<u>[Signature]</u> 12/02/25

CONTROLLED COPY

UNCONTROLLED COPY

Safetab Life Science					
BATCH PACKING RECORD					
Page 104 of 109					
Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR No.:	ST/MPR/005/R3	Supersedes MPR No.:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

21.0 TOTAL BATCH RECONCILIATION:

A	Batch Size	1250000	No's
B	Quantity received for packing	1227181	
C	Total quantity packed	1224160	
D	Control samples + Physical Sample	240 + 60	
E	Analytical Samples	52	
F	In process samples	1480	
G	Others samples if any (Stability, Regulatory Samples, etc.)	472	
H	Packing Loss = B - (C + D + E + F + G)	717	
I	Total Loss = A - (C + D + E + F + G)	23536	
J	Packing yield % (C + D + E + F + G) / B X 100 (Limit: NLT 99.00%)	99.94%	
K	Yield % (C + D + E + F + G) / A X 100 (Limit: NLT 97.00%)	98.12%	

21.1 TRANSFER OF FINISHED GOODS:

Date	FGR No	Pack Size of per shipper	Quantity transferred	No. of Shippers	Quantity sent to FGR	Issued by Packing	Checked by IPQA
14/02/25	1554	75	1224160	204+1	15302x805	✓	✓

21.2 RELEASE FOR DISPATCH:

Total batch Quantity released by (packing):

Checked by QA Head:


Job title:	Prepared by:	Checked by:	Approved by:
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee
Signature with			
UNCONTROLLED COPY CONTROLLED COPY			

		Safetab Life Science BATCH PACKING RECORD		Page 105 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No.:	078	Batch Size:	100000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

22.0 PROCESS HISTORY

Bar code label print given by 
12/02/25

Total - 2052100.

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY				



Safetab Life Science
BATCH PACKING RECORD

Page 106 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	9/8	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR No.:	ST/MPR/005/R3	Supersedes MPR No.:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

23.0 Following the BMR, BPR records in additional attached.		Checked by Production Head	Reviewed by QA Head
1	Filling stage analytical report	QA	QA
2	Stereo Impression approval sheet	✓	✓
3	Approved Specimen of Aluminium foil with the batch stereo overprinting	✓	✓
4	Approved Specimen of overprinting carton	✓	✓
5	Checked in Packing label paste	✓	✓
6	FGR transferred slip attached	✓	✓
7	Finished product Certificate of Analysis attached	✓	✓
8	Any other items attached	✓	✓
Note: Tick mark (✓) if Complies and (X) if Does not comply.			
This Batch Manufacturing Record submitted to QA Head on date <u>17/02/25</u>			
Production Head sign/date <u>M 17/02/25</u>			
This Batch Manufacturing Record received, reviewed and batch is released for distribution on date <u>17/02/2025</u>			
Reviewed by QA Head <u>[Signature]</u> <u>17/02/2025</u>			

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with	<u>[Signature]</u>	<u>[Signature]</u> <u>11/05/24</u>	<u>[Signature]</u>	

UNCONTROLLED COPY


Safetab Life Science
BATCH PACKING RECORD

Page 107 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size:	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg. Date:	01/2025	Exp. Date:	12/2027

24.0 ABBREVIATIONS

S.No.	Short Form	Abbreviation	S.No.	Short Form	Abbreviation
1	BPR	Batch Packing Record	16	Tab	Tablet
2	MFG	Manufacturing	17	Std	Standard
3	EXP	Expiry	18	Avg	Average
4	MRP	Maximum Retail Price	19	rpm	Rotation Per Minute
5	QA	Quality Assurance	20	No's	Numbers
6	Kg	Kilogram	21	Sign	Signature
7	g	Gram	22	IPQA	In Process Quality Assurance
8	Min	Minutes	23	SOP	Standard Operating Procedure
9	Hrs	Hours	24	mm	Millimeter
10	M/C	Machine	25	NLT	Not Less Than
11	Prod	Production	26	NMT	Not More Than
12	B.No	Batch number	27	QC	Quality Control
13	FGR	Finished Good Receipt	28	A.R.No.	Analytical Report Number
14	mg	Milligram	29	L	Lakhs
15	Qty	Quantity			


Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				
UNCONTROLLED COPY				

 Safetab Life Science BATCH PACKING RECORD		Page 108 of 109	
Product:	LITACOLD FLU TABLETS	Product ID No:	978
Batch Size:	1250000 Tablets		
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorphenamine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3
Supersedes MPR No:	ST/MPR/005/R2		
Batch No.:	G17250103	Mfg.Date:	01/2025
Exp.Date:	12/2027		

25.0 BATCH PACKING RECORD REVIEW CHECK LIST		
S.No.	CHECK POINTS	Status
1	The entire Page no: <u>1</u> to <u>109</u> are available, verified and found correct.	Yes / No
2	Batch numbers available in all pages verified and are correct.	Yes / No
3	Manufacturing and expiry is correctly allotted verified.	Yes / No
4	Dispensing carried out on calibrated balance verified.	Yes / No
5	Dispensing Cards attached.	Yes / No
6	Line clearance is taken prior to all packing activities.	Yes / No
7	Dispensing is carried out as per indent of packing materials and are correct.	Yes / No
8	Environmental conditions comply during all the packing stages.	Yes / No
9	All wrong entries if any have been strike out and signed and verified.	Yes / No
10	Actual equipment used is as specified and verified.	Yes / No
11	Individual specimen of Printed foil, Printed carton and shipper label, Inspected tablets container label is attached with signature, verified and are correct.	Yes / No
12	Destruction record of rejected packing material is filled and destruction is verified by IPQA.	Yes / No
13	Surplus packing materials return request is filled and surplus quantity of packing material is returned and verified.	Yes / No / NA
14	Packing material reconciliation is filled and variance is within acceptance limit.	Yes / No
15	In shipper weighing record none of the shippers are out of limit. Where any shipper weight is out of limit, has been opened and quantity verified and signed against that particular shipper number and correct.	Yes / No
16	If any deviation, change control or OOS initiated, related to this batch verified and corrective action implemented.	Yes / No / NA
Note:	If any abnormalities noticed during batch Packing record review it should be brought to the notice of QA Head.	

Batch Packing Record reviewed by IPQA Sign/Date: 17/02/25

QA Head Sign/Date: 17/02/2025

Job title:	Prepared by:	Checked by:	Approved by:	CONTROLLED COPY
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with Date				
UNCONTROLLED COPY				



Safetab Life Science
BATCH PACKING RECORD

Page 109 of 109

Product:	LITACOLD FLU TABLETS	Product ID No.:	978	Batch Size	1250000 Tablets
Generic:	Paracetamol, Phenylephrine Hydrochloride, Chlorpheniramine Maleate and Caffeine Anhydrous Tablets	MPR NO:	ST/MPR/005/R3	Supersedes MPR No:	ST/MPR/005/R2
Batch No.:	G17250103	Mfg.Date:	01/2025	Exp.Date:	12/2027

26.0 REVISION HISTORY:

Date of Revision	Revision No.	Changes Made	Effective Batch No.:	Effective Date
	R0	1. BPR is Prepared For 12.50 L and Plan 5-ply shipper change to Pre-Printed "Litacold" 5-ply shipper. Refer change control number: ST/CC/22/242	G1722115	20/12/2022
05/02/2024	R1	1. New shipper Litacold brand name printed on 5 ply shipper	G17240211	07/02/2024
11/05/2024	R2	1. Manufacturin from Safetab Life Science to GHPL Loan Licence. Refer change control No: ST/CC/24/039.	G17240415	11/05/2024

Job title:	Prepared by:	Checked by:	Approved by:	
Designation	QA Executive / Designee	Production Head / Designee	QA Head / Designee	
Signature with				
UNCONTROLLED COPY				CONTROLLED COPY